



TECHNICAL ASSOCIATES INSTRUMENTATION FOR NUCLEAR RESEARCH 140 WEST PROVIDENCIA AVENUE • BURBANK • CALIFORNIA

J	Son Mee Co	
	Juno Inspection Tag	
	Serial No. 2830	
-	1. Subchassis assembly and test with source MB	
	2. Install and check screens and slidesEV	
	3. Case finish and engraving <u>E.V.</u>	
	4. Install case and recheck slides <u>E.V.</u>	
	5. Install battery pack A and set zero	·
	6. Calibration: X1 X10 X100	
	7. Install back plate E.V.	
	8. Final Inspection	
	9. Completion Date	

•

MANUAL OF OPERATION JUNO RADIATION SURVEY METER

MODELS SRJ-7 AND HRJ-7

(REV. 10/60)



JUNO RADIATION SURVEY METER

MODELS SRJ-7 AND HRJ-7

I. GENERAL DESCRIPTION

The JUNO SURVEY METER IS A PORTABLE, BATTERY OPERATED INSTRUMENT FOR MEASURING THE INTENSITY OF, AND DISTINGUISHING BETWEEN ALPHA, BETA AND GAMMA RADIATION. THIS INSTRUMENT IS NORMALLY SUPPLIED IN TWO MODELS: SRJ-7 (standard range) which has a sensitivity range suitable for all routine applications; and HRJ-7 (high range) for use where exceptionally high radiation levels are likely to be encountered. To easily identify the High Range Juno (HRJ-7), the range selector switch knob and the meter face are red in color. Both models meet A.E.C. specifications.

The instrument comprises an ionization chamber, an electrometer circuit, an indicating meter calibrated directly in MR/HR, a removable plug-in battery pack, and two absorption filters for the rejection of alpha and beta particles. The Juno is self-contained in a hard-chrome plated aluminum case which can be easily decontaminated when necessary.

The Model 7 Juno is designed to permit operation under high humidity conditions. The switch box, which contains the Hi-meg grid resistors (R-1, R-2, and R-3) and the electrometer tube (V-1), is sealed by means of gaskets. A replaceable desiccant cartridge is used to take care of any leakage of moisture into the box.

II. SPECIFICATIONS

RANGES:

Model SRJ-7 Juno (standard range): 50, 500, 5000 MR/HR full-scale (Improved A.E.C. Model S1C-17B). White meter face.

Model HRJ-7 Juno (High Range): 250, 2500, 25,000 MR/HR full-scale. (Improved A.E.C. Model S1C-17D). Red meter face.

IONIZATION CHAMBER: VOLUME: 27 CUBIC INCHES

WINDOW OPENING: 3" x 4 5/8"

Alpha Screen: 0.3 mil Rubber hydrochloride (0.45 mg/cm^2) .

Alpha Rejection Absorber: $0.01^{\prime\prime}$ cellulose acetate sheet (36 mg/cm²).

BETA REJECTION: 0.102" ALUMINUM (720 MG/CM²)

DISTANCE FROM BOTTOM OF FEET TO THE ALPHA SCREEN 1.7 CM.

ACCURACY:

+10% TO -5% OF GAMMA RADIATION INTENSITY WHEN AIR DENSITY CORRECTIONS ARE MADE.

SRJ-7 HRJ-7

OPERATING TEMPERATURE Range: 350

35⁰F. то 135⁰F.

CASE:

HARD-CHROME PLATED ALUMINUM WITH ENGRAVED MARKINGS ON TOP OF CASE; GASKETED, DUST AND MOISTURE RESIST-ANT. DIMENSIONS: 11 1/2" L. x 5 3/4" W. x 6 1/2" H.

HANDLE: ALUMINUM ESPECIALLY CAST LOW-POROSITY, SMOOTHLY POLISHED.

WEIGHT:

6 LBS. 11 OZ. NET.

- III. THEORY OF OPERATION
 - A. IONIZATION CHAMBER.

THE IONIZATION CHAMBER IS A RECTANGULAR ALUMINUM CASE, (THE INSIDE OF WHICH IS COATED WITH AQUADAG) LOCATED IN THE FORWARD END, BOTTOM OF THE JUNO CASE. AIR AT ATMOSPHERIC PRESSURE IS THE IONIZATION MEDIUM. THE COLLECTOR, AN ALUMINUM ROD COATED WITH AQUADAG, IS EXTREMELY WELL INSULATED BY MEANS OF A TEFLON SUPPORT BUSHING. ALPHA, BETA AND/OR GAMMA RADIATION, WHEN PASSING THROUGH THE CHAMBER, COLLIDES WITH INDIVID-UAL MOLECULES OF AIR. A COLLISION CAUSES AN ELECTRON TO BE DISPLACED FROM ITS NORMAL ENERGY LEVEL. ONCE REMOVED SUFFICIENTLY FAR FROM ITS NORMAL ENERGY LEVEL IT IS SUBJECTED TO STATIC ELECTRIC FORCE. THIS FORCE PULLS THE FREE ELECTRON TO THE COLLECTOR WHICH IS 90 VOLTS POSI-TIVE WITH RESPECT TO THE CHAMBER WALLS. THE MOLECULE, OR ATOM, AFTER LOSING AN ELECTRON ASSUMES A POSITIVE POTENTIAL (SUCH & CHARGED ATOM IS CALLED AN ION) AND IS ATTRACTED BY THE NEGATIVELY CHARGED WALLS. SINCE THE COLLECTOR IS WELL INSULATED THE ELECTRONS WHICH LAND ON ITS SURFACE BUILD UP A CHARGE, RESULTING IN A MORE NEGATIVE POTENTIAL. THIS NEGA-TIVE SHIFT IN THE COLLECTOR IS WHAT IS DETECTED BY THE ELECTROMETER CIR-CUIT. THE REASON FOR CHOOSING ALUMINUM COATED WITH AQUADAG IS TO DUPLI-CATE THE PERFORMANCE OF A HYPOTHETICAL AIR WALL CHAMBER AS NEARLY AS POS-SIBLE. IN ORDER TO GIVE A RELATIVELY CORRECT READING ON ALL TYPES OF RADIATION AND ALL ENERGIES, THE SOLID MATERIAL OF WHICH A CHAMBER IS MADE MUST NEARLY DUPLICATE THE IONIZING PROPERTIES OF AIR AND THEREFORE HAVE APPROXIMATELY THE SAME ATOMIC WEIGHT. AQUADAG APPROACHES THIS VERY CLOSELY AND ALUMINUM IS THE MOST ACCEPTABLE OF THOSE METALS WHICH HAVE SATISFACTORY STRUCTURAL PROPERTIES.

B. ELECTROMETER TUBE AND GRID RESISTORS.

The electrometer tube is a specially constructed, high vacuum pentode which is operated as a triode. It requires a very low power signal on its grid to control its plate current. The low power input permits the use of an extremely high grid resistor, (R-1, R-2, or R-3). Hence, a very small current through the grid resistor will make a readily observed change in the plate current. In turn, the plate current is conducted through a 20 microampere meter which provides the visual indication of the level of radiation in the chamber. The CK 5886 electrometer tube filament operates on 10 ma., 1.25 volts. The plate voltage is 6.5V and the plate current (with no radiation) is 80 microamps. The tube is

SRJ-7 HRJ-7

SUPPORTED BY AN ANTI-VIBRATION TUBE CLAMP WITHIN THE SWITCH BOX. THE BASE IS SPECIALLY TREATED FOR HIGH INSULATION VALUE IN HIGH HUMIDITIES AND WILL PROBABLY BE PERMANENTLY DAMAGED IF TOUCHED WITH FINGERS. THE THREE GRID RESISTORS ARE ALSO SPECIALLY TREATED. THE HOUSING FOR THE SWITCH AND ELECTROMETER IS AT BIAS POTENTIAL SO THAT THE IONIZATION OF AIR IN THE TUBE COMPARTMENT WILL NOT AFFECT THE VOLTAGE OF THE ELECTRO-METER GRID.

This switch box has been sealed, by means of gaskets, to prevent moisture from affecting the operation of the instrument under high humidity conditions. A replaceable desiccant cartridge is used to take care of any leakage of moisture into the box. This cartridge is easily reached, for inspection and replacement, by removing the instrument back-plate (which also gives access to the battery pack). The desiccant material is effective when it is blue in color. When the material becomes pink, the cartridge should be replaced. A special wrench is supplied with each Juno to facilitate the removal of the cartridge if this ever becomes necessary.

CAUTION

IT IS NEITHER NECESSARY NOR DESIRABLE TO USE THE WRENCH WHEN INSTALLING THE CARTRIDGE. FINGER TIGHTNESS IS ADEQUATE AND POSSIBLE THREAD DAMAGE IS AVOIDED.

C. JUNO CIRCUIT.

B-3 is in series with 82,000 ohms and supplies 80 microamps to the METER. THE ELECTROMETER TUBE WITH A PLATE VOLTAGE OF 6.5 VOLTS, AND FILAMENT VOLTAGE OF APPROXIMATELY 1.25 VOLTS, HAS A PLATE CURRENT OF 80 MICROAMPS ALSO. THESE TWO CURRENTS EXACTLY CANCEL IN THE METER, PRO-VIDING THE INSTRUMENT IS PROPERLY ZEROED. ON INCREASING THE RADIATION SUFFICIENTLY TO GO FROM ZERO TO FULL SCALE, THE CURRENT THROUGH THE ELECTROMETER TUBE CHANGES FROM 80 MICROAMPS TO 60 MICROAMPS. IT IS THIS DIFFERENCE THAT IS INDICATED ON THE METER. THE METER IS SHORTED IN THE OFF AND ON POSITIONS SO AS TO LESSEN THE CHANCE OF DAMAGE FROM SHOCK AND VIBRATION. WITH THE SELECTOR SWITCH IN THE ON POSITION THE FILAMENT OF THE ELECTROMETER IS ALLOWED TO WARM UP. THIS INCREASES THE LIFE OF THE TUBE AND REDUCES ZERO DRIFT. WITH THE SELECTOR SWITCH IN THE SET POSI-TION, THE BIAS VOLTAGE, FILAMENT VOLTAGE, PLATE VOLTAGE, AND BUCKING VOLTAGE ARE APPLIED. IN ADDITION THE GRID RESISTANCE IS SHORTED OUT. THIS PERMITS ACCURATE ZEROING OF THE INSTRUMENT IN THE PRESENCE OF RADIA-TION. WHEN THE SELECTOR SWITCH IS TURNED TO THE X1, X10 OR X100 POSITION, THE IONIZATION CURRENT GENERATED IN THE CHAMBER CAUSES VOLTAGE DROP ACROSS R-1, R-2, OR R-3, RESPECTIVELY. SLIGHT VARIATIONS IN THESE RESISTANCES ARE COMPENSATED FOR BY THE CALIBRATION ADJUSTMENTS R-10, R-11 AND R-12, RESPECTIVELY. ZEROING IS ACCOMPLISHED BY CHANGING THE FILAMENT VOLTAGE AND THE TUBE BIAS VOLTAGE.

IV. OPERATION

- A. TO MEASURE RADIATION INTENSITIES:
 - 1. TURN THE SELECTOR SWITCH TO THE ON POSITION AND WAIT FIVE SECONDS FOR THE ELECTROMETER TUBE FILAMENT TO WARM UP.
 - 2. TURN THE SELECTOR SWITCH TO THE SET POSITION AND ADJUST THE ZERO CONTROL UNTIL THE METER READS EXACTLY ZERO.
 - 3. TURN THE SELECTOR SWITCH TO THE "X1" POSITION. IF THE METER READS OFF SCALE, THE RADIATION PRESENT IS GREATER THAN 50 MR/HR IN THE CASE OF THE SRJ-7 (STANDARD RANGE), AND THE SELECTOR SWITCH SHOULD BE TURNED TO EITHER THE "X10" OR "X100" POSITION.
 - 4. THE INSTRUMENT IS NOW READY FOR READING GAMMA RADIATION INTENSITY. TO READ GAMMA AND BETA SIMULTANEOUSLY, SLIDE BACK THE ALUMINUM SCREEN BY PULLING UP RIGHT HAND TAB "G" IN THE HANDLE. TO READ ALPHA, BETA AND GAMMA SIMULTANEOUSLY, SLIDE BACK BOTH THE ALUMINUM AND ACETATE FILTERS BY PULLING UP BOTH TABS IN THE HANDLE.

CAUTION

WITH BOTH FILTERS BACK, THE ALPHA SCREEN IS EXPOSED AND CAN BE VERY EASILY DAMAGED. A VERY DELICATE TOUCH MAY CAUSE IT TO RIP.

5. These instruments are calibrated in International Roentgens corrected to Standard Temperature and Pressure (32°F. and 29.92 inches Hg.) when used at 71.6°F. (22°C) and 29.92 inches Hg. (760 mm. Hg.) For temperature other than 71.6°F. and pressures other than 29.92" Hg., a correction factor must be applied to the meter reading to compensate for ionization chamber air density changes.

CORRECTION FACTOR = $\frac{459.7 + T^{\circ}F}{531.3} \times \frac{29.92}{P. \text{ Inches Hg.}}$

EXAMPLE: A RADIATION MEASUREMENT IS TAKEN AND THE METER READING IS 25 MR/HR ON THE X10 RANGE = 250 MR/HR. THE AIR TEMPERATURE IS 90°F. AND THE BAROMETRIC PRESSURE IS 29.13" Hg.

$$\frac{459.7 + 90^{\circ}}{531.3} = \frac{549.7}{531.3} = 1.035$$

 $\frac{29.92}{29.13} = 1.027$

CORRECTION FACTOR =1.035 x 1.027 = 1.063

CORRECTED READING = 250 MR/HR x 1.063 = 266 MR/HR.

THE EFFECTS OF AIR TEMPERATURE AND PRESSURE MUST BE TAKEN INTO

CONSIDERATION WHEN THESE INSTRUMENTS ARE RECALIBRATED OR WHEN THE CALIBRATION IS CHECKED.

V. CALIBRATION

A. INTERNAL ADJUSTMENTS.

Calibration. There are three internal calibration controls; one for each range. These controls are accessible, through holes in the top of case, by removing the cover plugs labeled X1, X10, and X100. These controls should not be changed unless there are calibrated radium or $\rm Co^{60}$ sources available for accurate recalibration. The proper procedure for calibrating the Model SRJ-7 (Standard Range) is as follows:

- 1. TURN THE INSTRUMENT ON AND ALLOW SEVERAL MINUTES WARM-UP.
- 2. TURN TO SET POSITION AND ZERO THE METER WITH THE ZERO CONTROL.
- 3. PLACE THE INSTRUMENT IN A RADIATION FIELD OF 15 MR/HR AND TURN THE SWITCH TO THE X1 POSITION.
- 4. MEASURE AIR TEMPERATURE AND BAROMETRIC PRESSURE. ADJUST METER READ-INGS SHOWN BELOW TO COMPENSATE FOR AIR DENSITY IN IONIZATION CHAMBER. (REFER TO SECTION IV., PARAGRAPH 5).
- 5. REMOVE THE COVER PLUG LABELED X1 ON THE TOP OF THE CASE AND ADJUST THE INTERNAL CONTROL UNTIL THE METER READS 15 MR/HR.
- 6. Place the instrument in a radiation field of 45 MR/HR and check the meter reading. Readjust the calibration control if necessary. Accuracy of the calibration on all three ranges should be such that indications on the meter will not be more than 5% lower nor more than 10% higher than the radiation intensity to which the chamber is exposed if Air density corrections are made.
- 7. Place the instrument in radiation fields of 150, 450, 1500 and 4500 MR/HR. With the range selector switch in the appropriate positions, and using the corresponding calibration controls, calibrate the X10 and X100 ranges using the procedure given above.
- 8. To calibrate the Model HRJ-7 (High Range), use the above procedures and increase the radiation intensities by a factor of five.
- B. INTERNAL ZERO ADJUSTMENT.

IF THE EXTERNAL ZERO CONTROL WILL NOT ZERO THE METER, IT IS AN INDI-CATION THAT THE 1.35 VOLT FILAMENT CELLS (B-2) HAVE DROPPED SLIGHTLY IN VOLTAGE. TO COMPENSATE FOR THIS, AN INTERNAL ZERO ADJUSTMENT IS PROVIDED WHICH IS ACCESSIBLE THROUGH A HOLE IN THE SIDE OF THE CASE AS MENTIONED IN SECTION V1, PARAGRAPH A. THIS CONTROL SHOULD BE ADJUSTED IN THE FOL-LOWING MANNER:

1. TURN THE SELECTOR SWITCH TO THE SET POSITION.

- 2. TURN THE EXTERNAL ZERO CONTROL TO APPROXIMATELY 3/4 OF A REVOLUTION FROM THE EXTREME COUNTER-CLOCKWISE POSITION.
- 3. Adjust the internal Zero adjustment until the meter reads zero.

VI. MAINTENANCE

A. BATTERY REPLACEMENT.

ALL BATTERIES ARE CONTAINED IN A REMOVABLE PLUG-IN BATTERY PACK. TO GAIN ACCESS TO THE PACK, REMOVE THE FOUR SCREWS FROM THE BACK-PLATE (NOT THE BOTTOM PLATE). MOVE THE FILTER CONTROL TABS TO THE UP POSITION AND REMOVE THE END PLATE. SPECIAL ATTENTION SHOULD BE GIVEN TO PLACING THE BATTERIES IN THEIR PROPER POSITIONS. THE INSTRUMENT WILL NOT FUNCTION PROPERLY IF POLARITY IS REVERSED. REFER TO BATTERY PACK SKETCH FOR BATTERY TYPES, LOCATION AND POLARITY.

The 1.35 volt mercury cells (B-2) which supply the tube filament current, will need to be replaced more frequently than the other batteries. These cells must be replaced when it becomes impossible to zero the instrument by means of both the external zero control and the Internal Zero control (accessible through a hole in the side of the case).

B. ALPHA SCREEN REPLACEMENT.

The Alpha screen, which covers the bottom of the chamber, is very delicate. If this screen is touched, or allowed to touch a surface being measured, it may be perforated or torn.

REPLACEMENT SCREENS ARE AVAILABLE FROM TECHNICAL ASSOCIATES AND MAY BE EASILY REPLACED BY FOLLOWING THE PROCEDURE GIVEN BELOW:

- 1. REMOVE THE FOUR NUTS WHICH FORM THE FEET OF THE INSTRUMENT.
- 2. REMOVE THE BOTTOM PLATE AND SLIDE THE TWO FILTERS BACK TO EXPOSE THE ALPHA SCREEN.
- 3. Remove the two small screws on either side of the metal frame of the Alpha screen.
- 4. SLIDE THE FRAME TOWARDS THE FILTERS, THEN TILT IT UP AND REMOVE BY SLIDING TOWARDS THE FRONT OF THE INSTRUMENT.
- 5. THE NEW ALPHA SCREEN CAN NOW BE INSTALLED BY REVERSING THE ABOVE PROCEDURE.

PARTS LIST

	Standard Range Juno SRJ-7	High Range Juno HRJ-7
R 1	4 x 10 ¹¹ Онмз	8 x 10 ¹⁰ Онмз
R 2	4 x 10 ¹⁰ Онмs	8 x 10 ⁹ Онмз
R 3	4 x 10 ⁹ Онмз	8 x 10 ⁸ Он м s
R 4 RESISTOR	470K, 1/2 watt, 10% - Alle	en Bradley Type EB 4741
R 5 RESISTOR	470K, 1/2 watt, 5% - Stema	G TYPE SLAD
R 6 ZERO CONTROL	50 Ohms, 2 watts - Clarost	TAT TYPE CM 10042
R 7 INTERNAL ZERO	250K, 2 watts - Allen Brad	DLEY TYPE JU 2541
R 8 RESISTOR	270K, 1/2 watt, 5% - Stema	G TYPE SLAD
R 9 RESISTOR	82K, 1/2 watt, 5% - Stemag	TYPE SLAD
R 10, R 11, R 12	Calibration Controls - 25K Type CM 10040	COHMS, 2 WATTS - CLAROSTAT
C 1 CAPACITOR	0.1 MFD, 400 volts - Sprag	UE TYPE 4 TM-P1, OR EQUAL
J 1 RECEPTACLE	Amphenol Type 26-183	
P 1 PLUG	Amphenol Type 26-182	
V 1 TUBE	RAYTHEON ELECTROMETER TYPE	CK 5886
M 1 METER	0-20 microamperes - Genera (SRJ-7 Scale calibration 0 (HRJ-7 Scale calibration 0	l Electric Type DO-91 (special) -50 MR/HR, White face) -250 MR/HR, Red face)
S 1A,B,C Switch	GRIGSBY ALLISON TYPE 12260	-4MLW-1 (SPECIAL)
S 1D Switch	TECHNICAL ASSOCIATES TYPE	SW-5
B 1 BATTERY	22 1/2 volts (4 required) VS084, or equal	- Eveready No. 412, or RCA No.
B 2 BATTERY	1.35 volts (2 required) - No. RM-12R, or equal	Eveready No. E-12, or Mallory
B3, B4 BATTERY	6.5 volts (2 required) - M	allory No. TR-115, or equal.
Alpha Screen	TECHNICAL ASSOCIATES TYPE	A-19204-B
HANDLE ASSEMBLY	TECHNICAL ASSOCIATES TYPE	C-19222
BETA FILTER	TECHNICAL ASSOCIATES TYPE	A-1111-A

- 7 -

PARTS LIST (CONTINUED)

Desiccant Cartridge	TECHNICAL ASSOCIATES TYPE 6248		
Desiccant Cartridge Wrench	Technical Associates Type A-23227		
0-Ring	Linear, Inc. Part No. 1	Linear, Inc. Part No. 1820-11	
	Recommended Spare Parts (Based on Yearly Per	s List iods)	
		QUANTITY	
B 1 BATTERIES		4 EACH	
B 2 BATTERIES		QUANTITY DEPENDS ON INSTRU- ment usage - normal life of these batteries is approxi- mately 800 hours.	
B 3, B 4 BATTERIES		Depends on instrument usage – normal life is approximately 3000 hours.	
V 1 TUBE		1 EACH	
Alpha Screen		2 EACH	
Desiccant Cartridge] EACH	



BATTERY PACK



TECHNICAL ASSOCIATES

BATTERY PACK SHAFT SEAL DESICCANT 00 CLAMP 0 Ē Ş ()) ŧ 5) 9)) Ð Ο ଚ ଢ SWITCH BOX > SWITCH BOX BASE ~ SWITCH ~ GASKET INSULATOR ANTENNA ANTENNA CONNECTOR

TECHNICAL ASSOCIATES

JUNO MODEL 7





JUND SURVEY METER, Models SRJ-7 and HRJ-7

Incorporating design features which permit effective operation under high humidity conditions

APPLICATION. The Juno Survey Meter is a portable instrument for measuring the intensity of, and discriminating between alpha, beta, and gamma radiation. It is used to protect personnel from the danger of over-exposure to radiation from radioactive materials or X-rays. While primarily intended for inspection of flat surfaces, the instrument is suitable for most uses where a high degree of accuracy is desired.

The T/A Juno is available in two models: SRJ-7 (standard range) for all normal applications; and HRJ-7 (high range) for use where exceptionally high intensity radiation is likely to be encountered. Both models meet A.E.C. specifications.

Models SRJ-7 and HRJ-7 are improved versions of the original Hanford instrument. The high impedance circuit switch box includes a desiccant cartridge and is sealed with gaskets. These design improvements insure high efficiency performance under adverse humidity conditions.

DESCRIPTION: The instrument comprises an ionization chamber, an electrometer circuit, absorption filters for the rejection of either alpha or beta particles, suitable batteries mounted in a removable power pack, and an indicating meter. The unit is battery operated and is self-contained. To easily identify the High Range Juno, its knobs and meter dials are finished in brilliant red.

ALPHA WINDOW

The ionization chamber has a volume of approximately

27 cubic inches. All surfaces within the chamber are coated with aquadag. The chamber is covered by a screen of .0003" (approximately 0.45 mg/cm²) rubber hydrochloride film. The alpha screen is within 7/16" of any flat surface on which the instrument may be placed, and is easily replaced by simply removing the bottom plate and two retaining screws.

Two absorbers are provided to reject either alpha or beta radiation. These are readily moved in and out of position by means of sliding tabs fitted in rails which form part of the handle. The tab marked "G," with a square end, operates the absorber which rejects alpha and beta, thus permitting a reading of gamma only. The tab marked "B," with a rounded end, operates the absorber which rejects alpha, permitting a reading of beta and gamma. The total of all three types of radiation is read, when both tabs are in "open" position.

The high quality microammeter, which is calibrated in milliroentgens per hour for gamma radiation, has a large easy-to-read face and is mounted in position to permit excellent visibility. Battery life is approximately 800 hours in normal intermittent use. An easily removable battery pack, with simple positive contacts, assures trouble-free operation over long periods of service.



TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA

Models SRJ-7 and HRJ-7

have similar circuits, the principal difference being the value of grid resistances used.

Radiation incident upon the ionization chamber produces a minute current which flows through a very high resistance in the grid circuit of the electrometer tube. The voltage thus produced at the grid causes a corresponding change in plate current which is indicated by the panel meter. A bucking current is provided through the meter in order that the no-signal plate current of the electrometer may be balanced out and readings of radiation intensity may start from the meter zero reading. Sensitivity is varied by switching appropriate values of grid resistance in the electrometer circuit.

The instrument is calibrated by adjusting a resistance in series with the meter. An individual adjustment is provided for each range. Zero setting is accomplished by means of a rheostat in the filament circuit of the electrometer tube. By varying the filament voltage, the plate current may be varied and thus adjusted to a value equal to the bucking current flowing through the meter. All high resistance points in the circuit are insulated with Teflon to insure minimum leakage.

SENSITIVITY. Both models are calibrated in three separate full-scale ranges in easily read increments of the meter scale, covering the total range of which the instrument is capable.

Ranges are based on radium gamma radiation intensity. Accuracy of calibration is such that indications on the meter will not be more than 5% lower nor more than 10% higher than the radiation intensity to which the chamber is exposed. For use in abnormal environmental conditions, air temperature and density correction data and curves are included with each instrument. Sensitivity dependency upon battery aging is limited to a 10% variation while the unit can be zeroed by means of the panel zero control.

Illustration below shows the T/A Juno in use at the Nuclear Radiation Laboratory of Admiral Corporation, monitoring for radio-activity as sample is transferred from "coffin" to shipping container.



SPECIFICATIONS:

IONIZATION CHAMBER:

Volume: 27 cubic inches.

Window Opening: 3" x 45%"

Alpha Window: 0.3 mil (0.45 mg/cm²) rubber hydrochloride.

Alpha Absorber: 0.01" cellulose acetate sheet.

Beta Absorber: 0.102" aluminum.

BATTERIES:

4 Eveready No. 412 22¹/₂ Volt "B" Batteries. 2 Eveready No. E12 1.35 Volt "A" Batteries. 2 Mallory No. TR-115 6.5 Volt "B" Batteries.

TUBE:

1 Sub-Miniature Electrometer Type CK 5886.

RANGES:

Model SRJ-7 Juno (standard range): 50, 500, 5000 MR/HR full-scale. (Improved A.E.C. Model SIC-17B). Model HRJ-7 Juno (high range): 250, 2500, 25,000 MR/HR full-scale. (Improved A.E.C. Model SIC-17D).

TIME CONSTANTS:

50 MR/HR - 18 Seconds500 MR/HR - 4 Seconds5000 MR/HR - 2.5 Seconds

OPERATING TEMPERATURE RANGE: 35°F to 135°F

CASE:

Hard-chrome plated aluminum, with engraved markings on top of case; gasketed, dust and moisture-resistant. Dimensions: $9\frac{1}{2}$ " x $5\frac{3}{4}$ " x 4".

HANDLE:

Aluminum, especially cast low-porosity, smoothly polished.

WEIGHT: Net 6 lbs. 11 oz. Shipping: 10 lbs.



The Juno can be used for X-ray detection and measurement by reference to the above curves.



Jech

JUNO RADIATION SURVEY METER

Models SRJ-6 and HRJ-6



TECHNICAL ASSOCIATES. INSTRUMENTATION FOR NUCLEAR RESEARCH 140 WEST PROVIDENCIA AVENUE . BURBANK . CALIFORNIA

JUNO RADIATION SURVEY METER

Models SRJ-6 and HRJ-6

8

۲

Technical Associates JUNO RADIATION SURVEY METER

Models SRJ-6 and HRJ-6

1. GENERAL DESCRIPTION

The Juno Survey Meter is a portable, battery operated instrument for measuring the intensity of, and distinguishing between alpha, beta and gamma radiation. This instrument is normally supplied in two models: SRJ-6 (standard range) which has a sensitivity range suitable for all routine applications; and HRJ-6 (high range) for use where exceptionally high radiation levels are likely to be encountered. Both models meet A.E.C. specifications.

The instrument comprises an ionization chamber, an electrometer circuit, an indicating meter calibrated directly in MR/HR, a removable plug-in battery pack, and two absorption filters for the rejection of either alpha or beta particles. The Juno is self-contained in a hard-chrome plated aluminum case which can be easily decontaminated when necessary.

To easily identify the High Range Juno (HRJ-6), the range selector switch knob and the meter face are red in color.

11. SPECIFICATIONS

Batteries:

Ionization Chamber: Volume: 27 cubic inches

Window Opening: 3" X 4 5/8"

- Alpha Screen: 0.4 mil (0.68 mg/cm²) rubber hydrochloride, or 0.25 mil (0.89 mg/cm²) mylar film
- Alpha Rejection Absorber: 0.01" cellulose acetate sheet

Beta Rejection: 0.102" aluminum

4 each 22를 Volt "B" Batteries

- 2 each 1.35 Volt "A" Batteries
- 2 each 6.5 Volt "B" Batteries

Tube:	l each Sub-Miniature Electrometer Type CK 5886
Ranges:	Model SRJ-6 Juno (standard range): 50, 500, 5000 MR/HR full-scale (Improved A.E.C. Model SlC-17B). White meter face.
	Model HRJ-6 Juno (high range): 250, 2500, 25,000 MR/HR full-scale. (Improved A.E.C. Model SIC-17D). Red meter face.
Case:	Hard-chrome plated aluminum with engraved markings on top of case; gasketed, dust and moisture resis- tant. Dimensions: $9\frac{1}{2}$ " X 5 3/4" X 4".
Handle:	Aluminum, especially cast low-porosity, smoothly polished.
Weight:	6 lbs. ll oz. net.

III. THEORY OF OPERATION

Refer to Fig. No. 1 for the circuit diagram.

A. Ionization Chamber

The ionization chamber is a rectangular aluminum case. (the inside of which is coated with aquadag) located in the forward end, bottom of the Juno case. Air at atmospheric pressure is the ionization medium. The collector, an aluminum rod coated with aquadag, is extremely well insulated by means of a teflon support bushing. Alpha, beta and/or gamma radiation, when passing through the chamber, collides with individual molecules of air. A collision causes an electron to be displaced from its normal energy level. Once removed sufficiently far from its normal energy level it is subjected to static electric force. This force pulls the free electron to the collector which is 90 volts positive with respect to the chamber walls. The molecule, or atom, after losing an electron assumes a positive potential (such a charged atom is called an ion) and is attracted by the negatively charged walls. Since the collector is well insulated the electrons which land on its surface build up a charge, resulting in a more negative potential. This negative shift in the collector is what is detected by the electrometer circuit.

- 2 -

The reason for choosing aluminum coated with aquadag is to duplicate the performance of a hypothetical air wall chamber as nearly as possible. In order to give a relatively correct reading on all types of radiation and all energies, the solid material of which a chamber is made must nearly duplicate the ionizing properties of air and therefore have approximately the same atomic weight. Aquadag approaches this very closely and aluminum is the most acceptable of those metals which have satisfactory structural properties.

B. Electrometer Tube and Grid Resistors

The electrometer tube is a specially constructed, high vacuum triode. It requires a very low power signal on its grid to control its plate current. The low power input permits the use of an extremely high grid resistor, (R-1, R-2, or R-3). Hence, a very small current through the grid resistor will make a readily observed change in the plate current. In turn, the plate current is conducted through a 20 microampere meter which provides the visual indication of the level of radiation in the chamber. The CK 5886 electrometer tube filament operates on 10 ma., 1.25 volts. The plate voltage is 6.5V and the plate current (with no radiation) is 75 micro amps. The tube is supported on its leads. The base is specially treated for high insulation value in high humidities and will probably be permanently damaged if touched with fingers. The three grid resistors are also specially treated. The housing for the switch and electrometer is at bias potential so that the ionization of air in the tube compartment will not affect the voltage of the electrometer grid.

C. Juno Circuit

B-3 is in series with 82,000 ohms and supplies 75 micro amps to the meter. The electrometer tube with a plate voltage of 6.5volts, and filament voltage of approximately 1.25 volts, has a plate current of

- 3 -

75 micro amps also. These two currents exactly cancel in the meter, providing the instrument is properly zeroed. On increasing the radiation sufficiently to go from zero to full scale, the current through the electrometer tube changes from 75 micro amps to 55 micro amps. It is this difference that is indicated on the meter. The meter is shorted in the OFF and ON positions so as to lessen the chance of damage from shock and vibration. With the Selector Switch in the ON position the filament of the electrometer is allowed to warm up. This increases the life of the tube and reduces Zero drift. With the Selector Switch in the SET position, the bias voltage, filament voltage, plate voltage, and bucking voltage are applied. In addition the grid resistance is shorted out. This permits accurate zeroing of the instrument in the presence of radiation. When the Selector Switch is turned to the xl, xl0, or xl00 position, the ionization current generated in the chamber causes voltage drop across R-1, R-2, or R-3, respectively. Slight variations in these resistances are compensated for by the calibration adjustments R-10, R-11 and R-12, respectively. Zeroing is accomplished by changing the filament voltage. R-6 and R-7 should never be turned fully counter-clockwise as this would place the full potential of B-2 on the sensitive electrometer tube filament and might damage the tube. Care must also be taken not to let either meter lead touch the case of the instrument as this would

place 90 volts across the electrometer filament.

IV. OPERATION

- A. To measure radiation intensities:
 - 1. Turn the Selector Switch to the ON position and wait five seconds for the electrometer tube filament to warm up.
 - 2. Turn the Selector Switch to the SET position and adjust the ZERO control until the meter reads exactly zero.

- 4 -

- 3. Turn the Selector Switch to the "xl" position. If the meter reads off scale the radiation present is greater than 50 mr/hr and the Selector Switch should be turned to either the "xl0" or "xl00" position.
- 4. The instrument is now ready for reading gamma radiation intensity. To read gamma and beta simultaneously, slide back the aluminum screen by pulling up right hand tab "G" in the handle. To read alpha, beta and gamma simultaneously, slide back both the aluminum and acetate filters by pulling up both tabs in the handle.

CAUTION

WITH BOTH FILTERS BACK, THE ALPHA SCREEN CAN BE VERY EASILY DAMAGED. A VERY DELICATE TOUCH MAY CAUSE IT TO RIP.

V. CALIBRATION

A. Internal Adjustments

Calibration. There are three internal calibration controls; one for each range. These controls are accessible, through holes in the top of the case, by removing the cover plugs labeled x1, x10, and x100. These controls should not be changed unless there are calibrated radium or co^{60} sources available for accurate recalibration. The proper procedure for calibrating the Model SRJ-6 (Standard Range) is as follows: 1. Turn the instrument ON and allow several minutes warm-up.

- 2. Turn to the SET position and zero the meter with the Zero Control.
- 3. Place the instrument in a radiation field of 15 MR/HR and turn the switch to the xl position.
- 4. Remove the cover plug labled xl on the top of the case and adjust the internal control until the meter reads 15 MR/HR.
- 5. Place the instrument in a radiation field of 45 MR/HR and check the meter reading. Readjust the calibration control if necessary. Accuracy of the calibration on all three ranges should be such that

- 5 -

indications on the meter will not be more than 5% lower nor more than 10% higher than the radiation intensity to which the chamber is exposed.

- 6. Place the instrument in radiation fields of 150, 450, 1500 and 4500 MR/HR. With the range selector switch in the appropriate positions, and using the corresponding calibration controls, calibrate the x10 and x100 ranges using the procedure given above.
- 7. To calibrate the Model HRJ-6 (High Range), use the above procedures and increase the radiation intensities by a factor of five.
- B. Internal Zero Adjustment

If the external Zero Control will not zero the meter, it is an indication that the 1.35 volt filament cells (B-2) have dropped slightly in voltage. To compensate for this, an internal Zero adjustment is provided which is accessible through a hole in the side of the case as mentioned in Section VI, paragraph A. This control should be adjusted in the following manner:

- 1. Turn the Selector Switch to the SET position.
- 2. Turn the external Zero control to approximately 3/4 of a revolution from the extreme counter-clockwise position.
- 3. Adjust the internal Zero adjustment until the meter reads zero.

VI. MAINTENANCE

A. Battery Replacement

All batteries are contained in a removable plug-in battery pack. To gain access to the pack, remove the four screws from the back, end plate (not the bottom plate). Move the filter control tabs to the up position and remove the end plate. Special attention should be given to placing the batteries in their proper positions. The instrument will not function properly if polarity is reversed. Refer to Figure 2 for battery types, location, and polarity.

- 6 -

The 1.35 volt mercury cells (B-2) which supply the tube filament current, will need to be replaced more frequently than the other batteries. These cells must be replaced when it becomes impossible to zero the instrument by means of both the external Zero control and the Internal Zero control (accessible through a hole in the side of the case).

B. Alpha Screen Replacement

The Alpha screen, which covers the bottom of the chamber, is very delicate. If this screen is touched, or allowed to touch a surface being measured, it may be perforated or torn.

Replacement screens are available from Technical Associates and may be easily replaced by following the procedure given below:

- 1. Remove the four nuts which form the feet of the instrument.
- 2. Remove the bottom plate and slide the two filters back to expose the Alpha screen.
- 3. Remove the two small screws on either side of the metal frame of the Alpha screen. The frame can now be slid towards the filters, and out of the slot which retains one end of the frame.
- 4. Slide the frame until it is out of the slot, then tilt it up and remove by sliding towards the front of the instrument.
- 5. The new Alpha screen can now be installed by reversing the above procedure.

- 7 -

TECHNICAL		ASSOCIATES		
JUNO	MODEL 6	PARTS	LIST	

Standard Range Juno SRJ-6 High Range Juno HRJ-6		
Rl	1.4 x 10 ¹¹ Ohms	4×10^{10} Ohms
R 2	1.4×10^{10} Ohms	4 x 10 ⁹ Ohms
R 3	1.4×10^9 Ohms	4 x 10 ⁸ Ohms
R 4 Resistor	470K, ½ watt, 10% - Allen Bradley Type	e EB 4741
R 5 Resistor	Selected on test (nominal value 680K,	<u>ਤ</u> watt, 10%)
R 6 Zero Control	50 Ohms, 2 watts - Clarostat Type CM 2	10042
R 7 Internal Zero	50 Ohms, 2 watts - Clarostat Type CM 1	10042
R 8 Resistor	Selected on test (nominal value 390K,	$\frac{1}{2}$ watt, 10%)
R 9 Resistor	82K, ½ watt, 5% - Allen Bradley Type H	EB 8235
R 10, R 11, R 12	Calibration Controls - 25 K Ohms, 2 wa	atts - Clarostat Type CM 10040
C 1 Capacitor	0.1 MFD., 400 volts - Sprague Type 4 7	IM-Pl, or equal
J l Receptacle	Amphenol Type 78-545	
J 2 Receptacle	Amphenol Type 26-183	
P l Plug	Amphenol Type 86-CP-48	
P 2 Plug	Amphenol Type 26-182	
V l Tube	Raytheon Electrometer Type CK 5886	
M l Meter	0-20 microamperes - Marion Type MM3 or (SRJ-6 Scale calibration 0-50 MR/HR, W1 (HRJ-6 Scale calibration 0-250 MR/HR, F	r Simpson Type 1327 nite face) Red face)
S 1A,B,C Switch	Grigsby Allison Type 12260-4MLW-1 (spe	ecial)
S 1D Switch	Technical Associates Type SW-3	
B 1 Battery	22 ¹ / ₂ volts (4 required) - Eveready No. or equal.	412, or RCA No. VS084,
B 2 Battery	1.35 volts (2 required) - Eveready No. RM-12R, or e	E-12, or Mallory No. equal.
B 3, B 4 Battery	6.5 volts (2 required) - Mallory No. 7	R-115, or equal.
Alpha Screen	Technical Associates Type A-19204	
Handle Assembly	Technical Associates Type B-1000-B	
Beta Filter	Technical Associates Type A-1111	





FIGURE 2

UNIVERSITY 4-4958

(HICAGO -IN 3-3538

BEN Z. RUBIN CO.

Electronic Manufacturers Representatives

MEYERS RU 16130

BEN Z. RUBIN

16282-046XENNE-AVENUE Detroit 35, Michigan





QUOTATION FROM



TECHNICAL ASSOCIATES

INSTRUMENTATION FOR NUCLEAR RESEARCH

140 WEST PROVIDENCIA AVENUE • BURBANK • CALIFORNIA

TELEPHONES: Victoria 9-1994 • Thornwall 8-8133

To: General Electric Company 4855 Electric Avenue Milwaukee 1, Wisconsin Date: January 23, 1959

Your Inquiry No.: Letter 1/19/59

Our No.: 394

Attention	: J.J. X-Rav	Jech, Coolidge Lab. Department		
ltem	Quantity	Description	Unit Price	Total
1.	l each	Technical Associates' Model SRJ-6 or HRJ-6 Juno Radiation Survey Meter as described in Bulletin No. 159 attached.	295.00	295. 00

Terms: Net 30 days

F.O.B.: Burbank, California

Delivery Schedule: 15 days after receipt of order

- Alexandra

Shipment via: Railway Express

Enclosures: Bulletin No. 159

Copies: (1)

TECHNICAL ASSOCIATES

by

Baron Getman, Asst. Sales Mo



TECHNICAL ASSOCIATES INSTRUMENTATION FOR NUCLEAR RESEARCH 140 WEST PROVIDENCIA AVENUE • BURBANK • CALIFORNIA TELEPHONES: Victoria 9-5838 • Thornwall 8-6649

December 1, 1959

To All Technical Associates' Catalog Holders:

Enclosed herewith are five (5) new pages describing the latest products added to the Technical Associates' line. Included also is our latest price list revision. Please insert these pages and price list in the catalog which, our records show, was previously issued to you.

Bulletin No.	Model No.	Product
167	HSM-10	New Beta/Gamma Hand and Shoe Monitor
	* * *	· ×
168	Р РМ- 8	Beta/Gamma Portal Type Personnel Monitor
	LIM-18	Beta/Gamma Laundry Inspection Monitor
	* * *	* * * * *
169	SA-20	Single Channel Pulse Height Analyzer
	la-6	Non-overloading Linear Amplifier
	* * *	* * * *
170	SM-1 0	Gamma Ray Spectrometer
	SS-3 0	Motor-Driven Spectrum Scanner
	* * *	* * * * *
171		Star Series of Gamma Ray Spectrometry Systems

Please notify us if your catalog has been misplaced and we will send another copy to you.

Very truly yours,

TECHNICAL ASSOCIATES



Beta-Gamma HAND and SHOE MONITOR

Model HSM-10

Completely automatic monitoring. Safe, error-proof, and easy-to-use. Wide range of warning levels. Positive warning of incomplete check. Includes external clothing probe.





No skill is required in the operation of the Model HSM-10 Monitor. The user simply follows the instructions on the illuminated multi-colored panels at the top of the instrument.

When the monitor is ready for testing, the blue "READY FOR USE" panel will be illuminated. To start the operation, the user merely steps on the shoe deck and

inserts hands into the two waist-high probe openings. The counting process is started and maintained by pressure of the finger-tips at the rear of the probe openings, and by the weight of the user on the shoe deck. As soon as the counting starts, the "READY FOR USE" panel darkens, and the yellow "COUNTER IN OPERATION" panel lights up.

NOTE: If at any time during the counting cycle, the user should remove either hand, or step off the shoe deck, the cycle will automatically stop, and the orange panel, reading "CHECK INCOMPLETE - RESET AND RE-PEAT," will be illuminated. This positive warning prevents erroneous readings.

After a short length of time - as preset by the mechanical timer-the "COUNTER IN OPERATION" panel darkens and either the green "CHECK O.K." panel or the red "DECONTAMINATION REQUIRED" panel lights up, depending upon whether the radioactivity present is above or below the preset maximum allowable level. The user then pushes the large red "RESET" pushbutton on the front panel, and the instrument is ready for the next user.

If the "DECONTAMINATION REQUIRED" panel lights up, one or more of the small red panels marked "LEFT HAND," "SHOES," or "RIGHT HAND" will also be illuminated, depending upon the location of the contamination. The degree of contamination can be checked by noting the register readings.

The clothing probe, which is located on the right side of the instrument cabinet, is actuated by a switch located just above the probe holder. When this is done, the white "CLOTHING PROBE IN USE" panel lights up. The output of the clothing probe is read on a count rate meter which is calibrated 0 to 10,000 counts per minute.



TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA
5-DIGIT SODECO REGISTERS

MULTI-COLORED ILLUMINATED INSTRUCTION PANELS

HAND PROBES

GLOW DECADE TUBES FOR SCALING CHANNELS

MECHANICAL TIMER

POWER ON/OFF

SLOTS FOR DISPOSABLE PROTECTIVE PAPER



CLOTHING PROBE RATE METER

SWITCH FOR CLOTHING PROBE

CLOTHING PROBE

AUTOMATIC RESET PUSH-BUTTON

SWITCHES TO VARY THE SCALING FACTORS

SWITCH FOR BACKGROUND TEST AND 60-CYCLE TESTING

TIMER ON/OFF

HIGH VOLTAGE METER

SHOE DECK

FLEXIBILITY

The Model HSM-10 offers an unusually wide range of operating limits and warning levels, and thus can serve any health physics requirement. Use of three switches in each channel provides as many as 15 different scaling factors. Time intervals on the mechanical timer are variable from 1 to 120 seconds.

CONSTRUCTION: The Model HSM-10 is housed in a metal cabinet with a smooth grey hammertone finish. A door, with lock, covers all controls, protecting the instrument from misadjustment by unauthorized personnel. Adjustments, tube replacements, and routine servicing can be done without removing the chassis from the cabinet.

SPECIFICATIONS:

SENSITIVITY: Minimum Beta energy 0.2 MEV. TUBES: T/A Type T1100 halogen-quenched.

ALARM SETTINGS: Hand Channels - 100 to 10,000 counts in

15 steps. Shoe Channel – 1000 to 100,000 counts in 15 steps.

COUNTING PERIOD: Variable from 1 to 120 seconds.

REGISTER RANGE: 0 to 99999, five digits on each of three Sodeco Registers.

WARNING LEVEL INDICATOR: Adjustable by the Health Physicist.

POWER SUPPLY REQUIRED: 300 watts, 95-125 volts 60 cycle A.C. **DIMENSIONS:** Cabinet: 65" high x 27" wide x 18" deep.

Foot Deck: 25" long x 18" wide x 4" high.

SHIPPING WEIGHT: 650 lbs.



TECHNICAL ASSOCIATES

140 West Providencia Ave. • Burbank, California



BETA-GAMMA PERSONNEL and LAUNDRY MONITORS Model PPM-8 Portal Type • Model LIM-18 Laundry Type



- Alarm lights on portal and meter indicators on console panel provide double check on location of contamination.
- Lead shielding around individual detectors minimizes background.
- Watertight threshold with detachable plastic foot mat for easy decontamination.
 - Audible (buzzer) and visual (red light) alarms.
- Eight separate counting rate circuit channels with individual alarm settings.
- Single button to reset all channels.
- Single H. V. power supply for all channels.
- Easy-to-service console cabinet.

APPLICATION: Model PPM-8 Portal Monitor provides a quick, efficient, economical method of "head to toe" monitoring of personnel entering or leaving an area. To prevent the spread of contamination, the Portal Monitor is placed at a control point so that personnel entering or leaving the controlled area pass through the portal. Should the contamination on any person passing through the portal exceed a preset radiation level, immediate detection and alarm will result.

OPERATION: No skill is required by personnel being monitored. As person passes through the portal, any detected change in the level of radiation at 8 different body areas is indicated on a corresponding meter. An audible



SEVEN DETECTORS are located around the sides and top of portal to check head, shoulders, waist, and legs; and four detectors are grouped in the threshold to check the soles and heels of shoes.

> The insert at right shows the lead shielding behind each individual detector (T/A #1120 Tube) in portal. This shielding minimizes the background count.



alarm sounds when the preset radiation level is exceeded, and indicator lights on *both* portal.and console panel glow "RED" to indicate the exact spots of contamination. The inside dimension of the portal is purposely narrow to prevent personnel from moving through too rapidly to complete a good check; however, passage can usually be made at 2-second intervals.



TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE • BURBANK, CALIFORNIA

PECIFICATIONS

of PPM-8 PORTAL MONITOR

RANGE: Portal Channels (7), 0-1000 CPM full scale. Threshold Channel (1), 0-2000 CPM full scale.

RESPONSE TIME: 2 to 4 seconds continuously variable.

ALARM SETTING: Minimum 1.5 x background reading.

SENSITIVITY: Less than 0.15 μ c (Beta/Gamma) on body surfaces. Less than 2 μ c (Beta) on soles of shoes.

DETECTORS: Portal Frame - 7 T/A No. T-1120 Halogen-quenched tubes. Threshold -- 4 T/A No. T-1100 Halogen-quenched tubes.

REMOTE ALARM RESET: Back panel mounted jack for remote reset accessory.

POWER SUPPLY REQUIRED: 110 volt 60 cycle 50 watts.

DIMENSIONS: Portal Frame (Outside) – 20" wide x 15" deep x 84" high. (Inside) – 16" wide x 15" deep x 77" high. Console Cabinet – 22" wide x 15" deep x 11" high.



FINISH: Portal and Console Cabinet - Hammertone grey. Console Panel - Flat dark grey.

WEIGHT: Portal, 98 lbs. Console, 54 lbs. Total packed for shipment, 200 lbs.

T/A Easy-to-Service CONSOLE CABINET

Top section raised to service tubes and upper chassis.

Complete chassis and panel raised to service circuit boards and under-chassis components.

LIM-18 LAUNDRY INSPECTION MONITOR

This is a simplified, low-cost instrument for monitoring clothing and equipment worn in Beta/Gamma contaminated areas. After being washed and prior to re-issue, the material is passed through the monitoring frame where it is scanned by six lead-shielded detectors. Laundered material having residual contamination in excess of permissible preset level will cause an immediate alarm at the Single Channel Radiation Console. A T/A Model P-7 Hand Probe is standard equipment with the Model LIM-18 Monitor for the purpose of localizing detected contamination.

SPECIFICATIONS

RANGE: 0-2000 CPM full scale. **RESPONSE TIME:** 2 to 4 seconds. ALARM SETTING: Minimum 1.5 x background reading. DETECTORS: 6 T/A No. T-1120 Halogen-quenched tubes in frame. 1 Victoreen 1885 G/M tube in hand probe. SENSITIVITY: Less than 0.15 micro-curies (Beta/Gamma). POWER SUPPLY REQUIRED: 110 volts 60 cycle 35 watts. DIMENSIONS: Frame approximately 44" wide x 4" high inside. (Can be varied within reasonable limits.) WEIGHT: Total shipping - 175 lbs.



TECHNICAL ASSOCIATES

140 West Providencia Avenue • Burbank, California



TECHNICAL ASSOCIATES INSTRUMENTATION FOR NUCLEAR RESEARCH 140 WEST PROVIDENCIA AVENUE • BURBANK • CALIFORNIA TELEPHONES: Victoria 9-5838 • Thornwall 8-6649

March 1, 1959

To All Technical Associates' Catalog Holders:

Herewith four (4) new pages and our latest price list revision for insertion in the catalog which, our records show, was previously issued to you:

<u>Bulletin</u>	Model	Product	Replaces
164	GS-7	New Binary Scaler	New Product
166	DS - 5B	New Decade Scaler	Bulletin 140
144 (revised)	IC-1	Ionization Chamber and Counter Tubes	Bulletin 144
165	CP-3	New Cutie Pie Survey Meter	Bulletin 147

Our latest price list revision includes changes in all personnel monitors and prices for new instruments listed above.

Please notify us if your catalog has been misplaced and we will send another copy to you.

Very truly yours,

TECHNICAL ASSOCIATES

PRICE LIST

TECHNICAL ASSOCIATES • 140 W. Providencia Ave. • Burbank, Calif.

Revised March 1, 1959 Subject to Change Without Notice

PERSONNEL MONITORS

ALM-1X	Alpha Hand Monitor (Argonne Count Rate Meter Design) with external Alpha Probe	\$2,875.00
ALM-2X	Beta/Gamma Hand and Foot Monitor (Argonne Count Rate Meter Design) with External Beta/Gamma Probe	2,975.00
AHC-2 HSC-1	Alpha Hand Counter (Hanford Binary Scaler Type) Beta/Gamma Hand and Shoe Counter (Hanford Binary Scaler Type)	2,675.00 2,650.00
	Note: Add \$75.00 per instrument for 50 cycle, 115 volt.	

ANALYTICAL INSTRUMENTS

* DS-5B	Decade Scaling Unit (5 microseconds)	835.0 0
* DS-58A	Decade Scaling Unit (1 microsecond)	885.00
DU-1	One microsecond plug-in decade unit	107.50
DU-5	Five microsecond plug-in decade unit	57.50
* *CR-2	Count Register used with DS-5B or DS-5BA	65.00
PT-1	Pre-determined Timer (Long Counting Intervals—999 Minutes)	95.00
**PT-2	Pre-determined Timer (Short Counting Intervals—120 Seconds)	95.00
* * PT-3	Pre-determined Timer, Liebel-Flarsheim (Range 1 second to 60 minutes)	110.00
* GS-7	Binary Scaling Unit (Scale of 256)	655.00
* RM-7	Dual Linear and Log Count Rate Meter	595.00
RM-7C	Dual Linear and Log Count Rate Meter with special time constants	610. 0 0
* SM-10	Medical Spectrometer (Single Channel Pulse Height Analyzer, with	
	built-in Linear Amplifier and Ratemeter)	675.00
* SA-20	Analytical Spectrometer (Single Channel Pulse Height Analyzer)	595.00
AP-20	Alpha Poppy Laboratory Monitoring Instrument—rack mounted 83/4" x 19"	550.00
* LA-5a	Linear Amplifier (specify delay line 1 or 2 microseconds) cable lengths to 10 ft	535.00
* LA-5b	Linear Amplifier (specify delay line 1 or 2 microseconds) cable lengths 10' to 100'	565.00
PA-5	Pre-Amplifier for LA-5 with cables	115.00
RHV-1	High Voltage Power Supply — rack mounted 51/4" x 19"	325.00

*Prices on all analytical instruments are for cabinet mounted models, unless otherwise noted. If a rack mounted instrument is desired, deduct \$10.00 from the price and add the letter "R" to the instrument designation.

**These accessories will be supplied in a cabinet or panel mounted at customer's option.

DETECTORS, PROBES AND ACCESSORIES

SD-1	Scintillation Detector with 1" x 1" Nal crystal, complete with cables	455.00
SD-2	Scintillation Detector with 2" x 2" Nal crystal, complete with cables	765.00
SD-2W	Well-type Scintillation Detector with $13/4$ x 2" Nal crystal (5%" x $11/2$ " well) complete with cables for use with LS-8 or LS-8X Shields	695.00
SD-3	Scintillation Detector with 11/2" Alpha Phosphor, complete with cables	495.00
SD-4	Scintillation Detector with 11/2" Beta Anthracene Crystal, complete with cables	545.00
\$D-12	Scintillation Detector for Analytical Spectrometry with 2" x 2" Nal crystal, complete with cables	825.00
SD-12₩	Well-type Scintillation Detector for Analytical Spectrometry, with $13/4$ x 2" Nal crystal (5%" x $11/2$ " well), complete with cables for use with LS-8 or LS-8X Shields	755.00
CH-1	Crystal Housing for SD-1 Scintillation Detector	25.00
CH-2	Crystal Housing for SD-2 or SD-12 Scintillation Detector	30.00
CH-2W	Crystal Housing for SD-2W and SD-12W Scintillation Detectors	30.00
NS-1	Nose Shield for Model SD-1 Detector	95.00
NS-2	Nose Shield for Model SD-2 Detector	115.00
C-1a	Collimator Type A (20° Flat Field) for SD-1 Detector	29.50
C-20	Collimator Type A (20° Flat Field) for SD-2 Detector	39.50
С-1Ь	Collimator Type B (Straight Bore) for SD-1 Detector	35.00

PRICE LIST (Continued)

r

	C-2b	Collimator Type B (Straight Bore) for SD-2 Detector	45.00
	C-1c	Collimator Type C (Focusing) for SD-1 Detector	75.00
	C-2c	Collimator Type C (Honeycomb) for SD-2 Detector	95.00
	AS-11	Alpha Detector (1") for SD-1, complete with mount	79.50
	AS-12	Alpha Detector $(1\frac{1}{2})$ for SD-2, complete with mount	85.00
	BS-11	Beta Detector (1") for SD-1, complete with mount	79.50
	BS-12	Beta Detector $(1\frac{1}{2})$ for SD-2, complete with mount	95.00
	MS-1	Mobile Stand for Scintillation Detector, counter weighted and with shelf for Scaler or Rate Meter	295.00
	DL-1	Distance Locator (Plastic disc with 10 c.m. rubber pointer) for use with NS-1 Nose Shield	10.00
	DL-2	Distance Locator (Plastic disc with 10 c.m. rubber pointer) for use with NS-2 Nose Shield	10.00
	P-7	Beta/Gamma Probe Assembly with 1885 Tube, cable and connector	39.50
	P-8	End Window Probe (less tube) complete with cable and connector	45.00
	P-8X	End Window Probe Assembly with T/A T-1180 Tube, cable and connector	75.00
	P-9	Alpha Probe for use with AP-20	135.00
IONIZ	ATION CHA	MBERS	
	IC-1	Ionization Chamber with windows (Beta/Gamma)	115.00
	IC-2	Ionization Chamber—Solid Wall (Gamma)	105.00
SURV	EY METERS	· · · ·	
50101	CP-3	New Cutie Pie Portable Survey Meter	265.00
	CP-3A	New Cutie Pie Portable Survey Meter with special ranges	
		(0-25, 0-250, 0-2500 mr/hr)	295.00
	SRJ-6	New Juno Survey Meter with plug-in battery pack and replaceable alpha window (Standard Range)	295.00
	HRJ-6	New Juno Survey Meter with plug-in battery pack and replaceable alpha window (High Range)	295.00
LEAD	SHIELDS		
	LS-1	Type TM-1 Tube Mount	185.00
	LS-2A	Standard with Type TM-6 Tube Mount and Sample Tray Holder, Type ST-2A Sample Tray and 10 Type PL-2 Aluminum Planchets	200.00
	LS-2B	Micrometric with 10 Type PL-1 Planchets	285.00
	LS-4A	For liquids complete with Type TM-1 Tube Mount, Type SK-1 Socket, and one Type ST-5 Marinelli Beaker	200.00
	LS-4B	Modified to accommodate special tubes with aluminum inner liner	235.00
	LS-4C	For solids, complete with Type TM-1 Tube Mount, Type SK-1 Socket, one Type ST-6 Ore Container and 100 Type ST-6A Paper Sleeves	225.00
· .	LS-5	Complete with Type TM-3 Tube Mount, Type SK-2 Socket, Type STH-1 Sample Tray Holder, Type ST-4A Sample Tray, and 10 Type PL-2 Aluminum Planchets	260.00
	LS-6	Complete with Type TM-7 Tube Mount, one Type ST-3B Sample Tray and 10 Type PL-3 Aluminum Planchets	275.00
	LS-6X	Eight-inch height extension for LS-6 Shield	98.50
• .:	LS-6L	Lid with 3" opening and retaining collar for Models SD-1, SD-2, SD-12 Scintillation Detectors	72.50
	LS-6P	Plug for 3" opening in LS-6L	12.50
	LS-7A	Multi-Purpose Type complete with SD-1 Scintillation Detector, three ST-7 Sample Travs and 10 PL-3 Planchets	720.00
. •	LS-7B	Multi-Purpose Type complete with SD-2 Scintillation Detector, three ST-7 Sample Trays and 10 PL-3 Planchets	1,030.00
	LS-7C	Multi-Purpose Type complete with TM-8 Tube Mount, three ST-7 Sample Trays and 10 PL-3 Planchets	325.00
	LS-7M	Multi-Purpose Type—Shield only, with adapter ring, three ST-7 Sample Trays and 10 PL-3 Planchets	265.00
· · .	LS-8	For well counting (shield only, less detector)	225.00

.

PRICE LIST (Continued) TECHNICAL ASSOCIATES • 140 W. Providencia Ave. • Burbank, Calif.

	LS-8W	For well counting, complete with SD-2W Scintillation		920.00
	10.07	For well equating with 2" chielding in all directions (Shield only less	Detector)	290.00
	LS-8X	For well counting with 2" chielding in all directions complete with		270.00
	L2-8WX	For well counting with 2 shielding in an unections complete with SD_2W scintillation Detector and $134'' \times 2''$ well crystal		985.00
COUN	TER TUBES			
	T-1100	Thin wall type tube 3" long used in Model ALM-2, ALM-2X and HSC	-1	20.00
	T-1120	Thin wall type tube 7" long used in Models ALM-2, ALM-2X and HS	C-1	35.00
	T-1140	Mica End Window, used in TM-7, TM-8 Tube Mounts and P-8 Probe		47.50
	T-1160	Mica End Window (Organic Quenched) used in TM-7, TM-8 Tube Mount and P-8 Probe	ts	47.50
	T-1180	Mica End Window Alpha/Beta/Gamma sensitive used in P-8X Probe		50.00
LAROI	RATORY AC	CESSORIES		
	AB-2	Absorber Set		85.00
	ST-1	Sample Tray	· · · · · · · · · · · · · · · · · · ·	6.50
	ST-2 ST-3 ST	T-4 Sample Trays		3.50
	CT 5	Marinelli Beaker		4.25
	ST-5	Are Container		12.50
	51-0 ST (A	Paras Slaves for Ore Container	Der C	2.00
• • • •	ST-0A	Cample Trace for 1 5 7	per c	2 50
1	51-/	Sample (ray for LS-/		13.50
	STH-T	Sample Tray Holder with tray and ten planchels		5.00
	SK-1, SK-2, S	SK-3 SOCKETS		15.00
	TM-1, 1M-2,	TM-3, TM-4 Tube mount		100.00
	TM-5	Tube Mount (micrometric for LS-2 Lead Shield)		22 50
	TM-6	Tube Mount and Sample Tray Holder		22.30
	TM-7	Tube Mount and Sample Tray Holder		27.50
	TM-8	Tube Mount and Adapting Ring for LS-7 Multi-Purpose Shield		60.00
	LC-1	Lead Source Container for four radium needles		45.00
	LC-3	Lead Source Container 1" lead walls, $1\frac{1}{2}$ " x 3" inside dimensions		40.00
	LC-3A	Lead Source Container 13%" lead walls, $34'' \times 214''$ inside dimension	15	47.50
	LC-4	Lead Source Container $1\frac{1}{2}$ " lead walls, $1\frac{1}{2}$ " x 3" inside dimensions		50.00
	LC-4A	Lead Source Container $134''$ lead walls, $34'' \times 2\frac{1}{2}''$ inside dimensions	• • • • • • • • • • • • • • • • • • • •	57.50
	LC-9	Mobile Source Unit — compartment 3" x 6" in shield of 2" lead.		
		Mounted on heavy duty hand cart		225.00
Lead Br	ick, machined al	Il sides 2'' x 3'' x 6''		8.50
Lead Br	ick, machined al	II sides 2'' × 4'' × 8''		12.50
Lead Br	ick, extruded sh	ape	Prices upon	application
Planche	ts, Type PL-1, S	heet steel, tinned	per C	3.25
			per M	30.00
			5M	27.00/M
Planche	ets, Type PL-2 &	PL-3, aluminum	per C	3.25
			per M 5M	27.00/M
<i>.</i>		P. DI 20	Der C	4 50
Planche	ets, Type PL-2C a	& PL-3C, copper	Der M	40.00
			per 5M	36.00/M
Planche	ate PL-25 stainle	ess steel	per C	9.00
(fanch			per M	65.00
			per 5M	58.50/M
Planche	ets, Type PL-3S,	stainless steel	per C	12.00
	-		per M	85.00
			per 5M	76.50/M
Planche	ets, Type PL-4, c	aluminum	per C	5.00
			per M	45.00
.		· · · · · · · · · · · · · · · · · · ·	tions quated .	
Planch	ets, special mate	Priors and sizesPriors and sizes	nces quotea L	ipon request

IN GENERAL

Unit prices are given in this price list, unless otherwise noted.

All prices quoted herein are f.o.b. Burbank, California. Payment terms, net 30 days upon approved credit.

No extra charge is made for packing for domestic shipments. Export charges will be quoted upon application.

We reserve the right to alter specifications at any time without incurring the obligation of incorporating new features in previously manufactured equipment.

Prices are subject to change without notice. Quotations remain firm for 30 days.

Shipping instructions: Please include shipping instructions when ordering; in the event that shipping instructions are not given, we will use our best judgment in the matter.

SERVICE CHARGES

Service charges for repairs are usually billed on an hourly basis, however, if desired, an estimate can be given before the work is undertaken. The unit to be repaired should be shipped prepaid to us.

WARRANTY

Technical Associates warrants instruments and equipment (except tubes, fuses, batteries and crystals), manufactured

by it to be free from defects in workmanship or materials under normal use for a period of one year from the date of shipment from the factory to the buyer. Tubes, fuses, batteries and crystals are subject to the guarantee established by the manufacturer of them, however, Technical Associates will assist the customer to obtain full benefits of these guarantees.

If, within the one year warranty period, any Technical Associates' instrumentation or equipment requires service as a result of a defect, the buyer may return it to the factory of Technical Associates at Burbank, California or to a service station designated by Technical Associates, transportation charges prepaid, for service at no charge under the warranty. The buyer is urged to communicate with Technical Associates when warranty service is required, stating the nature of the difficulty and giving model and serial number of instrument. It may be possible to diagnose the trouble and send a replacement part or assembly, thereby avoiding the expense of shipment.

Technical Associates will return the instrument to the buyer, transportation charges prepaid, after repairs or replacement under warranty are completed. The liability to Technical Associates under this warranty is limited to the cost of replacement of defective parts upon prompt notification of such defect.



CUTIE PIE, MODEL CP-3

PORTABLE ALPHA, BETA, and GAMMA SURVEY METER



- 3 ranges: 0-50, 500, or 5000 mr/hr.
- Selector switch positions permit testing of all batteries before use.
- New battery pack provides over 800 hours of operating life.
- Thin end window permits alpha detection.

APPLICATION: The Model CP-3 Cutie Pie Survey Meter is an improved version of the Model CP-2. This portable instrument measures and distinguishes between alpha, beta, and gamma radiation. Designed for laboratory, reactor, and industrial use, it has found wide acceptance wherever radioisotopes are handled. The Model CP-3 is especially useful in health physics work to monitor such inaccessible spots as corners, behind pipes, beams, etc.

Three full-scale ranges of 50, 500, and 5000 milliroentgens per hour provide excellent coverage of radiation levels encountered in most laboratories. The tolerance rate of 7.5 mr/hr for a 40-hour week is easily read on the 50 mr/hr range. To assure maximum reliability, the range selector switch includes 3 test positions for checking batteries prior to use.

DESCRIPTION: The Model CP-3 Cutie Pie is a gun-type survey meter comprising an ionization chamber, an electrometer circuit, absorption filters for the rejection of alpha or beta particles, an indicating meter, and a battery complement identical to that of the Juno Radiation Survey Meter, Models SRJ-6 or HRJ-6. The undesired radiations are easily rejected by swinging the proper absorber into place in the absorber bracket mounted on the front of the ionization chamber.

The Model CP-3 has a bright, chromeplated all-aluminum ionization chamber and case. The end of the chamber is closed with a rubber hydrochloride alpha screen .0003" thick (0.45 mg/ cm²), which is frame-mounted and held in place by clips, thus permitting its easy replacement. The case is chrome-plated aluminum for easy cleaning and decontamination. All edges are rounded. Rubber gaskets provide protection against high humidity conditions.

The T/A Model CP-3 incorporates circuit improvements and battery pack which extends the battery life to more than 800 operating hours and permits the testing of all batteries prior to use. This is done by turning the selector switch to each of the 3 battery test

positions. Batteries of proper voltage produce a reading in the green sector on the meter. A set position permits the meter to be adjusted to read zero even in a radiation field. The remaining 3 switch positions permit the selection of radiation measuring ranges.





CUTIE PIE

SERIAL NO

TECHNICAL ASSOCIATES

Cutie Pie is shown testing radioactivity of vacuum melted steel.

TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE • BURBANK, CALIFORNIA

PECIFICATIONS

RADIATION RANGES: 50, 500, and 5000 mr/hr full scale.

- CALIBRATION: Factory calibrated, using gamma standard calibrated by the National Bureau of Standards. Calibration accuracy $\pm 10\%$. Internal individual calibration controls for each range, screwdriver adjusted from outside of case. Access is permitted by means of snap plugs.
- CIRCUIT: Reliable single tube electrometer circuit. Ranges of sensitivity are obtained by switching input grid resistors. This is accomplished by a special Teflon insulated switch. All high resistance points in the circuit are insulated with Teflon or Kel-F to insure minimum leakage. Ranges and corresponding input resistance are:

a oonooponani,	
50 mr/hr	1.4 x 10 ¹¹ ohms
500 mr/hr	1.4 x 10 ¹⁰ ohms
5000 mr/hr	1.4 x 10 ⁹ ohms
anntaine five	individual notantiamators

The circuit contains five individual potentiometers for purposes of zeroing and calibration.

- METER: High quality 31/2" meter, scale length 2.36", with 50 scale divisions. Appropriately calibrated to read in milliroentgens per hour for gamma radiation. Mounted in position to permit excellent visibility.
- CONTROLS: Single control switches meter to zero position, battery test points, and 3 operating ranges. Meter Zero Control is located directly below the meter and zeros the meter by adjusting the filament voltage. Coarse Meter Zero Control and Calibration Control are screwdriver adjustments located beneath snap plugs on sides of instrument case and are provided for use only after battery voltages have dropped considerably.
- TIME CONSTANTS: Range 50 mr/hr 6 seconds Range 500 mr/hr less than 1 second Range 5000 mr/hr less than 1 second

ZERO DRIFT: Negligible after 15 minutes warm-up period.

BATTERY LIFE: Over 800 operating hours.

IONIZATION CHAMBER:

Aluminum Cylinder: 27/8" inside diameter, 63/8" long. Volume: Approximately 36 cubic inches. Window opening: 23/4" in diameter. Alpha Screen: Easily removable, ring-mounted rubber hydrochloride .0003" thick (0.45 mg/cm²). Alpha Rejection Absorber: Hinge held, ring-mounted cellulose acetate

0.01" thick. Beta Rejection Absorber: Hinge held, aluminum disc 0.102" thick.

- CASE: Formed and welded sheet aluminum with chrome finish and large, clearly engraved markings.
- HANDLE: Cast aluminum, free from porosity and highly polished. Hollow construction to reduce weight.

BATTERY COMPLEMENT:

4 Eveready No. 412 221/2 Volt "B" Batteries 2 Eveready No. E12 1.35 Volt "A" Batteries 2 Mallory No. TR-115 6.5 Volt "B" Batteries

VACUUM TUBE: Raytheon Type 5886.

WEIGHT: 4 lbs. 12 oz. net. Shipping Weight: 8 lbs.

A special Model CP-3A is available at slightly higher cost, for applications where higher sensitivity is required. This instrument is equipped with ranges of 0-25, 0-250, 0-2500 mr/hr.



TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE **BURBANK, CALIFORNIA**





MODEL CP-3 CUTIE PIE. with front tripod feet in position for bench or table use.



MODEL GS-7 BINARY SCALER

- with automatic electric reset
- predetermined count
- elasped time indicator
- binary scale of 256
- 300-3000 volt detector supply



APPLICATION and SCOPE OF USE...

The Model GS-7 Binary Scaler is a general purpose, precision counting instrument which features extreme simplicity in operation, as well as unusual flexibility in counting procedures. It is designed for use with Geiger, scintillation, or proportional detectors, and has many applications in diagnostic and therapeutic procedures in nuclear medicine, radiochemical studies, and other laboratory counting problems. Exceptional reliability is assured through the use of newly developed components of the highest quality. The GS-7 has a resolving time of 5 microseconds. For counting problems requiring shorter resolving time, the GS-7A, with 1 microsecond resolution, is available.

OPERATION: The GS-7 uses 8 binary scaling stages, providing an electronic scaling factor up to 256. A scale selector switch allows a wide choice of scaling speeds. The total recorded count is determined by multiplying the reading on the five-digit register by the scaling factor, and adding the values of the



with scanner for thyroid "mapping."

interpolation lights which are illuminated. Built-in features include an elapsed time indicator, which gives the exact amount of time elapsed during any counting run. The GS-7 is available with either a timer reading in seconds, or in hundredths of a minute. Either timer indicates a total elapsed time of 60 minutes.

An outstanding feature of the GS-7 is the automatic electric reset button, which simultaneously resets the register, timer, and binary scale indicators in one single operation. Laboratory technicians will quickly recognize the tremendous value of this feature in simplifying operation and assuring faster, more efficient counting procedures.

The GS-7 incorporates an accurately regulated 300-3000 volts D.C. detector supply, with extremely accurate settings accomplished by "coarse" and "fine" adjustments. A panel switch turns the high voltage on and off, without changing the voltage setting.

Predetermined count selection is provided by a switch which enables the user to automatically stop the scaling action after 1, 10, or 100 register counts. When the preset count is reached, the scaler stops, and the time required to accumulate the count is read on the timer. The GS-7 may also be coupled with such external accessories as the Predetermined Timers PT-1, PT-2, or PT-3, to gain the added advantage of preselected counting time.

TECHNICAL ASSOCIATES

for use with

Technical Associates Scalers



PREDETERMINED TIMER, Model PT-1

Provides preset time control for GS-7 Binary Scaler from 1 minute to 999 minutes, with 1 minute increments. For the convenience of the operator, an indicator light goes off when counting period ends. ACCURACY of any selected time

interval is $\pm 0.2\%$. SIZE: 5" x $3\frac{1}{2}$ " x 5". WEIGHT: 2 lbs.

Pre

PREDETERMINED TIMER, Model PT-2

For use where short counting intervals prevail. Provides pre-set time control from 0 to 120 seconds, with 1 second increments. Available in other ranges and dial divisions upon request. ACCURACY: Overail ± ½

of 1% of full scale. Repeat ± ¼ of 1% of full scale. SIZE: 8" x 8" x 8". WEIGHT: 5 lbs.

PREDETERMINED TIMER, Model PT-3



Contains Liebel-Flarsheim Model 2-D Nuclear Timer. Combines pre-set and elapsed time with timing range from 1 second minimum to 60 minutes maximum, in one second divisions. ACCURACY: $\pm 1/10$ th second

> SIZE: 8" High 8" Wide 8" Deep WEIGHT:

> > 6 lbs.

pecifications

of GS-7 Binary Scaler

- **RANGE:** Diode coupled binary scale of 256, with scale selection of 4, 8, 16, 32, 64, 128, and 256, followed by a 5-digit Sodeco Register providing an aggregate count capability of 6400 cps.
- **ELAPSED TIMER:** New Cramer #691, clock face type, with 60 minute total elapsed time and accuracy of .01 second. Choice of increments of one second or one hundredth of a minute, as ordered.
- **RESOLVING TIME:** 5 microseconds on standard GS-7. 1 microsecond resolution time on Model GS-7A.
- HIGH VOLTAGE POWER SUPPLY: Variable from 300 to 3000 volts by means of "Coarse" and "Fine" Controls.
- HIGH VOLTAGE METER RANGES: 300 to 1500 volts full-scale; 600 to 3000 volts full-scale.
- **REGULATION:** Less than .002% change in high voltage for 1% line voltage change between 95 and 130 volts A.C.
- **ELECTRIC RESET:** High speed electrical reset instantly resets register, timer, and binary scale indicators by simply depressing a single push-button.
- **PREDETERMINED COUNT:** Stops count register at 1, 10 or 100, and also stops elapsed time indicator.

DN REAR PANEL: Latch-lock connectors for PT-1, PT-2 or PT-3 Predetermined Timers, and SD-1 or SD-2 Scintillation Detectors. Signal input and high voltage output connector. Test signal switch, permitting check of scaling stages. Oscilloscope connector.

INPUT SENSITIVITY: Variable from 0.15 volts to 3.0 volts negative pulse, with chassis-mounted control. Factory adjusted for 0.25 volts.

POWER REQUIREMENTS: 100 to 130 volts, 60 cycle 180 watts. Available for 50 cycle, 110 volts or 220 volts, on request.

CONNECTORS:

Signal Input (front and rear)	UG-931/U
Preamplifier	Cannon XL-3-13
External Register	Cannon XL-4-13
External Timer	Cannon P5-13
Scope Output	Switchcraft #11
Power Input	Amphenol 61-M10

DIMENSIONS: In cabinet: 201/2"W x 103/4"H x 133/4"D (Model GS-7). For rack mounting: 19"W x 83/4"H x 13"D (Model GS-7R).

WEIGHT: 35 lbs.

FINISH: Smooth grey hammertone.



GS-7 AND TIMER AVAILABLE AS A COMPLETE UNIT

All of the timers shown can be supplied on a 19" panel for mounting with the GS-7 Binary Scaler in a single cabinet, at customer's option.

TECHNICAL ASSOCIATES



DECADE SCALER Model DS-5B and Model DS-5BA



Built-in elapsed time indicator. Choice of 1 or 5 microsecond resolving time.

APPLICATION. Technical Associates Models DS-5B and DS-5BA Decade Scalers are precision counting instruments for use with Geiger, Scintillation, or Proportional Detectors, and have been carefully designed to combine exceptional versatility with extreme reliability and simplicity of operation. Model DS-5B provides 5 microsecond resolving time; Model DS-5BA provides 1 microsecond resolving time for higher counting rates. In addition to operation as a conventional scaler, they may be used with a printing medical scanner or with printing timers or similar external accessories.

DECADE UNITS. Five Model DU-5 decades provide a count capability of 99,999. By use of the Model CR-2 External Register, this may be extended to 1 x 10⁸ counts. A special Model DU-1 one-microsecond decade is used as the first decade in Model DS-5BA. HIGH VOLTAGE DETECTOR SUPPLY. Built-in high voltage detector supply provides 300 to 3000 volts of exceptionally stable direct current. Desired voltages are obtained by a dual voltage range selector switch in conjunction with coarse and fine controls. PREDETERMINED COUNT. A single control provides choice of ten predetermined counts (4-10-40-100-400-1K-4K-10K-40K-100K) or manual operation for counting against time.

ELAPSED TIME INDICATOR. The time elapsed during any counting operation is registered on a five digit (odometer type) indicator which will indicate up to 999.99 minutes of elapsed time. Elapsed time may be accurately read to hundredths of a minute. PRESET TIME. A connector is provided at rear of instrument to accept T/A Models PT-1, PT-2, or PT-3 Predetermined Timers, for applications where the total number of counts for a specific period of time is desired.

PUSH-BUTTON RESET. Resetting the entire instrument is accomplished by the operation of a single push-button, which electrically resets the 5 decades and the elapsed time indicator, as well as external accessories. A two-position stop-count lever type switch is provided for manual control of the instrument.

INSTRUMENT TEST. A 60 cycle test switch permits checking the instrument's operation with the 3600 cpm line frequency.

AUTO-NORMAL CONTROL. When the Auto-Normal switch is in the normal position, the Model DS-5B or DS-5BA operates as a standard scaler and offers predetermined count and manual operation. When the switch is in auto position, the instrument may be used with printing medical scanners in making body scans, or with print-out timers in half-life decay studies. In auto position, the scaler will count to any of the 10 predetermined counts to which the predetermined count control has been set, then automatically reset the entire instrument to zero, and resume counting and recycling. Each time the predetermined count is reached, a positive pulse is provided for operating external accessories.

TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA

PECIFICATIONS

MODELS DS-5B and DS-5BA

INPUT SENSITIVITY: Chassis mounted pulse height control variable from 0.1 to 6 volts. (Factory adjusted to 0.25V.)

INPUT SIGNAL POLARITY: Negative pulses.

- RESOLVING TIME: Model DS-5B: 5 microseconds. Model DS-5BA: 1 microsecond.
- RANGE: Count capability of 99,999. This may be extended to 1 x 10⁸ counts by using a Model CR-2 External Count Register.
- ELAPSED TIME INDICATOR: 5 digit odometer type, indicating in hundredths of a minute with a maximum capability of 999.99 minutes.

HIGH VOLTAGE POWER SUPPLY:

 ${\bf Range:+300}$ volts to +3000 volts D.C. continuously variable by adjustment of "course" and "fine" panel mounted potentiometer controls.

Polarity: Positive output (negative ground).

Stability: 0.01% drift per day maximum.

Regulation: Less than 0.002% change in high voltage for a 1% line voltage change between 95 to 130 volts A.C.

Ripple: Less than 0.01% of voltage

Current: Maximum output 1 ma on 1500 volt range; 0.3 ma on 3000 volt range.

Automatic Time Delay: Prevents high voltage from being turned on until regulator tubes are in operation, thereby eliminating excessive voltage output.

High Voltage Meter: Large 4" meter face calibrated to within $\pm 1.5\%$ full scale. Scales are 300V to 1500V, and 600V to 3000V full scale.

LOW VOLTAGE SUPPLY: Used for Scintillation Detector Preamplifier. Provides 250V D.C. plate voltage and 6.3V A.C. filament voltage.

PANEL MOUNTED CONTROLS:

Combined Power and H. V. Switch: Turns unit on and selects the high voltage power supply range.

Coarse and Fine Controls: Provide continuous variation of high voltage.

Auto-Normal Control: Normal position for conventional scaler use; auto position for use with external printing devices.

Stop-Count Control: Starts and stops counting operations. Reset Button: Electrically resets all 5 decades, elapsed time indicator, and external accessories (if in use).

Predetermined Count: Stops count and elapsed time indicator after 4-10-40-100-400-1K-10K-40K-100K counts. Has position for manual control by operator or predetermined timer unit. 60-Cycle Test Switch: Permits checking the instrument's operation with the 3600 cpm line frequency.

ON REAR PANEL:

Latch Lock Connectors: Connect T/A Model CR-2 Count Register, T/A Models PT-1, PT-2, or PT-3 Predetermined Timers, and T/A Models SD-1, SD-2, or SD-2W Scintillation Detectors to the DS-5B or DS-5BA.

Decade Scaler shown with

Model DS-5B

Model SD-1 Scintillation Detector and Model MS-1 Stand for thyroid uptake work



CONNECTORS:

High Voltage and Signal Input	UG-931/U
Separate High Voltage Output	Amphenol 83-798
Detector Preamplifier Power Supply	Cannon XL-3-13
External Count Register	Cannon XL-4-13
External Predetermined Timer	Cannon P5-13
Power Input.	Bryant 7486
Scope or Ratemeter Jack	Switchcraft No. 11

POWER REQUIREMENTS: 95 to 130 volts, 60 cycle, 200 watts. Available for 50 cycle and/or 220 volt operation on special order.

DIMENSIONS: Cabinet 20" wide x 11" high x 13" deep. Panel 19" wide x 8³/4" high x 13" deep.

WEIGHT: 40 lbs. net. Shipping weight: 46 lbs.

CABINET: Smooth grey hammertone finish.



DS-5B or DS-5BA available in Single Cabinet with Timer and/or External Register

T/A Model PT-2 or Model PT-3 Predetermined Timer, as well as Model CR-2 External Register, can be supplied on a 19" panel for mounting with Models DS-5B or DS-5BA Decade Scaler in a single cabinet, at customer's option.

TECHNICAL ASSOCIATES 140 West Providencia Avenue Burbank, California



IONIZATION CHAMBERS (MODELS IC-1 and IC-2)



APPLICATION: These ionization chambers are used for detection and measurement of Beta and Gamma radiation or Gamma radiation only. They can be installed in "hot" locations in a laboratory or reactor installation where it is desirable to monitor for control, health, or safety purposes. Model IC-1 has approximately one-half of the wall area cut away in the form of four windows which are fitted with .005" thick cellulose acetate. It is sensitive to soft radiation and is used for the detection and measurement of Beta and Gamma radiation. Model IC-2 is of solid wall construction and is used for the detection and measurement of Gamma radiation.

CONSTRUCTION: These chambers consist of a bakelite tube 5" in diameter x $18\frac{3}{4}$ " long, with bakelite end plates. One end plate holds the high voltage input and signal output connectors. At the other end of the chamber, a space is provided to hold a dessicant to assure reliable operation under conditions of high humidity. Silica jell is used as the dessicant and is supplied in a cloth bag suitable for insertion in the cavity provided.

The collector is an aluminum tube 1/2" in diameter x 17 %" long, and is supported by Teflon insulation at the back end, and by polystyrene insulation at the connector end. Guard rings are provided at either end to further reduce leakage. Leakage resistance is 10⁹ ohms or higher to either ground or high voltage.

The inside of the 5" cylinder and the inside of the end plates are coated with aquadag to form a conducting surface. The aquadag coating is connected to the high voltage input connector. The electrical continuity resistance of the aquadag coating is held to less than 5000 ohms.

Careful craftsmanship throughout assembly and testing insures that rigid insulation specifications are maintained.

SPECIFICATIONS:

- OPERATING VOLTAGE RANGE: 275 to 325 volts D.C.
- · SENSITIVITY: Less than 1 mr/hr with micro-micro ammeter on 0.1 volt 10-11 ampere range.
- · ENERGY DEPENDENCE: Small corrective factor below 0.3 mev; flat above 0.3 mev.
- DIMENSIONS: 5" dia. x 20" long. WEIGHT: Net 3 lbs. Shipping 6 lbs.

INQUIRIES INVITED

In addition to the types described, Technical Associates is prepared to supply ionization chambers for special purposes to customer specifications.



140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA

4M-10/58

ALOGEN QUENCHED permanently sensitive **Radiation Counter Tubes used in T/A Equipment**

Beta-Gamma sensitive stainless steel wall type tubes

... for greater sensitivity and pulse height

Used in hand probe and foot deck of monitoring instruments. (T/A Model ALM-2X Hand & Foot Monitor hand probe only, and T/A Model HSC-1 Hand & Shoe Counter.) Cathode size: length 3", dia. %", wall thickness 30-40 mg/cm². Tube base standard small 4 pin.

T1120

T1100

Used in foot deck of monitoring instruments. (T/A Model ALM-2X Hand & Foot Monitor.) Cathode size: length 7", dia. %", wall thickness 30-40 mg/cm². Tube base standard small 4 pin.



T1180



Alpha-beta sensitive. Used in general radio-assay work and radiochemical analysis. (T/A Model LS-6 Lead Shield, T/A Model

TM-8 Tube Mount in conjunction with LS-7 Lead Shield, and T/A Model P-8 Probe.)

End window diameter 13/32"; thickness 1.4-2 mg/cm².

Tube base standard medium 4 pin.

Alpha-beta-gamma sensitive. Used in general laboratory survey instruments where areater alpha and beta sensitivity are required. (T/A Model P-8X Probe in conjunction with a count rate meter or scaler.) End window diameter ²⁵/₃₂"; thickness 1.4-2 mg/cm². Tube base connector standard single

pin.

PECIFICATIONS

	T1100	T1120	T1140	T1160 *	T1180
OPERATING VOLTAGE:	900V DC	900V DC	700V DC	1300V DC	700V DC
PLATEAU LENGTH IN EXCESS OF:	200 Volts	200 Volts	200 Volts	250 Volts	180 Volts
SLOPE OF PLATEAU:	10% per 100V	10% per 100V	5% to 10% per 100V	1.5% per 100V	10% per 100V
STARTING VOLTAGE (0.3V Pulses):	825V max.	825V max.	625V max.	11 80V m ax.	620V max.
MAXIMUM COUNTING RATE:	1700 cps	1700 сря	830 cps	1100 cps	1100 cps
BACKGROUND (Shielded with 2"	•				
lead and 1/8" aluminum):	50 cpm max.	75 cpm max.	50 cpm max.	50 cpm max.	75 cpm max.
DEAD TIME (Approximate):	100 µ sec.	100 µ sec.	200 д. зес.	150 да зес.	150 д. sec.
OPERATING TEMPERATURE:	—55° to +75°C	-55° to $+75^\circ$ C	—55° to + 75°C	\pm 15° to \pm 50°C	-55° to $+$ 75°C
OVERALL DIMENSIONS:	5‰″ long, 1⅔2″ dia.	11 ² 5 ₃₂ " long, 15 ₃₂ " dia.	4 ¹ / ₃₂ " long 1%" dia.	4 ¹¹ ⁄ ₃₂ " long, 1 %" dia.	6″ long, 1″ dia.
LIFE EXPECTANCY IN COUNTS:	Unlimited	Unlimited	Unlimited	Approx. 1.5 x 10 ⁸	Unlimited

* T1160 is identical to T1140, except that it contains an organic quenching agent.

TECHNICAL ASSOCIATES

PRICE LIST

•

TECHNICAL ASSOCIATES

140 W. Providencia Ave. • Burbank, Calif.

Revised February 1, 1961

Subject to Change Without Notice

Bulletin No.		Model No.	Description of Product F.O.B.	Price . Burbank
			AUTOMATED SAMPLE CHANCING SYSTEMS	
174		ASC-1	Gas Flow Detection of Alpha/Beta activity (proportional or Geiger, window or win- dowless) complete with 200 planchet holders and planchets (less cylinder	\$5375.00
174		ASC-2	Gamma Scintillation Detection complete with 200 planchet holders and planchets	5575.00
174		ASC-3	Both Gas Flow (Alpha/Beta) and Scintillation (Gamma) Detection, simultaneously or separately, complete with 200 planchet holders and planchets (less cylinder	7005.00
173 173 173 173 173 177 166 166 166 166 166 166 166 164 164 164	* * * * * * * * * * * *	LLB-40A LLB-40B PT-1R ACS-42 GTS-6 DS-5B DS-5BA DU-1 DU-5 CR-2 PT-1 PT-2 PT-3 PT-4 GS-7S GS-7D SM-10 SA-20 LA-6 SS-30 PA-6	of gas) ANALYTICAL AND COUNTING INSTRUMENTS Low Level Beta System with Amperex Type 18516/18518 Tubes complete with stand Low Level Beta System with Amperex Type 18515/18517 Tubes complete with stand Predetermined Timer for LLB-40 Anti-Coincidence Scaler for use with Amperex Beta Counters, mounted in T/A Model C-1225 cabinet High Speed Decade Scaler (1 microsecond) with built-in preset time and preset count (less H.V. Supply) Decade Scaler (1 microsecond) One microsecond plug-in decade unit Five microsecond plug-in decade unit Count Register used with DS-5B or DS-5BA Pre-determined Timer (Long Counting Intervals — 999 minutes) Pre-determined Timer, Liebel-Flarsheim (1 second/60 minutes) Pre-determined Timer, Liebel-Flarsheim (1 second/60 minutes) Binary Scaler (Elapsed time in seconds) Binary Scaler (Elapsed time in seconds) Spectrometer (Single Channel Pulse Height Analyzer Linear Amplifier (Non-overloading type) Spectrum Scanner Pre-Amplifier for LA-6	7985.00 3750.00 3690.00 110.00 1030.50 750.00 825.00 915.00 147.50 57.50 65.00 95.00 100.00 100.00 725.00 725.00 765.00 685.00 595.00 475.00 115.00
163 178 178	* * *	RHV-1B RHV-2 RHV-3	 HIGH VOLTAGE SUPPLIES Regulated High Voltage Supply (1500-3000 volts) Two Ranges, Continuously Variable Reversible Polarity High Voltage Supply (500V to 2500V) Continuously Variable Regulated High Voltage Supply (500V to 5000V) Two Ranges for Geiger and Proportional * Prices shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Eary to Sonice" Correle Cabinet add \$20.00 to prices shown See runnes See r	375.00 450.00 425.00
171 171 171 171		Polaris Lyra (Decade) Lyra (Binary) Arcturus	 of Bulletin No. 168 for description of "Easy-to-Service" cabinet. ** These accessories will be supplied in a cabinet or panel mounted at customer's option. STAR SERIES SPECTROMETRY SYSTEMS Includes SM-10 and RHV-1 mounted in C-1400 cabinet, complete with interconnecting cables and system test Includes SM-10 and DS-5B mounted in C-1750 cabinet, complete with interconnecting cables and system test Includes SM-10 and GS-7 mounted in C-1750 cabinet, complete with interconnecting cables and system test Includes SM-10 and GS-7 mounted in C-1750 cabinet, complete with interconnecting cables and system test Includes SM-10, RHV-1, SS-30 and Texas Instruments Rectilinear Recorder Model RR1M-A16 in P-1 Panel Mount installed in C-3150 cabinet, with interconnecting cables and system test Includes SM-10, C-3150 cabinet, with interconnecting cables and system test Includes SM-10, RHV-1, SS-30 and Texas Instruments Rectilinear Recorder Model RR1M-A16 in P-1 Panel Mount installed in C-3150 cabinet, with interconnecting cables and system test Includes IA-6 SA-20, DS-5B, PT-3 or PT-4 and CR-2 mounted in C-3150 cabinet. 	1180.00 1635.00 1535.00 2250.00
		(Decade)	with interconnecting cables and system test	2330.00

PAGE 1

LIST (Continued) PRICE

Price F.O.B. Burbank

Bulletin	Model
No.	No.
171	Capella (Binary)

Description of Product

Includes LA-6, SA-20, GS-7, CR-2 and PT-3 or PT-4 mounted in C-3150 cabinet, with interconnecting cables and system test 2230.00

NOTE: For special applications requiring other combinations of T/A instruments, please request a quotation.

SCINTILLATION DETECTORS AND ACCESSORIES

157 157 157	SD-1 SD-2 SD-2W	Scintillation Detector with $1'' \times 1''$ Nal crystal, complete with cables Scintillation Detector with $2'' \times 2''$ Nal crystal, complete with cables Well-type Scintillation Detector with $13/4'' \times 2''$ Nal crystal ($5/6'' \times 11/2''$ well)	455.00 765.00
157	50.3	Scintillation Detector with 116" Alpha Phoephor, complete with cables	495.00
157	SD-3	Scintilation Detector with 1/2 Apha Phosphol, complete with cables	545.00
137	50.11	Branium Desolution Scillation Detector (Cathode Follower Type) with 1" x 1"	343.00
179	30-11	Nal crystal, complete with cables	515.00
179	SD-12	Nal crystal, complete with cables	825.00
179	SD-12W	Premium Resolution Well-type Scintillation Detector (Cathode Follower Type) with $13/4$ x 2" Nal crystal ($5/8$ " x $11/2$ " well) complete with cables for use with LS-8 or LS-8X Shields	755.00
157	CH-1	Crystal Housing for SD-1 or SD-11 Scintillation Detector	25.00
157	CH-2	Crystal Housing for SD-2 or SD-12 Scintillation Detector	30.00
157	CH-2W	Crystal Housing for SD-2W and SD-12W Scintillation Detectors	30.00
157	NS-1	Nose Shield for Model SD-1 or SD-11 Detector	95.00
157	NS-2	Nose Shield for Model SD-2 or SD-12 Detector	115.00
157	C-1a	Collimator Type A (20° Flat Field) for SD-1 or SD-11 Detector	29.50
157	C-2a	Collimator Type A (20° Flat Field) for SD-2 or SD-12 Detector	39.50
157	C-1b	Collimator Type B (Straight Bore) for SD-1 or SD-11 Detector	35.00
157	C-2b	Collimator Type B (Straight Bore) for SD-2 or SD-12 Detector	45.00
157	C-1c	Collimator Type C (Eccusing) for SD-1 or SD-11 Detector	75.00
157	C-2c	Collimator Type C (Honeycomb) for SD-2 or SD-12 Detector	95.00
157	AS_11	Alpha Detector (1") for SD-1 complete with mount	79.50
157	AC 12	Alpha Detector (11/4') for SD-2 complete with mount	85.00
157	PC 11	Apple Detector $(1^{\prime})_{2}$, or SD-1, complete with mount	79 50
157	21 28	Beta Detector (1)/6") for SD-2 complete with mount	95.00
166	NAS 1	Mobile Stand for Sciptillation Detector counter weighted and with shelf for Scaler	20.00
		or Rate Meter	2 95 .00
176	MSC-1	Shielded Gas Flow Counter (window and windowless) complete with GD-1 Detector	750.00
159	D 7	Bota Gamma Probe Assembly with T-1090 Tube cable and connector	39 50
150	F-/	End Window Probe (lass tube) complete with cable and connector	45.00
100		End Window Probe Assembly with T/A T 1180 Tube cable and connector	87 50
100	P-0A P-10	Beta/Gamma Probe Assembly with T-1110 G.M. Tube, cable and connector	75.00
		COUNTER TUBES	
	T-1090	Metal Thin wall, 3-11/16" long, used in Model F-6 Geiger Counter, P-7 Probe or other Beta/Gamma applications	15.00
144	T-1100	Thin wall type 3" long used in Models ALM-2. ALM-2X. PPM-8 and HSC-1	27.00
144	T-1110	Thin wall 3" long used in P-10 Probe	27.00
144	T-1120	Thin wall type 7" long used in Model PPM-8	53.00
144	T-1140	Mica End Window used in TM-7, TM-8 Tube Mounts and P-8 Probe	58.50
144	T-1160	Mica End Window (Organic Quenched) used in TM-7, TM-8 Tube Mounts and P-8 Probe	67.00
144	T-1180	Mica End Window, Alpha/Beta/Gamma sensitive, used in P-8X Probe	61.00
172	A 18515	Beta Counter (Amereca) for low activity low energy	186 50
173	A-18514	Reta Counter (Amperex) for low activity higher energy	206.25
173	A 19517	Guard Tube (Amperex) for use with A-18515	373.00
173	A 19519	Guard Tube (Amperex) for use with A-18516	412 50
175	N-10010	Guard Fube (Amperex) for use with A-10310	412.30
144		IONIZATION CHAMBERS	115.00
144 111	10-1	Ionization Chamber with windows (beta/Gamma)	105.00
177	10-2		, 05.00

PERSONNEL MONITORS Bata/Gamma Hand and Shoe Monito

167	HSM-10	Beta/Gamma Hand and Shoe Monitor (Dresden Automatic Design) with external
		Beta/Gamma Probe
161	ALM-2X	Beta/Gamma Hand and Foot Monitor (Argonne Count Rate Meter Design) with
		external Beta/Gamma Probes

TST (Continued) DRICF

Bulletin		Model	F	Price
No.	· · · ·	No.	Description of Product F.O.B.	Burbank
175		AHM-11	Alpha Hand Monitor (Scintillation Type - Dresden Design)	2975.00
1/6		HSC.1	Beta/Gamma Hand and Shoe Counter (Hanford Binary Scaler Type)	2750.00
140			Beta/Gamma Dartal Type Monitor	2150.00
100		FF/V1-0	Beta/Gamma Fortal Type Monitor	1150.00
108		L1/VI-10	Deta/Gamma Laundry Monitor	1130.00
			SURVEY METERS	
165		CP-3	New Cutie Pie Portable Survey Meter	265.00
165		CP-3A	New Cutie Pie Portable Survey Meter with special ranges (0-25, 0-250, 0-2500 mr/hr)	295.00
180	н. н. н. Н	SRJ-7	New Juno Survey Meter with plug-in battery pack and replaceable Alpha window (Standard Range)	325.00
180		HRJ-7	New Juno Survey Meter with plug-in battery pack and replaceable Alpha window (High Range)	325.00
136		F-6	Geiger Counter for Uranium Prospecting	159.50
137	14	FS_11	Scintillation Counter for Uranium Prospecting	449.50
137	÷	13-11		
			LEAD SHIELDS	
152		LS-1	Complete with Type ST-1 Sample Tray, Type SK-1 Socket and Type TM-1 Tube	105.00
			Mount	185.00
152	г. Тъ	LS-2A	Standard with Type TM-6 Tube Mount and Sample Tray Holder, Type ST-27	200.00
		LC OR	Alignments with 10 Time PI 1 Displate	285.00
152			For Unrite with Top TA 1 Tube Mount Type SK-1 Socket and one Type	200.00
153		L3-4A	ST-5 Marinelli Beaker	200.00
153	,	LS-4B	Modified to accommodate special tubes with aluminum inner liner	235.00
153	24. 24.	LS-4C	For solids, complete with Type TM-1 Tube Mount, Type SK-1 Socket, one Type ST-6 Ore Container and 100 Type ST-6A Paper Sleeves	225.00
151	e de l'Ara L	LS-5	Complete with Type TM-3 Tube Mount, Type SK-2 Socket, Type STH-1 Sample	260.00
150		15.4	Complete with Type TH-7 type Mount one Type ST-38 Sample Tray and 10	
152	,	63-0	Type PL-3 Aluminum Planchets	295.00
	2	1.5.6X	Fight-inch beight extension for LS-6 Shield	98.50
		LS-6L	Lid with 3" opening and retaining collar for Models SD-1, SD-2, SD-12 Scintilla-	
			tion Detectors	72.50
		LS-6P	Plug for 3" opening in LS-6L	12.50
160		LS-7A	Multi-Purpose Type complete with SD-1 Scintillation Detector, three ST-7 Sample Trave and 10 PL-3 Planchets	780.00
160		LS-7B	Multi-Purpose Type complete with SD-2 Scintillation Detector, three ST-7 Sample	
			Trays and 10 PL-3 Planchets	1090.00
160		LS-7C	Multi-Purpose Type complete with TM-8 Tube Mount, three S1-7 Sample Trays and 10 PL-3 Planchets	385.00
160		LS-7M	Multi-Purpose Type — Shield only, with adapter ring, three ST-7 Sample Trays and 10 PL-3 Planchets	325.00
148		LS-8	For well counting (shield only, less detector)	265.00
148		LS-8W	For well counting, complete with SD-2W Scintillation Detector and 13/4" x 2"	
		10.1014/	well crystal	960.00
	1 - 1 - 1	L3-10W	well crystal	1020.00
148		LS-8X	For well counting with 2" shielding in all directions (Shield only less Detector)	350.00
148		LS-8WX	For well counting with 2" shielding in all directions complete with SD-2W Scintilla-	1035.00
		10 104/94	For well counting with 2" chielding in all directions complete with SD 12W	1000.00
	. :	L2-18WX	Scintillation Detector and 134" x 2" well crystal	1095.00
173		LS-41	Lead, Steel and Mercury Shield for Amperex Beta Counters (less tubes), but com-	2100.00

A REAL PROPERTY OF A REAL PROPER

LABORATORY ACCESSORIES

162	AB-2	Absorber Set	85.00
153	ST-1	Sample Tray	6.50
154	ST-2, ST-3,	ST-4 Sample Trays	3.50
153	ST-5	Marinelli Beaker	4.25
153	ST-6	Ore Container	12.50
153	ST-6A	Paper Sleeves for Ore Containerper C	2.00
160	ST-7	Sample Tray for LS-7	2.50
154	STH-1	Sample Tray Holder with tray and ten planchets	13.50
154	 SK-1, SK-2,	SK-3 Sockets	5.00
154	TM-1, TM-2	2, TM-3, TM-4 Tube Mount	15.00
154	TM-5	Tube Mount (micrometric for LS-2 Lead Shield)	100.00
154	TM-6	Tube Mount and Sample Tray Holder	22.50
154	TM-7	Tube Mount and Sample Tray Holder	27.50
160	TM-8	Tube Mount and Adapting Ring for LS-7 Multi-Purpose Shield	60.00

plete with stand _____ 2100.00

PAGE 3

PRICE LIST (Continued)

Bulletin No.	Model No.	Description of Product		F.O.	Price B. Burbank
162	LC-1	Lead Source Container for four radium needles			45.00
162	LC-3	Lead Source Container 1" lead walls, 11/2" x 3" inside dimensions	- • • - • • • • • •		40.00
162	LC-3A	Lead Source Container 13% " lead walls $34'' \times 214''$ inside dimensions			47.50
162	LC-4	Lead Source Container $1\frac{1}{2}$ " lead walls, $1\frac{1}{2}$ " x 3" inside dimensions			50.00
162	LC-4A	Lead Source Container $1\frac{3}{4}$ " lead walls $\frac{3}{4}$ " x $2\frac{1}{2}$ " inside dimension	s		57.50
162	LC-9	Mobile Source Unit — compartment 3" x 6" in shield of 2" lead. heavy duty hand cart	Moun	ted oi	n 250.00
130		Lead Brick, machined all sides 2" x 3" x 6"			8.50
130		Lead Brick, machined all sides 2" x 4" x 8"			12.50
130		Lead Brick, extruded shape	Prices	upon	application
153		Planchets, Type PL-1, Sheet steel, tinned	. per	Ċ.	3.25
			per	м	30.00
			per	5M	27.00/M
153		Planchets, Type PL-2 & PL-3, aluminum	. per	С	3.25
			per	Μ	30.00
			per	5M	27.00/M
153		Planchets, Type PL-2C & PL-3C, copper	per	С	4.50
			per	м	40.00
			per	5M	36.00/M
153		Planchets, PL-2S, stainless steel	. per	С	9.00
		T = 0 , $T = 0$, $T =$	per	M	65.00
			per	5M	58.50/M
153		Planchets, Type PL-3S, stainless steel	. per	С	12.00
		$(x,y) \in \{x,y\}$, $(x,y) \in \{x$	per	м	85.00
			per	5M	76.50/M
153		Planchets, Type PL-4, aluminum	. per	С	5.00
			per	M	45.00
			per	5M	40.50/M
153		Planchets, special materials and sizes Price	s quoi	ted up	oon request
174		Planchet Holder for ASC-1, ASC-2, ASC-3	. per	C	6.00
			per	M	55.00
174		Sample Tray for ASC-1, ASC-2, ASC-3	. per	C	8.00
			per	M	75.00

CABINETS

C-875	8 ³ /4" height "Easy-to-Service" console cabinet for standard 19" width panel	30.00
C-1050	10/2 neight console cabinet for standard 19 width panels	37.50
C-1400	14" height console cabinet for standard 19" width panels	40.00
C-1750	$17\frac{1}{2}$ " height console cabinet for standard 19" width panels	45.00
C-2100	21" height console cabinet for standard 19" width panels	50.00
C-3150	$31\frac{1}{2}$ " height console cabinet for standard 19" width panels	60.00
C-4375	433/4" height console cabinet for standard 19" width panels	70.00

IN GENERAL

Unit prices are given in this price list unless otherwise noted.

All prices quoted herein are f.o.b. Burbank, California. Payment terms, net 30 days upon approved credit.

No extra charge is made for packing for domestic shipments. Export charges will be quoted upon application. We reserve the right to alter specifications at any time without incurring the obligation of incorporating new features in previously manufactured equipment. Prices are subject to change without notice. Quotations remain firm for 30 days.

Shipping instruments: Please include shipping instructions when ordering; in the event that shipping instructions are not given, we will use our best judgment in the matter.

SERVICE CHARGES

Services charges for repairs are usually billed on an hourly basis, however, if desired, an estimate can be given before the work is undertaken. The unit to be repaired should be shipped prepaid to us.

WARRANTY

Technical Associates warrants instruments and equipment (except tubes, fuses, batteries and crystals), manufactured by them to be free from defects in workmanship or materials under normal use for a period of one year from the date of shipment from the factory to the buyer. Tubes, fuses, batteries and crystals are subject to the guarantee established by the manufacturer of them, however, Technical Associates will assist the customer to obtain full benefits of these guarantees.

If, within the one year warranty period, any Technical Associates instrumentation or equipment requires service as a result of a defect, the buyer may return it to the factory of Technical Associates at Burbank, California or to a service station designated by Technical Associates, transportation charges prepaid, for service at no charge under the warranty. The buyer is urged to communicate with Technical Associates when warranty service is required, stating the nature of the difficulty and giving model and serial number of instrument. It may be possible to diagnose the trouble and send a replacement part or assembly, thereby avoiding the expense of shipment.

Technical Associates will return the instrument to the buyer, transportation charges prepaid, after repairs or replacement under warranty are completed. The liability to Technical Associates under this warranty is limited to the cost of replacement of defective parts upon prompt notification of such defect.



TECHNICAL ASSOCIATES INSTRUMENTATION FOR NUCLEAR RESEARCH 140 WEST PROVIDENCIA AVENUE • BURBANK • CALIFORNIA

TELEPHONES: Victoria 9-1994 • Thornwall 8-8133

January 23, 1959

General Electric Company 4855 Electric Avenue Milwaukee 1, Wisconsin

Attention: Mr. J. J. Jech, Coolidge Lab. X-Ray Department

Gentlemen:

Your letter of January 19 was received and read with considerable interest. In response we are enclosing our descriptive Bulletin No. 159 which explains the Juno Survey Meter in complete detail. You will also find enclosed our quotation and a X-Ray Calibration Chart which was made by the National Eureau of Standards on our Juno.

The calibration chart provides the user of the Juno with a correction factor so that he may use the instrument to determine X-Ray dosage rates.

To use the calibration chart it is necessary to know the KEV level at which the x-ray unit is operating. We have red lined an example on the chart. In the example, the Juno is being used with the alpha rejection absorber in place and is therefore positioned for detection of beta. Further, the KEV level of the x-ray unit being monitored is 30 as marked on the base line. In this case the correction factor would be .8, therefore the reading of the Juno should be multiplied by .8 to obtain the correct x-ray dosage rate.

Should you find the Juno acceptable for your particular application, it is advisable to check the calibration with some known dosage rates.

We trust you will feel free to contact us should you have any further questions.

Very truly yours,

TECHNICAL ASSOCIATES

Baron Mimar

BG:rt Encl.

Baron Getman Assistant Sales Manager

P.S. If you would like to have a copy of our complete catalog, please fill in and return the enclosed postcard.

TECHNICAL ASSOCIATES 140 West Providencia Avenue

. . .

Burbank, California REQUEST FOR CATALOG

Name	
Organization	
Department	
Street Address	
City	Zone

TECHNICAL ASSOCIATES

.

-ign.

140 West Providencia Avenue Burbank, California



JUNO RADIATION SURVEY METER MODELS SRJ-6 and HRJ-6



APPLICATION. The Juno Survey Meter is a portable instrument for measuring the intensity of, and discriminating between alpha, beta, and gamma radiation. It is used to protect personnel from the danger of over-exposure to radiation from radioactive materials or X-rays. While primarily intended for inspection of flat surfaces, the instrument is suitable for most uses where a high degree of accuracy is desired.

The T/A Juno is an improved version of the original Hanford instrument and is available in two models: SRJ-6 (standard range) for all normal applications; and HRJ-6 (high range) for use where exceptionally high intensity radiation is likely to be encountered. Both models meet A.E.C. specifications.

DESCRIPTION: The instrument comprises an ionization chamber, an electrometer circuit, absorption filters for the rejection of either alpha or beta particles, suitable batteries mounted in a removable power pack, and an indicating meter. The unit is battery operated and is self-contained. To easily identify the High Range Juno, its knobs and meter dials are finished in brilliant red.

The ionization chamber has a volume of approximately 27 cubic inches. All surfaces within the chamber are coated with aquadag. The chamber is covered by a screen of .0003'' (approximately 0.45 mg/cm^2) rubber hydrochloride film. The alpha screen is within 7/16'' of any flat surface on which the instrument may be placed, and is easily replaced by simply removing the bottom plate and two retaining screws.

Two absorbers are provided to reject either alpha or beta radiation. These are readily moved in and out of position by means of sliding tabs fitted in rails which form part of the handle. The tab marked "G," with a square end, operates the absorber which rejects alpha and beta, thus permitting a reading of gamma only. The tab marked "B," with a rounded end, operates the absorber which rejects alpha, permitting a reading of beta and gamma. The total of all three types of radiation is read, when both tabs are in "open" position.

The high quality microammeter, which is calibrated in milliroentgens per hour for gamma radiation, has a large easy-to-read face and is mounted in position to permit excellent visibility. Battery life is approximately 300 hours in normal intermittent use. An easily removable battery pack, with simple positive contacts, assures trouble-free operation over long periods of service.



TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA

Models SRJ-6 and HRJ-6 have similar circuits, the principal difference being the value of grid resistances used.

Radiation incident upon the ionization chamber produces a minute current which flows through a very high resistance in the grid circuit of the electrometer tube. The voltage thus produced at the grid causes a corresponding change in plate current which is indicated by the panel meter. A bucking current is provided through the meter in order that the no-signal plate current of the electrometer may be balanced out and readings of radiation intensity may start from the meter zero reading. Sensitivity is varied by switching appropriate values of grid resistance in the electrometer circuit.

The instrument is calibrated by adjusting a resistance in series with the meter. An individual adjustment is provided for each range. Zero setting is accomplished by means of a rheostat in the filament circuit of the electrometer tube. By varying the filament voltage, the plate current may be varied and thus adjusted to a value equal to the bucking current flowing through the meter. All high resistance points in the circuit are insulated with Teflon to insure minimum leakage.

SENSITIVITY. Both models are calibrated in three separate full-scale ranges in easily read increments of the meter scale, covering the total range of which the instrument is capable.

Ranges are based on radium gamma radiation intensity. Accuracy of calibration is such that indications on the meter will not be more than 5% lower nor more than 10% higher than the radiation intensity to which the chamber is exposed. Sensitivity dependency upon battery aging is limited to a 10% variation while the unit can be zeroed by means of the panel zero control.

Illustration below shows the T/A Juno in use at

SPECIFICATIONS:

IONIZATION CHAMBER:

Volume: 27 cubic inches. Window Opening: 3" x 4%". Alpha Screen: 0.3 mil (0.45 mg/cm²) rubber

hydrochloride. Alpha Rejection Absorber: 0.01" cellulose acetate

sheet.

Beta Rejection Absorber: 0.102" aluminum.

BATTERIES:

4 Eveready No. 412 22¹/₂ Volt "B" Batteries. 2 Eveready No. E12 1.35 Volt "A" Batteries.

2 Mallory No. TR-115 6.5 Volt "B" Batteries.

TUBE:

1 Sub-Miniature Electrometer Type CK 5886.

RANGES:

Model SRJ-6 Juno (standard range): 50, 500, 5000 MR/HR full-scale. (Improved A.E.C. Model SIC-17B). Model HRJ-6 Juno (high range): 250, 2500, 25,000 MR/HR full-scale. (Improved A.E.C. Model SIC-17D).

CASE:

Hard-chrome plated aluminum, with engraved markings on top of case; gasketed, dust and moistureresistant. Dimensions: 91/2" x 53/4" x 4".

HANDLE:

Aluminum, especially cast low-porosity, smoothly polished.

WEIGHT:

Net 6 lbs. 11 oz. Shipping: 10 lbs.





The Juno can be used for X-ray detection and measurement by reference to the above curves.

140 W. PROVIDENCIA AVE. . BURBANK, CALIF.

TECHNICAL ASSOCIATES



JUNO RADIATION SURVEY METER



MODELS SRJ-6 and HRJ-6 Cost 295:00

APPLICATION. The Juno Survey Meter is a portable instrument for measuring the intensity of, and discriminating between alpha, beta, and gamma radiation. It is used to protect personnel from the danger of over-exposure to radiation from radioactive materials or X-rays. While primarily intended for inspection of flat surfaces, the instrument is suitable for most uses where a high degree of accuracy is desired.

The T/A Juno is an improved version of the original Hanford instrument and is available in two models: SRJ-6 (standard range) for all normal applications; and HRJ-6 (high range) for use where exceptionally high intensity radiation is likely to be encountered. Both models meet A.E.C. specifications.

DESCRIPTION: The instrument comprises an ionization chamber, an electrometer circuit, absorption filters for the rejection of either alpha or beta particles, suitable batteries mounted in a removable power pack, and an indicating meter. The unit is battery operated and is self-contained. To easily identify the High Range Juno, its knobs and meter dials are finished in brilliant red.

The ionization chamber has a volume of approximately 27 cubic inches. All surfaces within the chamber are coated with aquadag. The chamber is covered by a screen of .0003'' (approximately 0.45 mg/cm^2) rubber hydrochloride film. The alpha screen is within 7/16'' of any flat surface on which the instrument may be placed, and is easily replaced by simply removing the bottom plate and two retaining screws.

Two absorbers are provided to reject either alpha or beta radiation. These are readily moved in and out of position by means

of sliding tabs fitted in rails which form part of the handle. The tab marked "G," with a square end, operates the absorber which rejects alpha and beta, thus permitting a reading of gamma only. The tab marked "B," with a rounded end, operates the absorber which rejects alpha, permitting a reading of beta and gamma. The total of all three types of radiation is read, when both tabs are in "open" position.

The high quality microammeter, which is calibrated in milliroentgens per hour for gamma radiation, has a large easy-to-read face and is mounted in position to permit excellent visibility. Battery life is approximately 300 hours in normal intermittent use. An easily removable battery pack, with simple positive contacts, assures trouble-free operation over long periods of service.



TECHNICAL ASSOCIATES



ALPHA SCREEN

ALPHA REJECTION

140 WEST PROVIDENCIA AVENUE • BURBANK, CALIFORNIA

Models SRJ-6 and HRJ-6

have similar circuits, the principal difference being the value of grid resistances used.

Radiation incident upon the ionization chamber produces a minute current which flows through a very high resistance in the grid circuit of the electrometer tube. The voltage thus produced at the grid causes a corresponding change in plate current which is indicated by the panel meter. A bucking current is provided through the meter in order that the no-signal plate current of the electrometer may be balanced out and readings of radiation intensity may start from the meter zero reading. Sensitivity is varied by switching appropriate values of grid resistance in the electrometer circuit.

The instrument is calibrated by adjusting a resistance in series with the meter. An individual adjustment is provided for each range. Zero setting is accomplished by means of a rheostat in the filament circuit of the electrometer tube. By varying the filament voltage, the plate current may be varied and thus adjusted to a value equal to the bucking current flowing through the meter. All high resistance points in the circuit are insulated with Teflon to insure minimum leakage.

SENSITIVITY. Both models are calibrated in three separate full-scale ranges in easily read increments of the meter scale, covering the total range of which the instrument is capable.

Ranges are based on radium gamma radiation intensity. Accuracy of calibration is such that indications on the meter will not be more than 5% lower nor more than 10% higher than the radiation intensity to which the chamber is exposed. Sensitivity dependency upon battery aging is limited to a 10% variation while the unit can be zeroed by means of the panel zero control.

Illustration below shows the T/A Juno in use at the Nuclear Radiation Laboratory of Admiral Corporation, monitoring for radio-activity as sample is transferred from "coffin" to shipping container.

SPECIFICATIONS:

IONIZATION CHAMBER:

Volume: 27 cubic inches.

Window Opening: 3" x 4%". Alpha Screen: 0.3 mil (0.45 mg/cm²) rubber

hydrochloride.

Alpha Rejection Absorber: 0.01" cellulose acetate sheet.

Beta Rejection Absorber: 0.102" aluminum.

BATTERIES:

4 Eveready No. 412 22¹/₂ Volt "B" Batteries. 2 Eveready No. E12 1.35 Volt "A" Batteries. 2 Mallory No. TR-115 6.5 Volt "B" Batteries.

TUBE:

1 Sub-Miniature Electrometer Type CK 5886.

RANGES:

Model SRJ-6 Juno (standard range): 50, 500, 5000 MR/HR full-scale. (Improved A.E.C. Model SIC-17B). Model HRJ-6 Juno (high range): 250, 2500, 25,000 MR/HR full-scale. (Improved A.E.C. Model SIC-17D).

CASE:

Hard-chrome plated aluminum, with engraved markings on top of case; gasketed, dust and moistureresistant. Dimensions: $9\frac{1}{2}$ " x $5\frac{3}{4}$ " x 4".

HANDLE:

Aluminum, especially cast low-porosity, smoothly polished.

WEIGHT:

Net 6 lbs. 11 oz. Shipping: 10 lbs.



The Juno can be used for X-ray detection and measurement by reference to the above curves.





140 W. PROVIDENCIA AVE. . BURBANK, CALIF.

TECHNICAL ASSOCIATES



HIGH RANGE SURVEY METER MODELS CP-TP-1A and CP-TP-1B

- Full-scale ranges to 5000 R/HR. CP-TP-1A: 0 to .5, 5, or 50 R/HR. CP-TP-1B: 0 to 50, 500, or 5000 R/HR.
- Interchangeable ionization chambers.
- 3-position rotating chamber.

APPLICATION. The T/A Model CP-TP is the only portable, battery operated survey meter especially designed to enable the health physicist to measure high intensity beta and gamma radiation fields. The instrument features a 40" extension rod with swivel bracket for rotating the chamber, thus allowing the user to stand behind a shield or out of the direct line of a beam.

DESCRIPTION. Technical Associates' CP-TP High Range Survey Meter is available in two models, each having three linear ranges: CP-TP-1A provides 0 to .5, 5, and 50 R/HR (beta and gamma); CP-TP-1B provides 0 to 50, 500, and 5000 R/HR (gamma only).

The basic instrument is the same for both models, the only difference being the type of ionization chamber selected. The lower and higher ranges are determined by the size and configuration of the ionization chambers and their electrometer circuits. Each chamber contains its own electrometer circuit. The chambers are interchangeable on the instrument. The CP-TP-1A chamber has a removable beta absorber to allow measurement of beta radiation through an acetate end-window.

The instrument case is chrome-plated aluminum for easy decontamination. It contains the range selector switch, zeroing and calibration controls, batteries, and a high quality wide-face meter. The unit has a pistol grip handle for easy manipulation. Each model comes complete with a 40" extension rod, which is detachable from the instrument case. Either ionization chamber can be connected directly to the case, if desired. The swivel bracket at the end of the extension rod has two connectors for use in attaching the chamber for the desired ranges. One connector is at the forward end and the other is on the side. The swivel bracket rotates, thus moving the chamber connected to the side of bracket into 3 separate positions: right, left, and downward.

SPECIFICATIONS:

RANGES: Model CP-TP-1A: 0.5, 5, 50 R/HR. Model CP-TP-1B: 50, 500, 5000 R/HR.

- CONTROLS: Combined "Off-On" and range selector switch. Individual calibration controls for each range. Zero control can be set in radiation field.
- **CIRCUIT:** Reliable single tube electrometer circuit with all high resistance points insulated with Teflon or Kel-F to insure minimum leakage.
- METER: Wide-face 3¹/₂" meter with dual scale calibrated in R/HR (black scale for 1A chamber, red scale for 1B chamber). Divided into 50 divisions. Mounted to provide excellent visibility.
- ZERO DRIFT: Negligible after 15 minute warm-up. IONIZATION CHAMBERS:
- Model CP-TP-1A: Cylindrical bakelite internally coated with aquadag. Volume approx. 42 cubic inches. Chamber has a removable beta absorber

(432 mg/cm²) to allow beta measurements through an acetate end-window (6 mg/cm²). **Model CP-TP-1B:** Cylindrical aluminum internally coated with aquadag. Volume approx. 3 cubic inches. Sensitive to gamma only.

- NOTE: Both chambers include their own electrometer circuit for specific ranges.
- **EXTENSION:** Made of polished aluminum tubing $40^{\prime\prime}$ long, with connector on one end for mounting on instrument case and swivel connector at the other end for mounting chamber.
- **INSTRUMENT CASE:** Chrome-plated aluminum with all markings engraved for easy decontamination.
- **BATTERIES:** 4 Eveready No. 412, 22.5 volts. 2 Burgess No. W5BP, 7.5 volts. 2 Eveready "D" Cells, 1.5 volts.
- WEIGHT: 6 lbs.





EXTENSION ROD enables user to stay at safe distance from beam.

afe distance from beam.

MODEL CP-TP Shown with 1A Chamber

TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE • BURBANK, CALIFORNIA



READY FOR US COUNTER IN OPERATION

CHECK INCOMPLETE - RESET AND SEREN

CHECKON

Ta

Beta-Gamma HAND and SHOE MONITOR **Model HSM-10A**

- Completely automatic monitoring.
- Error-proof, illuminated read-out panel.
- Alarm settings as low as background.
- Positive warning of incomplete check.
- External clothing probe with separate circuitry, meter, and speaker.

APPLICATION. The Model HSM-10A Hand and Shoe Monitor is a health physics instrument designed to provide

completely automatic detection and measurement of beta-gamma contamination on the palm and back of each hand, the bottom surfaces of both shoes, and also the hair and clothing. Its simple, fool-proof operation allows an individual to perform this check without direct supervision of the health physicist.

DESCRIPTION. The Model HSM-10A employs the decade scaling principle and features illuminated multi-colored panels that enable the user to tell - in one look - whether he is contaminated, and where. Separate register read-outs are provided on the monitor panel, showing the total count on each hand and the combined count for both shoes.

An outstanding feature of the instrument is the new lower warning levels, which permit alarm settings as low as the background count. Count ranges are adjustable from 10 to 190, in increments of 10, with variable time settings from 1 to 120 seconds.

A separate rate meter circuit is utilized for the external clothing probe, including a speaker for audible check and a meter for visual check. All the detecting elements and sections exposed to radiation are designed for easy decontamination. A roll of kraft paper is mounted in the base for use over the shoe deck to prevent contaminated material from dropping into the detector chambers.





OPERATION. No skill is required in the operation of the Model HSM-10A Monitor. The user simply checks the illuminated multi-colored panels on the front of the instrument.

When the monitor is turned on, the blue "Ready for Use" panel will be illuminated. The user merely steps on the shoe deck and inserts hands into the two waist-high probe openings. The counting process is started and continued by pressure of the fingertips at the rear of the probe openings,

and by the weight of the person on the shoe deck. As soon as the counting starts, the "Ready for Use" panel darkens, and the yellow "Counter in Operation" panel lights up.

NOTE: If at any time during the counting cycle, the user removes either hand from the hand probe or steps off the shoe deck, the cycle will automatically stop and the orange panel reading "Check Incomplete - Reset and

Repeat" will light up. This positive warning prevents erroneous readings.

As soon as the preset time has elapsed, the "Counter in Operation" panel darkens and either the green "Check O.K." panel or the bold red "Decontamination Required" panel lights up, depending upon whether the radioactivity present is above or below the alarm setting. If decontamination is required, bold red panels marked "Left Hand," "Right Hand," or "Shoes" will light up, indicating the location of contamination. The degree of contamination is indicated by the registers on the monitor panel.

The clothing probe, mounted on the right side of the instrument, is used to detect contamination on the hair and clothing. The output of the clothing probe is read visually on the count rate meter and aurally by the speaker, both of which are mounted on the monitor panel.

After the check is completed, the user pushes the large red "Reset" button on the front panel, and the instrument is ready for the next user.



TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE • BURBANK, CALIFORNIA



Alarm Settings: Preset Count Ranges: 10 to 190 counts, in increments of 10. Preset Time Periods: Variable from 1 to 120 seconds.

Sensitivity: Less than 0.15 μ c of 0.2 MEV or higher energy.

- Detectors: Hand Probes: 4 T/A Type T-1100 halogen quenched Geiger tubes in each hand probe.
- Shoe Deck: 2 T/A Type T-1100 halogen quenched Geiger tubes in each shoe probe.

Clothing Probe: 1 T/A Type T-1110 halogen quenched Geiger tube. Shielding: Hand and shoe detectors are shielded with $\frac{1}{2}$ " of lead.

Register Range: 0 to 9999, four digits on each of three Sodeco Registers.

Speaker for Clothing Probe: 21/2" adjustable volume.

- Meter for Clothing Probe: 5" Assembly Products meter, registering to 10,000 cpm.
- **Power Requirements:** 105 to 125 volts, 60 cycle A.C., 300 watts. **Dimensions:** Cabinet: $67\frac{3}{4}$ " high x $27\frac{1}{6}$ " wide x 24" deep.

Shoe Deck: 23" long x 18" wide x $4\frac{1}{4}$ " high.

Finish: Grey baked enamel. Shipping Weight: 650 lbs.

The complete unit is housed in a sturdy metal cabinet with grey baked enamel finish. Locked doors cover all controls, protecting the instrument from misadjustment by

unauthorized personnel.



TECHNICAL ASSOCIATES

140 West Providencia Avenue • Burbank, California



REGULATED HIGH VOLTAGE POWER SUPPLY MODEL RHV-1B





APPLICATION. The Technical Associates Model RHV-1B High Voltage Power Supply is specifically designed for nuclear applications. It supplies operating current for gas flow counters, scintillation detectors, Geiger tubes, ionization chambers, and other applications where high voltage and low current are required. Because of its outstanding stability and reliability, the instrument is ideal for use as a component in spectrometry systems.

DESCRIPTION. The Model RHV-1B is highly stable and free from jitter, bounce, and corona. The residual ripple is less than 10 millivolts under all conditions. The instrument employs precision deposited carbon resistors which are unaffected by temperature changes. The output voltage is essentially unchanged by line or load fluctuations. An automatic time delay prevents the high voltage from being turned on until regulator tubes are in operation, thereby eliminating damage to the detectors.

Front panel controls are provided for power, range selection, and variable coarse and fine adjustments of output voltage. An "off" position is provided on the range switch to allow the high voltage to be turned off independently of the main power switch and without disturbing the settings of the "Coarse" and "Fine" controls. When the range switch is returned to range position, no warm-up time is required.

Separate indicator lights are provided for power switch and range switch. The 41/2" meter has a tilted dial and was selected for easier and more accurate reading. For maximum convenience, output receptacles are provided on both front panel and rear chassis.

- 2 RANGES: +300V to +1500 V, and +600V to +3000V D.C. continuously variable by adjust-ment of "coarse" and "fine" potentiometer controls.
- POLARITY: Positive output (negative ground). STABILITY: Less than .01% drift per day after warm-up
- LINE REGULATION: Less than .01% change in high voltage for 1 volt line change between 105 to 125 volts.
- LOAD REGULATION: Output voltage decrease is less than .02% from no load to full load.
- CURRENT OUTPUT: Maximum output 1 ma on 1500 volt range: 0.3 ma on 3000 volt range. RIPPLE: Less than 10 millivolts at any voltage.
- POWER REQUIREMENTS: 105 to 125 Volts A.C., 50/ 60 cycles, 120 watts.

PANEL CONTROLS: Power Switch: "On" and "Off?" 3-Position High Voltage Switch: "Off?" "1.5KV?" and "3KV?"

High Voltage Adjustments: "Coarse" and "Fine!" HIGH VOLTAGE METER: 41/2" meter with 2 scales.

- HIGH VOLTAGE RECEPTACLES: One each Amphe-nol 83-798 and one each UG-931/U located on rear of chassis. Identical set of receptacles on
 - front panel.
- INPUT CONNECTOR: 3-prong safety plug Bryant 7486 with 6 foot cord.
- DIMENSIONS: 19" wide x 51/4" high x 13" deep. For rack mounting.

WEIGHT: 18 pounds net; 26 pounds packed for shipment.

ECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA

HIGH VOLTAGE POWER SUPPLY **MODEL RHV-2**



APPLICATION. The T/A Model RHV-2 is designed expressly for applications in analytical spectrometry and general systems where great stability and minimum ripple content are needed and a choice of positive or negative output polarity is desirable. It will supply operating current for Geiger tubes, scintillation detectors, gas flow counters, ionization chambers, and other sensing components requiring stable high voltage and low current.

DESCRIPTION. The Model RHV-2 High Voltage Power Supply is engineered for high performance and reliability. The circuit design is conservative and the components selected operate well below their rated limits. The instrument is highly stable and free from noise and corona. An automatic time delay prevents the high voltage from being turned on until regulator tubes are in operation, thus eliminating damage to detectors.

Front panel controls are provided for power, positive or negative range selection, and variable coarse and fine adjustments of output voltage. An "off" position is provided on the range switch to allow the high voltage to be turned off independently of the main power switch and without disturbing the coarse and fine adjustments. When the range switch is returned to range position, no warm-up time is required. The instrument has separate indicator lights for power switch and range selector switch. The 41/2" meter has a tilted dial and was selected for easier, more accurate reading.

The continuously variable type of voltage control featured in the Model RHV-2 permits the precise selection of the optimum voltage for a particular scintillation detector to insure maximum gain and, at the same time, minimum noise level. The high voltage transformer is hermetically sealed for long life.

이 제 이 제 그 야 과 지 않는 것을 들었다.

2 RANGES: +500 volts to +2500 volts D.C. -500 volts to --2500 volts D.C.

- POLARITY: Positive or negative, selected by front panel switch.
- LONG TERM STABILITY: Less than .02% drift per day after warm-up.
- LINE REGULATION: Less than 1% for variations from 105 volts to 125 volts.
- LOAD REGULATION: Output voltage decrease is less than .01% from no load to full load.
- RIPPLE: Less than 10 millivolts on either polarity. CURRENT OUTPUT: Maximum output 3 ma at 500 volts, 1 ma at 1500 volts, 0.6 ma at 2500 volts.
- POWER REQUIREMENTS: 105 to 125 volts A.C., 50/ 60 cycles, 130 watts,

- PANEL CONTROLS: Power Switch: "On" and "Off!" 3-Position High Voltage Switch: "Off," "-2.5KV," and "+2.5KV." High Voltage Adjustments: "Coarse" 4-posi-tion switch with overlapping ranges; "Fine" continuously variable potentiometer.
- HIGH VOLTAGE METER: 41/2" meter with 0 to 2.5KV
- scale. RECEPTACLE: High voltage output UG-931/U.
- INPUT CONNECTOR: 3-prong safety plug Bryant 7486 with 6-foot cord.
- DIMENSIONS: 19" wide x $5\frac{1}{4}$ " high x $14\frac{1}{2}$ " deep. For rack mounting.
- WEIGHT: 20 pounds net; 26 pounds packed for shipment.



140 West Providencia Avenue • Burbank, California

TECHNICAL ABSOCIATES



IONIZATION CHAMBERS (MODELS IC-1 and IC-2)



APPLICATION: These ionization chambers are used for detection and measurement of Beta and Gamma radiation or Gamma radiation only. They can be installed in "hot" locations in a laboratory or reactor installation where it is desirable to monitor for control, health, or safety purposes. Model IC-1 has approximately one-half of the wall area cut away in the form of four windows which are fitted with .005" thick cellulose acetate. It is sensitive to soft radiation and is used for the detection and measurement of Beta and Gamma radiation. Model IC-2 is of solid wall construction and is used for the detection and measurement of Gamma radiation.

CONSTRUCTION: These chambers consist of a bakelite tube 5" in diameter x 18¾" long, with bakelite end plates. One end plate holds the high voltage input and signal output connectors. At the other end of the chamber, a space is provided to hold a dessicant to assure reliable operation under conditions of high humidity. Silica jell is used as the dessicant and is supplied in a cloth bag suitable for insertion in the cavity provided.

The collector is an aluminum tube $\frac{1}{2}$ " in diameter x 17%" long, and is supported by Teflon insulation at the back end, and by polystyrene insulation at the connector end. Guard rings are provided at either end to further reduce leakage. Leakage resistance is 109 ohms or higher to either ground or high voltage.

The inside of the 5" cylinder and the inside of the end plates are coated with aquadag to form a conducting surface. The aquadag coating is connected to the high voltage input connector. The electrical continuity resistance of the aquadag coating is held to less than 5000 ohms.

Careful craftsmanship throughout assembly and testing insures that rigid insulation specifications are maintained.

SPECIFICATIONS:

- OPERATING VOLTAGE RANGE: 275 to 325 volts D.C.
- SENSITIVITY: Less than 1 mr/hr with micro-micro ammeter on 0.1 volt 10-11 ampere range.
- ENERGY DEPENDENCE: Small corrective factor below 0.3 mev; flat above 0.3 mey.
- DIMENSIONS: 5" dia. x 20" long. WEIGHT: Net 3 lbs. Shipping 6 lbs.

INQUIRIES INVITED

In addition to the types described, Technical Associates is prepared to supply ionization chambers for special purposes to customer specifications.



140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA

TECHNICAL ASSOCIATES

12.53

ALOGEN QUENCHED permanently sensitive

Radiation Counter Tubes used in T/A Equipment

Beta-Gamma sensitive stainless steel wall type tubes



T1100

Used in hand probes and foot deck of monitoring instruments. (T/A Models HSM-10A, HSC-1, and ALM-2X Hand & Shoe Monitors, and Model PPM-8 Portal Monitor).

Cathode size: length 3", diameter %", wall thickness 30-40 mg/cm².

Tube base: standard small 4 pin.

T1110

- Used in clothing probe of monitoring instruments & Model P-10 Probe Assembly.
- Cathode size: Same as T1100, except small 3 pin base.

T1120

Used in portal frame detectors of PPM-8 Portal Monitor and frame detectors of LIM-18 Laundry Monitor. Cathode size: length 7", diameter %", wall thickness 30-40 mg/cm². Tube base: standard small 4 pin.

T1140

Mica end-window type tubes

T1180



Alpha-beta sensitive.

Used in general radio-assay work and radiochemical analysis. (T/A Model LS-6 Lead Shield, T/A Model TM-8 Tube Mount in conjunction with LS-7 Lead Shield, and T/A Model P-8 Probe.)

End window diameter 13/32"; thickness 1.4-2 mg/cm².

Tube base standard medium 4 pin.

Alpha-beta-gamma sensitive. Used in general laboratory survey instruments where greater alpha and beta sensitivity are required. (T/A Model P-8X Probe in conjunction with a count rate meter or scaler.) End window diameter ²⁵/₃₂"; thickness

 $1.4-2 \text{ mg/cm}^2$.

Tube base connector standard single pin.

CIFICATIONS S P E

	T1090*	T1100	T1120	T1140	T1160**	T1180
OPERATING VOLTAGE: PLATEAU LENGTH IN EXCESS OF: SLOPE OF PLATEAU: STARTING VOLTAGE (0.3V Pulses): MAXIMUM COUNTING RATE:	900V DC 200 Volts 3 % per 100V 825V max. 1700 cps	900V DC 200 Volts 10% per 100V 825V max. 1700 cps	900V DC 200 Volts 10% per 100V 825V max. 1700 cps	900V DC 200 Volts 5% to 10% per 100V 825V max. 830 cps	1300V DC 250 Volts 1.5% per 100V 1180V max. 1100 cps	700V DC 180 Volts 10% per 100V 620V max. 1100 cps
BACKGROUND (Shielded with 2" lead and 1/s" aluminum): DEAD TIME (Approximate): OPERATING TEMPERATURE: OVERALL DIMENSIONS: LIFE EXPECTANCY IN COUNTS:	50 cpm max. 100 μ sec. 	50 cpm max. 100 µ sec. —55° to +75°C 5¾″ long, 1⅔″ dia. Unlimited	75 cpm max. 100 µ sec. 	50 cpm max. 200 µ sec. —55° to +75°C 4 ¹¹ ⁄ ₂ ″ long 1 ‰″ dia. Unlimited	50 cpm max. 150 µ sec. +15° to +50°C 4 ¹¹ ⁄ ₃₂ " long, 1¾" dia. Approx. 1.5 x 10 ⁸	75 cpm max. 150 µ sec. —55° to +75°C 6″ long, 1″ dia. Unlimited

*T1090 has organic quenching agent; base connector is single pin as in T1180.

**T1160 is identical to T1140, except it contains an organic quenching agent.

TECHNICAL ASSOCIATES

THE NU-TEC POCKET MONITOR ON MODEL PM-100

Ramin Longtor



NUCOR'S new, inexpensive, wide range fully transistorized pocket monitor is a Geiger-Mueller tube type radiation detector of the most modern design.

Packaged in a lightweight nonbreakable plastic case, it is the ideal instrument for carrying in a pocket or belt holder for health physics purposes.

SPECIFICATIONS

RANGES: Four scales: 0-10 mr/hr; 0-100 mr/hr; 0-1000 mr/hr; 0-100 r/hr.

ENERGY DEPENDENCE:

 \neq 15% for gamma or X-ray energies between 80KEV and 1.2 MEV.

SATURATION CHARACTERISTIC: Will not jam or saturate in any field.

OPERATION CHECK:

All electronic circuits can be checked by a test position on the function selector switch. Complete radiological operation can be checked with a small radioactive source, such as a luftimus ratch diat.

READOUT: 1-100" mater provides immediate and visual reador Rhine tack provided for earphone.

Hi-impac plastic. Waterproof and sock resistant.

SIZE: 2 - 5/8" x 4 - 1/6" x 1 - 3/8" deep.

WEIGHT: 10 ounces

POWER SUPPLY: One 4 voit battery provides more than 40 hours of continued operation.

MANUAL: Complete operating and instruction manual provided.

ACCESSORY KIT: Earphone and Radioactive Test Source Included.

PRICE COMPLETE: \$52.95 Includes Accessory Kit



CORPO

FROM: NUCLEAR CORPORATION OF AMERICA Instrument and Control Division 2 Richwood Place Denville, New Jersey

FOR RELEASE:

POCKET SIZE DETECTOR MONITORS RADIATION

The NU-TEC, a pocket size radiation detector, designed for industrial applications and as a personal survival aid in atomic war, has been marketed by Nuclear Corporation of America, Denville, New Jersey.

Fully transistorized, and inexpensive (\$52.95), the battery operated NU-TEC is a Geiger-Mueller tube type radiation detector sensitive to a wide range of gamma radiation energies.

Four ranges of 0 to 10/100/1000 mr/hr and 0 to 100 r/hr provide a direct and instanteous reading of the gamma or x-ray field present.

Packaged in a lightweight non-breakable plastic case, the "NU-TEC" is slightly larger than a pack of cigarettes and is the ideal instrument for carrying in a pocket or belt holder for health physics purposes.

An accessory kit containing an earphone and a radiation check source is included with the instrument. Used in situations where it is inconvenient to observe readings, the NU-TEC earphone can indicate changes in radiation levels. Radiation fields will produce a series of clicks in the earphone; the higher the radiation field, the greater the number of clicks.
A radiation check source can be used to obtain positive indications on the instrument. The check source provided with the NU-TEC utilizes a small harmless radiation field for checking the complete radiological response of the pocket monitor.

Further details and specification sheets may be obtained by writing to the Nuclear Corporation of America, Denville, New Jersey.

000

PRICE LIST

-

TECHNICAL ASSOCIATES

ς.

140 W. Providencia Ave. • Burbank, Calif.

Revised May 1, 1963

Subject to Change Without Notice

No. Description of Product Freston NEL Monitor Complete Statistics 168 PPM-82 Console for Model PPM-8 Portal Monitor Statistics Statistics <t< th=""><th>Bulletin</th><th></th><th>Model</th><th>Description of Deschart FO</th><th>Price B Burbank</th></t<>	Bulletin		Model	Description of Deschart FO	Price B Burbank
PERSONNEL MONITORS PERSONNEL MONITORS \$1650.00 168 PPM-8 Console for Model PPM-8 Portal Monitor 1295.00 168 LIM-18 Beta/Camma Laundry Monitor 1295.00 168 LIM-18 Beta/Camma Laundry Monitor 1150.00 164 HSM-10A Beta/Camma Laundry Monitor 365.00 207 FM-1 Frisker - Monitor (basic instrument with Rustrak Recorder and cabinet) 355.00 207 RAL-3 Remote Alarm for Frisker - Monitor 80.00 207 RAL-3 Remote Alarm for Frisker - Monitor 80.00 207 RAL-3 Remote Alarm for Frisker - Monitor 80.00 208 CP-TP-1A Cutie Pic Totem Pole Portable Survey 3 Range Meter, 0 to 5, 50 and 500 R/HR 440.00 185 18 Chamber Cutie Pic Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 225.00 186 CP-3 Cutie Pic Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 255.00 186 CP-3A Cutie Pic Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 225.00 186 CP-4A Cutie Pic Portable Survey	No.		No.	Description of Product F.O	.D. Durbank
168 PPM-8 Beta/Camma Portal Monitor complete				PERSONNEL MONITORS	
168 PPM-8C Console for Model PPM-8 Portal Monitor 1293.00 168 LIM-18 Beta/Camma Landry Monitor 1150.00 164 HSM-10A Beta/Camma Landry Monitor 1150.00 167 Finisker - Monitor (Basic instrument with Rustrak Recorder and cabinet) 355.00 168 LIM-18 Recorder and cabinet) 355.00 169 Cutie Pic Totem for Frisker - Monitor 80.00 160 Cutie Pic Totem Pole Portable Survey 3 Range Meter, 0 to 5, 5, and 50 R/HR 460.00 161 Cutie Pic Totem Pole Portable Survey 3 Range Meter, 0 to 50, 500 and 5000 R/HR 440.00 163 LA Chamber Cutie Pic Totem Pole Portable Survey 3 Range 0 to 50, 500 and 5000 R/HR 225.00 164 CP-3 Cutie Pic Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 295.00 166 CP-4 Cutie Pic Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 275.00 166 CP-4 Cutie Pic Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 225.00 166 CP-4 Cutie Pic Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 225.00 166 DS-58 <td>168</td> <td></td> <td>PPM-8</td> <td>Beta/Gamma Portal Monitor complete</td> <td>\$1650.00</td>	168		PPM-8	Beta/Gamma Portal Monitor complete	\$1650.00
168 LIA-18 Beta/Gamma Laundry Monitor 1190.00 184 HSA-10A Beta/Gamma Hand and Shce Monitor with external Beta/Gamma clothing probe. 3375.00 207 FM-1A Frisker - Monitor (basic instrument with Rustrak Recorder and cabinet) 355.00 207 FM-1A Frisker - Monitor (basic instrument with Rustrak Recorder and cabinet) 80.00 207 RAL-3 Remote Alarm for Frisker - Monitor 80.00 207 RAL-3 Remote Alarm for Frisker - Monitor 80.00 207 RAL-3 Remote Alarm for Frisker - Monitor 80.00 208 CP-TP-1A Cutie Pie Totem Pole Portable Survey 3 Range Meter, 0 to 5, 5, and 50 R/HR 460.00 185 IA Chamber Chamber only for CP-TP, range 0 to 5, 500 and 5000 R/HR 225.00 186 CP-3A Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 275.00 186 CP-4A Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 225.00 180 CP-4A Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 500 and 5000 MR/HR 225.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges	168		PPM-8C	Console for Model PPM-8 Portal Monitor	1295.00
184 HSM-10A Beta/Camma Hand and Shoe Monitor with external Beta/Camma Chiling probe 3575.00 207 FM-1 Frisker - Monitor (basic instrument exprobes and alarms) mounted in cabinet 365.00 207 RAL-3.2 Plug-in Alarm for Frisker - Monitor 80.00 207 RAL-3 Remote Alarm for Frisker - Monitor 80.00 207 RAL-3 Remote Alarm for Frisker - Monitor 80.00 208 CP-TP-1A Cutie Pie Totem Pole Portable Survey 3 Range Meter, 0 to 50, 500 and 500 R/HR 440.00 185 I A Chamber Chamber only for CP-TP, range 0 to 50, 500 and 500 MR/HR 225.00 186 CP-3 Cutie Pie Totable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 295.00 186 CP-3 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 275.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 275.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 225.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 225.00 180 Spl.1-7	168		LIM-18	Beta/Gamma Laundry Monitor	1150.00
207 FM-1 Frisker - Monitor (basic instrument less probes and alarms) mounted in cabinet	184		HSM-10A	Beta/Gamma Hand and Shoe Monitor with external Beta/Gamma clothing probe	3575.00
207 FM-1A Frisker - Monitor (basic instrument with Rustrak Recorder and cabinet) 515.00 207 RL3-2 Remote Alarm for Frisker - Monitor 80.00 207 RL3-3 Remote Alarm for Frisker - Monitor 80.00 208 SURVEY METERS 80.00 80.00 185 CP-TP-1A Cuite Pie Totem Pole Portable Survey 3 Range Meter, 0 to 50, 50.00 and 5000 R/HR 440.00 185 IB Chamber only for CP-TP, range 0 to 5, 5 and 50 R/H 225.00 186 CP-3 Cuite Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 225.00 186 CP-3A Cuite Pie Portable Survey Meter 3 Ranges 0 to 52, 250 and 2500 MR/HR 275.00 186 CP-3A Cuite Pie Portable Survey Meter 3 Ranges 0 to 52, 250 and 2500 MR/HR 275.00 186 CP-4A Cuite Pie Portable Survey Meter 3 Ranges 0 to 52, 250 and 2500 MR/HR 225.00 180 SRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 52, 250 and 2500 MR/HR 225.00 180 SRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 52, 250 and 2500 MR/HR 225.00 180 SRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 52, 250 and 2500 MR/	207		FM-1	Frisker - Monitor (basic instrument less probes and alarms) mounted in cabinet .	365.00
207 ALS-2 Plug-in Alarm for Frisker-Monitor 80.00 207 RAL-3 Remote Alarm for Frisker-Monitor 55.00 307 RAL-3 Remote Alarm for Frisker-Monitor 55.00 308 Cutie Pic Totem Pole Portable Survey 3 Range Meter, 0 to 5, 5, and 50 R/HR 460.00 185 C-TP-1-B Cutie Pic Totem Pole Portable Survey 3 Range Meter, 0 to 50, 500 and 5000 R/HR 225.00 185 1B Chamber Chamber only for CP-TP, range 0 to 50, 500, and 5000 R/HR 225.00 186 CP-3 Cutie Pic Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 225.00 186 CP-4 Cutie Pic Portable Survey Meter 3 Ranges 0 to 52, 250 and 2500 MR/HR 275.00 186 CP-4A Cutie Pic Portable Survey Meter 3 Ranges 0 to 52, 250 and 2500 MR/HR 275.00 180 RRL-7 Juno Portable Survey Meter 3 Ranges 0 to 52, 500 and 5000 MR/HR 225.00 180 HRL-7 Juno Portable Survey Meter 3 Ranges 0 to 252, 250 and 2500 MR/HR 225.00 180 HRL-7 Juno Portable Survey Meter 3 Ranges 0 to 252, 250 and 2500 MR/HR 225.00 180 HRL-7 Juno Portable Survey Meter 3 Ranges 0	207		FM-1A	Frisker - Monitor (basic instrument with Rustrak Recorder and cabinet)	515.00
207 RAL-3 Remote Alarm for Frisker - Monitor 55.00 208 SURVEY METERS 55.00 185 CP-TP-1A Cutie Pie Totem Pole Portable Survey 3 Range Meter, 0 to 5, 5, and 500 R/HR 460.00 185 CP-TP-1B Cutie Pie Totem Pole Portable Survey 3 Range Meter, 0 to 50, 500 and 5000 R/HR 225.00 185 1B Chamber Chamber only for CP-TP, range 0 to 50, 500 and 5000 RR/HR 225.00 186 CP-3 Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 295.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 275.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 225.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 02 and 2500 MR/HR 225.00 180 SRJ.7 Juno Portable Survey Meter 3 Ranges 0 to 25, 02 and 25,000 MR/HR 225.00 180 F6 Geiger Counter — Ranges 0-500, 0-5000 CPM 325.00 180 SRJ.7 Juno Portable Survey Meter 3 Ranges 0 20.00 OR/HR 205 F5-11 Super Sensitive Gamma Scintillation Survey Meter = Ranges 0.10 ov	207		ALS-2	Plug-in Alarm for Frisker - Monitor	80.00
SURVEY METERS 185 CP-TP-1A Cutie Pie Totem Pole Portable Survey 3 Range Meter, 0 to 5, 5, and 50 R/H 460.00 185 CP-TP-1B Cutie Pie Totem Pole Portable Survey 3 Range Meter, 0 to 50, 500 and 5000 R/HR 225.00 185 18 Chamber Chamber only for CP-TP, range 0 to 5, 5 and 500 R/H 225.00 186 CP-3 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 225.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 275.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 275.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 250, 250 and 2500 MR/HR 225.00 180 SRJ.7 Juno Portable Survey Meter 3 Ranges 0 to 250, 250 and 25000 MR/HR 225.00 180 F-6 Geiger Counter — Ranges 0-500, 0-5000 CPM 250 180 SRJ.7 Juno Portable Survey Meter 3 Ranges 0 250.00 MR/HR 250.00 205 FS-11 Super Sensitive Gamma Scintillation Survey Meter — Ranges .01 over background .1, 2, 5, 1.0, 2.5 MR/HR 250.00 260 206 Decade Scale	207		RAL-3	Remote Alarm for Frisker - Monitor	55.00
185 CP.TP-1A Cutie Pie Totem Pole Portable Survey 3 Range Meter, 0 to 5, 5, and 50 R/HR. 460.00 185 IA Chamber Chamber only for CP-TP, range 0 to 5, 500, and 5000 R/HR. 225.00 186 CP-3 Cutie Pie Totam Pole Portable Survey 3 Range Meter, 0 to 5, 500 and 5000 R/HR. 225.00 186 CP-3 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR. 295.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR. 275.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR. 275.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR. 275.00 180 SRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 25, 250 and 25000 MR/HR. 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 25, 250 and 25000 MR/HR. 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 25, 250 and 25000 MR/HR. 325.00 205 FS-11 Super Sensitive Gamma Scintillation Survey Meter 3 Ranges 0 to 25, 250 and 25000 MR/HR. 325.00 206 DS-58 Decade Scaler (1 microsecond) 965.00				SURVEY METERS	
185 CP-TP-1B Cutie Pie Totem Pole Portable Survey 3 Range Meter, 0 to 50, 500 and 5000 R/HR 245.00 185 1B Chamber Chamber only for CP-TP, range 0 to 50, 500, and 5000 R/HR 225.00 186 CP-3 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 295.00 186 CP-3 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 295.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 275.00 180 SRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 2500 MR/HR 275.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 2500 MR/HR 2325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 2500 MR/HR 2325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 2500 MR/HR 2325.00 180 F-6 Geiger Counter — Ranges 0-500, 0-5000 CPM 160.00 205 FS-11 Super Sensitive Gamma Scintillation Survey Meter — Ranges 0.1 over background 1, 2, 5, 1, 0, 2, 5 MR/HR 250.00 166 DS-58 Decade Scaler (1 microsecond) 965.00	185		CP-TP-1A	Cutie Pie Totem Pole Portable Survey 3 Range Meter, 0 to .5, 5, and 50 R/HR	460.00
185 1.A Chamber Chamber only for CP-TP, range 0 to 5, 5 and 50 R/H 225.00 185 11B Chamber Chamber only for CP-TP, ranges 0 to 50, 500, and 5000 R/HR 205.00 186 CP-3 Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 5000 MR/HR 295.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 5000 MR/HR 275.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 5000 MR/HR 275.00 180 SRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 25, 0200 and 5000 MR/HR 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 25000 MR/HR 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 250.00 MR/HR 325.00 205 FS-11 Super Sensitive Gamma Scintillation Survey Meter — Ranges .01 over background, 1, 25, 1.0, 2.5 MR/HR 450.00 206 DS-58 Decade Scaler (1 microsecond) 965.00 181 RM-8 Log-Linear Ratemeter 810.00 190 SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 765.00 <	185		CP-TP-1B	Cutie Pie Totem Pole Portable Survey 3 Range Meter, 0 to 50, 500 and 5000 R/H	R 440.00
185 18 Chamber Chamber only for CP-TP, ranges 0 to 50, 500, and 5000 R/HR 205.00 186 CP-3 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 295.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 275.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 275.00 180 SRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 250,00 MR/HR 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 25,000 MR/HR 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 25,000 MR/HR 325.00 205 FS-11 Super Sensitive Gamma Scintillation Survey Meter — Ranges .01 over background .1, .25, 1.0, 2.5 MR/HR 450.00 206 DS-5B Decade Scaler (1 microsecond) 915.00 206 DS-9 Decade Scaler (1 microsecond) 915.00 207 SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 810.00 181 RM-8 Log-Linear Ratemeter) 120.00 375.00	185		1A Chamber	Chamber only for CP-TP, range 0 to .5, 5 and 50 R/H	225.00
186 CP-3 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 295.00 186 CP-3A Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 295.00 186 CP-4A Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 275.00 186 CP-4A Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 275.00 180 SRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 325.00 205 FS-11 Super Sensitive Gamma Scintillation Survey Meter — Ranges .01 over background .1, .25, 1.0, 2.5 MR/HR 450.00 206 DS-5B Decade Scaler (1 microsecond) 965.00 915.00 190 SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 765.00 190 SA-20 Single Channel Pulse Height Analyzer (1500 and 3000 volts) 375.00 190 SA-30 Spectrum Scintillation Detector, counte	185		1B Chamber	Chamber only for CP-TP, ranges 0 to 50, 500, and 5000 R/HR	205.00
186 CP-3A Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 295.00 186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 5000 MR/HR 275.00 180 SRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 2500 MR/HR 325.00 205 F5-11 Super Sensitive Gamma Scintillation Survey Meter — Ranges .01 over background .1, .25, 1.0, 2.5 MR/HR 450.00 206 ANALYTICAL AND COUNTING INSTRUMENTS 825.00 166 DS-58 Decade Scaler (5 microsecond) 915.00 200 DS-9 Decade Scaler (1 microsecond) 915.00 190 SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 810.00 190 SA-20 Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 355.00 190 SA-20 Single Channel Pulse Height Analyzer with built-in Linear Amplifier scaler or Ratemeter) 765.00 1	186		CP-3	Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR	295.00
186 CP-4 Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR 275.00 180 SRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 275.00 180 SRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 25000 MR/HR 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 25.000 MR/HR 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 25.000 MR/HR 325.00 205 FS-11 Super Sensitive Gamma Scintillation Survey Meter — Ranges .01 over background .1, .25, 1.0, 2.5 MR/HR 450.00 206 DS-5B Decade Scaler (5 microsecond) 915.00 200 DS-9 Decade Scaler (1 microsecond) 965.00 181 RM-8 Log-Linear Ratemeter 810.00 190 SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 765.00 190 SA-20 Single Channel Pulse Height Analyzer with built-in Linear Amplifier or Ratemeter 755.00 190 MS-1 Mobil Stand for Scintillation Ductor, counter-weighted and with shelf for scaler or Ratemeter 325.00	186		CP-3A	Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR	295.00
186 CP-4A Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR 275.00 180 SRJ.7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 25000 MR/HR 325.00 180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 25,000 MR/HR 325.00 205 F5-11 Super Sensitive Gamma Scintillation Survey Meter — Ranges .01 over background .1, .25, 1.0, 2.5 MR/HR 450.00 206 PS-58 Decade Scaler (5 microseconds) 825.00 200 DS-59 Decade Scaler (1 microsecond) 915.00 200 DS-59 Decade Scaler (1 microsecond) 965.00 200 DS-9 Decade Scaler (1 microsecond) 965.00 200 SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 765.00 169 SA-20 Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 765.00 178 RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Resetter Social instrument supplied in 19" panel for rack mounting unters or Ratemeter or Ratemeter 350.00 <t< td=""><td>186</td><td></td><td>CP-4</td><td>Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR</td><td> 275.00</td></t<>	186		CP-4	Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR	275.00
180 SRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR	186		CP-4A	Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR	275.00
180 HRJ-7 Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 25,000 MR/HR 325.00 205 F-6 Geiger Counter — Ranges 0-500, 0-5000 CPM 160.00 205 FS-11 Super Sensitive Gamma Scintillation Survey Meter — Ranges .01 over background 450.00 206 FS-58 Decade Scaler (5 microseconds) 825.00 206 DS-58 Decade Scaler (1 microsecond) 915.00 200 DS-9 Decade Scaler (1 microsecond) 965.00 181 RM-8 Log-Unrear Ratemeter 810.00 190 SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier 765.00 169 SA-20 Single Channel Pulse Height Analyzer 685.00 190 PS-31 Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-11 Regulated High Voltage Power Supply 2 Ranges (1500 and with shelf for scaler or ratemeter 350.00 190 MS-11 Regulated High Voltage Power Supply 2 Ranges (1500 and with shelf for scaler or ratemeter 350.00 190 MS-11 Regulated High Voltage Power Supply 2 Ranges (1500 and with shelf for scaler or ratemeter 350.00 190	180		SRJ-7	Juno Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR	325.00
F-6 Geiger Counter — Ranges 0-500, 0-5000 CPM	180		HRJ-7	Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 25,000 MR/HR	325.00
205 FS-11 Super Sensitive Gamma Scintillation Survey Meter — Ranges .01 over background 1, .25, 1.0, 2.5 MR/HR 450.00 ANALYTICAL AND COUNTING INSTRUMENTS ANALYTICAL AND COUNTING INSTRUMENTS 825.00 166 DS-58 Decade Scaler (5 microseconds) 915.00 200 DS-9 Decade Scaler (1 microsecond) 915.00 200 DS-9 Decade Scaler (1 microsecond) 915.00 181 RM-8 Log-Linear Ratemeter 810.00 190 SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 745.00 190 SA-20 Single Channel Pulse Height Analyzer 685.00 200 PA-6B Pre-Amplifier 120.00 178 RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Mobile Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 190 MS-1 Prices shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-to- Service" Console cabinet add \$30.00 to prices shown. 95.00 203 ** PT-1 Predetermined Timer (Long counting intervals			F-6	Geiger Counter — Ranges 0-500, 0-5000 CPM	160.00
ANALYTICAL AND COUNTING INSTRUMENTS 166 DS-5B Decade Scaler (5 microseconds) 825.00 166 DS-5BA Decade Scaler (1 microsecond) 915.00 200 DS-9 Decade Scaler (1 microsecond) 965.00 181 RM-8 Log-Linear Ratemeter 810.00 190 SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 765.00 169 SA-20 Single Channel Pulse Height Analyzer 685.00 190 SS-30 Spectrum Scanner 475.00 200 PA-6B Pre-Amplifier 120.00 178 RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 275.00 190 MS-1 Mobile Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 * Prices shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-tor Service" Console cabinet add \$30.00 to prices shown. 95.00 203 ** PT-1 Predetermined Timer (Long counting intervals 999 minutes) 95.00 203 ** PT-2 Predetermined Timer (Lobe	205		FS-11	Super Sensitive Gamma Scintillation Survey Meter — Ranges .01 over backgrou .1, .25, 1.0, 2.5 MR/HR	nd 450.00
166 DS-58 Decade Scaler (5 microsecond) 825.00 166 DS-58A Decade Scaler (1 microsecond) 915.00 200 DS-9 Decade Scaler (1 microsecond) 965.00 181 RM-8 Log-Linear Ratemeter 810.00 190 SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 765.00 166 SA-20 Single Channel Pulse Height Analyzer 685.00 190 SS-30 Spectrum Scanner 475.00 200 PA-68 Pre-Amplifier 120.00 178 RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Redulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 355.00 190 MS-1 Preces shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-to-Service" Console cabinet add \$30.00 to prices shown. 203 ** PT-1 Predetermined Timer (Long counting intervals 120 seconds) 95.00 203 ** PT-2 Predet				ANALYTICAL AND COUNTING INSTRUMENTS	
166 DS-5B Decade Scaler (1 microseconds) 022.50 166 DS-5BA Decade Scaler (1 microsecond) 915.00 200 DS-9 Decade Scaler (1 microsecond) 965.00 181 RM-8 Log-Linear Ratemeter 810.00 190 SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 765.00 169 SA-20 Single Channel Pulse Height Analyzer 685.00 200 PA-6B Pre-Amplifier 775.00 200 PA-6B Pre-Amplifier 120.00 178 RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Mobile Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 * Prices shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-to-Service" Console cabinet add \$30.00 to prices shown. 95.00 203 ** PT-1 Predetermined Timer (Short counting intervals 120 seconds) 95.00 203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 203 ** PT		•	DC 50	ANALITICAL AND COUNTING INSTRUMENTS	825 00
166 DS-SDA Decade Scaler (1 microsecond) 713.00 200 DS-9 Decade Scaler (1 microsecond) 965.00 181 RM-8 Log-Linear Ratemeter 810.00 190 SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 765.00 169 SA-20 Single Channel Pulse Height Analyzer 685.00 190 SS-30 Spectrom Scanner 475.00 200 PA-6B Pre-Amplifier 120.00 178 RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Mobile Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 * Prices shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-to-Service" Console cabinet add \$30.00 to prices shown. 203 ** PT-1 Predetermined Timer (Long counting intervals 999 minutes) 95.00 203 ** PT-2 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 </td <td>100</td> <td>•</td> <td>DS-DB</td> <td>Decade Scaler (5 microseconds)</td> <td>915.00</td>	100	•	DS-DB	Decade Scaler (5 microseconds)	915.00
200 - Ds-9 Decade Scaler (1 microsecond) 903.00 181 RM-8 Log-Linear Ratemeter 810.00 190 * SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 765.00 169 * SA-20 Single Channel Pulse Height Analyzer 685.00 190 * SS-30 Spectrum Scanner 475.00 200 PA-68 Pre-Amplifier 120.00 178 * RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Mobile Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 190 MS-1 Prices shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-to-Service" Console cabinet add \$30.00 to prices shown. 203 ** PT-1 Predetermined Timer (Long counting intervals 299 minutes) 95.00 203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 203 ** PT-6 Electronic Timer 500.00 203 **	100	÷	DS-SBA	Decade Scaler (1 microsecond)	915.00
181 - KM-3 Log-Linear Katemeter 610.00 190 * SM-10 Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Amplifier and Ratemeter) 765.00 169 * SA-20 Single Channel Pulse Height Analyzer 685.00 190 * SS-30 Spectrum Scanner 475.00 200 PA-6B Pre-Amplifier 120.00 178 * RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Mobile Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 190 MS-1 Prices shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-to-Service" Console cabinet add \$30.00 to prices shown. 203 ** PT-1 Predetermined Timer (Long counting intervals 999 minutes) 95.00 203 ** PT-2 Predetermined Timer (Short counting intervals 120 seconds) 95.00 203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 203 ** PT-6 Electronic Timer 500.00	200	•	05-9	Decade Scaler (1 microsecond)	905.00 810.00
190 * SM-10 Spectrometer (Single Channel Pulse Height Analyzer with Dulit-in Linear Amplifier and Ratemeter) 765.00 169 * SA-20 Single Channel Pulse Height Analyzer 685.00 190 * SS-30 Spectrum Scanner 475.00 200 PA-6B Pre-Amplifier 120.00 178 * RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Mobife Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 190 MS-1 Mobife Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 * Prices shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-to- Service" Console cabinet add \$30.00 to prices shown. 95.00 203 ** PT-1 Predetermined Timer (Long counting intervals 999 minutes) 95.00 203 ** PT-2 Predetermined Timer (Long counting intervals 120 seconds) 95.00 203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 203 ** PT-6 <td< td=""><td>181</td><td>-</td><td>RM-8</td><td>Log-Linear Katemeter</td><td> 010.00</td></td<>	181	-	RM-8	Log-Linear Katemeter	010.00
169 • SA-20 Single Channel Pulse Height Analyzer 685.00 190 * SS-30 Spectrum Scanner 475.00 200 PA-6B Pre-Amplifier 120.00 178 * RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Mobife Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 190 MS-1 Mobife Stand for Scintillation Detector, counter-weighted in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-to-Service" Console cabinet add \$30.00 to prices shown. 203 ** PT-1 Predetermined Timer (Long counting intervals 999 minutes) 95.00 203 ** PT-2 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 200 PT-6 Electronic Timer 500.00 203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 203 ** PT-4 Extended Count Register (4 digit) 65.00 203 ** CR-2 Extended Count Register (4 digit) 147.50 204 ** CR-2 Extended Count Register (4 digit) 65.00 203 ** C	190	<u>^</u>	SM-10	and Ratemeter)	765.00
190 * SS-30 Spectrum Scanner 475.00 200 PA-6B Pre-Amplifier 120.00 178 * RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Mobile Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 190 MS-1 Mobile Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 190 MS-1 Prices shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-to-Service" Console cabinet add \$30.00 to prices shown. 203 ** PT-1 Predetermined Timer (Long counting intervals 999 minutes) 95.00 203 ** PT-2 Predetermined Timer (Short counting intervals 999 minutes) 95.00 203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 203 ** PT-6 Electronic Timer 500.00 203 ** CR-2 Extended Count Register (4 digit) 65.00 203 ** CR-2 Extended Count Register (4 digit) 65.00 204 PT-6 Electronic Timer 500.00 2	169	٠	SA-20	Single Channel Pulse Height Analyzer	685.00
200 PA-6B Pre-Amplifier 120.00 178 * RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Mobife Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 * Prices shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-to-Service" Console cabinet add \$30.00 to prices shown. \$50.00 203 ** PT-1 Predetermined Timer (Long counting intervals 999 minutes) 95.00 203 ** PT-2 Predetermined Timer (Short counting intervals 120 seconds) 95.00 203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 203 ** CR-2 Extended Count Register (4 digit) 65.00 203 ** PT-3 Predetermined Timer 500.00 203 ** PT-3 Predetermined Timer 500.00 203 ** CR-2 Extended Count Register (4 digit) 65.00 204 PT-6 Electronic Timer 500.00 50.00 205 Five microsecond plug-in de	190	*	SS-30	Spectrum Scanner	475.00
178 * RHV-1B Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) 375.00 190 MS-1 Mobile Stand for Scintillation Detector, counter-weighted and with shelf for scaler or Ratemeter 350.00 * Prices shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-to-Service" Console cabinet add \$30.00 to prices shown. \$95.00 203 ** PT-1 Predetermined Timer (Long counting intervals 999 minutes) 95.00 203 ** PT-2 Predetermined Timer (Short counting intervals 120 seconds) 95.00 203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 203 ** CR-2 Extended Count Register (4 digit) 65.00 166 DU-1 One microsecond plug-in decade unit 147.50	200		PA-6B	Pre-Amplifier	120.00
 * Prices shown are for instruments supplied in 19" panel for rack mounting unless otherwise noted. For single instrument mounted in T/A Model C-875 "Easy-to-Service" Console cabinet add \$30.00 to prices shown. COUNTING INSTRUMENT ACCESSORIES 203 ** PT-1 Predetermined Timer (Long counting intervals 999 minutes)	178 190	*	RHV-1B MS-1	Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts) Mobile Stand for Scintillation Detector, counter-weighted and with shelf for sca or Ratemeter	375.00 ller 350.00
COUNTING INSTRUMENT ACCESSORIES 203 ** PT-1 Predetermined Timer (Long counting intervals 999 minutes) 95.00 203 ** PT-2 Predetermined Timer (Short counting intervals 120 seconds) 95.00 203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 200 PT-6 Electronic Timer 500.00 203 ** CR-2 Extended Count Register (4 digit) 65.00 166 DU-1 One microsecond plug-in decade unit 147.50 166 DU-5 Five microsecond plug-in decade unit 57.50				* Prices shown are for instruments supplied in 19" panel for rack mounting unl otherwise noted. For single instrument mounted in T/A Model C-875 "Easy- Service" Console cabinet add \$30.00 to prices shown.	ess to-
203**PT-1Predetermined Timer (Long counting intervals 999 minutes)95.00203**PT-2Predetermined Timer (Short counting intervals 120 seconds)95.00203**PT-3Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes)100.00200PT-6Electronic Timer500.00203**CR-2Extended Count Register (4 digit)65.00166DU-1One microsecond plug-in decade unit147.50166DU-5Five microsecond plug-in decade unit57.50				COUNTING INSTRUMENT ACCESSORIES	
203 ** PT-2 Predetermined Timer (Short counting intervals 120 seconds) 95.00 203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 200 PT-6 Electronic Timer 500.00 203 ** CR-2 Extended Count Register (4 digit) 65.00 166 DU-1 One microsecond plug-in decade unit 147.50 166 DU-5 Five microsecond plug-in decade unit 57.50	203	**	PT-1	Predetermined Timer (Long counting intervals 999 minutes)	95.00
203 ** PT-3 Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes) 100.00 200 PT-6 Electronic Timer 500.00 203 ** CR-2 Extended Count Register (4 digit) 65.00 166 DU-1 One microsecond plug-in decade unit 147.50 166 DU-5 Five microsecond plug-in decade unit 57.50	203	* *	PT-2	Predetermined Timer (Short counting intervals 120 seconds)	95.00
200 PT-6 Electronic Timer 500.00 203 ** CR-2 Extended Count Register (4 digit) 65.00 166 DU-1 One microsecond plug-in decade unit 147.50 166 DU-5 Five microsecond plug-in decade unit 57.50	203	* *	PT-3	Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes)	100.00
203 ** CR-2 Extended Count Register (4 digit) 65.00 166 DU-1 One microsecond plug-in decade unit 147.50 166 DU-5 Five microsecond plug-in decade unit 57.50	200		PT-6	Electronic Timer	500.00
166 DU-1 One microsecond plug-in decade unit	203	* *	CR-2	Extended Count Register (4 digit)	65.00
166 DU-5 Five microsecond plug-in decade unit	166		DU-1	One microsecond plug-in decade unit	147.50
	166		DU-5	Five microsecond plug-in decade unit	57.50

** These accessories will be supplied in a cabinet or panel mounted at customer's option.

PAGE ONE

dideen and the second second

-

Bulletin	Model	P	rice
No.	No.	Description of Product F.O.B.	Burbank
		AUTOMATED COUNTING SYSTEMS	
189	MST-201	Multi-Sample Tabulator, 2" sample diameter, single channel complete with one DS-6 Scaler and Victor digital printer (less detector)	6420.00
189	MST-202	Multi-Sample Tabulator, 2" Sample diameter, dual channel complete with two DS-6 Scalers and Victor digital printer (less detector)	7250.00
189	MST-301	Multi-Sample Tabulator, 3" sample diameter, single channel complete with one DS-6 Scaler and Victor digital printer (less detector)	7670.00
189	MST-302	Multi-Sample Tabulator, 3" sample diameter, dual channel complete with two DS-6 Scalers and Victor digital printer (less detector)	8500.00
189	DS-6	Discriminating Decade Scaler for use with all model Multi-Sample Tabulator series	1050.00
189	MST-GF-2	Gas Flow Detector for all models of Multi-Sample Tabulator including lead shielding	440.00
189	MST-GM-2	Geiger Tube Detector with 1.860 dia. Window x 1.4 to 2 mg/cm ² thick for all models of Multi-Sample Tabulator including lead shielding	380.00
189	MST-AS-2	Alpha Scintillation Detector with 3" diameter zinc sulfide phosphor for all models of Multi-Sample Tabulator	778.00
189	MST-BS-2	Beta Scintillation Detector with 21/2" diameter Stilbene crystal for all models of Multi-Sample Tabulator	1075.00
189	MST-GS-2	Gamma Scintillation Detector with 2" diameter x 1" thick Nal crystal for all models of Multi-Sample Tabulator including lead shielding	685.00
201	GR-1	Gas Regulator for gas flow detectors used with Multi-Sample Tabulator and model GD-6 Flow Counter	47.50
		CTAD CEDIC CRECTROMETRY SYSTEMS	
171	Polaris	Includes SM-10 and RHV-1B mounted in C-1400 cabinet, complete with inter-	1100.00
171	Lyra	connecting cables and system test Includes SM-10 and DS-5B mounted in C-1750 cabinet, complete with inter-	1180.00
171	Arcturus	connecting cables and system test Includes SM-10_RHV-1B_SS-30 and Texas Instruments Rectilinear Recorder Model	1635.00
	,	FRRIM-A16T in P-1 Panel Mount installed in C-3150 cabinet, with intercon- necting cables and system test	2250.00
171	Orion	Includes SM-10, DS-9 and PT-6 mounted in C-2100 cabinet complete with cables	2280.00
		NOTE: For special applications requiring other combinations of T/A instruments please request a quotation.	
		SCINTILLATION DETECTORS AND ACCESSORIES	
202	SD-1	Scintillation Detector with $1'' \times 1''$ Nal crystal, complete with cables	455.00
202	SD-2	Scintillation Detector with 2" x 2" Nal crystal, complete with cables	765.00
202	SD-2W	Well-type Scintillation Detector with 13/4" x 2" Nal crystal (5/8" x 11/2" well) complete with cables for use with LS-8 or LS-8X Shields	695.00
202	SD-3	Scintillation Detector with 11/2" Alpha Phosphor, complete with cables	495.00
202	SD-4	Scintillation Detector with 11/2" Beta Anthracene crystal, complete with cables	545.00
183	SD-6-1	Scintillation Detector with transistorized pre-amplifier, 2" photomultiplier tube and 1" x 1" Nal crystal, complete with cables. (For use in LS-6 Lead Shield with extension or LS-6) Universal Shield	445 00
183	SD-6-1-5	Scintillation Detector with transistorized pre-amplifier 2" photomultiplier tube and	
165	30-0-1.3	$1/2'' \times 1''$ Nal crystal, complete with cables. (For use in LS-6 Lead Shield with extension or LS-66) Universal Shield	575.00
183	SD-6-2	Scintillation Detector with transistorized pre-amplifier, 2" photomultiplier tube and	
		2" x 2" Nal crystal, complete with cables. (For use in LS-6 Lead Shield with extension or LS-66) Universal Shield	765.00
183	SD-6W	Well Scintillation Detector with transistorized pre-amplifier, 2" photomultiplier tube and $134'' \times 2"$ Nal crystal ($56'' \times 112''$ well) complete with cables. (For use	<u> 695 00</u>
202		m worde Universal Smith	30.00
202		Crystal Housing 50-1 and 50-0-1 Scintillation Detector	30.00
183		Crystal mousing for SD-0-1.5 scintiliation Detector	30.00
202	CH-2	Crystal mousing for SD-22 and SD-0-2 scintillation Detector	30.00
202	CH-2W	Crystal Housing for SU-2W and SU-6W Scintillation Detector	30.00
202	NS-1	Nose Shield for Model SD-1 Detector	95.00
202	NS-2	Nose Shield for Model SD-2 Detector	115.00
202	C-1a	Collimator Type A (20° Flat Field) for SD-1 Detector	35.00
202	C-2a	Collimator Type A (20° Flat Field) for SD-2 Detector	45.00
202	C-16	Collimator Type B (Straight Bore) for SD-1 Detector	40.00
202	С-2ь	Collimator Type B (Straight Bore) for SD-2 Detector	50.00

PAGE TWO

L

-

L

Bulletin	Model		Price
No.	No.	Description of Product F.O.	B. Burbank
202	C-1c	Collimator Type C (Forcusing) for SD-1 Detector	90.00
202	C-2c	Collimator Type C (Honeycomb) for SD-2 Detector	
202	AS-11	Alpha Phosphor (1") for SD-1 and SD-6-1, complete with mount	79.50
202	AS-12	Alpha Phosphor $(1\%')$ for SD-2 and SD-6-2 complete with mount	
202	BS-11	Beta Stillege crystal $(1'')$ for SD-1 and SD-6-1 complete with mount	79.50
202	03-11	Bota Stillbong crystal (11/611) for SD-2 and SD 6-2, complete with mount	95.00
202	D3-12	beta Stilbene crystal (172 / 101 50-2 and 50-6-2, complete with modification	75.00
		GAS FLOW DETECTORS AND G.M. PROBES	
182	GD-6	Gas Flow Counter with built-in transistorized pre-amplifier (for use with LS-6 Lea Shield) includes extra window	id 395.00
204	P-7	Beta/Gamma Probe Assembly with T-1090 Tube, cable and connector	39.50
204	P-8	End Window Probe complete with T-1140 Tube, cable and connector	100.00
204	P-10	Beta/Gamma Probe Assembly with T-1110 G.M. Tube, cable and connector	
204	P-AS-2	Alpha Scintillation Probe	315.00
201	GR-1	Gas Regulator for Gas Flow Detectors	47.50
204	AMS-2	Beta/Gamma Area Monitoring Station	
		COUNTED THREE	
	7 1000	COUNTER TUBES	~ -
144	1-1090	other Beta/Gamma applications	or 9.00
144	T 1100	This wall type 3" long used in Models PPM-8 and HSM-10A	27.00
144	T-1100	Thin wall 3" long used in P-10 Probe	27.00
144	T 1120	This wall tree 7" long used in Model DDM 9	53.00
144	T-1120	thin wait type 7 long used in Model Fride Maunta and D.9 Broba	59.50
144	1-1140	Mica End Window used in TM-7, TM-8 Tube Mounts and P-6 Probe	50.50
144	1-1160	Mica End Window (Organic Quenched) used in TM-7, TM-8 Tube Mounts ar P-8 Probe	67.00
144	T-1180	Mica End Window, Alpha/Beta/Gamma sensitive	61.00
		IONIZATION CHAMBERS	
144	IC-1	Ionization Chamber with windows (Beta/Gamma)	115.00
144	IC-2	Ionization Chamber — Solid Wall (Gamma)	105.00
		LEAD SHIELDS	
191	LS-1	Complete with Type ST-1 Sample Tray, Type SK-1 Socket and Type TM-1 Tul Mount	oe 185.00
191	LS-2A	Standard with Type TM-6 Tube Mount and Sample Tray Holder, Type ST-2A Samp Tray and 10 Type PL-2 Aluminum Planchets	le 200.00
191	LS-2B	Micrometric with 10 Type PL-1 Planchets	285.00
192	LS-4A	For liquids complete with Type TM-1 Tube Mount, Type SK-1 Socket, and or Type ST-5 Marinelli Beaker	ne 200.00
102		Nodified to accommodate special tubes with aluminum inner liner	235.00
192	L3-4D	For selide several to accommodate special tubes with administration inter the several to accommodate special tubes with administration inter the several tubes are tubes with administration inter the several tubes are tubes and tubes with administration inter tubes are tubes are tubes and tubes are tubes and tubes are tubes a	233.00
192	LS-4C	ST-6 Ore Container and 100 Type ST-6A Paper Sleeves	225.00
192	LS-5	Complete with Type TM-3 Tube Mount, Type SK-2 Socket, Type STH-1 Samp Tray Holder, Type ST-4A Sample Tray, and 10 Type PL-2 Aluminum Planche	le ts 275.00
193	LS-6	Complete with Type TM-7 Tube Mount, one Type ST-3B Sample Tray and 10 Ty PL-3 Aluminum Planchets	pe 295.00
193	LS-6X8	Eight-inch height extension for LS-6 shield used with SD-1, SD-2 and SD-6 Scintill tion Detector	a-
193	LS-6L	Lid with 3" opening and retaining collar for Models SD-1 and SD-2 Scintill tion Detectors	a-
103	1 S-6P	Plug for 3" opening in 15-61	12.50
204	[S_KA	Ton section of 15.66 for the nurness of converting a model 15.6 to universal	265.00
200	10 44	Universal Land Shield includes Detector Mount land ring for well detector as	200.00 nd
200	L2-00	splash pan	540.00
194	LS-7A	Multi-Purpose Type complete with SD-1 Scintillation Detector, three ST-7 Samp Trays and 10 PL-3 Planchets	le 780.00
194	LS-7B	Multi-Purpose Type complete with SD-2 Scintillation Detector, three ST-7 Samp Trays and 10 PL-3 Planchets	ole 1090.00
194	LS-7C	Multi-Purpose Type complete with TM-8 Tube Mount, three ST-7 Sample Tra and 10 PL-3 Planchets	ys 385.00

Construction of the second second

Bulletin No.	Model No.	Pr Description of Product F.O.B.	ice Burbank
194	LS-7M	Multi-Purpose Type — Shield only, with adapter ring, three ST-7 Sample Trays and 10 PL-3 Planchets	325.00
195	LS-8	For well counting (shield only, less detector)	265.00
195	LS-8W	For well counting, complete with SD-2W Scintillation Detector and $1\frac{3}{4}$ x 2" well crystal	960.00
195	LS-8X	For well counting with $2^{\prime\prime}$ shielding in all directions (Shield only less Detector)	350.00
195	LS-8WX	For well counting with 2" shielding in all directions complete with SD-2W Scintillation Detector and $13/4$ " x 2" well crystal	1035.00

LABORATORY ACCESSORIES

and the second sec

198	AB-2	Absorber Set		• • • • • • • • • • • • •	85.00
199	ST-1	Sample Tray		·····	6.50
199	ST-2, \ST-3,	ST-3X, ST-3A, ST-3B, ST-4, ST-4A Sample Trays		•••••	3.50
1 99	ST-5	Marinelli Beaker		· · · · · · · · · · · · · · · · · · ·	4.25
199	ST-6	Ore Container			12.50
199	ST-6A	Paper Sleeves for Ore Container		per	C 2.00
199	ST-7	Sample Tray for LS-7		·	2.50
199	STH-1	Sample Tray Holder with tray and ten planchets	•••••		13.50
199	SK-1, SK-2,	SK-3 Sockets		••••	12.00
199	тм-1, тм-2	, TM-3, TM-4 Tube Mount		••••••	15.00
199	тм-5	Tube Mount (micrometric for LS-2 Lead Shield)			100.00
199	ТМ-6	Tube Mount and Sample Tray Holder	••••••	····	22.50
199	TM-7	Tube Mount and Sample Tray Holder		••••••	27.50
199	тм-8	Tube Mount and Adapting Ring for LS-7 Multi-Purpose Shield		.	70.00
206	SDM-9	Mount and Support Ring when SD-6 is used in LS-66 Shield			27.50
197	LC-1	Lead Source Container for four radium needles		•••••	45.00
197	LC-3	Lead Source Container 1" lead walls, $1\frac{1}{2}$ " x 3" inside dimensions			40.00
197	LC-3A	Lead Source Container $13/8''$ lead walls $3/4'' \times 21/4''$ inside dimensions			47.50
197	LC-4	Lead Source Container $1\frac{1}{2}$ " lead walls $1\frac{1}{2}$ " x 3" inside dimensions			60.00
197	LC-4A	Lead Source Container $1\frac{3}{4}$ " lead walls $\frac{3}{4}$ " x $2\frac{1}{2}$ " inside dimensions			67.50
197	LC-9	Mobile Source Unit — compartment 3" x 6" in shield of 2" lead.	Мони	nted c	m
		Lead Brick $2'' \times 3'' \times 6''$			9.00
		Lead Brick, $2 \times 6 \times 6$			12 50
					14.00
		Lond Details and advantage of the	n		
100		Lead Brick, extruded shapes	Prices	upon	application
199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned	Prices per per	upon C M	application 3.25 30.00
199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned	Prices . per per per	upon C M 5M	application 3.25 30.00 27.00/M
199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum	Prices per per per	upon C M 5M C	application 3.25 30.00 27.00/M 3.25 30.00
199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum	Prices per per per per per per	upon C M 5M C M 5M	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M
199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum	Prices - per per per per per per per	upon C 5M C M 5M C	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75
199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum	Prices per per per per per per per per	ироп С 5М С 5М С 5М С М	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00
199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum	Prices per per per per per per per per per	ироп С 5 М С М 5 М С М 5 М 5 М	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M
199 199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper	Prices per per per per per per per per per	ироп С 5 М 5 М С М 5 М С М 5 М С М 5 М С	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25
199 199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper	Prices - per per per - per per - per per per per per	ирол С М 5 М 5 С М 5 М 5 М 5 М	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 47.27
199 199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper	Prices - per per per per per per per per	upon C M 5M C M 5M C M 5M C M 5M C M 5M	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 43.50/M
199 199 199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper Planchets, Type PL-3C, Copper	Prices - per per per per per per per per per per	900 90 90 90 90 90 90 90 90 90 90 90 90	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 43.50/M 6.00
199 199 199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper Planchets, Type PL-3C, Copper	Prices per per per per per per per per per per per	900 50 50 50 50 50 50 50 50 50 50 50 50 5	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 43.50/M 6.00 54.00 54.00 50.00/M
199 199 199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper Planchets, Type PL-3C, Copper Planchets, Type PL-3C, Stainless Steel	Prices per per per per per per per per per per	upon しんか しんか しんか しんか しんか しんか	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 43.50/M 6.00 54.00 50.00/M 9.00
199 199 199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper Planchets, Type PL-3C, Copper Planchets, Type PL-3C, Stainless Steel	Prices per per per per per per per per	upon しんか しんか しんか しんか しんか しん	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 43.50/M 6.00 54.00 50.00/M 9.00 81.00
199 199 199 199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper Planchets, Type PL-3C, Copper Planchets, Type PL-3C, Stainless Steel	 Prices per 	upon し んち し んち し んち し んち し んち し んち	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 43.50/M 6.00 54.00 50.00/M 9.00 81.00 73.00/M
199 199 199 199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper Planchets, Type PL-3C, Copper Planchets, Type PL-3S, Stainless Steel Planchets, Type PL-3S, Stainless Steel	Prices - per per per per per per per per	upon し へ 55 し ん 55 し	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 43.50/M 6.00 54.00 50.00/M 9.00 81.00 73.00/M 12.00
199 199 199 199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper Planchets, Type PL-3C, Copper Planchets, Type PL-2S, Stainless Steel Planchets, Type PL-3S, Stainless Steel	 Prices per per per per per per per per per per	u つ A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 43.50/M 6.00 54.00 50.00/M 9.00 81.00 73.00/M 12.00 108.00
199 199 199 199 199 199		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper Planchets, Type PL-3C, Copper Planchets, Type PL-3C, Stainless Steel Planchets, Type PL-2S, Stainless Steel	 Prices per per per per per per per per per per	u つ A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5 し A 5 U	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 43.50/M 6.00 54.00 50.00/M 9.00 81.00 73.00/M 12.00 108.00 98.00/M
 199 		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper Planchets, Type PL-3C, Copper Planchets, Type PL-3S, Stainless Steel Planchets, Type PL-3S, Stainless Steel Planchets, Type PL-4, Aluminum	Prices per per per per per per per per	90 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 43.50/M 6.00 54.00 50.00/M 9.00 81.00 73.00/M 12.00 108.00 98.00/M 5.00
 199 		Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper Planchets, Type PL-3C, Copper Planchets, Type PL-3C, Stainless Steel Planchets, Type PL-3S, Stainless Steel Planchets, Type PL-4, Aluminum	 Prices per per per per per per per per per per	90 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 43.50/M 6.00 54.00 50.00/M 9.00 81.00 73.00/M 12.00 108.00 98.00/M 5.00 45.00 45.00
 199 	DH 3	Lead Brick, extruded shapes Planchets, Type PL-1, Sheet steel, tinned Planchets, Type PL-2, Aluminum Planchets, Type PL-3, Aluminum Planchets, Type PL-2C, Copper Planchets, Type PL-3C, Copper Planchets, Type PL-3C, Copper Planchets, Type PL-2S, Stainless Steel Planchets, Type PL-3S, Stainless Steel Planchets, Type PL-4, Aluminum Planchets, Type PL-4, Aluminum	Prices per per per per per per per per per per	90 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5 5	application 3.25 30.00 27.00/M 3.25 30.00 27.00/M 3.75 34.00 31.00/M 5.25 47.27 43.50/M 6.00 54.00 50.00/M 9.00 81.00 73.00/M 12.00 108.00 98.00/M 5.00 40.50/M

Bulletin No.	Model No.	Description of Product	Price F.O.B. Burbank
189	PH-2	Plastic Holder for 2" diameter filter paper and wipe samples used with MST-200 series per	C 30.00
189	PLH-30	Plastic Holder for 3" diameter planchets used with MST-300 series per	C 36.00
189	PLH-20	Plastic Holder for 2" diameter planchets used with MST-200 series per	C 30.00
189	PL-5	Aluminum Planchet 2" diameter x 3/16" deep used with MST-200 series	C 8.00
189	PL-6	Aluminum Planchet 3'' diameter x 3/16'' deep used with MST-300 series	C 11.00

CABINETS

C-875	8¾" height "Easy-to-Service" console cabinet for standard 19" width panel	30.00
C-1050	$10\frac{1}{2}$ " height console cabinet for standard 19" width panels	35.00
C-1225	121/4" height console cabinet for standard 19" width panels	37.50
C-1400	14" height console cabinet for standard 19" width panels	40.00
C-1750	$17\frac{1}{2}$ " height console cabinet for standard 19" width panels	45.00
C-2100	21" height console cabinet for standard 19" width panels	50.00
C-2275	223/4" height console cabinet for standard 19" width panels	55.00
C-2625	261/4" height console cabinet for standard 19" width panels	55.00
C-3150	31 V_2 " height console cabinet for standard 19" width panels	60.00
C-4375	433/4'' height console cabinet for standard 19'' width panels	70.00

IN GENERAL

Unit prices are given in this price list unless otherwise noted.

All prices quoted herein are f.o.b. Burbank, California. Payment terms, net 30 days upon approved credit.

No extra charge is made for packing for domestic shipments. Export charges will be quoted upon application. We reserve the right to alter specifications at any time without incurring the obligation of incorporating new features in previously manufactured equipment. Prices are subject to change without notice. Quotations remain firm for 30 days.

Shipping instruments: Please include shipping instructions when ordering; in the event that shipping instructions are not given, we will use our best judgment in the matter.

SERVICE CHARGES

Service charges for repairs are usually billed on an hourly basis, however, if desired, an estimate can be given before the work is undertaken. The unit to be repaired should be shipped prepaid to us.

WARRANTY

Technical Associates warrants instruments and equipment (except tubes, fuses, batteries and crystals), manufactured by them to be free from defects in workmanship or materials under normal use for a period of one year from the date of shipment from the factory to the buyer. Tubes, fuses, batteries and crystals are subject to the guarantee established by the manufacturer of them, however, Technical Associates will assist the customer to obtain full benefits of these guarantees.

If, within the one year warranty period, any Technical Associates instrumentation or equipment requires service as a result of a defect, the buyer may return it to the factory of Technical Associates at Burbank, California or to a service station designated by Technical Associates, transportation charges prepaid, for service at no change under the warranty. The buyer is urged to communicate with Technical Associates when warranty service is required, stating the nature of the difficulty and giving model and serial number of instrument. It may be possible to diagnose the trouble and send a replacement part or assembly, thereby avoiding the expense of shipment.

Technical Associates will return the instrument to the buyer, transportation charges prepaid, after repairs or replacement under warranty are completed. The liability to Technical Associates under this warranty is limited to the cost of replacements of defective parts upon prompt notification of such defect.



TECHNICAL ASSOCIATES INSTRUMENTATION FOR NUCLEAR RESEARCH 140 WEST PROVIDENCIA AVENUE • BURBANK • CALIFORNIA TELEPHONES: Victoria 9-5838 • THornwall 8-6649

Dear Sir:

Technical Associates is pleased to announce the addition of several new products to our growing catalog.

In keeping with our practice to up-date catalogs in the field we are enclosing Bulletins No. 182, 183, 200, 205, 206 and 207 describing our new products and a current price list.

For the Counting Laboratory, we offer a new family of instruments. The Series SD-6 Scintillation Detectors, Gas Flow Counter Model GD-6, the Universal Lead Shield Model LS-66 and Decade Scaler Model DS-9, all have been engineered to meet the requirements of most sample counting applications.

Technical Associates' Model LS-6 Lead Shield, which has been the standard of the nuclear industry, was used as a basis for this unique, interchangeable package for sample measurement. Any Model LS-6 can be easily converted to the Universal Lead Shield Model LS-66 by adding the top section designated as LS-60.

For Monitoring, we offer the new Frisker-Monitor Model FM-1, a completely transistorized instrument with interchangeable probes for use when monitoring personnel or equipment for Alpha or Beta/Gamma activity, or for use as a Beta/Gamma laboratory monitor. Also offered is the new Super Sensitive Gamma Scintillation Survey Meter Model FS-11, which has a full scale range as low as .01 MR/HR over background.

Inquiries for further information or additional catalogs will be handled promptly.

Sincerely,

TECHNICAL ASSOCIATES

John E. Borg Sales Manager

Improved Model PPM-8 PORTAL MONITOR (Bulletin No./168) with eight P.S. individually adjustable channels is available from stock. We offer this instrument for \$1650.00 complete.

JEB: rt Encls.



UNIVERSAL LEAD SHIELD MODEL LS-66

APPLICATION. The Model LS-66 Universal Shield provides unequalled versatility. Its unique design and 2-section construction enable the researcher to adapt the shield to any type of counting application. This new, all-purpose model reflects Technical Associates' 22 years of experience in manufacturing lead shields for the nuclear industry.

The Universal Shield is provided with mountings and sample tray holders for Geiger tubes, window or windowless flow counters, as well as scintillation detectors with well or solid crystals. Inside dimensions accommodate well scintillation assemblies with crystals up to 3''x3'' in size. All detector assemblies can be installed without dismantling the shield.

The bottom section of the Universal Shield is identical to the T/A Model LS-6, the most popular and widely used shield in the industry. Addition of the top section to the Model LS-6 easily converts it to the Universal Model LS-66.

a



MODEL SD-6 SCINTILLATION DETECTOR with SOLID CRYSTAL

may be used in the Model LS-66 Universal Shield for counting samples in planchets. Mountings in bottom section of shield hold sample tray holders.

DESCRIPTION. The Model LS-66 is constructed of low porosity virgin lead meeting Federal Specification QQ-L-171. The upper plane of the door in the lower section is cut at a 30° angle and provides an efficient wedge-type seal for blocking stray radiation and light. The top section is designated as Model LS-60 and provides the well counting capability. The lid is counterbalanced for smooth easy operation and contains removable plugs to accept test tubes and a 50 ML beaker. A lead ring and splash pan are provided for installation around the well crystal, thus adding 2" to the shielding. The unit is lined with aluminum to minimize backscatter.



MODEL SD-6W SCINTILLATION DETECTOR with WELL CRYSTAL may be inserted in Model LS-66 Shield for counting liquid samples. A lead ring to further reduce background is provided for use around well crystal.



TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA

pecifications of MODEL LS-66 UNIVERSAL SHIELD

INSIDE DIMENSIONS: 5%" diameter x 15" high.
DOOR OPENING AREA: 4%" high (mean) x 5" wide.
GEOMETRIC REPRODUCIBILITY: 0.1%
SHIELDING: Well Counter — lead 3⁵/32", brass ¼". Solid Crystal — lead 1%", brass ¼".
WEIGHT: 355 lbs. Shipping weight, 375 lbs.
FINISH: Grey enamel. Chrome plated handles, hinges, and plugs.



Į 🕡

may be installed in the Universal Shield for window or windowless counting of alpha or beta emitting samples. A gas flow meter can be mounted on the shield.







Addition of the Model LS-60 Top Section to the LS-6 Lead Shield easily converts it to the Model LS-66 Universal Shield.





Model LS-66 Universal Shield shown with the Orion Spectrometry System, one of the many versatile T/A systems available for spectrometry applications.

ECHNICAL ASSOCIATES 140 West Providencia Avenue • Burbank, California



DESCRIPTION. The scintillation detector assembly contains a 2" photomultiplier tube with a mu metal shield and a transistorized dual preamplifier built as an integral part of the probe. Short leads from the photomultiplier tube assure high counting efficiency, even when the signals are very weak. A switch is provided to select the output polarity. The preamplifier for the negative output is extremely linear and has a gain of 10. The negative output of the preamplifier gives the SD-6 Series Detectors the versatility of operating into standard counting instruments with ¼ volt negative input. The positive output for spectrometry utilizes a double emitter follower circuit with unity gain. High voltage and signal cables are electrostatically shielded.

The crystals are hermetically sealed and optically coupled to the photomultipler tube by a rigid crystal mount. The crystal mount and probe casing are light tight and ruggedly constructed of brass. The assembly is hard chrome plated for durability and easy decontamination. **MODEL SD-6-1** employs a 1''x1'' Nal crystal with a 2'' photomultiplier tube for gamma ray detection.

MODEL SD-6-1.5 is identical to Model SD-6-1 except that it uses a $1^{1}\!\!/_{2}''$ dia. x 1'' thick Nal crystal.

MODEL SD-6-2 employs a 2"x2" Nal crystal with a 2" photomultiplier tube. This larger crystal significantly increases the counting efficiencies and is recommended for spectrometry applications. **MODEL SD-6-2W** uses a well type Nal crystal (134" dia v 2" thick

MODEL SD-6-2W uses a well type Nal crystal $(1\frac{3}{4})^{"}$ dia, x 2" thick with a well 21/32" dia, x $1\frac{3}{4}$ " deep). The well crystal is a highly efficient detector for spectrometry and measurement of low level gamma-emitting samples in liquid or powder form.

MODEL SD-6-B1.5 has a sealed light-tight $1\frac{1}{2}n''$ dia. Beta Stilbene crystal $\frac{1}{2}mm$ thick mounted behind an end window of polished aluminum foil .0005" thick. This detector permits the counting of beta particles with energy greater than 40KV and is comparatively insensitive to gamma rays.

MODEL SD-6-A1.5 has a sealed light-tight $1\frac{1}{2}$ " dia. Alpha Phosphor mounted behind an aluminum foil and protected by a perforated screen. This detector is not sensitive to gamma rays and yields a pulse height large enough to eliminate background count due to tube noise.

	SD-6-1 SD-6-1.5 SD-6-2 SD-6-2W SD-6-B1.5 SD-6-A1.5
CRYSTAL:	$\begin{array}{cccccccccccccccccccccccccccccccccccc$
RESOLUTION:	Does not apply Pos. output only. Does not apply Better than 10% for Cesium 137.
PHOTOMULTIPLIER TUBE:	RCA 6655-A (Same for all Series SD-6 Models)
PRE-AMPLIFIER:	All transistor dual polarity: Negative output has a gain of 10. Positive output is double emitter follower with unity gain.
PLATEAU LENGTH:	200 volts minimum (negative output only)
PLATEAU SLOPE:	Less than 5% per 100 volts using Cobalt 60 (negative output only)
HIGH VOLTAGE REQUIREMENTS:	Stable positive polarity supply continuously variable from 500 to 1500 volts.
LOW VOLTAGE REQUIREMENTS:	Probes are designed for use with instruments having 250 volt low voltage supply. Adaptable for instruments having 7 volt supply.
CONNECTORS:	H.V. and Sig. Output UG-932/U Pre-Amp. Power Cannon XLR-3-12
PROBE SIZE:	3½" dia. x 9¾" 3½" dia. x 10½" 3½" dia. x 9"
WEIGHT:	6½ pounds (Same for all Series SD-6 Models)

AN IMPORTANT RESEARCH TOOL

To secure a permanent record of the spectrum of the sample being studied, a Series SD-6 Probe may be installed in the LS-66 Lead

Shield, and the signal fed into the Arcturus Spectrometry System. (See Bulletin #171)





TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA

GAS FLOW COUNTER, MODEL GD-6

Exceptional plateau and high efficiency. Proportional or Geiger operation. Built-in transistorized preamplifier. Ultra-thin nylon window. Fits inside standard LS-6 Shield. 2" diameter sample capacity.

APPLICATION. The Model GD-6 Flow Counter is engineered to produce high counting efficiencies where extremely low levels of alpha or beta radioactivity are to be measured. It is ideal for applications in biological, physical, and chemical research, or in environmental sampling for health physics and public health services. The instrument includes a sample tray holder which accepts samples up to 2" in diameter and provides 5 reproducible geometries. The Model GD-6 is designed to fit inside the standard Model LS-6 Shield for background reduction during beta counting.

DESCRIPTION. The combination of the full 2 pi hemisphere, the precision collector wire, and the ultra-thin nylon window is the key to the high counting efficiency of the Model GD-6 Flow Counter. The instrument provides counting efficiencies which are 85% of internal counters. No purging time is required between counting cycles, thus allowing faster sample counting and reducing gas consumption.

The built-in transistorized preamplifier is coupled directly to the collector wire with a 1" lead to eliminate loss of low amplitude signals. The preamplifier circuit is designed to operate with ± 250 Volts D.C., which is supplied by most scalers and ratemeters. A switch is provided to bypass the preamplifier when the Model GD-6 is used as a Geiger counter.

SPECIFICATIONS:

SAMPLE PANS: 2" diameter x $1\!\!\!/_4$ " deep, 2" diameter x $1\!\!/_4$ " deep, and 2" diameter x flat (aluminum or stainless).

HEMISPHERE: 2³/₈" diameter stainless steel. COLLECTOR WIRE: 1 mil stainless steel.

WINDOW: $2\frac{1}{4}$ " diameter, 90 to 125 µg/cm² thick. (One extra window is supplied.)

COUNT YIELD: Alpha 38%, Beta 43% with C¹⁴. PLATEAU: Length, 300V; slope, 1% per 100V. BACKGROUND: Alpha less than 0.1 cpm, Beta

less than 30 cpm in LS-6 Shield. COUNTING GAS: P-10 proportional, Q gas Geiger.

(Incl. flow meter & low pressure regulator.)



GAS FLOW COUNTING SYSTEM

OPERATING VOLTAGE: Proportional – Alpha 1200V, Beta 1800V, Geiger 1300V.

PREAMPLIFIER POWER REQUIREMENT:

PREAMPLIFIER OUTPUT: 2 volts nega-

tive. (Triggers standard instruments

CABLE CONNECTORS: Signal and High

Voitage-UG-932/U Preamplifier Power

+250 Volts D.C. @ 10 ma.

having .25V negative input.)

SIZE: 51/2" diameter x 8" high.

- Cannon XLR-3-12C

FINISH: Chrome.

WEIGHT: 5 pounds.

offers extremely high counting efficiency. The system combines the Model GD-6 Flow Counter in the Model LS-6 Shield and the Model DS-9 Decade Scaler with electromechanical timer.





This is an actual plateau from the Model GD-6 Flow Counter using a C¹⁴ sample.



DECADE SCALER, MODEL DS-9

- 1 microsecond resolving time.
- Total count capacity to 100,000,000.
- 5000 volt power supply.
- Manual, remote, and preset count operation.
- Push-button controls.



APPLICATION. The Model DS-9 Decade Scaler is a precision counting instrument designed for use with any type of scintillation detector, Geiger tube, or proportional flow counter. This versatile scaler, combined with the appropriate detector and shield, provides an ideal radioisotope counting system for applications in research, health physics, and medicine. The instrument may be included as a component in analytical spectrometry systems as well. (See Bulletin No. 171.) **DESCRIPTION.** The Model DS-9 Decade Scaler features 1 microsecond resolving time and provides count storage capacity to 100,000,000 with 4 direct reading decades and a four digit electro-mechanical register. The instrument has manual and preset count modes of operation and is adaptable to predetermined or elapsed timing with an electronic or electro-mechanical timer. Input sensitivity is variable from 0.2 to 2.0 negative volts. A stable regulated high voltage power supply features dual ranges 500 to 2500 and 1000 to 5000 volts, together with coarse and fine adjustments. Push-button controls are provided on the front panel for easy operation of "Start," "Stop," and "Reset."



RADIOISOTOPE COUNTING SYSTEM utilizes the Model DS-9 Decade Scaler and Model PT-6 Electronic Timer. System is shown with LS-66 Universal Shield which accommodates any type of detector.

SPECIFICATIONS

INPUT SENSITIVITY: 0.2 to 2.0 negative volts. (Set at factory 0.25V) INPUT SIGNAL POLARITY: Negative. **RESOLVING TIME:** 1 microsecond to pulse pair. MAXIMUM COUNTING RATE: 250 KC. COUNT CAPACITY: 100,000,000. PRESET COUNT RANGES: 10, 40, 100, 400, 1000, 4000, 10,000, 40,000. PUSH-BUTTON CONTROLS: Start, stop, reset. HIGH VOLTAGE POWER SUPPLY: Range: 500 to 2500 volts and 1000 to 5000 volts continuously variable with "coarse" and "fine" controls. Polarity: Positive output (negative ground). Stability: Less than .01% drift per day after warm-up. Regulation: Less than .01% change in high voltage for a 1 volt line change between 105 to 125 volts A.C. Ripple: Less than 0.005% of output voltage. Automatic Time Delay: Prevents high voltage from being turned on until regulator tubes are in operation, thereby eliminating damage to detector. High Voltage Meter: Large 5" meter face calibrated to within $\pm 2\%$ full scale. POWER REQUIREMENTS: 105 to 125 Volts A.C., 250 watts. FINISH: Grey baked enamel. SIZE: 8¾" high x 19" long x 15" deep. WEIGHT: 40 lbs. Shipping weight: 46 lbs.

TA

ELECTRONIC TIMER, MODEL PT-6



The T/A Model PT-6 Electronic Timer is designed for use with counting systems where short-time error must be minimized and longtime counting is required. The instrument incorporates fast start and stop circuitry, as well as push-button reset. This timer may be used with any electronically gated scalers or any scalers that stop and start with mechanical contact closures.

SPECIFICATIONS

TIME INDICATION: By 4 glow transfer tubes.

PRESET TIME: 10, 20, 50, 100, 200, 500, 1000, 5000, 10,000 seconds. TIME RESOLUTION: To the nearest 1 second on times longer than 1000 seconds and 0.1 second on times less than 1000 seconds.

- (The X1.X10 switch changes the frequency of the timer drive by a factor of 10, thus moving the decimal point one place to the left or to the right as desired.)
- ELAPSED TIME: 0 to 9999 seconds.
- ACCURACY: Limited entirely by line frequency variations.

- **OUTPUT:** Contact closure on accumulation of preset time as well as electronic gate output. A gate output is also provided on start and stop.
- CONTROLS: Power, stop-start, reset, preset time, and elapsed time selector switches.
- **DIMENSIONS:** $3\frac{1}{2}'' \times 19'' \times 9''$ deep for rack or cabinet mounting. A line cord is included with the instrument and cabinets are available at slight extra cost.

POWER REQUIREMENTS: 115 volts, single phase, 60 cycle, 75 watts. CONNECTORS: Output-AN-3102A-2029S • AC Power-Amphenol 160-3

PREAMPLIFIER, MODEL PA-6B. This instrument is designed to have low noise and good non-overloading characteristics for use with proportional flow counters and scintillation detectors where the signal output may be too low in amplitude to feed directly into a scaler, ratemeter, or other counting instrument. The Model PA-6B Preamplifier has an adjustable input sensitivity from 1 millivolt to 1000 millivolts and a gain of 250, thus allowing it to be used with any standard scaler or ratemeter. Cables are provided for power input and signal-high voltage input. The instrument is offered mounted on a standard 19" panel or self-contained in a 4" x 5" x 8" cabinet. The cable connectors and receptacles are as follows: power input (cable), Cannon XLR-3-12; high voltage input and signal input(cable), UG-932 A/U; and detector input(receptacle), UG-931/U.





GAS FLOW COUNTING SYSTEM

offers extremely high counting efficiency. The system combines the Model GD-6 Flow Counter in the Model LS-6 Lead Shield and the Model DS-9 Decade Scaler with Electro-Mechanical Timer.



FRISKER-MONITOR, MODEL FM-1





MODEL FM-1A FRISKER MONITOR is shown mounted with chart recorder on a standard 19" panel.

- RANGES: 0 to 100, 300, 1000, 3000, 10,000, 30,000.
- TIME CONSTANTS: Fast, medium, slow (5 sec., 25 sec., 50 sec.)

ACCURACY: Better than 5% on all ranges.

INPUT SENSITIVITY: Negative, variable from .2 to 2 volts. (Factory set at .25 volts)

- **RESOLUTION:** Less than 10 microseconds.
- VISUAL READ-OUT: 4 1/2" dual purpose meter for count rate read-out and high voltage.
- AURAL READ-OUT: Audio amplifier and 3" speaker with volume control on panel.
- FAST METER RETURN: Button on front panel.
- **RECORDER OUTPUTS:** Connections to drive 10 MV and 1 MA strip chart recorder.
- ZERO STABILITY: Grounded output, no drift.
- TEST SIGNAL: 3600 cpm calibration check provided on front panel.
- HIGH VOLTAGE: Regulated 500 to 1600 volt supply with coarse and fine
- adjustments on front panel.
- POWER REQUIREMENTS: 115 volts, 60 cycle, 6 watts.
- CONNECTORS: UG-931/U Signal & H.V. jack mates with UG-932/U plug. AN 3102A-14S-5S Alarm receptacle mates with AN 3106A-14S-5P. Switchcraft 13B recorder jack.

 - Amphenol 160-2 auxiliary power receptacle. Cinch-Jones Type 3-141Y barrier terminal strip for 10 MV recorder output and external test signal input.
- SIZE: 81/4" high x 17" wide x 10" deep.
- WEIGHT: 20 pounds. Shipping weight: 28 pounds



Chart recorder outputs to provide permanent records.

APPLICATION. The T/A Model FM-1 Frisker Monitor is a multi-purpose ratemeter designed primarily to detect radiation contamination on the clothes, hair, hands, and shoes of personnel moving about from active areas to clean areas. The basic unit is also used in checking laboratory equipment during decontamination processes.

With the inclusion of the Model ALS-2 Plug-In Alarm, the Frisker Monitor can be used as a continuous area monitor for laboratories and source rooms. The alarm will alert personnel when the radioactivity reaches the preset warning level.

In addition to personnel and laboratory monitoring, this flexible instrument is ideal for educational demonstrations and medical isotope studies.

DESCRIPTION. The Model FM-1 Frisker Monitor features completely transistorized circuits on plug-in boards to assure long dependable service. Its 6 linear ranges and 3 time constants were selected to offer efficient capability for counting alpha, as well as beta/gamma, with a single instrument using interchangeable probe assemblies. The count rate is audible through a built-in 31/2" loudspeaker and easily visible on a large 41/2" meter face, which displays the full scale ranges. Included on the front panel are a fast zero control, calibration test signal, and high voltage controls with coarse and fine adjustments.



PLUG-IN ALARM To increase the versatility of the Frisker Monitor, a Model ALS-2 Plug-In Alarm is offered on an optional basis. This transistorized accessory is designed to alarm at a preset level and includes a red light indicator and reset button. The alarm settings are variable from 10% to 100% of full count range. The meter continues to read the exact count rate.



REMOTE ALARM Also available for use with the Frisker Monitor is the Model RAL-3 Remote Alarm Indicator. This accessory contains a clearly visible flashing red light and buzzer type audio signal which are activated when the radioactivity reaches the warning level set on the alarm. A front panel switch is provided to switch off the audio signal. The alarm is furnished with 25 feet of connecting cable, but can be used at greater distances.





PROBE ASSEMBLIES and AREA MONITORING STATION



MODEL P-AS-2 ALPHA PROBE

MODEL P-AS-2 PROBE is ideal for surveying contaminated areas or for "frisking" personnel, clothing, and equipment. It can be used with ratemeters, scalers, portable survey meters, or monitoring stations. The unit consists of a $4\frac{1}{2}$ " $x4\frac{1}{2}$ " alpha phosphor covered with 1.0 mg/cm² aluminized mylar and optically coupled to a 2" photomultiplier tube. The probe is specially designed to assure high efficiency over the entire window surface. Power requirements are +900 to +1000 volts at approximately 25 micro-amps.



BETA/GAMMA PROBE

is primarily designed for use as a "frisker" with a ratemeter, scaler, portable survey meter, or as an accessory with monitoring stations. The unit includes a T/A T-1090 GM Tube (with organic quenching agent) and features a 180° rotating beta shield. The probe is ruggedly constructed of brass with chrome plating.



MODEL P-8 BETA/GAMMA PROBE

This end window probe includes the T/A Model T-1140 GM Tube and is ideal for detecting soft beta radiation. It can be used with counting instruments such as ratemeters, scalers, and portable survey meters. The unit is chrome-plated and features a snap action end cap which protects window and acts as a beta absorber when surveying for gamma only.



MODEL P-10

This probe is similar in application to the Model P-7 Probe, except that it utilizes a T/A T-1110 GM Tube (stainless steel, with halogen quenching agent). The unit features a 180° rotating beta shield and is constructed of durable chrome-plated brass.

NOTE: All probes are supplied with 8' cable and a UG-932/U connector so that they can be connected directly to T/A equipment. Other connectors will be supplied on request.





Super Sensitive GAMMA SCINTILLATION SURVEY METER



- Full-scale range as low as .01 mr/hr over background.
- Large clear meter for easy reading of weak emitters.
- Built-in battery checker indicates remaining "life"in batteries.
- Light, gun-type probe can easily be placed in the best position for maximum sensitivity.

APPLICATION.

The T/A Model FS-11 Scintillation Survey Meter is an extremely sensitive instrument with high statistical accuracy. It has many features that make it suitable for use in health physics, research, medicine, and ore prospecting. The design includes a top quality Sodium Iodide (Thalium activated) crystal to assure optimum efficiency for detecting and measuring very small amounts of gamma ray emitting radioactive material.

DESCRIPTION. The Model FS-11 Scintillation Survey Meter consists of a compact, lightweight, gun-type probe with a hermetically sealed 1"x1" Sodium Iodide crystal and photomultiplier tube. The instrument has five sensitivity ranges which are displayed on a large clear meter face. Each range is individually calibrated and holds calibrations for hundreds of hours. There is no zero drift because of the inherent stability of the circuit. Meter response is adjusted by a 3-position time constant control. A single control knob provides both range selection and battery check. The electronic circuit, meter, controls, and batteries are contained in a ruggedly constructed weatherproof case. Each instrument is supplied with a calibration source.

SPECIFICATIONS:

- 5 FULL-SCALE RANGES: .01 over back-ground, .1, .25, 1.0, 2.5 MR/HR abso-lute.
- Iute. **TIME CONSTANTS:** 2, 5, 15 seconds. **CONTROLS:** Range Selector, "A" and "B" Battery Tester, "B" Adjust, Fast Meter Return, 01 MR/HR Range Background Adjust, Time Constant. **BRORE:** Observe olated steel convenient
- PROBE: Chrome plated steel, convenient
- pistol grip, sealed construction. CRYSTAL: 1"x1" Nal hermetically sealed, shock-mounted.

PHOTOMULTIPLIER TUBE: Selected RCA 6199, shock-mounted.

- BATTERIES: (2) Eveready #467 67½ volt, (4) Eveready #D-99 1½ volt flash-light cells, (1) Eveready #412 22½ volt, (4) Eveready #912 penlight cells.
- CASE DIMENSIONS: 101/4 " long x 43/4" wide x 5⁷/₈" high.

FINISH: Hammertone grey. WEIGHT 10 lbs. Shipping weight: 13 lbs.





HIGH RANGE SURVEY METER MODELS CP-TP-1A and CP-TP-1B

- Full-scale ranges to 5000 R/HR. CP-TP-1A: 0 to .5, 5, or 50 R/HR. CP-TP-1B: 0 to 50, 500, or 5000 R/HR.
- Interchangeable ionization chambers.
- 3-position rotating chamber.

APPLICATION. The T/A Model CP-TP is the only portable, battery operated survey meter especially designed to enable the health physicist to measure high intensity beta and gamma radiation fields. The instrument features a 40" extension rod with swivel bracket for rotating the chamber, thus allowing the user to stand behind a shield or out of the direct line of a beam.

DESCRIPTION, Technical Associates' CP-TP High Range Survey Meter is available in two models, each having three linear ranges: CP-TP-1A provides 0 to .5, 5, and 50 R/HR (beta and gamma); CP-TP-1B provides 0 to 50, 500, and 5000 R/HR (gamma only).

The basic instrument is the same for both models, the only difference being the type of ionization chamber selected. The lower and higher ranges are determined by the size and configuration of the ionization chambers and their electrometer circuits. Each chamber contains its own electrometer circuit. The chambers are interchangeable on the instrument. The CP-TP-1A chamber has a removable beta absorber to allow measurement of beta radiation through an acetate end-window.

The instrument case is chrome-plated aluminum for easy decontamination. It contains the range selector switch, zeroing and calibration controls, batteries, and a high quality wide-face meter. The unit has a pistol grip handle for easy manipulation. Each model comes complete with a 40" extension rod, which is detachable from the instrument case. Either ionization chamber can be connected directly to the case, if desired. The swivel bracket at the end of the extension rod has two connectors for use in attaching the chamber for the desired ranges. One connector is at the forward end and the other is on the side. The swivel bracket rotates, thus moving the chamber connected to the side of bracket into 3 separate positions: right, left, and downward.

SPECIFICATIONS:

RANGES: Model CP-TP-1A: 0.5, 5, 50 R/HR. Model CP-TP-1B: 50, 500, 5000 R/HR.

- CONTROLS: Combined "Off-On" and range selector switch. Individual calibration controls for each range. Zero control can be set in radiation field.
- CIRCUIT: Reliable single tube electrometer circuit with all high resistance points insulated with Teflon or Kel-F to insure minimum leakage.
- METER: Wide-face 31/2" meter with dual scale calibrated in R/HR (black scale for 1A chamber, red scale for 1B chamber). Divided into 50 divisions. Mounted to provide excellent visibility.
- ZERO DRIFT: Negligible after 15 minute warm-up. **IONIZATION CHAMBERS:**

Model CP-TP-1A: Cylindrical bakelite internally coated with aquadag. Volume approx. 42 cubic inches. Chamber has a removable beta absorber (432 mg/cm²) to allow beta measurements through an acetate end-window (6 mg/cm²). Model CP-TP-1B: Cylindrical aluminum internally coated with aquadag, Volume approx. 3 cubic inches. Sensitive to gamma only.

- NOTE: Both chambers include their own electrometer circuit for specific ranges.
- EXTENSION: Made of polished aluminum tubing 40" long, with connector on one end for mounting on instrument case and swivel connector at the other end for mounting chamber.
- INSTRUMENT CASE: Chrome-plated aluminum with all markings engraved for easy decontamination.
- BATTERIES: 4 Eveready No. 412, 22.5 volts. 2 Burgess No. W5BP, 7.5 volts. 2 Eveready "D" Cells, 1.5 volts.
- WEIGHT: 6 lbs.



is shown connected to forward end of bracket. **Dotted lines** indicate the 3 positions to which the chamber can be rotated.



MODEL CP-TP SHOWN WITH **1A CHAMBER**



TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA



JUND SURVEY METER, Models SRJ-7 and HRJ-7

Incorporating design features which permit effective operation under high humidity conditions

APPLICATION. The Juno Survey Meter is a portable instrument for measuring the intensity of, and discriminating between alpha, beta, and gamma radiation. It is used to protect personnel from the danger of over-exposure to radiation from radioactive materials or X-rays. While primarily intended for inspection of flat surfaces, the instrument is suitable for most uses where a high degree of accuracy is desired.

The T/A Juno is available in two models: SRJ-7 (standard range) for all normal applications; and HRJ-7 (high range) for use where exceptionally high intensity radiation is likely to be encountered. Both models meet A.E.C. specifications.

Models SRJ-7 and HRJ-7 are improved versions of the original Hanford instrument. The high impedance circuit switch box includes a desiccant cartridge and is sealed with gaskets. These design improvements insure high efficiency performance under adverse humidity conditions.

DESCRIPTION: The instrument comprises an ionization chamber, an electrometer circuit, absorption filters for the rejection of either alpha or beta particles, suitable batteries mounted in a removable power pack, and an indicating meter. The unit is battery operated and is self-contained. To easily identify the High Range Juno, its knobs and meter dials are finished in brilliant red.

ALPHA WINDOW

The ionization chamber has a volume of approximately 27 cubic inches. All surfaces within the chamber are coated with aquadag. The chamber is covered by a screen of .0003'' (approximately 0.45 mg/cm^2) rubber hydrochloride film. The alpha screen is within 7/16'' of any flat surface on which the instrument may be placed, and is easily replaced by simply removing the bottom plate and two retaining screws.

Two absorbers are provided to reject either alpha or beta radiation. These are readily moved in and out of position by means of sliding tabs fitted in rails which form part of the handle. The tab marked "G," with a square end, operates the absorber which rejects alpha and beta, thus permitting a reading of gamma only. The tab marked "B," with a rounded end, operates the absorber which rejects alpha, permitting a reading of beta and gamma. The total of all three types of radiation is read, when both tabs are in "open" position.

The high quality microammeter, which is calibrated in milliroentgens per hour for gamma radiation, has a large easy-to-read face and is mounted in position to permit excellent visibility. Battery life is approximately 800 hours in normal intermittent use. An easily removable battery pack, with simple positive contacts, assures trouble-free operation over long periods of service.



Ą

TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA

Models SRJ-7 and HRJ-7

have similar circuits, the principal difference being the value of grid resistances used.

Radiation incident upon the ionization chamber produces a minute current which flows through a very high resistance in the grid circuit of the electrometer tube. The voltage thus produced at the grid causes a corresponding change in plate current which is indicated by the panel meter. A bucking current is provided through the meter in order that the no-signal plate current of the electrometer may be balanced out and readings of radiation intensity may start from the meter zero reading. Sensitivity is varied by switching appropriate values of grid resistance in the electrometer circuit.

The instrument is calibrated by adjusting a resistance in series with the meter. An individual adjustment is provided for each range. Zero setting is accomplished by means of a rheostat in the filament circuit of the electrometer tube. By varying the filament voltage, the plate current may be varied and thus adjusted to a value equal to the bucking current flowing through the meter. All high resistance points in the circuit are insulated with Teflon to insure minimum leakage.

SENSITIVITY. Both models are calibrated in three separate full-scale ranges in easily read increments of the meter scale, covering the total range of which the instrument is capable.

Ranges are based on radium gamma radiation intensity. Accuracy of calibration is such that indications on the meter will not be more than 5% lower nor more than 10% higher than the radiation intensity to which the chamber is exposed. For use in abnormal environmental conditions, air temperature and density correction data and curves are included with each instrument. Sensitivity dependency upon battery aging is limited to a 10% variation while the unit can be zeroed by means of the panel zero control.

Hustration below the Nuclear Red

ration, monitorin

SPECIFICATIONS:

IONIZATION CHAMBER:

Volume: 27 cubic inches.
Window Opening: 3" x 4%".
Alpha Window: 0.3 mil (0.45 mg/cm²) rubber hydrochloride.
Alpha Absorber: 0.01" cellulose acetate sheet.

Beta Absorber: 0.102" aluminum.

BATTERIES:

4 Eveready No. 412 22¹/₂ Volt "B" Batteries. 2 Eveready No. E12 1.35 Volt "A" Batteries. 2 Mallory No. TR-115 6.5 Volt "B" Batteries.

TUBE:

1 Sub-Miniature Electrometer Type CK 5886.

RANGES:

Model SRJ-7 Juno (standard range): 50, 500, 5000 MR/HR full-scale. (Improved A.E.C. Model SIC-17B). Model HRJ-7 Juno (high range): 250, 2500, 25,000 MR/HR full-scale. (Improved A.E.C. Model SIC-17D).

TIME CONSTANTS:

50 MR/HR – 18 Seconds 500 MR/HR – 4 Seconds 5000 MR/HR – 2.5 Seconds

OPERATING TEMPERATURE RANGE: 35°F to 135°F

CASE:

HANDLE:

Aluminum, especially cast low-porosity, smoothly polished.

WEIGHT: Net 6 lbs. 11 oz. Shipping: 10 lbs.



The Juno can be used for X-ray detection and measurement by reference to the above curves.



ry of J



Hard-chrome plated aluminum, with engraved markings on top of case; gasketed, dust and moistureresistant. Dimensions: 91/2" x 53/4" x 4".



CUTIE PIE

MODELS CP-3 and CP-3A

PORTABLE ALPHA, BETA, and GAMMA SURVEY METER

- 3 linear ranges to 5000 mr/hr. CP-3: 0 to 50, 500, or 5000 mr/hr. CP-3A: 0 to 25, 250, or 2500 mr/hr.
- Thin window permits alpha detection.
- Built-in alpha and beta absorbers.
- Selector switch positions permit checking all batteries.
- Zero adjustment in radiation field.
- New battery pack provides over 800-hour operating life.

APPLICATION The T/A Models CP-3 and CP-3A Cutie Pie Survey Meters measure alpha, beta, and gamma radiation with excellent energy independence. These lightweight, portable instruments are improved versions of the original Cutie Pie design, which was developed under the Manhattan Project. They have found wide acceptance by health physicists for surveys at reactor sites and research and industrial laboratories, and are extensively used for determining shielding effectiveness, checking source shipping containers, monitoring areas, as well as for decontaminating and cleanup purposes. They are especially useful for monitoring such inaccessible spots as corners, behind pipes, beams, etc.

DESCRIPTION. The Models CP-3 and CP-3A comprise an ionization chamber, an electrometer circuit, alpha and beta absorbers, battery pack, and an indicating meter. Both units are battery operated and self-contained.

The ionization chamber is made of aluminum and coated internally with aquadag. The end of the chamber is closed with a rubber hydrochloride alpha window (0.45 mg/cm^2) , which is frame-mounted and held in place with clips, thus permitting easy replacement. Alpha and beta absorbers are mounted with hinges on the front of the chamber. Undesired radiations are easily rejected by swinging the proper absorber into place.

Each instrument has three linear ranges: the CP-3 provides 0 to 50, 500, or 5000 mr/hr; the CP-3A provides 0 to 25, 250, 2500 mr/hr. These ranges provide excellent coverage of radiation levels normally encountered. To assure maximum reliability, the range selector switch includes three test positions for checking batteries prior to use. Batteries of proper voltage produce a reading in the green sector on the meter. A "Set" position permits the meter to be adjusted to read zero even in radiation fields. The remaining three switch positions permit the selection of ranges.

The instrument case is made of chrome-plated aluminum, with engraved lettering for easy decontamination. Rubber gaskets are used for protection against high humidity. The case contains the range selector switch, the zeroing control, and a large-face meter calibrated to read milliroentgens per hour, with 50 scale divisions. The meter is mounted in position for excellent visibility. A reliable high impedance electrometer circuit with an improved 800-hour battery complement assures long uninterrupted service. Both Models CP-3 and CP-3A provide front tripod feet for bench or table use.





140 WEST PROVIDENCIA AVENUE .

BURBANK, CALIFORNIA

ASSOCIATES

Cutie Pie is shown testing radioactivity at reactor site.

Printed in U.S.A

Buletin No. 186

CUTIE

TECHNICAL ASSOCIATI

PIE



- RADIATION RANGES: CP-3: 50, 500, and 5000 mr/hr full scale. CP-3A: 25, 250, and 2500 mr/hr full scale.
- **CALIBRATION:** Factory calibrated, using gamma standard calibrated by Nat'l Bureau of Standards. Calibration accuracy \pm 10%. Individual calibration control for each range.
- **CIRCUIT:** Reliable single tube electrometer circuit. All high resistance points are insulated with Teflon or Kel-F to insure minimum leakage.
- METER: High quality 31/2" meter, with 50 scale divisions. Appropriately calibrated to read in milliroentgens per hour for gamma radiation.
- **CONTROLS:** Single control switches meter to battery test points, zero position, and 3 operating ranges. Meter Zero Control is located directly below the meter.

ZERO DRIFT: Negligible after 15 minute warm-up.

TIME CONSTANTS:

- CP-3: Range 50 mr/hr, 6 seconds Range 500 mr/hr, less than 1 second Range 5000 mr/hr, less than 1 second.
- CP-3A: Range 25 mr/hr, 12 seconds Range 250 mr/hr, 2 seconds Range 2500 mr/hr, 1 second.

IONIZATION CHAMBER:

Aluminum Cylinder: 27%" inside diameter, 6%" long. Volume: Approximately 36 cubic inches. Window Opening: 23%" diameter. Alpha Window: Removable, ring-mounted rubber hydrochloride (0.45 mg/cm²). Alpha Absorber: Hinge-held, ring-mounted cellulose acetate (36 mg/cm²). Beta Absorber: Hinge-held, aluminum disc (720 mg/cm²).

CASE: Chrome-plated aluminum, with clearly engraved markings.

BATTERY LIFE: Over 800 operating hours.

BATTERY COMPLEMENT:

4 Eveready No. 412 22¹/₂ Volt "B" Batteries 2 Eveready No. E12 1.35 Volt "A" Batteries 2 Mallory No. TR-115 6.5 Volt "B" Batteries (Battery complement is identical to the Model 7 Juno and CP-4.)

WEIGHT: 4 lbs. 12 oz. net. Shipping Weight: 8 lbs.

CUTIE PIE, MODELS CP-4 and CP-4A PORTABLE BETA and GAMMA SURVEY METER

For those customers who prefer a Cutie Pie with a chamber made of bakelite, instead of aluminum, and who require measurements of beta and gamma only, T/A offers Models CP-4 and CP-4A. These instruments utilize the same electrometer circuitry and include the outstanding features and specifications of Models CP-3 and CP-3A. The Model CP-4 has the same ranges as the Model CP-3 (0 to 50, 500, or 5000 mr/hr); while the Model CP-4A has the same ranges as the Model CP-3A (0 to 25, 250, or 2500 mr/hr).

The bakelite chamber of the Models CP-4 and CP-4A has a beta end-window (6 mg/cm²), permitting detection of low energy beta particles. They have a bakelite beta absorber (432 mg/cm²), in the form of a cap, which is easily placed in position over the end, thus permitting

measurement of gamma only.



TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE 🔗 BURBANK, CALIFORNIA

with

RETA ARSORRER

ALPHA ABSORBE

s locked at th

sides of

ALPHA WINDOW



ASSOCIATE

MODEL CP-4 CUTIE PIE with front tripod feet in position for bench or table use.



LOG-LINEAR COUNT RATEMETER MODEL RM-8



- Linear and log scales.
- 9 overlapping linear ranges.
- Choice of 7 time constants.
- 2% accuracy on all linear ranges.
- Outputs for galvenometer or potentiometer recorders.
- Speaker with "howler" or "clicker."
- Stable 2500 volt power supply.

APPLICATION. The Technical Associates Model RM-8 Log-Linear Count Ratemeter is an extremely precise and versatile radiation counting instrument. It converts randomly spaced pulses from an external detector into an average count rate per minute. The Model RM-8 presents this average on a wide panel meter in a logarithmic scale or in one of 9 different linear ranges. The count rate information may be fed to a chart recorder for permanent record. A front panel mounted speaker provides an aural indication of the count rate by means of a "clicker" or "howler."

The Model RM-8 is recommended for applications in laboratory research, reactor studies, medical radioisotope diagnostic studies, analytical procedures, etc., that require accurate measurement of counting rates, and is particularly suitable for studies concerned with changes in counting rates. The instrument is designed for use with Geiger tubes, scintillation detectors, proportional counters, and gas chromatograph scanners.

DESCRIPTION. The Model RM-8 Log-Linear Count Ratemeter features a 5-decade logarithmic scale which permits counts to be made over a full range of 10 to 1,000,000 cpm. The log scale makes it easy to chart rapid changes in count rates or to determine quickly the appropriate linear scale. 9 separate linear scales are avail-



able, ranging from 0-100 to 0-1,000,000 cpm. This wide choice of ranges permits the selection of an optimum range for a specific application.

The instrument provides six time constants, with one additional position for a special time constant, if desired. The availability of extremely short time constants on the upper ranges allows rapid response with no significant increase in probable error. Other front panel controls provide fast meter return, a meter zero control, a 60-cycle test, and a speaker for aural indication of the count rate.

The Model RM-8 has an exceptionally stable 600 to 2500 volt power supply with its own meter and with separate coarse and fine voltage adjustment controls. An outstanding feature of the instrument is its monostable circuit utilizing the EFP-60 secondary-emission tube. This circuit minimizes drift and assures excellent accuracy. The Model RM-8 will drive standard potentiometer or galvenometer chart recorders from either the logarithmic or linear ranges.



TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA

PECIFICATIONS of Model RM-8 Log-Linear Count Ratemeter

RANGES:

Log: 10 to 1,000,000 counts per minute over 5 decades. Linear: 9 full scale ranges - 0 to 100, 0 to 300,

0 to 1000, 0 to 3000, 0 to 10,000, 0 to 30,000, 0 to 100,000, 0 to 300,000, and 0 to 1,000,000.

Range Meter: Easy-to-read 41/2" meter provides direct indications of either log or linear count rate. Color coding of meter scales (black and green) identifies the scale in use.

ACCURACY: Linear Ranges: 2% Log Ranges: 5%

TIME CONSTANTS:

Log: Variable with counting rate.

- Linear: 6 individual constants may be selected for each range -0.3, 1, 3.0, 10, 30, 100 seconds. An additional position is provided for a special time constant, if desired.
- INPUT: Negative. Pulse adjustable between 0.1 and 5.0 volts. Factory preset at 0.25 volts.

RESOLUTION: 2 microseconds to pulse pairs.

AURAL SIGNAL: Panel mounted 3" speaker with internal volume control. Variable pitch control for "Howler" or "Clicker." **RECORDER OUTPUT:** Cathode follower output - 0 to 1.0 MA, 0 to 10 MV, 0 to 100 MV for standard potentiometer or galvenometer type graphic recorders.

ZERO CONTROL: Front pariel knob.

FAST RETURN: Front panel button.

TEST SIGNAL: 3600 counts per minute.

HIGH VOLTAGE SUPPLY:

Polarity: Positive

Range: 600 to 2500 volts

Current Output: 1 MA up to 1500 volts. 0.5 MA up to 2500 volts.

Meter: 41/2" meter with 50 division scale.

Ripple: Less than 10 millivolts.

Line Regulation: Less than .01% output shift for 1% change in line voltage between 105 and 125 volts.

POWER REQUIREMENTS: 105 to 125 volts, 150 watts.

DIMENSIONS: 19" x 834" x 14"

WEIGHT: 32 lbs. Shipping weight: 38 lbs.



GRAPHIC RECORDER Model GRP-1 (Potentiometer)

MODEL RM-8 RATEMETER shown with Model GRG-2F Graphic Recorder (Galvenometer) and Model LS-6 Lead Shield. This highly versatile ratemeter and lead shield combination meets numerous radiation counting applications in all phases of nuclear work.

TECHNICAL ASSOCIATES

140 W. PROVIDENCIA AVE. . BURBANK, CALIF.



MULTI-SAMPLE TABULATOR

MODELS MST-201, MST-202, MST-301, and MST-302

- Handles samples up to 3" in diameter.
- Counts 200 samples and provides full data on every sample.
- Features dual or single channels with discriminating scalers.
- Accommodates any 2 of five different types of detectors.
- Provides binary-coded outputs to operate printers or converters.





TABULATOR TAPE shows sample number and count from each detector (color-coded for easy identification).

APPLICATION. The T/A Multi-Sample Tabulator is an

automatic sample changing system designed to offer unusual capacity and speed in detecting and measuring radiation contamination contained in large-size filter paper samples, wipe samples, solid samples, and evaporated samples. The exclusive dual channel feature allows the selection of any 2 of five types of detectors and thus the instrument can measure two different kinds of radiation or two different gamma energies of each sample. Complete data on each sample is obtained by a single programming.

DESCRIPTION. The Multi-Sample Tabulator consists of an automatic sample changing mechanism with input and output storage towers that accommodate 200 samples and positions for mounting two separate detectors. The electronic console of the tabulator includes an automatic control unit and two Model DS-6 Decade Scalers (one for each detector). The control unit contains an interval timer, three decades for sample number count, and circuits to automatically start and stop the sample changer and reset the scaler decades. The scalers feature pulse height discrimination by means of a ten turn helipot mounted on the front panel and incorporate exceptionally stable 500V to 3000V power supplies.

OPERATION. The changer mechanism moves each sample from the input storage tower to the first detector position, then to the second detector position, and finally to the ouput storage tower. Accuracy of the lift mechanism assures reproducible counting of the sample in either detector position. The sample change time is 15 to 20 seconds. The samples are programmed through the system by the automatic control unit. They are counted to preset time by the Model DS-6 Scalers which provide outputs to feed sample numbers and binary-coded count information to digital printers or data converters to drive electric typewriters, teletypewriters, and card or tape punches.

5 TYPES OF DETECTORS AVAILABLE:





GAMMA SCINTILLATION



TA

TECHNICAL ASSOCIATES

140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA

MULTI-SAMPLE TABULATOR

CAPACITY: 200 samples

SAMPLE SIZE:

- MST-300 SERIES: 3" diam. holders for filter papers or wipe samples. 3" diam. x 3/16" deep planchets for solid samples.
- MST-200 SERIES: 2" diam. holders for filter papers or wipe samples. 2" diam. x 3/16" deep planchets for solid samples.
- (Unit comes with 200 filter paper holders or planchets)

SAMPLE CHANGE TIME: 15 to 20 seconds. PRE-SET TIME: 0 to 20 minutes in 1 second increments. **PRINT-OUT:** System includes Victor Digital Printer with 2-color ribbon. (Data Converters available.)

DIMENSIONS:

- MST-300 SERIES: Sample Changer: 32" x 32" x 73" high. Console: 19" wide x 27½" high. MST-200 SERIES: Sample Changer: 29" x 73" high. Console: 19" wide x 27½" high. WEIGHT: 450 lbs.

POWER REQUIREMENTS: 115 volts, 60 cycle, 450 watts.

MODEL DS-6 SCALER

INPUT POLARITY: Negative

INPUT SENSITIVITY: 0.2 volts negative. AMPLIFIER: Non-overloading with X25 gain.

PULSE HEIGHT DISCRIMINATOR: Front panel control variable from +5 volts to +55 volts.

RESOLVING TIME: 5 microseconds

COUNT CAPACITY: 100,000 counts.

- **OUTPUT:** Provides four line binary coded output to operate digital printers or converters to drive electric typewriters, card or tape punches.
- HIGH VOLTAGE: 500 to 3000 volts in two ranges (1.5KV and 3.0KV with coarse and fine controls). STABILITY: Less than .01% drift per day after warm-up.

REGULATION: Less than .01% change for 1 volt variation in line voltage. RIPPLE: Less than 10 millivolts.

INPUT SELECTOR SWITCH: Selects Test, G.M., or Amp. input. COUNT SWITCH: Up position permits instrument to count. Down position prevents counting.

MANUAL REMOTE SWITCH: Selects mode of operation.

DETECTORS

ALPHA SCINTILLATION DETECTOR (Model MST-AS-3)

Uses a 3" diameter photomultiplier tube and a 3" diameter Zinc Sulphide phosphor with transistorized preamplifier. BETA SCINTILLATION DETECTOR (Model MST-BS-3)

- Identical to Alpha Scintillation Detector, except that it uses a $2\frac{1}{2}$ " diameter Stilbene phosphor.
- MICA END WINDOW G-M TUBE (Model MST-GM-2)

Uses Anton 1001 G-M tube with $1\frac{1}{2}''$ diameter mica window 1.4 to 2 mg/cm² thick. Shielding: $\frac{7}{8}''$ lead, $\frac{1}{4}''$ brass. GAMMA SCINTILLATION DETECTOR (Model MST-GS-2)

Uses 2" diameter x 1" thick solid NaI crystal with 2" photomultiplier tube and transistorized preamplifier. Shielding: 1" lead, $\frac{1}{4}$ " brass. GAS FLOW COUNTER (Model MST-GF-2)

Consists of a 2 pi hemispherical chamber with 2" diameter window and transistorized preamplifier. For Geiger or proportional counting. Shielding: $1\frac{1}{6}$ " lead, $\frac{1}{4}$ " brass.

NOTE: 'Light-tight' design of sample changer permits use of windowless scintillation detectors.



SEPARATE MODELS AVAILABLE

MODEL PH-3

3" FILTER PAPER HOLDER

for use with

MST-300 Series

MODEL PH-2

PAPER HOLDER

2" FILTER

for use with MST-200 Series

The T/A Multi-Sample Tabulator is available with single or dual channel system, and with capacity to hold samples up to 2" or up to 3" in diameter. (NOTE: The tabulator can be furnished with a single channel as a starter system. It can be converted to a dual channel system at any time by the addition of a detector and DS-6 Scaler.) Molded plastic holders are provided for paper samples or for 3/16" deep metal planchets.

MST-201 - Single channel system. Has capacity of 200 samples up to 2" in diameter. MST-202 - Dual channel system. Has capacity of 200 samples up to 2" in diameter. MST-301 — Single channel system. Has capacity of 200 samples up to 3" in diameter. MST-302 --- Dual channel system. Has capacity of 200 samples up to 3" in diameter.

MODEL PLH-20 2" PLANCHET HOLDER* for use with MST-200 Series

MODEL PL-5 **PLANCHET*** (2" diam. x 3/16" deep) for use with MST-200 Series

*3" Planchets and Holders available for MST-300 Series



ECHNICAL ASSOCIATES 140 WEST PROVIDENCIA AVENUE . BURBANK, CALIFORNIA



PRICE LIST

ullet

TECHNICAL ASSOCIATES

140 W. Providencia Ave. • Burbank, Calif.

Revised August 1, 1962

Subject to Change Without Notice

Bulletin No.		Model No.	Description of Product F	Price O.B. Burbank
			PERSONNEL MONITORS	
168		PPM-8	Beta/Gamma Portal Type Monitor	\$1975.00
168		LIM-18	Beta/Gamma Laundry Monitor	1150.00
184		HSM-10A	Beta/Gamma Hand and Shoe Monitor with external Beta/Gamma clothing probe	3575.00
			SURVEY METERS	
185		CP-TP-1A	Cutie Pie Totem Pole Portable Survey 3 Range Meter, 0 to .5, 5, and 50 R/HF	460.00
185		CP-TP-1B	Cutie Pie Totem Pole Portable Survey 3 Range Meter, 0 to 50, 500 and 5000 R/ $$	HR 440.00
185		1A Chamber	Chamber only for CP-TP, ranges 0 to .5, 5 and 50 R/HR	225.00
185		1B Chamber	Chamber only for CP-TP, ranges 0 to 50, 500, and 5000 R/HR	205.00
186		CP-3	Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR	295.00
186		CP-3A	Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR	295.00
186		CP-4	Cutie Pie Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR	275.00
186		CP-4A	Cutie Pie Portable Survey Meter 3 Ranges 0 to 25, 250 and 2500 MR/HR	275.00
180		SRJ-7	Juno Portable Survey Meter 3 Ranges 0 to 50, 500 and 5000 MR/HR	325.00 ;
180		HRJ-7	Juno Portable Survey Meter 3 Ranges 0 to 250, 2500 and 25,000 MR/HR	325.00
136 137		F-6 FS-11	Scintillation Counter Ranges .01 over background .1, .25, 1.0, 2.5 MR/HR	159.50
			ANALYTICAL AND COUNTING INSTRUMENTS	
166	*	DS-5B	Decade Scaler (5 microseconds)	825.00
166	*	DS-5BA	Decade Scaler (1 microsecond)	915.00
200	*	DS-9	Decade Scaler (1 microsecond)	825.00
181	*	RM-8	Log-Linear Ratemeter	810.00
190	*	SM-10	Spectrometer (Single Channel Pulse Height Analyzer with built-in Linear Ampl and Ratemeter)	ifier 765.00
169	*	SA-20	Single Channel Pulse Height Analyzer	685.00
190	*	SS-30	Spectrum Scanner	475.00
		PA-6	Pre-Amplifier	115.00
178	*	RHV-1B	Regulated High Voltage Power Supply 2 Ranges (1500 and 3000 volts)	375.00
190		MS-1	Mobile Stand for Scintillation Detector, counter-weighted and with shelf for so or Ratemeter	aler 295.00
			* Prices shown are for instruments supplied in 19" panel for rack mounting ur otherwise noted. For single instrument mounted in T/A Model C-875 "Eas Service" Console cabinet add \$30.00 to prices shown. See reverse side of Bul No. 168 for description of "Easy-to-Service" cabinet.	nless 7-to- letin
			COUNTING INSTRUMENT ACCESSORIES	
203	* *	PT-1	Predetermined Timer (Long counting intervals 999 minutes)	95.00
203	* *	PT-2	Predetermined Timer (Short counting intervals 120 seconds)	95.00
203	* *	PT-3	Predetermined Timer (Liebel-Flarsheim, 1 second-60 minutes)	100.00
203	* *	CR-2	Extended Count Register (4 digit)	65.00
166		DU-1	One microsecond plug-in decade unit	147.50
166		DU-5	Five microsecond plug-in decade unit	57.50
			** These accessories will be supplied in a cabinet or panel mounted at custom option.	ner's
			AUTOMATED COUNTING SYSTEMS	
189		MST-201	Multi-Sample Tabulator, 2" sample diameter, single channel complete with DS-6 Scaler and Victor Digital printer (less detector)	one 6420.00
189		MST-202	Multi-Sample Tabulator, 2" Sample diameter, dual channel complete with DS-6 Scalers and Victor digital printer (less detector)	two 7250.00
189		MST-301	Multi-Sample Tabulator, 3" sample diameter, single channel complete with DS-6 Scaler and Victor digital printer (less detector)	one 7670.00
189		MST-302	Multi-Sample Tabulator, 3" sample diameter, dual channel complete with DS-6 Scalers and Victor digital printer (less detector)	two 8500.00
189		DS-6	Discriminating Decade Scaler for use with all model Multi-Sample Tabulator series	es 830.00
			PAGE ONE	

Ĵ

1

ノ

.

Bulletin No.	Model No.	Description of Product F.O.B.	Price Burbank
189	MST-GF-2	Gas Flow Detector for all models of Multi-Sample Tabulator including lead shielding	440.00
189	MST-GM-2	Geiger Tube Detector with Anton 1001 tube for all models of Multi-Sample Tabu- lator including lead shielding	380.00
189	MST-AS-2	Alpha Scintillation Detector with 3" diameter zinc sulfide phosphor for all models of Multi-Sample Tabulator	778.00
189	MST-BS-2	Beta Scintillation Detector with 21⁄2" diameter Stilbene crystal for all models of Multi-Sample Tabulator	1075.00
189	MST-GS-2	Gamma Scintillation Detector with 2" diameter x 1" thick Nal crystal for all models of Multi-Sample Tabulator including lead shielding	685.00
201	ASC-21	Automatic Sample Changer, Single Channel System equipped with Gas Flow Detec- tor, one power supply, one scaler and designed for time print-out	6975.00
201	ASC-22	Automatic Sample Changer, Single Channel System equipped with Gamma Scintilla- tion Detector, one power supply, one scaler and designed for time print-out	7275.00
201	ASC-23	Automatic Sample Changer, Dual Channel System equipped with Gas Flow Detector and Gamma Scintillation Detector, two power supplies, two scalers and designed for time print-out	9675.00
201	GR-1	Gas Regulator for gas flow detectors used with Multi-Sample Tabulator and Auto- matic Sample Changer.	42.50
		STAR SERIES SPECTROMETRY SYSTEMS	
171	Polaris	Includes SM-10 and RHV-1B mounted in C-1400 cabinet, complete with inter- connecting cables and system test	1180.00
171	Lyra	Includes SM-10 and DS-5B mounted in C-1750 cabinet, complete with inter- connecting cables and system test	1635.00
171	Arcturus	Includes SM-10, RHV-1B, SS-30 and Texas Instruments Rectilinear Recorder Model RR1M-A16 in P-1 Panel Mount installed in C-3150 cabinet, with interconnect- ing cables and system test	2250.00
		NOTE: For special applications requiring other combinations of T/A instruments, please request a quotation.	
		SCINTILLATION DETECTORS AND ACCESSORIES	
188	SD-1	Scintillation Detector with 1" x 1" Nal crystal, complete with cables	455.00
188	SD-2	Scintillation Detector with 2" x 2" Nal crystal, complete with cables	765.00
188	SD-2W	Well-type Scintillation Detector with $13/4$ x 2" Nal crystal ($5/6$ " x $1/2$ " well) complete with cables for use with LS-8 or LS-8X Shields	695.00
179	SD-3	Scintillation Detector with $1\frac{1}{2}$ " Alpha Phosphor, complete with cables	495.00
179	SD-4	Scintillation Detector with 11/2" Beta Anthracene crystal, complete with cables	545.00
183	SD-6-1	Scintillation Detector with transistorized pre-amplifier, 2" photomultiplier tube and 1" x 1" Nal crystal, complete wth cables. (For use in LS-6 Lead Shield with extension.)	445.00
183	SD-6-1.5	Scintillation Detector with transistorized pre-amplifier, 2" photomultiplier tube and 11/2" x 1" Nal crystal, complete with cables. (For use in LS-6 Lead Shield with extension.)	575.00
183	SD-6-2	Scintillation Detector with transistorized pre-amplifier, 2" photomultiplier tube and 2" x 2" Nal crystal, complete with cables. (For use in LS-6 Lead Shield with extension.)	765.00
183	SD-6W	Well Scintillation Detector with transistorized pre-amplifier, 2" photomultiplier tube and $134'' \times 2''$ Nal crystal ($5/8'' \times 1\frac{1}{2}$ " well) complete with cables. (For use in LS-6 Lead Shield with extension and lid for well counting.)	695.00
202	CH-1	Crystal Housing for SD-1 and SD-6-1 Scintillation Detector	25.00
202	CH-2	Crystal Housing for SD-2 and SD-6-2 Scintillation Detector	30.00
202	CH-2W	Crystal Housing for SD-2W and SD-6W Scintillation Detector	30.00
202	NS-1	Nose Shield for Model SD-1 Detector	95.00
202	NS-2	Nose Shield for Model SD-2 Detector	115.00
202	C-la	Collimator Type A (20° Flat Field) for SD-1 Detector	29.50
202	C-2a	Collimator Type A (20° Flat Field) for SD-2 Detector	39.50
202	C-15	Collimator Type B (Straight Bore) for SD-1 Detector	35.00
202	C-2b	Collimator Type B (Straight Bore) for SD-2 Detector	45.00
202	C-lc	Collimator Type C (Focusing) for SD-1 Detector	75.00
202	C-2c	Collimator Type C (Honevcomb) for SD-2 Detector	95.00
202	AS-11	Alpha Phosphor (1") for SD-1 and SD-6-1 complete with mount	79.50
202	AS-12	Alpha Phosphor (11/2") for SD-2 and SD-6-2 complete with mount	85.00
202	BS-11	Reta Stillene crystal (1") for SD-1 and SD-6-1 complete with mount	70 50
202	BS-12	Beta Stilbene crystal $(1\frac{1}{2}'')$ for SD-2 and SD-6-2, complete with mount	95.00

Bulletin No.	Model No.	P Description of Product F.O.B.	rice Burbank
		GAS FLOW DETECTORS AND G.M. PROBES	
182	GD-6	Gas Flow Counter with built-in transistorized pre-amplifier (for use with LS-6 Lead Shield with extension) includes extra window.	295.00
204	P-7	Beta/Gamma Probe Assembly with T-1090 Tube, cable and connector	39.50
204	P-8	End Window Probe (less tube) complete with cable and connector	45.00
204	P-8X	End Window Probe Assembly with T/A T-1180 Tube, cable and connector	87 50
204	P-10	Beta/Gamma Probe Assembly with T-1110 G.M. Tube, cable and connector	75.00
204	P_AS_2	Alpha Scintillation Probe	275.00
201	GR-1	Gas Regulator for Gas Flow Detectors	42.50
		COUNTER TUBES	
144	T-1090	Metal Thin wall, 3 11/16" long, used in Model F-6 Geiger Counter, P-7 Probe or other Beta/Gamma applications.	9.00
144	T-1100	Thin wall type 3" long used in Models PPM-8 and HSM-10A	27.00
144	T-1110	Thin wall 3" long used in P-10 Probe	27.00
144	T-1120	Thin wall type 7" long used in Model PPM-8	53.00
144	T-1140	Mica End Window used in TM-7, TM-8 Tube Mounts and P-8 Probe	58.50
144	T-1160	Mica End Window (Organic Quenched) used in TM-7, TM-8 Tube Mounts and	67.00
144	T-1180	Mica End Window, Alpha/Beta/Gamma sensitive, used in P-8X Probe	61.00
		IONIZATION CHAMBERS	
144	IC-1	Ionization Chamber with windows (Beta/Gamma)	115.00
144	10-1	Ionization Chamber — Solid Wall (Gamma)	105.00
	10-2		100.00
191	LS-1	Complete with Type ST-1 Sample Tray, Type SK-1 Socket and Type TM-1 Tube Mount	185.00
191	LS-2A	Standard with Type TM-6 Tube Mount and Sample Tray Holder, Type ST-2A Sample Tray and 10 Type PL-2 Aluminum Planchets	200.00
191	LS-2B	Micrometric with 10 Type PL-1 Planchets	285.00
192	LS-4A	For liquids complete with Type TM-1 Tube Mount, Type SK-1 Socket, and one Type ST-5 Marinelli Beaker	200.00
192	LS-4B	Modified to accommodate special tubes with aluminum inner liner	235.00
192	LS-4C	For solids, complete with Type TM-1 Tube Mount, Type SK-1 Socket, one Type ST-6 Ore Container and 100 Type ST-6A Paper Sleeves	225.00
192	LS-5	Complete with Type TM-3 Tube Mount, Type SK-2 Socket, Type STH-1 Sample	240.00
193	LS-6	Complete with Type TM-7 Tube Mount, one Type ST-3B Sample Tray and 10 Type	200.00
100		PL-3 Aluminum Planchets	295.00
193	LS-6XT	One-inch height extension for LS-6 shield when used with GD-6 Gas Flow Detector	48.00
193	LS-6X8	Eight-inch height extension for LS-6 shield used with SD-1, SD-2 and SD-6 Scintilla- tion Detector	98.50
193	LS-6L	Lid with 3" opening and retaining collar for Models SD-1 and SD-2 Scintillation Detectors	72.50
193	LS-6P	Plug for 3" opening in LS-6L	12.50
194	LS-7A	Multi-Purpose Type complete with SD-1 Scintillation Detector, three ST-7 Sample Trays and 10 PL-3 Planchets	780.00
194	LS-7B	Multi-Purpose Type complete with SD-2 Scintillation Detector, three ST-7 Sample Trays and 10 PL-3 Planchets	1090.00
194	LS-7C	Multi-Purpose Type complete with TM-8 Tube Mount, three ST-7 Sample Trays and 10 PL-3 Planchets	385.00
194	LS-7M	Multi-Purpose Type — Shield only, with adapter ring, three ST-7 Sample Trays and 10 PL-3 Planchets	325.00
195	LS-8	For well counting (shield only, less detector)	265.00
195	LS-8W	For well counting, complete with SD-2W Scintillation Detector and $13/4$ x 2"	
-		well crystal	960.00
195	LS-8X	For well counting with 2" shielding in all directions (Shield only less Detector)	350.00
195	LS-8WX	For well counting with 2" shielding in all directions complete with SD-2W Scintillation Detector and $134'' \times 2''$ well crystal	1035.00

LABORATORY ACCESSORIES

198

AB-2

C

L

Absorber Set

85.00

Bulletin No.	Model No.	Description of Product		F.0	Price .B. Burbank
199	ST-1	Sample Tray			6.50
199	ST-2, ST-3	3, ST-3X, ST-3A, ST-3B, ST-4, ST-4A Sample Trays			3.50
199	ST-5	Marinelli Beaker			4.25
199	ST-6	Ore Container			12.50
199	ST-6A	Paper Sleeves for Ore Container		per	C 2.00
199	ST-7	Sample Tray for LS-7	•••••		2.50
199	STH-1	Sample Tray Holder with tray and ten planchets			13.50
199	SK-1, SK-	2, SK-3 Sockets			5.00
19 9	TM-1, TM	1-2, TM-3, TM-4 Tube Mount			15.00
199	TM-5	Tube Mount (micrometric for LS-2 Lead Shield)			100.00
199	TM-6	Tube Mount and Sample Tray Holder			22.50
199	TM-7	Tube Mount and Sample Tray Holder	••••••		27.50
199	TM-8	Tube Mount and Adapting Ring for LS-7 Multi-Purpose Shield			60.00
	SDM-9	Mount and Support Ring when SD-6 is used in LS-6 Shield	••••••••••••••••••••••••••••••••••••••		27.50
197	LC-1	Lead Source Container for four radium needles			45.00
197	LC-3	Lead Source Container 1" lead walls, $1\frac{1}{2}$ " x 3" inside dimensions			40.00
197	LC-3A	Lead Source Container $13\%''$ lead walls $34'' \times 214''$ inside dimensions	;		47.50
197	LC-4	Lead Source Container $1\frac{1}{2}$ " lead walls $1\frac{1}{2}$ " x 3" inside dimensions			60.00
197	LC-4A	Lead Source Container $1\frac{3}{4}$ " lead walls $\frac{3}{4}$ " x $2\frac{1}{2}$ " inside dimensions	• • • • •		67.50
197	LC-9	Mobile Source Unit — compartment 3" x 6" in shield of 2" lead.	Moun	ted o	250.00
130		Lead Brick $2'' \times 3'' \times 6''$		•••••	9 00
130		Lead Brick, 2" x 4" x 8"	•••••		12 50
130		Lead Brick, extruded shapes	Prices	upon	application
199		Planchets, Type PL-1, Sheet steel, tinned	per	c	3.25
			per	Μ	30.00
			per	5M	27.00/M
199		Planchets, Type PL-2, Aluminum	per	C M	3.25
			per	5M	27.00/M
199		Planchets, Type PL-3, Aluminum	per	с	3.75
			per	M	35.00
100			per	5M	31.50/M
199		Planchets, Type PL-2C, Copper	. per	м	5.25 47.00
			per	5M	42.30/M
199		Planchets, Type PL-3C, Copper	. per	С	6.00
			per	M 5M	55.00
100		Planchets Type PL-25 Stainless Steel	Der	5/m	9 00
177			per	м	65.00
			per	5M	58.50/M
199		Planchets, Type PL-3S, Stainless Steel	per	C	12.00
			per	5M	76.50/M
199		Planchets, Type PL-4, Aluminum	per	C	5.00
			per	Μ	45.00
			per	5M	40.50/M
		Planchet Holder for ASC-21, ASC-22, ASC-23	. per	С м	6.00 55.00
		Sample Trav for ASC-21 ASC-22 ASC-23	. Der	c	8 00
			per	M	75.00
189	PH-3	Plastic Holder for 3" diameter filter paper and wipe samples used wit MST-300 series	h per	с	30.00
189	PH-2	Plastic Holder for 2" diameter filter paper and wipe samples used wit	h 	c	20.00
180		Mol-200 series	. per	č	30.00
189	Р∟П-30 ріц эл	Plastic Holder for 2" diameter planchets used with MST-200 series	. per	c	30.00
189	PI_5	Aluminum Planchet 2" diameter x 3/16" deen used with MST-200		0	50.00
100	. E.J	series	. per	С	8.00
189	PL-6	Aluminum Planchet 3'' diameter x 3/16'' deep used with MST-30 series	J per	с	11.00
		Planchets, special materials and sizes Price	s quot	ed u	pon request

PAGE FOUR

/

Ĵ

Bulletin	Model	Description of Product F(
NO.	140.		
168	C-875	8 ³ / ₄ " height "Easy-to-Service" console cabinet for standard 19" width panel	30.00
168	C-1050	101/2" height console cabinet for standard 19" width panels	35.00
168	C-1225	121/4" height console cabinet for standard 19" width panels	37.50
168	C-1400	14" height console cabinet for standard 19" width panels	40.00
168	C-1750	171/2" height console cabinet for standard 19" width panels	45.00
168	C-2100	21" height console cabinet for standard 19" width panels	50.00
168	C-2625	261/4" height console cabinet for standard 19" width panels	55.00
168	C-3150	311/2" height console cabinet for standard 19" width panels	60.00
168	C-4375	433/4" height console cabinet for standard 19" width panels	70.00

IN GENERAL

Unit prices are given in this price list unless otherwise noted.

All prices quoted here in are f.o.b. Burbank, California. Payment terms, net 30 days upon approved credit.

No extra charge is made for packing for domestic shipments. Export charges will be quoted upon application. We reserve the right to alter specifications at any time without incurring the obligation of incorporating new features in previously manufactured equipment. Prices are subject to change without notice. Quotations remain firm for 30 days.

Shipping instruments: Please include shipping instructions when ordering; in the event that shipping instructions are not given, we will use our best judgment in the matter.

SERVICE CHARGES

Services charges for repairs are usually billed on an hourly basis, however, if desired, an estimate can be given before the work is undertaken. The unit to be repaired should be shipped prepaid to us.

WARRANTY

Technical Associates warrants instruments and equipment (except tubes, fuses, batteries and crystals), manufactured by them to be free from defects in workmanship or materials under normal use for a period of one year from the date of shipment from the factory to the buyer. Tubes, fuses, batteries and crystals are subject to the guarantee established by the manufacturer of them, however, Technical Associates will assist the customer to obtain full benefits of these guarantees.

If, within the one year warranty period, any Technical Associates instrumentation or equipment requires service as a result of a defect, the buyer may return it to the factory of Technical Associates at Burbank, California or to a service station designated by Technical Associates, transportation charges prepaid, for service at no charge under the warranty. The buyer is urged to communicate with Technical Associates when warranty service is required, stating the nature of the difficulty and giving model and serial number of instrument. It may be possible to diagnose the trouble and send a replacement part or assembly, thereby avoiding the expense of shipment.

Technical Associates will return the instrument to the buyer, transportation charges prepaid, after repairs or replacement under warranty are completed. The liability to Technical Associates under this warranty is limited to the cost of replacement of defective parts upon prompt notification of such defect.

Manual of Operation and Service Instructions



NUCLEAR INSTRUMENT & CHEMICAL CORPORATION

INSTRUCTION BOOK

MODEL 2610A ≠ /747

SECTION I

GENERAL DESCRIPTION

PURPOSE

Model 2610A count-rate meter is a lightweight battery operated instrument designed for general survey work and for the location of small amounts of beta and gamma radiation in rooms, laboratories, on desks, laboratory coats, and similar areas. It has also been found useful for X-ray monitoring and geological surveying for radioactive ore.

GENERAL SPECIFICATIONS

- 1. The instrument consists of a single unit containing all circuit components, electron tubes, batteries, and a count-rate meter.
- 2. The detector is a thin-wall Geiger counter mounted in a metal probe at the end of a three foot cable. The probe has a rotary shield which, when covering the sensitive area of the counter, effectively cuts out beta radiation. The instrument then reads only the gamma ray component of the incident radiation, excluding all but the highest energy beta rays commonly encountered.
- 3. Three full scale ranges of 20, 2, and 0.2 milliroentgens per hour (mr/hr) have been provided to permit quick and convenient measurements of radiation dosage rates.
- 4. Each instrument is calibrated at the factory with a known ionization intensity produced by gamma rays from radium in equilibrium with its short-lived decay products.
- 5. When it is inconvenient to read the meter, the earphones may be plugged into the jack on the instrument. One "click" is produced in the phones for each ionizing event occuring in the counter.

MECHANICAL SPECIFICATIONS

- 1. The count-rate meter is housed in an aluminum case with a smooth enamel finish to make surface decontamination easy.
- 2. The instrument is ruggedly constructed for long operating life with both case and probe sealed against moisture.
- 3. The probe is conveniently mounted in the handle of the unit and is easily removed.
- 4. The range switch is located below the handle and can be controlled with a fingertip while the probe is carried in the other hand.

~~]___

- 5. The plug-in Geiger tube (model D50) has a three pin "pee wee" base, requires no soldering, and the probe is quickly disassembled by a simple twist when necessary.
- 6. Cabinet dimensions: 10" x 4-3/4" x 5-3/4" deep.

7. Weight: 9-1/2 pounds.

ELECTRICAL SPECIFICATIONS

All the batteries necessary for proper operation of the instrument are shipped with the equipment. They consist of:

- 1. Three 300 wolt battery packs which are used for GM tube potential. The life of these batteries is equal to shelf life.
- 2. One 67.5 volt battery which is used for B-plus potential.
- 3. One 1.5 volt battery which is used for filament potential.

The life of the 67.5 volt and 1.5 volt batteries is adequate to maintain sensitivities within 10% when operated continuously for 100 hours. When operated four hours a day, the life of these batteries is over 250 hours.

SECTION II

OPERATION

FACTORY CALIBRATION

The instrument has been calibrated using gamma rays from a radium source with the rotary shield on the metal probe closed. Therefore, the shield should be closed when reading gamma radiation in order to preserve the calibration. With the shield closed, no beta rays will be able to penetrate into the chamber and affect the reading. Gamma rays can thus be monitored in the presence of a beta background. When it is desired to admit beta rays into the chamber, the shield should be set in the open position.

The calibration of the instrument may be checked by use of the radium calibration source supplied with the unit.

RADIUM CALIBRATION SOURCE

1.4

The radium calibration source provides a convenient means of checking the calibration of the instrument. It is contained in a small plastic cylinder mounted at the rear of the instrument Case.

To use the calibrating source, <u>open</u> the probe shield and hold a flat side of the source against the unshielded section of the probe. Move the source until the maximum reading of the meter is observed. The numerals stamped on that side of the source next to the Geiger counter represent the reading in mr/hr. which should be indicated on the meter. If the error is greater than 5% of full scale, readjust the instrument by means of the CALIBRATE control located at the top of the instrument.

The background reading of the meter will be increased if the calibrating source is mounted in its holder while the instrument is in use. If the source is not moved some distance from the instrument, this increase in meter reading should be taken into account.

TAKING A COUNT

When taking a count the Geiger counter may be used in its attached position if desired, or removed for probing in confined spaces or small areas.

The range switch has four positions: OFF and three full scale sensitivity selections of 20, 2, and 0.2 milliroentgens per hour. The maximum total dose to which any part of the body of a person shall be exposed continuously or intermittently in a given time is now 300 milliroentgens per week. On the basis of 48 hours per week of uniform exposure, the permissible dosage rate is 6.25 mr per hour. This dosage rate is indicated on the 20 mr/hr. scale.

The range switch should be turned first to the least sensitive scale (20 mr/hr) to prevent harming the meter if a strong radiation field is present. No warm-up period is necessary. The countrate meter is ready to operate almost instantly. If no reading appears on the 20 mr/hr scale, the switch should be turned to the next more sensitive scale, etc. With the radium calibrating source removed from the vicinity of the instrument, cosmic ray background will indicate about 1/10 of full scale on the 0.2 mr/hr. range.

SECTION III

CIRCUIT THEORY

GEIGER TUBE

The Geiger tube used with the instrument detects beta and gamma radiation from a radioactive source. The three 300 volt batteries in series produce a 900 volt electric field between the anode and cathode of the Geiger tube. Entering beta or gamma rays

collide with the atoms of gas in the tube, causing ionization of the gas and resulting in a current flow between the electrodes. The random height, random width voltage pulses produced by this minute current flow are then fed to a two tube selfrelaxing trigger circuit. Refer to the schematic diagram at the rear of this manual.

TRIGGER CIRCUIT

The output of the GM tube is differentiated (sharpened) by Cl and R3 and fed to a one-shot multivibrator, or trigger circuit. The first tube of the trigger circuit (VT-1) is normally conducting, developing a DC voltage with respect to ground on the filaments of both tubes and biasing the second tube of the trigger circuit to cut-off.

When a negative pulse from the Geiger tube is impressed on the grid of VT=1, the conduction of the tube decreases and its plate goes in a positive direction. This positive surge in plate voltage is transmitted to VT=2 through C4, C5, or C6 and VT=2 then leaves cut-off. The plate current through VT=2 will further bias the first tube. This process will be self-maintaining and will terminate when the first tube is fully cut-off and the second tube is conducting heavily.

As C4, C5, or C6 becomes charged, the grid of VT=2 approaches the cathode potential. At this point, grid-cathode conduction ceases and the charging proceeds at a slower rate through R9, R10, and R5. The grid of VT=2 gradually goes negative decreasing the current flow through the tube and decreasing the bias on VT=1. When this bias decreases to the point where plate current will once again flow in the first tube, the circuit will rapidly trigger back to the initial condition with VT=1 conducting heavily and VT=2 cut=off. The time during which the second tube is conducting is varied by the range selector switch which changes the time constant of the coupling circuit. Since the range selector switch varies the "on" time of the second tube, it changes the number of input pulses necessary to give the same average current flow through VT=2.

METER-INTEGRATING CIRCUIT

Each time VT_{-2} conducts, a pulse of current is put into the meter integrating circuit made up of R8, C7, and the meter in the plate circuit of VT_{-2} . These pulses of current are averaged by the charging and gradual discharging of the condenser through the resistor and meter. The time constant of the meter-integrating circuit, i.e., the time it takes the meter to reach 63% of its final reading is about 5 seconds. It will be noted that for any given amount of radiation, the meter needle fluctuates around an average reading; the largest fluctuations being on the most sensitive scale (O_{2} mr/hr). This is due to the randomness of the radiation flux. The observer's final reading should be the average position, as seen by the eye, of the fluctuating needle.
LINEARITY

The linearity of the electronic count rate circuit corresponds very favorably (about 2%) with the linearity of the 0-20 microamp meter. However, the resolution time of the Geiger counter causes coincidence losses at the higher counting rates. The counter is insensitive for a time equal to:

N
$$-\frac{1}{60}$$
 Minutes during each minute.

Where N = observed count rate in counts per minute T = resolution time in seconds

Two particles or quanta entering the counter within this time interval produce only one pulse since the second particle enters the counter before most of the positive ions from the previous discharge are collected. Because of the random nature of emissions from a radioactive source, there is as much chance for pulses to occur during this time as during any other of equal duration.

Thus, the higher the counting rate, the higher is the number of counts "lost". Coincidence loss has been minimized on the 2 and 20 mr/hr. ranges by the introduction of slightly larger capacitors (C4 and C5) in the coupling circuit. Figure 1 shows the deviation of an average instrument from true linear calibration.



CALIBRATION CORRECTIONS *

A plot of correction factor versus photon energy for a Geiger-Müller counter is shown below. For energies above 750 kv (.75 mev), the correction factor is equal to or near 1.0. For less energetic radiations the plot exhibits maximum and minimum values of correction factors caused by a change in efficiency of the Compton and photoelectric processes.

NOTE

This graph represents average readings taken on a number of 2610A count-rate meters. The graph serves merely to indicate the approximate error which may be expected from low energy radiation and is not to be taken as an absolute guide for each instrument.

This instrument was carefully calibrated at the factory using gamma rays from a radium source in equilibrium with its shortlived decay products. The calibration of the instrument will be correct for radium radiation since the gamma energies of RaC (which is the principle source of gammas in radium decay) fall for the most part on the flat portion of the curve below.

The relative energy and intensity of RaC gammas is given in the



* Excerpt from: National Bureau of Standards Circular No. 507(7/25/51)

-6-

following chart:

ENERGY	RELATIVE INTENSITY
°6 mev	8,9
1.12 mev	2.9
1.76 mev	2.8
2 . 19 me v	1

From the graph on the opposite page it may be seen that as the energy of radiation is reduced the correction factor rises and reaches a maximum value in the region of .25 to .3 mev. Upon further decreases in energy the correction factor drops rapidly and passes through a minimum around 0.1 mev, after which it again rises sharply.

The response of a Geiger-Müller counter depends primarily upon the number of electrons that traverse the sensitive volume of the counter. This number is proportional to the number of electrons released in the vicinity of the sensitive volume and to the range of these electrons. The correction factor of a GM counter varies primarily in proportion to the ratio

Energy of electrons generated in the counter vicinity

Range of electrons generated in the counter vicinity

The range of an electron increases faster than its energy increases at moderately high energies and in direct proportion to its energy in the multimillion-volt region. Therefore, for energies within the plot, the energy-to-range ratio (the correction factor) must increase as the energy decreases. This assumption is approximately valid in the high energy region where most of the electrons arise from the Compton effect, because this effect is almost independent of the mode of binding of the electrons within matter.

At lower energies, electrons are released primarily because of the photoelectric effect. This effect is the more intense the more tightly bound are the electrons within matter, that is, the higher is the atomic number of the material. Therefore, the presence of high-atomic number materials in the proximity of the sensitive volume of a GM counter increases the response to low energy particles sharply. Accordingly, the correction factor must be expected to drop as one proceeds from high to low energies, as soon as the photoelectric effect becomes important. This drop may be seen in the plot and represents increased sensitivity of the GM counter to radiation energies in the 100 kv region.

From the chart it can be seen that the meter reading for low energy gamma emitters is essentially erroneous. To obtain the correct gamma dose rate, it is necessary to multiply the meter reading by the correction factor.



SECTION IV

MAINTENANCE

CALIBRATION WITHOUT A STANDARD SOURCE

If a standard source is not available, an approximate calibration may be obtained by the use of an electronic pulse generator. The pulses should be negative, no more than 75 microseconds wide (preferably much narrower) with a fast rising time.

Feed the pulses into the input grid circuit. The screwdriver adjustment should be made so that full scale on the second range (2 mr/hr.) is indicated when 110 pulses per second are fed into the input. Full scale on the 0.2 mr/hr. range should then be 11.0 pulses per second.

This method of calibration was determined experimentally by averaging the pulses per minute at 1 mr/hr. for a large number of Geiger counters. Individual counters may differ by as much as a factor of two, so this method of calibration should only be used if a radium standard is not available.

GEIGER TUBE REPLACEMENT

- 1. With the tube socket housing held firmly in one hand and pointed to the operator, grasp the tube housing in the other hand and turn clockwise.
- 2. Slide the probe tube housing off the Geiger counter, being careful not to break the counter tube.
- 3. Grasp the Geiger counter by the base and remove it from the socket. CAUTION The Geiger tube is very fragile and should be handled carefully.
- 4. See that the "O" ring is in the tube socket housing.
- 5. Insert a new Geiger counter in the tube socket. "O" ring now forms a seal between the housing and the base of the counter.
- 6. Carefully slide the probe tube housing over the Geiger counter and turn counter-clockwise on the tube socket housing.

HOW TO TEAR DOWN

The instrument should be "torn down" in the following manner for servicing.

- 1. Remove the two screws from the bottom of the cabinet.
- 2. Holding the cabinet with one hand and pulling on the handle with the other will remove the top section and chassis (which is attached to the top).

3. The top piece may be folded to one side to facilitate working on the circuit. This is done by removing the two screws that hold the circuit shelf. The top can then be rotated on the rivets to a position which exposes the circuit for examination.

In reassembling the instrument, care must be taken in lining up the top of the instrument with the groove in the top of the cabinet. Uniform pressure should be applied to the top to bring the bottom section evenly to the rubber seal. If the seal is not completely contacted, the cabinet is no longer waterproof.

BATTERY REPLACEMENT

- 1. BE CAREFUL the high voltage batteries can give a shock which may be lethal.
- 2. Turn the range selector switch to the OFF position, Remove the negative and positive leads from the high voltage batteries and disconnect the jumpers between the batteries.
- 3. If the low voltage batteries are to be replaced, their leads must be disconnected. The 67.5 volt battery has snap-type connections which are easily removed with a screwdriver or sharp-pointed instrument.
- 4. Remove the two screws that hold the center shelf. Tilt back the shelf and remove the batteries.
- 5. Replace the batteries as shown below. The low voltage batteries should be replaced first. Then insert the two end high voltage batteries and slide the center high voltage battery in position.
- 6. Replace the battery leads and jumper connections.



Battery Replacement and Wiring

-10-

SECTION V

TROUBLE SHOOTING

GENERAL

A visual inspection will often help in determining the cause of improper operation of the instrument. Listed below are some of the defects that may be found:

1. Loose screws and nuts.

2. Broken solder connections.

3. Damaged Insulation.

4. Pitted or dirty switch contacts.

5. Burned or broken resistors.

6. Exposed wiring, causing shorts with wiring, terminals, or ground.

7. Corroded or swollen batteries.

The following notes are given to aid in the diagnosis and repair of any failures not due to any of the above defects. The steps are given in the order in which they should be tried. If the first step does not locate the trouble, proceed to the second, etc.

NO METER INDICATION OF BACKGROUND OR OTHER RADIATION

- 1. Check high voltage batteries and plateau of Geiger counter to see if the counter is being operated in the plateau region. The ordinary life of the self-quenching organic vapor Geiger counter is about 10⁸ counts. The counter should have a threshold voltage of from 800 to 825 volts when new. As the batteries age, the voltage will decrease. The high voltage should be at least 25 volts and not more than 200 volts above the threshold voltage of the Geiger tube.
- 2. Check the operation of the circuit with a pulse generator fed into the input grid circuit (see page 9). The circuit should trigger and give a meter indication with a 0.5 to 1.0 volt negative pulse.
- 3. Check circuit batteries and voltages.

4. Replace tubes.

ABNORMALLY HIGH METER READING WHEN NOT IN RADIATION FLUX

1. Put probe with counter in complete darkness. If meter reading

falls to the background value, the Geiger tube is photosensitive. Paint the tube with black glyptol or cover with black paper.

- 2. Remove the GM tube from the probe. If the meter no longer indicates, the counter was probably in continuous discharge and must be replaced.
- 3. Circuit may be oscillating. This can be detected in the headphones. Replace the tubes.
- 4. VT-2 may not be properly cut-off and thus may be drawing plate current with no signal. Replace VT-2.
- 5. VT-1 may have an open filament and thus not draw plate current necessary to bias VT-2 to cut-off. Replace VT-1.
- 6. Check voltages and test for shorted or open circuits.

ERRATIC METER MOVEMENT

- 1. If the meter needle moves to full scale and then starts to come down to zero again when brought into a more intense radiation field, the high voltage may be too low. If the high voltage is not enough above the threshold voltage at high radiation intensities, the output pulses will become too small to trigger the circuit. Normally, an instrument will not exhibit such symptoms unless the radiation intensity exceeds 2 to 20 r/hr (100 to 1000 times full scale on the least sensitive range). At this intensity, the Geiger counter becomes paralyzed.
- 2. Check operation of the circuit with a pulse generator. The circuit should trigger and give a meter indication with a 0.5 to 1.0 volt negative pulse fed into the input grid circuit (see page 9).
- 3. Check low voltage batteries and circuit voltages.

4. Replace tubes.

VOLTAGE MEASUREMENTS

A voltmeter of at least 5 megohms input resistance should be used for measuring circuit voltages. The voltages given on the circuit diagram (rear of this manual) and on the terminal board diagram (opposite page) were measured with new batteries. The voltages, excluding high voltage for the Geiger counter tube, may vary by as much as 15% to 20% without affecting circuit operation (except for calibration).

-All voltages were measured with respect to ground.

-12-



1

Voltages to chassis measured on voltmeter with total resistance of 5 megohms.

SECTION VI

PARTS LIST

NUCLEAR PART NO.

BATTERIES

B1	1.5 volt	Burgess 2F	BA003
B2	67.5 volt	Eveready 467	BA005
B3	300 volt (3 required)	Eveready 493	BAOOO

RESISTORS

All resistors are Allen-Bradley, $\frac{1}{2}$ watt, 10% fixed composition type unless otherwise indicated.

RI	2.2 M		RC20AE225K
R2	2.7 M		RC20AE275K
R3	330 K		RC2OAE334K
RÁ	33 K	•	RC2OAE333K
R5	200 K		RC20AE204K
RŚ	1.2 M		RC20AE125K
R7	150 K		RC20AE154K
R8	33 K		RC20AE333K
R9	200 K		RC20AE204K
RÍO	500 K	2 Watt Potentiometer	RV030

CAPACITORS

C1	270 mmfd	1600 V	Mica	CM003
C2	01 mfd	300 V	Mica 20%	CM35A103M
C3	005 mfd	500 V	Mica 20%	CM35A502M
C4	005 mfd	500 V	Silver Mica 2%	CM35C502G
C5	450 mmfd	500 V	Silver Mica 2%	CM35C451G
C6	33 mmfd	500 V	Silver Mica 2%	CM20C330G
Č6	33 mmfd	500 V	Silver Mica 2%	CM2OC33OG
C7	50 mfd	6 V	Electrolytic	C E OO8

TUBES

VT1	CK522	Raytheon	VTCK522AX
VT2	CK522	Raytheon	VTCK522AX
GM	CK1020 (D50)	Raytheon	VT=010
GM	CKT020 (D20)	nayoneon	

MISCELLANEOUS

Meter 0-20 milliamps	(waterproof) Simpson 185	MP003
Selector Switch	Centralab 2CHW11486	SW-011
46 ¹ / ₂ -in. Probe Cable	Amphenol RG590	RA=001
Hadioactive Source	Nuclear	PR-005
Phones (double headset) Brush 200A	PH-001

NUCLEAR PART NO.

B1 Battery Plug B2 Battery Connectors United Carr 52383	PL-002 C0-023
United Carr 52384 3 Black Battery Plugs (for B3) Amphenol 71-1M	CO-024 PL-006
3 Red Battery Plugs (for B3) Amphenol 71-1L	PL-007
Midget Phone Jack Switchcraft 2J1065	JA-002
2 Phone Jack Covers Croname A23559	JA-005
Phone Plug Switchcraft 2P-1020	PL=004
Female Connector (GM Socket) Amphenol 91MPF-3S	CO-012
2 Socket Nuts (for above) Amphenol 91MC-3F	CG-002
2 Tube Sockets (5 pin min.) Cinch 54A11953	S0-006
2 Retainer Rings (for above) Cinch 20K12446	CL-016
Tube Holder Nuclear	BC-024
2 Rubber "O" Rings for #10 screw Lavelle	GR-007
"O" Ring-21/32 ID x 29/32 OD x 7/16 Lavelle	GR-008
2 1" Rubber Grommets Lavelle	GR-004
2 $\tilde{1}0$ " Rings 3/8 OD x 1/4 ID x 1/16 Lavelle	GR-009
Handle, with Probe Holder Nuclear	HA-005
Plastic Battery Insulator Nuclear	IN-040
Selector Switch Knob H. Davies 2110AZ	KN-015
Potentiometer Bracket Nuclear	BC-023
2 Self-tap screws 10-32 Slotted, 1" Cad. Plated	SC-010
3 ft. Fungus Resistant Waxed Cable Cord	TW-001

SUGGESTED MAINTENANCE KIT OF PARTS FOR FIVE INSTRUMENTS, FOR ONE YEAR

2	Vacuum Tubes	Raytheon CK522	VTCK522AX
5	Geiger Counter Tubes	Raytheon CK1020	VT-010
2	GM Tube Sockets	Amphenol 91MPF3S	CO-012
2	Socket Housing Nuts	Amphenol 91MC-3F	CG-002
2	Capacitors, 50 mfd, 6 V, Elect	rolytic	CE-008
2	Capacitors, 270 mmfd, 1600 V,	Mica	CM-003
2	Phone Jack Covers	Croname A23559	JA-005
1	Switch - wafer, 6 circuit, 4]	osition	SW-011
1	Meter, 0-20 ma (waterproof)	Simpson 185	MP-003
4	Screws, 10-32 Slotted, self-ta	pping	SC-010
10) rubber "O" Rings, $3/8$ OD x L	/4 ID x 1/16	GR-009
5	Rubber "O" Rings, 21/32 ID x	29/32 OD x 7/16	GR-008

9

NOTE

Since battery replacements depend upon rate of usage as well as shelf life, no general recommendations can be made for spares.



с, 2 ----

Beckman

Process Instruments Division

2500 Fullerton Road, Fullerton, California : LAmbert 5-8241 : from Los Angeles, OWen 7-1771 a division of Beckman Instruments, Inc.

October 31, 1958

General Electric Company 4855 Electric Avenue Milwaukee 1, Wisconsin

Attention: J. J. Jech X-Ray Department

Reference: Your letter dated October 13, 1958

Gentlemen:

In reply to your referenced inquiry, please be advised that we do not have a readily available parts list or kit to convert the Model V to a logarithmic type. (Model VL).

Due to the required special handling, the final cost for this conversion would exceed that of a new Model VL. In view of this fact, we obviously recommend the purchase of a new Model VL for \$425.00 net exclusive of freight. Shipment for this item is seven days.

We trust that the above meets with your approval, and if any further questions should arise, please do not hesitate to contact this office or our local sales engineer, Mr. Herbert Feitler, 4914 West Belmont Avenue, Chicago 41, Illinois.

Very truly yours,

B.B. Hoffini

B. B. Goffin Product Specialist Process Instruments

BBG:ch cc: H. Feitler Encl.

micro- microammeter

measures one millionth of a millionth of an ampere



d-c breaker amplifier

measures a-c and d-c voltages in the microvolt and fractional microvolt ranges



Beckman^{*}/model V "Vibrating Reed" micro-microammeters

Beckman/vibrode*

lope, changes a low-level d-c to a lowlevel a-c potential. The Vibrode holds the secret of the exceptional driftfree performance and low price of the Models V, VS and V2. Plates are perfectly parallel and are approximately .002 inches apart at rest position. The dry nitrogen gas in which they are sealed builds up a "cushion" and prevents any strain during operation. The Vibrode has a longer life than conventional vibrating reeds and is replaced as easily and inexpensively as a vacuum tube.

The Beckman Vibrode, a unique vibrating reed sealed in a metal enve-

model V specifications:

Input power. 60 cycles, 103 to 127 volts, 100 watts.

Detection range. $3 \ge 10^{-13}$ to $3 \ge 10^{-7}$ ampere.

Indicator. 200 microamperes, $4\frac{1}{2}$ inch scale.

Ambient Temperature requirements. 20°F to 115°F.

Zero drift. Less than 1 mv in 24 hours. Drift is not cumulative.

Response time. Time constant (4.0 to 0.12 second) varies with input capacity.

Warmup time. Amplifier, 3 to 5 minutes. Isothermal shield, 30 minutes.

Output connections. Recorder, 50 millivolt potentiometer type. Multiple switch, Beckman 1710 multiple switch.

Accuracy and Reproducibility. $\pm 1\%$ on all ranges.

Output noise. Less than 3% with no input capacitance and less than 8% with 5,000 µµf input capacitance on 3×10^{-13} ampere range.

When the huge Savannah River Atomic Energy Project was started, an instrument which would measure electrical currents less than one millionth of a millionth of an ampere was needed, to measure radiation involved in nuclear reactions. Two Beckman instruments were submitted to fill this need, and after months of tests, placed first and second against all other entries, in stability, accuracy, and reproducibility. One of these was the Beckman Model V Micro-Microammeter-a vibrating-reed electrometer.

The extreme sensitivity of this instrument is useful for measuring output of radioactive sensing elements in reactor control and health monitoring; for tracing output of ionization chambers while monitoring and controlling liquid levels, thickness of coating or impregnation applied to moving sheets; for measuring



Typical Model V application in density measuring system.

currents across capacitors and resistors, output of phototubes, static electricity generated in turbulent systems, and research work involving micro-currents.

The Beckman Micro-Microammeter is available in four models. Three of these models (V, VS, V2) employ an alternating current vibrating-capacitance modulator, the Beckman Vibrode. Thirteen measuring ranges from 3×10^{-7} (0.3 micromomperes) to 3×10^{-13} (0.3 micro-microamperes) are available. At the current sensitivity at which it is operated, zerodrift is negligible and an interruption of power does not affect the calibration.



model V multiple range model

Model V is an outstanding precision instrument for measuring, recording and controlling micro-currents in any of 13 ranges. Designed for easy relay-rack mounting or supplied with attractive cabinet for bench use. Range: 3 x 10^{-13} to 3 x 10^{-7} ampere. Weight 38 lbs.

model V2 limited range model

Model V2 is an inexpensive unit for industrial applications and original equipment manufacturers. The range is $10 \ge 10^{-13}$ to $3 \ge 10^{-11}$ ampere. Weight 38 lbs.

model VL logarithmic scale model

A companion laboratory instrument to the Model V but without the Vibrode, Model VL does not convert a signal to a-c. It has a single logarithmic scale for measuring, recording and controlling micro-currents within the range of 10^{-13} to 10^{-7} ampere. No data is lost by range switching or by the meter going "off scale." The Model VL gives almost instantaneous response to input signal fluctuations between 10^{-11} and 10^{-7} ampere. Weight 38 lbs.

available with all models

Regulated 200-volt d-c source available for polarizing ion chambers. Terminals provided for connecting external meter and for any 50 mv potentiometer-type recorder. Long input cables can be used without introducing large time constants. Source has maximum drift rate of 0.30 millivolts per second.







During World War II, a new method, superior in many respects to the suspension galvanometer, was developed for measuring d-c voltages and low-frequency a-c voltages in the microvolt and fractional microvolt regions. The Beckman d-c Breaker-Amplifier, which utilizes this new method, is insensitive to vibration, faster in response, and much more flexible than previous instruments. It has the lowest noise level and greatest zero stability of any d-c amplifier available commercially. It is capable of driving directly standard recorders, relays and meters. The Breaker-Amplifier is used extensively in both research and industry for replacing galvanometers in bridge circuits for measuring the Hall effect; for radiation measurements used in conjunction with barrier layer cells, phototubes, and thermocouples; in such diverse fields as optical pyrometry, precision tem perature measurements, colorimetry, high sensitivity photometry, spectroscopy, molecular-weight determinations, seismology, meteorology, petroleum exploration, and in physiology for the measurement of nerve-voltages in blood-flow meters, and in oximeters.

Beckman LB breaker-type ultrasensitive d-c amplifier



MODEL LB14 employs an 8-cycle breaker to provide discrimination against the effects of 60-cycle and higher frequency pickup of over 1000. It can be used as a null indicator on several ultraprecision Wheatstone and Mueller bridges without requiring additional shielding of the bridges.

Dimensions: 151/2" wide, 8" high, 10" deep. Weight: 40 lb.



MODEL LB10A makes use of an 80-cycle breaker in combination with a sharp cut-off filter to give a frequency response that is essentially flat to 10 cycles per second. This unit is more sensitive to the effects of 60-cycle pickup.

Dimensions: 22"wide, 101/2" high, 15" deep. Weight: 60 lb.

specifications:

Input characteristics. The input transformer can be designed for any given input impedance from 5 ohms to 100.000 ohms. Zero stability. Drift is less than .005 μ v per 8hour day. Gain stability. With a stable power line, the amplifier will vary less than 0.3% over an 8-hour day. Output characteristics. The unit is linear to approximately 2 volts when terminated into 500 ohms, to 4 volts when used with a meter of infinite impedance. Noise level. Approaches theoretical limitations. Immunity. Discrimination against 60-cycle pickup of induced a-c signals is over 1:1000. Controls. Coarse gain control of 20 positions, 4DB per position. 7 test signals are provided for calibration ranging from 0.1 microvolts to 100 microvolts. Power requirements. 110 volts, 50 or 60 cycles, 20 watts.

Beckman[,]

Scientific and Process

Instruments Division Beckman Instruments, Inc. 2500 Fullerton Road, Fullerton, California CLEVELAND PHILADELPHIA HOUSTON CHICAGO LOS ANGELES NEW YORK PITTSBURGH CINCINNATI

PRINTED IN U.S.A.





the Analyzer ®

January 1965

Volume VI, Number 1

Published quarterly by the

SCIENTIFIC AND PROCESS INSTRUMENTS DIVISION

OF BECKMAN INSTRUMENTS, INC.

Fullerton, California

Jim Notarian, Editor

Alan Chute, Art Director

Contents

- 5 How Johnson's Wax shines in quality control
- 11 Fire-gas analyzer used in training of New York firemen
- 15 Infrared spectrophotometers and resolution
- 18 Index to volumes 1, 2, 3, 4, and 5, The Analyzer
- 22 Beam and stream lines
- 23 Literature

The Analyzer is published quarterly to inform readers about the latest uses and applications of Beckman laboratory and process instrumentation. Information is obtained from actual users of Beckman instrumentation as well as from other sources. While the sources of this information are believed to be reliable, no representation is made by Beckman Instruments, Inc., that the information is accurate and correct, and no warranty, express or implied, in connection with such information and its use (including but not limited to claims for direct or contributory patent infringements) is made. Beckman Instruments, Inc., expressly disclaims any liability for special or consequential damages arising in any way from use of such information.

The staff of *The Analyzer* invites readers to submit material of general interest to instrument users. Articles may be either originals or abstracts of previously printed material. Unless otherwise stated in writing by the submitter, it will be assumed that such articles have not been copyrighted, published, or offered for publication elsewhere and that Beckman Instruments, Inc., is free to publish and disseminate accepted articles without incurring any liability to others, including the submitter. Manuscripts and all correspondence pertaining to *The Analyzer* should be addressed to: *Editor, The Analyzer, Scientific and Process Instruments Division, Beckman Instruments, Inc., California.*

AQUAMETER, BECKMAN, DESICOTE, DB, DK, DU, HYGROMITE, MEGACHROM, SPINCO, THERMOMATIC, THE ANALYZER, THERMOTRAC, VIBRODE, and ZEROMATIC are *registered* trademarks, and MICROZONE is a trademark of Beckman Instruments, Inc.

HIPPOCRATES: 460 B.C. - 377 B.C. Facts about the Father of Medicine, Hippocrates, are surprisingly few. Born on the Greek Island of Cos, in a medical environment, Hippocrates apparently became his father's apprentice at an early age.

The Hippocratic Oath embodying a code of medical ethics, generally taken by recipients of the M.D. Degree, is attributed to Hippocrates. There is frequent speculation that he did not write it. Also, careful analysis of the text has revealed that the oath was written for a medical family guild instead of the medical profession as a whole.

The essence of Hippocratic thinking is contained in the treatise "On Epidemics." Here Hippocrates viewed the medical arts as having three main components: the disease, the patient and the physician.

Hippocrates remains today a symbol of greatness in medicine.





How Johnson's Wax shines in quality control

William K. Miller, PhD: Analytical Supervisor, S. C. Johnson & Son, Inc.

Ever walk across the gleaming floors of a new office building and wonder how they got that way? Chances are, they are polished regularly with a very specially formulated mixture of carnauba wax and other key ingredients. Getting just the right "shine," wearing characteristics, non-slip surface, and other qualities is the continuing task of the Research & Development Division, S. C. Johnson and Son, Racine, Wisconsin.

Making such waxes for commercial and home use is the prime concern of the Johnson company. Founded in 1886 by Samuel Curtis Johnson, a parquet flooring salesman, the firm is now the world's leading manufacturer and marketer of floor waxes, furniture polishes, and a number of other home care products. Such trade names as Glo-Coat*, Klear*, Pledge*, Raid*, and Glade* have become household words. In commercial, industrial and institutional maintenance, Johnson markets such products as Step-Ahead*, Over and Under*, and Forward*.

Its line of more than 200 products, all manufactured in Racine, also includes coolants and lubricants for metalworking; finishing products for furniture manufacturers; agricultural waxes to protect fruits and vegetables, and chemical intermediates for printing inks, paints, and varnishes.

Johnson products are known throughout most of the free world. The Racine company formed its first international subsidiary in England in 1916 and now has 23 manufacturing and marketing associate companies operating in Europe, Africa, Latin America, and the Far East.

Behind this international reputation for high quality products is a unique balance of creative research, modern manufacturing techniques, and marketing skill.

The company's Research and Development Division, headquartered in the world-famous research tower designed by the late Frank Lloyd Wright, has more than 200 employees. In 1929, Johnson's lab consisted of only three people.

R & D activities at Johnson Wax are supported by a wide variety of analytical instruments. Many of these instruments are concentrated in the Analytical Section of the Quality Control Department. In addition to developing methods for quality control, this group performs analytical service for all R & D activity, the Production Division and the Household, Service Products, Chemical, and International Marketing Divisions.

Workhorse of the analytical instruments is a Beckman IR-4 Infrared Spectrophotometer. Johnson has been using infrared spectrophometry for 12 years. The present IR-4, more than five years old, turns out approximately 200 spectra a month, and the number is increasing. The Research Department recently added a Beckman IR-5A to its laboratory to lessen the load on the IR-4 and the Analytical Section.

One function of the Research Department is to study the composition of raw materials, especially waxes. Materials to be studied are separated into fractions by physical, chemical and chromatographic methods. The fractions may be subjected to chemical treatment for further separation. Each fraction is carefully analyzed, and infrared spectrophotometry plays a major role. The composition of carnauba wax and other natural products are much better known as a result of these studies.

Much of the confidence that customers place in Johnson products can be traced to rigid quality control, which extends from the analysis of raw materials, through production, packaging and finally to the finished products. For example, Glo-Coat Self-Polishing Floor Wax undergoes 75 quality control tests. Over and Under, a new floor sealer and undercoater marketed by the Service Products Division for commercial use, undergoes 57 tests on raw materials alone, and 27 additional tests on the finished product.

The Beckman IR-4 plays an important part in this quality control system, supplementing the chemical and physical methods used, to assure that raw materials meet specifications, and that the chemical composition and performance of finished products are consistent.

Adulteration in carnauba wax is detected by infrared spectrophotometry in the analysis of the incoming raw materials. By noting the relative depths of the absorption bands at 1460 and 1740 cm⁻¹, it is possible to detect adulteration of carnauba with paraffin wax. The effect of the paraffin in the carnauba spectrum can be seen in Figure 1-B. A quantitative determination of paraffin in carnauba can also be made by chromatography, but the detection and estimation by infrared is much faster.

Many production problems are also solved by the Beckman IR-4. For example, large particle size, instability, or other undesirable properties of polymer emulsions may be caused by incorrect ratio of monomers, insufficient emulsifier, wrong emulsifier system, or contamination with another material. The cause often can be determined simply by examining the infrared spectrum of the polymer. More often, it is necessary to separate the contaminant or emulsifier from the polymer, but the final step will be infrared examination of the spectra of the separated fractions.

Finished products can be analyzed in the same manner, if other quality control tests indicate the possibility of improper formulation.

Customer complaints are sometimes investigated by infra-



Ŗ

RANK

red spectrophotometry. The gellation of an emulsion floor polish can be caused by the housewife or custodian mixing two partially-used floor waxes that are not compatible. Infrared analysis quickly determines that contamination has occurred, and often identifies the contaminating product. The left part of Figure 2 shows a portion of the spectrum of Traffic Grade, an emulsion floor finish high in wax content, that is marketed by the Service Products Division. On the right is the spectrum of the same finish which has been contaminated by a small amount of a polymer floor finish. Such contamination, of course, would cause the product to perform differently on a floor.

Johnson's Service Products Division, which markets maintenance products for commercial, industrial, and institutional use, has made effective use of the infrared spectrophometer as a marketing tool.

Because labor costs take up to 95% of the commercial building maintenance budget, purchasing agents and maintenance supervisors are becoming convinced that performance is far more important than price, and that any floor finish should be selected on the basis of an on-the-floor test rather than chemical specifications. But, after selecting a product on its performance in a test, the purchaser wants proof that it will perform the same in all subsequent purchases.

The Service Products Division offers this proof with an infrared "fingerprint" of the product tested, plus an affidavit that all subsequent gallonage will have the same "fingerprint" and, hence, the same formulation and same performance.

Johnson's Analytical Section has developed a method for sample preparation of emulsion products that permits the running of an infrared spectrum of the entire non-volatile portion of the product. If the spectrum is identical to that of a previous batch, there can be no significant difference in the chemical composition. The same technique is used to

A. Dr. Miller makes a final check against the operating guide to prepare the Beckman IR-4 Infrared Spectrophotometer for an analysis. Figure 1. To maintain consistent quality in Johnson Wax products, purity checks are run against raw materials and finished goods on the infrared spectrophotometer. Graphs below show how adulteration of carnauba wax with paraffin can be detected.





detect formula changes of competitive products. Major differences indicate formula change, while minor differences may result from poor quality control. If a formula change is indicated, the product is separated into the individual components by solvent extraction techniques, and the fractions are examined by infrared spectrophotometry.

The film of an emulsion floor polish is too brittle to be removed from a surface for infrared analysis. It cannot be applied directly to a sodium chloride plate because the water in the product will cloud the plate. The product cannot be cast from a solvent because no solvent will dissolve all of the non-volatile components of an emulsion floor polish. However, a film can be cast directly on a silver chloride disc by the following technique. A disc of silver chloride, 25 mm in diameter and 1 mm in thickness, is placed in a die which has a cavity of 11/2" in diameter. A pressure of 10,000 to 12,000 psi is applied for two minutes, producing a thinner disc of greater transparency. Two holes are drilled through the disc near the edge and directly across from each other, and a thin steel or copper wire hook is inserted into one of the holes to allow the disc to be suspended by the wire. The disc is dipped into the test emulsion and suspended for two minutes in an oven at 110°C. After the disc is removed from the oven and allowed to cool, the wire hook is inserted in the opposite hole and the dipping and drying process repeated until a film of sufficient thickness is built up. This is indicated by testing the absorption on the spectrophotometer. An average of about six dips is required. The silver chloride disc can be taped over an opening cut in a piece of cardboard for insertion into the cell holder of the instrument.

The silver chloride disc can be re-used after cleaning with hot water and soap. In extreme cases, hot solvents or even scraping with a razor blade may be necessary. When the infrared transparency of the discs become low, they can generally be rejuvenated by immersion in concentrated ammonia for two minutes, rinsing, drying and pressing again in the die.

Polymer emulsion can be cast on the silver chloride discs in the same manner. Fewer dips are required because of the higher non-volatile content of most polymer emulsions. Alternatively, polymers can be cast from a solution in a suitable solvent. Most laboratories evaporate the water, and then dissolve the polymer. Johnson analytical chemists have found that a more efficient method of dissolving the polymer is to add the solvent directly to the emulsion, and heat until the water has been removed by azeotropic distillation. The polymer generally is in solution when all of the water has been removed. Dissolution is accelerated because of the extremely high surface area of the polymer in suspension. Chlorobenzene is an excellent solvent for most polymers. The final preparation of the sample for infrared analysis is accomplished by depositing the polymer solution on a sodium chloride plate, and evaporating the solvent by applying heat from an infrared lamp.

In addition to the IR-4 and IR-5A, Beckman analytical instruments at Johnson Wax include a DK[®]-2 Ultraviolet Spectrophotometer with flame attachment, two DU[®] Ultraviolet Spectrophotometers, a Beckman Hygromite[®] electrolytic hygrometer for determining moisture in aerosol propellents, and more than a dozen pH meters.



Figure 2. Contamination of finished products can be quickly detected by infrared spectrophotometry. The left spectrum is that of a water emulsion floor finish marketed by the Johnson Wax Service Products Division. At right is the same product that has been contaminated by a small amount of an ingredient used in a polymer floor finish. B. Dr. William K. Miller (right), Analytical Supervisor, and Chemist Lee R. Williamson examine an infrared spectrum in the Johnson Wax Quality Control Department.

c. Silver chloride disc is dipped into a water emulsion floor wax, top, to prepare a sample for infrared examination of the entire non-volatile content of a finished product. The disc is then oven-heated to evaporate the water. The dipping and drying process is repeated until a film of sufficient thickness is built up. The disc is then mounted, bottom, on a cardboard and inserted into the IR-4.



Fire-gas analyzer used in training of New York firemen

James T. Ward, Assistant Chief: Division of Training, New York Fire Department

Firemen entering burning buildings are faced with many more dangers than heat, flame, and falling debris. Among these additional hazards are the quick drop in oxygen content of the air and the rapid rise of carbon monoxide. These elements could render ordinary filter masks almost useless and leave the user at the mercy of the flames.

In order to study these hazards, the New York City Fire Department recently dedicated a new Fire-Gas Analysis System. Installed on the apparatus floor of the Department's Welfare Island training tower, the system reads out samplings of gases piped from adjoining fire-training rooms.

Initial findings of the system indicate that oxygen content drops as low as 10 per cent and carbon monoxide concentration rises as high as 2 per cent minutes after the start of an ordinary fire. These statistics, much higher that fire fighters' previous estimates, indicate that the fire service should take a much harder look at protective breathing. A filter mask could be potentially dangerous in a concentration of 2 per cent CO. And what happens to a fire fighter with a self-contained mask when his air supply is exhausted inside a vehicular tunnel?

Answers to such questions received high priority when plans for the recently completed training center were on the drawing boards. Live-fire training rooms were incorporated into the training tower, and doctors stressed that such training should be conducted under rigid control and supervision. Such rigid control called for a fire-gas analysis system. Collected samples in evacuated containers for later analysis were rejected as inadequate for planned operations. It was decided that continuous, industrial-type monitoring equipment would be the minimum acceptable for installation, and the contract was eventually let to Beckman Instruments, Inc.

To eliminate warm up time, the console, which uses ordinary house current, is never de-energized. In its present state of development, the fire-gas analysis unit consists of sampling pipes leading from monitoring stations in two live-fire rooms (the cellar and the first floor) to an instrument console. A thermocouple directly connected to an electrical pyrometer at the console is located at each monitoring station.

Equipment used in the unit is as follows:

1. An electrical pyrometer with a range from 0° to 2,500°F and a selective rotary-switching device to enable the operator to read temperature at any of the fixed stations.

2. A Beckman Model 115 Carbon Monoxide Analyzer (non-dispersive infrared type) with a range of 0 to 10 per cent carbon monoxide. This instrument compares the effect of two infrared beams of energy passing through gaseous atmospheres. One beam passes through a reference cell containing a gas that does not absorb infrared energy. The other beam passes through the sample cell with the gas under test. The difference in output energy is proportional to the amount of carbon monoxide present in the gas under test. The concentration of gas is recorded on a strip chart recorder, and is indicated by the meter on the instrument.

3. A Beckman Model F3 Oxygen Analyzer (magnetic susceptibility type) with a range of 0 to 25 per cent oxygen. This instrument measures the partial pressure of the oxygen present in the gas sample using the known magnetic properties of oxygen gas. Oxygen, when introduced into a magnetic field, is unique in that it is paramagnetic, and reacts in a similar manner to soft iron. Other common gases, with a few exceptions, are diamagnetic. The electrical energy, automatically expended in balancing the instrument under changing conditions, is read out directly on the instrument meter and actuates the pen on the strip-chart recorder.

4. In addition to the direct-reading meters furnished with each analyzer, a single-pen electronic recorder is also used with each.

5. A system of piping and individual control valves to enable the operator to sample at any fixed station.

6. A bank of compressed gas cylinders containing the zero and span gases used in the system.

Fire rooms from which samplings are taken are of Class I construction and are finished with a special concrete that can withstand repeated fires without spalling or cracking. With a coat of paint and some furniture, the rooms closely approximate the same type occupancy and construction found in high-rise buildings.

Both rooms are larger than normal but are suitable for testing the fire fighter and his equipment. They will be more useful when plans for compartmentation are completed. By shifting fire-resistive partitions, different room layouts can be arranged, which will closely simulate actual apartments and office suites.

No scientific conclusions can be drawn from the relatively small number of fire-gas samplings taken to date. But many highly practical findings have developed.

In Class A fires with a fuel load that included a foam rubber mattress (single-bed size), carbon monoxide concentrations reached more than 2 per cent by volume in a matter of minutes, and in one test reached 6 per cent. For the same fire, oxygen content gradually dropped off to below 10 per cent. Readouts from the analyzer indicated that the most



ţ

dangerous atmosphere was found above the fire with, of course, the greatest heat concentration. It is reasonable to assume that the same conditions, and probably worse, would meet fire fighters stretching into a high-rise apartment house.

Tests, made thus far, indicate a need for a prompt and adequate ventilation above the fire by men wearing selfcontained masks.

The gas analysis unit permits live-fire training under the strictest of controls. The training staff, under the direction of Battalion Chief Herbert Whyte and Captain Joseph Carroll, is planning to expand the monitoring capacity of the analyzer to include other fire gases. Eventually, gases generated by almost any type of fire will have their intensities registered by the fire-gas analysis unit.

The New York Fire Department plans to make test results, as well as the conclusions drawn from them, available to the entire fire service.

Portions of this article have been reprinted from *Fire Engineering*, June 1964. All rights reserved.

Copyright, the Reuben H. Donnelley Corporation, New York, New York. Reprinted with permission.

A. Fire-gas analysis system installed in the Welfare Island training tower, New York City Fire Department, includes a Beckman Model 115 carbon monoxide analyzer and a Beckman Model F3 oxygen analyzer. Figure 1. Diagram shows arrangement of vacuum pump and piping which draw samples from fire rooms. Note location of thermocouples, numbered one through eight.



Figure 1



- 11

Infrared spectrophotometers and resolution

Howard J. Sloane: Senior Applications Chemist, Beckman Instruments, Inc., Fullerton, California

Richard Cavenah: Associate Applications Chemist, Beckman Instruments, Inc., Fullerton, California

From time to time, we have been confronted with the question "What is the best resolution of the IR-9?" Sounds like a simple question, but the answer may be rather complex.

Perhaps the first thing that needs to be asked in return is "What do you mean by resolution?" and, second, "How shall we measure it?"

The most obvious aspect of this elusive word "resolution" is that it provides a means for measuring our instrument's capability for showing the separation of two closely-spaced lines or bands. Qualitatively, this idea is quite acceptable, but now we must ask what constitutes separation. Are bands resolved when it is barely possible to tell that there are two bands rather than one, or must our record indicate a degree of separation such that the valley between peaks is a certain given displacement above them? Certain criteria have been suggested for making this judgment of separation but, unfortunately, to date no standard is universally accepted. Moreover, in practice, it is rather impractical to attempt to measure resolution at all infrared wavelengths in this manner since this would require seeking out virtually an infinite number of absorption or emission line pairs of exactly the proper known separation.

To answer our original question, therefore, it would probably be best to resort to the instrument designer's view of resolution and to consider three terms in this regard. These are *theoretical resolution*, *limiting resolution*, and *practical resolution*.

Theoretical resolution is determined entirely from the number, type, and size of dispersing elements, how they are employed, and the wavelength under consideration without regard to the remainder of the real instrument. For example, the theoretical resolution of a single prism instrument at a given wavelength is determined from the size and angle of the prism, the number of passes through it, and the variation of refractive index with wavelength at the chosen wavelength. In practice, this term is useful only to the extent that it describes the ultimate theoretical resolving power of which this dispersing system is capable.

When our dispersing system is placed into a spectrophotometer the term *limiting resolution* becomes useful. The limiting resolution is the best resolution that would have been attainable if our system were not energy-limited. In this case, resolution may be limited by optical aberrations, misalignments and so-called "Rayleigh diffraction" which produces an optical widening of the actual slit image. Limiting resolution for many spectrophotometers in the infrared region, however, is frequently not attainable because of energy limitations. Limiting resolution is a particularly useful concept because, in practical operation, no improvement can be realized by narrowing the slits beyond this limiting width. A contributing factor to limiting resolution, especially in the infrared region, is the Rayleigh diffraction limit. Interestingly, the calculation of this limiting slit width is rather simple for a given spectrophotometer. To a good first approximation, the Rayleigh limit is simply the f-number of the system (a pure number) times the wavelength under consideration. Thus, in the IR-4, -7, -9, -11, and -12 series, the f-number is 10, and at 10 microns, for example, the diffraction limit is $10 \times 10 = 100$ microns, or 0.1 mm slit width. At 3 microns, limiting slits are 0.03 mm, etc. Thus, we now realize that at the wavelengths specified, there would never be any reason to narrow the slit widths to less than the above calculated limiting widths, for no further improvement in resolution could be realized.

We have still not, however, satisfactorily answered our initial question concerning the IR-9's "best resolution." In a real infrared spectrophotometer, we must first acknowledge the fact that we are normally dealing with an energy-limited system employing a noise-limited detector. In consideration of this, we may refer to the term *practical resolution*, or the resolution that can be obtained for a given specified signalto-noise ratio and instrument period (the latter determines the maximum scanning speed which may be used).

Dr. Wilbur Kaye, of Beckman's Research Department, has suggested that practical resolution may be denoted compactly and symbolically as $R_p^{s/n}$. Thus, if one were able to obtain 1 cm⁻¹ resolution at a 2% noise level (signal-to-noise is thus $\frac{100}{2}$ =50) and an 8-second period, one could express practical resolution as $R_8^{50} = 1$ cm⁻¹.

Here then, finally, is a suitable yardstick for measurements that will answer our question. Note that we must answer with qualifications, always specifying the conditions (i.e., s/n and p) under which the resolution is attainable.

The term commonly used to describe resolution in a quantitative fashion is spectral slit width*, which is the width in cm⁻¹ (or microns) measured at half peak intensity of the band of energy emerging from the instrument's exit slit. Technically speaking then, we wish to calculate the instrument's spectral slit width, and this can be done by reference to its dispersion curve.

The dispersion curve for the IR-9, for example, is shown in Figure 1 and it gives the dispersion in cm⁻¹ per mm of slit width for the entire wavenumber range covered by the instrument. These data are calculated from knowledge of the theoretical resolution provided by this particular monochromator and, thus depend on the prisms and gratings used and their manner of use.

Spectral slit width may now be estimated by multiplying the dispersion, J, at any particular wavenumber by the mechanical slit width taken from the instrument's slit width readout dial. Thus, for example, a mechanical slit width of 0.5 mm at 1600 cm⁻¹ where J=3 cm⁻¹/mm will produce a spectral slit width J=3(0.5)=1.5 cm⁻¹.

It should be realized, however, that the sample calculation shown above is really a first order approximation of the spectral slit width in that it does not include the effects of slit widening due to diffraction and slight imperfections in the optics and instrument alignment. At very narrow slits these factors, especially the diffraction term, can contribute appreciably to spectral slit width.

While this calculation does tell us approximately what bandwidth is emerging from the exit slit under a static (i.e., non-scanning) condition, it still falls short of satisfying our practical resolution requirement in that the parameters necessary to record at the calculated spectral slit width have not been specified.

In considering the dynamic or scanning conditions necessary to realize a given resolution, the operator must make a compromise in gain setting and period which determine noise level and permissable scanning speed. For example, the spectral slit width of 1.5 cm^{-1} might be used to record a spectrum with a 1% noise level and at 100 cm⁻¹/min. scanning speed, or at 0.5% noise and 25 cm⁻¹/min. scanning speed; however, the benefit of both the low noise level and fast scanning speed



cannot be realized simultaneously at the specified slit width.

In Figure 2 are plotted experimental data taken on an IR-9. Curve A illustrates the approximate spectral slit width that can be achieved for the particular conditions of a 1% noise level and two-second period; i.e., it is a plot of R_2^{100} as a function of wavenumber and is uncorrected for aberration and diffraction effects which are relatively minor contributors at wide slits.

Similarly, Curve B is a plot of R_{32}^{25} corrected for slit widening due to diffraction and it very nearly represents limiting resolution conditions for the IR-9.

The points for Figure 2 were taken by setting the gain and period switches to give the indicated noise levels and then experimentally determining at 50 or 100 cm⁻¹ intervals the mechanical slit setting necessary to produce the true indicated periods. The latter were measured as the time in seconds for the pen to travel 98% of full chart length in response to a sudden block or unblock signal in the sample beam.

The mechanical slit settings measured were then converted to spectral slit width values by means of the dispersion curve discussed above.

The results of Curve B, Figure 2, then go a long way toward answering the question of "best resolution." Note, however, that this is a qualified answer and corresponds only to the R_{25}^{25} conditions.

It should be pointed out again, however, that there is a limit to the practical resolution achievable under dynamic conditions; i.e., slit width cannot be decreased indefinitely at the expense of noise and/or period. The limiting factor for this type of instrument is, in fact, precisely the limiting resolution discussed earlier. In practice, the practical resolution achievable for an instrument like the IR-9 is normally within a factor of 2 of diffraction limits when operating under the most severe conditions.

That calculated diffraction limits cannot be reached exactly reflects the fact that we are dealing with a real instrument with imperfect optics. Nevertheless, the achievement of even this factor of 2 is a tribute to today's instrument design engineers.

It might be worthwhile now to discuss some of the features of the resolution curves shown in Figure 2.

Note first that under the conditions of noise and period specified, the curves are relatively flat through most of the instrument's range. This points out one of the truly important advantages of a properly-designed grating instrument; i.e., nearly constant resolution as a function of wavelength. This is in sharp contrast to the early simple sodium chloride prism instruments in which dispersion is a rather steep function of wavelength. Typically, such an instrument would yield 4 cm^{-1} spectral slit width at 1000 cm^{-1} and as much as $20-30 \text{ cm}^{-1}$.

A rather subtle point in the design of this instrument involves the choice of the scribing of gratings and their blaze wavelengths in such a manner that no marked change in resolution takes place over the regions of change in gratings or orders (at 671 cm⁻¹, 1200 cm and 2000 cm⁻¹ in the IR-9). This proper choice then avoids large discrepancies in the nature of spectra recorded on both sides of the grating or order changes because of radically different spectral slit widths.



We do note that the curve rises more sharply at the ends of

the instrument range. At the 400 cm⁻¹ end, this rise is due to the absorption of the KBr prism just before its transmission cut-off.

If necessary, considerable improvement of the practical resolution at the 4000 cm⁻¹ end can be accomplished by operating the source at a higher temperature, thus significantly increasing the energy input to the monochromator. This is done at some expense in source life, however, and this "trick" is only effective at high frequencies because of the nature of the blackbody emission curve.

Finally, it is observed that Curve B predicts the region of best achievable practical resolution is the 1200-1400 cm⁻¹ region. This is borne out in practice by the high resolution spectrum of the v_4 band of methane taken at 60 mm pressure in a 10 cm cell for the 1220-1390 cm⁻¹ region (Figure 3).

For this curve, \mathbf{R}_{64}^{25} is estimated at less than 0.3 cm⁻¹ as judged from the spacing of the doublets at 1237, 1260, and 1266 cm⁻¹. Diffraction limits at 1250 cm⁻¹ corresponds to

 0.12 cm^{-1} and, therefore, good agreement with the factor of 2 of achievable practical resolution for the above conditions is observed.

* In some of Beckman's previous literature and manuals, this quantity has been described as "effective band width" and the term spectral slit width was defined as two times this quantity. Most recent literature is in conformance with the above, i.e., half-intensity band width.

A. Application Engineers Howard J. Sloane and Richard Cavenah discuss resolution capabilities as they operate the Beckman IR-9 Infrared Spectrophotometer.

Figure 1. Dispersion curve for IR-9 showing dispersion in cm^{-1}/mm of slit width as a function of wavenumber. Data are calculated from knowledge of dispersing elements and their manner of use. Figure 2. Practical resolution of IR-9 for two given sets of conditions: (A) For 1% noise and 2 second period (R_{19}^{190}), slit widths uncorrected for slit widening due to diffraction. (B) For 4% noise and 32 second period (R_{19}^{150}), corrected for diffraction.

Figure 3. High resolution scan of v_4 band of methane in 1220-1390 cm⁻¹ region. \mathbf{R}_{44}^{ss} is less than 0.3 cm⁻¹ as judged from spacing of doublets at 1237, 1260, and 1266 cm⁻¹.





Volume	Number	Page	Title				
			AOUAMETER	1	4	1	Petrochemical analy
				1	1	8	Preparative gas chro
1	4	13	Jet fuel contamination	3	1	13	Process electrodes
				2	1	5	Sampling systems for
			CARBON DIOXIDE ANALYZERS	2	1	1	Textile technology
2	4	1	O, CO ₂ in diving studies	4	1	11	The trend to process
				5	2	9	Utilizing gas chroma
			CARBONACEOUS ANALYZER				organic chemistry
5	4	11	Use of the total carbon analyzer for pollution control	4	3	19	What is the moon m
			ELECTRODES				GENERAL
2	A	1	0 CO in diving studies	3	2	IBC*	Beckman-Toshiba, l
2	1	12	Process electrodes	1	3	16	Carbon dioxide sens
י ז	1	17	Process electrodes to Ag Li NH ⁴ Ions	2	2	15	Continuous trace ar
1	4	17	Sodium Ion activity measurements	5	2	11	Excellence at Emory
1	4	У	Solium fon activity measurements	2	2	8	Fisk University IR S
			EL COTROBUORESIS	1	1	12	Kruger acquisition
			ELECTROPHORESIS	2	4	7	Offner acquisition
5	3	9	Pioneering new medical instrumentation at	1	1	11	Pittsburgh papers
			Hermann Hospital	5	2	15	R for improvement
				1	1	12	Silica "U" Cells
			FLAME PHOTOMETRY	3	3	16	Space engineering g
1	1	2	Analysis for rorar elements	2	3	16	Spectrograms on m
T	1	3	Analysis for farer elements	3	1	4	USPHS air monitorin
			FLAME SPECTROSCOPY				HYGROMETERS
2	4	3	New horizons	2	2	1	Regenerable solid a
			FLUOROMETER	2	2	12	Transistor manufac
3	3	3	Double-beam ratio fluorometer				INFRARED ANALY
			CAS CURONATOCRABIN	1	2	20	Air sampling unit
			GAS CHROMATOGRAPHY	3	4	1	Environmental con
1	1	12	Abstracting service	5	1	17	Ford uses infrared a
1	1	12	Applications manual	1	2	12	Optical absorption
4	3	8	Canadian Petrofina moves toward optimum production	1	1	5	Petroleum prospect
4	2	1	Computer control-Sohio looks before it leaps	2	2	1	Regenerable solid a
2	4	15	Condensation and evaporation phenomena	5	3	17	USAF solves fuel sys
2	4	7	Conoco's use of megachrom				sac bombers
3	3	16	Detection of helium impurities				
2	3	4	Far UV spectroscopy and gas chromatography				INFRARED REFLE
2	1	14	Flame ionization detector background	•	2	14	Applications of spe
4	2	13	Gas chromatographic determination of blood gases	2	3	14	Applications of spe
2	3	1	Heart research				
3	2	11	Hydrocarbons analysis				INFRARED SPECT
1	2	3,8	Hydrogen flame detector	3	4	11	Analysis for oils in
4	1	1	LPG quality control at Lone Star Gas	4	2	5	Bristol-Myers relie
1	4	14	Model 320-C process gas chromatograph	1	4	16	Christiansen effect
3	4	8	New 520-D gas chromatograph	3	1	7	D₂O in body fluids
2	2	4	Optimized temperature programming	2	3	1	Heart research

4	1	Petrochemical analysis
1	8	Preparative gas chromatography
1	13	Process electrodes
1	5	Sampling systems for liquid injection
1	1	Textile technology
1	11	The trend to process control by gas chromatography
2	9	Utilizing gas chromatography in undergraduate
		organic chemistry
3	19	What is the moon made of?
		GENERAL
2	IBC*	Beckman-Toshiba, Ltd., formation
3	16	Carbon dioxide sensor
2	15	Continuous trace analyses
2	11	Excellence at Emory
2	8	Fisk University IR Spectroscopy Institute
1	12	Kruger acquisition
4	7	Offner acquisition
1	11	Pittsburgh papers
2	15	R for improvement: Short courses
1	12	Silica "U" Cells
3	16	Space engineering group
3	16	Spectrograms on microfilm
1	4	USPHS air monitoring systems
		HYGROMETERS
2	1	Regenerable solid adsorbent system
2	12	Transistor manufacturing
		INFRARED ANALYZERS
2	20	Air sampling unit
4	1	Environmental control systems
1	17	Ford uses infrared analyzers for exhaust studies
2	12	Optical absorption methods
1	5	Petroleum prospecting
2	1	Regenerable solid adsorbent system
3	17	USAF solves fuel system leak detection problem on
		sac bombers
		INFRARED REFLECTANCE
3	14	Applications of specular IR reflectance
		INFRARED SPECTROSCOPY
٨	11	Analysis for oils in solvents
7 7	5	Bristol-Myers relieves a quality control headache
ے م	16	Christiansen effect

1	3	9	High resolution	5	1	13	Experimental studies in gastric secretion
4	4	1	Infrared spectroscopy in medical diagnosis	3	3	6	Film processing
3	4	5	Inorganic studies	1	2	1	Heart-lung machine
2	4	13	Instrumentation in crime laboratory	2	4	13	Instrumentation in crime laboratory
3	2	8	IR-8 Spectrophotometer	4	1	5	The measurement of parts per billion quantities of
3	2	5	IR-9 Spectrophotometer				sodium ions in high purity water
1	2	18	Long wavelength interchange	4	4	4	New research pH Meter raises standards for accuracy
5	1	5	Maintaining quality control standards at	1	1	7	Pharmaceutical quality control
			U.S. Polymeric	2	3	10	Process control of color film
2	1	8	Medicinal values of plants studies	3	2	16	Redox potential standards
4	2	11	MIT-pioneer in infrared-continues its summer session	2	3	12	Simulated stomach
3	3	8	Oxygen purity control	2	1	1	Textile technology
4	4	12	Technology $+$ Instrumentation $=$ quality packaging				
3	2	1	Union Carbide technical center				PYCNOMETER
1	3	4	United States Pharmacopeia XVI	-			
-	-	-		2	1	11	Air comparison pychometer
			LABORATORY INSTRUMENTATION	3	3	14	Cell structure of foams
				2	1	10	Controlling specific gravity of barite
4	3	15	Air Products company-wide uniform quality program	2	4	11	Purgeable air comparison pycnometer
4	4	8	Analysis of packaging materials for	4	1	8	The purgeable air comparison pychometer:
		_	extractable substances				evaluation and application
3	1	9	Archer-Daniels-Midland Laboratory				
3	1	1	Coca-Cola quality control				SAMPLE HANDLING
3	1	11	Ecology studies	3	4	14	Modular sample conditioners
1	3	13	Food additives	5	4	17	Sample handling techniques in the far-infrared
1	3	6	Instrumental analysis in education	5	4	11	Sample handning considers in the fair finited
5	2	17	Instrumentation assures product quality at MC&B				LIT TRAMICRO
5	2	5	New high school program trains chemical technicians				ULIKAMICKO
4	3	11	The Prime Biological Instrument	5	3	5	Ultramicro: Lifesaver at Cook County Hospital
5	3	13	Undergraduate analytical chemistry courses at the				
			University of Colorado				ULTRAVIOLET SPECTROPHOTOMETRY
			OXYGEN ANALYZERS	3	4	16	Analysis of monoaromatics
1		-		2	3	8	Determination of iron in bismuth
1	4	1	Cement kiln oxygen control	1	2	16	Double-beam spectrophotometer
1	3	1	Citrus concentrate quality control	3	2	13	DU & DK fluorescence attachment
3	3	1	Copper ore recovery	1	2	9	DU Cell and cell compartment cover
3	4	1	Environmental control systems	1	4	4	Enzyme chemistry
4	3	3	Head space analysis at Pabst	3	2	3	Far UV absorption studies
1	4	16	Hillary's Himalayan expedition	2	3	4	Far UV spectroscopy and gas chromatography
5	4	19	Oxygen analyzers help retard food spoilage	1	2	10	FBI examinations
1	3	12	Packaged meat	1	2	14	Industrial spectroscopy
2	2	1	Regenerable solid adsorbent system	2	4	13	Instrumentation in crime laboratory
1	1	1	Therapy for premature babies	5	1	9	The lizard wizard
				2	1	8	Medicinal values of plants studied
			ph meters	1	3	IBC*	Nitrogen purging
1	4	16	Antarctica research laboratory	5	4	5	Organic productivity studies in Antarctic waters
3	3	11	Blood gas analysis system	2	2	15	Plotting device for DU spectrophotometer
1	1	10	Buffer solutions	2	3	10	Process control of color film
4	2	8	Chemithon cuts detergent costs	2	2	10	Quality control in aircraft industry
2	1	13	Cosmic ray effects on printing	2	1	1	Textile technology
the Analyzer Title Index Volumes I-V Numbers 1, 2, 3, 4

e	H		the mack volumes I v
lum	mpe	ŝ	
۰ ۲o	ñ	Pag	Title
		10	All the still a service of the shows to prophy
1	1	12	Abstracting service on gas chromatography
2	1	11	Air comparison pychometer
			by Dr. William F. Ulrich
4	3	15	Air Products company-wide uniform quality program
			by Willard L. Ent
3	4	1	AiResearch environmental control systems
1	2	20	Air sampling unit
3	3	14	Analysis of the cell structure of foam materials
			by Robert J. Joyce
4	4	8	Analyses of packaging materials for
			extractable substances
			by E. A. Garlock
2	3	4	The application of far ultraviolet spectroscopy
			to gas chromatography
			by Dr. Wilbur I. Kaye
2	3	14	Applications of specular infrared reflectance
			by Dr. William F. Ulrich
3	1	9	Beckman instrumentation at Archer-Daniels-Midland
			Company
			by W. E. Link, Ph.D.
3	3	16	Beckman space engineering group
3	2	IBC*	Beckman-Toshiba, Ltd.
3	3	11	A blood gas analysis system for the
			measurement of blood oxygen, carbon dioxide
			and pH on a single static sample
		-	by Allan L. Budd
4	2	5	Bristol-Myers relieves a quality control headache
•	•		by Norman Shane
3	2	11	Canadian Chemical uses name ionization
			detector to measure hydrocarbons in air
	2	0	<i>by W.A.Flook</i>
4	3	0	by I. E. Parkar and G. A. Irwin
2	1	13	Case of the cosmic cultrit
1	2	15	Characteristics of the hydrogen flame detector
1	2	0	for gas chromatography
			by Dr. William S. Gallaway and Maurice R. Burnell
1	1	Q	Chemical analysis by prenarative gas chromatography
1	1	0	by Dr William F Ulrich
4	2	1	Chemithon cuts detergent costs
1	4	16	The Christiansen effect in infrared spectroscopy
1	4	13	Compating jet fuel contamination
4	ว	1	Computer control—Sobio looks before it leaps
2	4	15	Condensation and evaporation phenomena in
-	•	10	gas chromatography
			by George S. Turner
2	2	15	Continuous trace analyses
	-		by Thomas J. Kehoe
2	1	10	Controlling specific gravity of barite
			by John Sloan
3	3	6	Data monitoring and display for automating
-			film processing—a system by Houston Fearless
1	2	16	The design and performance of a new double beam
			spectrophotometer
			by Donald G. Tipotsch
3	1	7	D_2O in body fluids studied with Beckman IR-5

1	4	14	Developments in process gas chromatography
	2	11	Excellence of Emory
ט ד	2	11	Excenence at Emory
3	1	15	by David A. Brodie
3	2	3	Far ultraviolet absorption studies on
			biologically active substances
			by Dr. William F. Ulrich
3	4	16	Far ultraviolet analysis of monoaromatics
1	2	10	The FBI laboratory and ultraviolet
-	_		spectrophotometric examinations
			by I Edgar Hoover
2	2	8	Fisk University Infrared Spectroscopy Institute
-	-	Ŭ	by Drs Nelson Fusion Ernest A. Jones, James R. Lawson
1	1	3	Elame analysis for the rarer elements
1	1	5	hy Paul T Gilbert Ir
2	1	14	Flame ionization detector background signal
2	1	14	There a Johns
1	3	1	Elavor control at Minute Maid
2	2	12	A fluorescence attachment for Beckman DI and DK
5	4	15	A nuclescence attachment for Deckman D'o and Dic
			by Milton Leikin and Lerome F. Thomas
5	1	17	Ford uses infrared analyzers for exhaust studies
1	1	12	Gas chromatography applications manual
1	2	13	Gas chromatographic determination of blood gases
7	~	15	by Theron Johns and Brian Thompson
1	4	9	Glass electrodes for the direct measurements of
1	т	,	sodium ion activity in aqueous solutions
			by Dr. John E. Leonard
4	3	3	Head space analysis at Pabst
·		_	by Martin E. Ward and William D. Ziegler
2	3	1	Heart disease cause and cure
1	3	9	High resolution infrared spectroscopy
			by James Stewart
1	4	16	Himalayan expedition
1	3	13	Industry's response to the 1958 food additive amendment
			by William H. Butz
3	4	11	Infrared analysis of residual hydrocarbon oils
			in missile systems wash solvent
			by Howard J. Sloane
1	3	4	Infrared spectrophotometry in U.S.P. XVI
			by Dr. Arthur Osol and Dr. A. R. Gennaro
4	4	1	Infrared spectroscopy in medical diagnosis
			by D. S. Erley and R. D. Stewart
1	3	6	Instrumental analysis at theUniversity of Maryland
5	2	17	Instrumentation assures product quality at MC&B
			by Floyd J. Green
3	3	8	Instrumentation at Southern Oxygen
1	1	12	Kruger acquisition
1	4	7	Lehigh controls cement kiln oxygen
5	1	9	The lizard wizard
1	2	18	The long wavelength interchange for the Beckman IR-7
			by Joseph Ashley
4	1	1	LPG quality control at Lone Star Gas
5	1	5	Maintaining quality control standards at
			U.S. Polymeric

	3	1	11	Martin Aircraft conducts closed ecology studies	2	3	10	Quality color film through process control
	4	1	5	The measurement of parts per billion quantities	3	1	1	Quality control for Coca-Cola
				of sodium ions in high purity water	3	2	16	Redox potential standards
				by Allan L. Budd and Robert H. Jones				by Allan L. Budd
	1	1	5	Methane seep surveyor	3	2	1	Research and development at Union Carbide's
	4	2	11	MIT-pioneer in infrared-continues its summer session				technical center
				by Dr. Richard C. Lord	2	4	17	Response of glass electrodes to silver, lithium,
	1	2	14	Modern industrial spectroscopy at Arizona State University				and ammonium ions
				by Dr. Jacob Fuchs				by Allan L. Budd
	3	1	13	New Beckman process electrodes	2	4	13	Role of instrumentation in Georgia crime laboratory
				by T. L. Kehoe	3	4	14	Sample handling: prerequisite to accurate
	5	2	5	New high school program trains chemical technicians				process analyses
				by Paul E. Malesky	5	4	17	Sample handling techniques in the far-infrared
	4	4	4	New research pH Meter raises standards for accuracy				by Howard J. Sloane and Kenneth E. Stine
				by Allan L. Budd	2	1	5	Sampling systems for liquid injection process
	3	3	3	A new double-beam ratio fluorometer				chromatographs
				by Milton Laikin				by E. A. Houser
	2	4	3	New horizons in flame spectrophotometry	1	1	12	Silica "U" Cells
	1	2	3	A new hydrogen flame detector	2	3	12	Simulated stomach
	-	~	5	by Maurice R Burnell	1	4	4	Some applications of spectrophotometry in
	3	2	8	The new IR-8 infrared spectrophotometer	•	•	•	enzyme chemistry
	5	2	U	by G T Keehl				by Bruce Mackler M D
	2	2	5	The new IR 0 infrared spectrophotometer	2	2	15	Special plotting device for DU spectrophotometer
	5	2	5	by Joseph Ashley and Norman Shiftin	2	2	15	Spectrograms on microfilm
	,		0	A new process and showed correction the Backman 520 D	2	2	10	Spectrophotometric analyses and identifications
	3	4	o	A new process gas chromatograph—The Beckman 520-D	2	2	10	in the circle of the dustry
				by Kaiph A. Dora				h Edward A. Calualli
	1	3	IBC*	Nitrogen purging in the near infrared				by Eawara A. Calvent
	2	4	17	Offner acquisition	1	2	1	Substitute heart
	1	4	16	Operation deep freeze	1	1	7	Synthetic drugs
	1	2	12	Optical absorption methods of continuous analysis	2	2	1	System for space survival
				by Robert L. Chapman	4	4	12	Technology + Instrumentation = quality packaging
	2	2	4	Optimized temperature programming in				by I. Lee Brazier
				gas chromatography	2	1	1	Textile technology at Deering Milliken
				by Maurice R. Burnell and Dr. Abdel S. Said	2	2	12	Transistor production at Hughes
	5	4	5	Organic productivity studies in Antarctic waters	4	1	11	The trend to process control by gas chromatography
				by Sayed Z. El-Sayed				by Robert F. Vogelsang
	5	4	19	Oxygen analyzers help retard food spoilage	2	4	1	The Trieste and inner space
	1	3	12	Oxygen measurement in packaged meats	3	4	5	The use of infrared spectroscopy in inorganic chemistry
	1	1	1	Oxygen therapy for infants				by John R. Ferraro
	5	3	9	Pioneering new medical instrumentation at	2	3	8	The use of spectrophotometry to determine iron
				Hermann Hospital				in Bismuth
	1	1	11	Pittsburgh papers	5	3	5	Ultramicro: Lifesaver at Cook County Hospital
	2	1	8	Potions, poultices, and prescriptions	5	3	13	Undergraduate analytical chemistry courses at the
	1	1	10	Precise pH standardization				University of Colorado
	2	4	7	Preparative gas chromatography at Conoco				by Harold F. Walton
	5	2	15	R for improvement: Short courses	5	3	17	USAF solves fuel system leak detection problem on
	4	3	11	The Prime Biological Instrument		-		SAC bombers
		-		hy Julian Pleasants	5	4	11	Use of the total carbon analyzer for pollution control
	3	3	16	Process gas chromatography applications of current interest		•	• •	by George Kramig and Robert B. Schaffer
	5	2	10	by Gary Nuss	3	1	4	LISPHS sets up network of air monitoring systems
	1	4	1	Product analysis at Dow Chemical	5	2	۰ ۵	Litilizing gas chromatography in undergraduate
	· ~	4	11	Purgeable air comparison pychometer	2	4	,	organic chemistry
-	4	-	11	hy Dohart I. Joyaa				by K N Carter and R B Huff
	٨	1	D	The purpose of comparison	2	2	1	Waste ore recovered by Bagdad Coppor
-	4	1	ō	and purgeable an comparison pychometer:	3 A	ა ი	10	What is the moon made of?
				evaluation and application	4	3	19	what is the moon made of :
				DV Kichard A. Landy and John F. McGahan				oy waiter Donner

ļ

ŀ

Beam and stream lines

Physiological measurements group formed

Beckman Instruments, Inc. has announced the formation of a new unit, the Physiological Measurements Group, to develop physiological instruments and systems for hospitals, clinics and biological research applications.

Joseph W. Lewis, Beckman vice president, said the new group's research director will be Dr. Alan R. Kahn, M.D., who also serves as corporate medical director. Dr. Kahn will coordinate physiological instrumentation programs of all Beckman divisions.

Frank J. Domingues will be projects manager for the group, and will supervise development and engineering programs for physiological instruments.

Lewis said the group will be concerned with the development of new concepts in patient monitoring systems for present and potential medical markets. He indicated that new areas of instrument research will include the development of fetal monitoring systems for use in difficult childbirth cases, and electrocardiogram units which permit heart rate readings under conditions of exercise or stress.

The group will utilize products and technologies of all Beckman divisions, and will be supported by the production and scientific capabilities of the Scientific and Process Instruments Division.

New medical research electrode monitors tiny signals

A new research electrode that can measure electrical current from body tissues in millionths of a volt readings has been announced by Beckman Instruments, Inc.

Dr. Alan R. Kahn, medical director, said the new device is designed for use during surgery and in advanced research projects to monitor and record tiny electrical signals generated by brain, heart and muscle tissue.

Dr. Kahn said the electrode, which provides unusual measurement sensitivity and stability, will enable the researcher to learn more about the vital organs of the body and how they work. The device is expected to be especially useful in evaluating the body's reaction to drugs over extended periods.

The wick-like tip of the inch-long glass electrode carries minute electrical pulses from its point of contact through an electrolyte solution to a silver-silver chloride sensing element, which then transmits the signal to a recorder for readout on chart paper. The device is manufactured by Beckman's Physiological Measurements Group in Fullerton, California.

Trade shows

January 25-28: Eighteenth Annual Symposium on Modern Methods of Analytical Chemistry, Louisiana State University, Department of Chemistry. (Baton Rouge, Louisiana) March 1-5: The 16th Pittsburgh Conference on Analytical Chemistry and Applied Spectroscopy, Penn-Sheraton Hotel. (Pittsburgh, Pennsylvania)

April 6-8: Midwest Chemical Exposition, Cobo Hall. (Detroit, Michigan)

April 10-14: Federation of American Societies for Experimental Biology (FASEB), Convention Hall. (Atlantic City, New Jersey)

Short courses

February 1-4: Third Annual Introductory Course in Gas Chromatography, Roosevelt University. (Chicago, Illinois)

April 19-23: Gas Chromatography Institute, Canisius College. (Buffalo, New York)





Literature

New literature

tory or pilot plant.

during processing.

772 TUBE FITTINGS OF TEFLON*

A STATE AND A STATE

ļ

This new four-page bulletin describes the recently-introduced TRENDOMATIC Series, a new "packaged system approach" to process gas chromatography. Designed for closed-loop control and trend recording, TRENDOMATIC provides a single stream, one or two component capability. The basic system consists

process gas chromatography. Designed for closed-loop control and trend recording, TRENDOMATIC provides a single stream, one or two component capability. The basic system consists of an E Analyzer, wall or panel mount 620 Programmer, and a miniature trend-recorder with Application Engineering. It is available in six system packages that will meet virtually any control requirement in gas chromatography.

This six-page bulletin describes Beckman's Tube Fittings of Teflon for joining all common diameters of glass, plastic and metal tubing. Tubing systems using these fittings are leakproof under pressures to 100 psig and temperatures to 100°C.

Tube Fittings of Teflon are ideal for assembling permanent,

semi-permanent or temporary tubing systems in the labora-

This newly revised four-page bulletin describes the Beckman

Chloride Ion Analyzer, which continuously and automatically

determines the chloride ion concentration in aqueous solu-

tions, over a range from 1.0 to 10,000 ppm. Applications include monitoring for chloride pollution in industrial wastes,

and for salt content of finished and intermediate products

GC-4058 TRENDOMATIC PROCESS CHROMATOGRAPH SYSTEMS

pH-4022-b beckman chloride ion analyzer

7056 beckman planchet counting systems

This illustrated bulletin describes the LOWBETA II and the new WIDEBETA II Radioactivity Counting Systems. Both systems are low level counting instruments designed specifically for measurement of alpha and beta emitting samples. The new WIDEBETA II permits accurate counting at rates from a few counts per hour up to 1,000,000 counts per minute with an effective dead time of less than 0.5 microseconds. Also, a newly developed absorber insert system allows automatic printout of alpha/beta ratio and beta/beta ratio.

Application data sheet

NI-8081 planchet counting of tritium

This four-page application data sheet describes the measurement of tritium with Beckman Planchet Counting Systems. The problems of window and air absorption, sample selfabsorption and windowless counting are described, as well as Beckman's screen window application which makes it possible to efficiently count tritium as a solid.

Technical reprints

reprint 6206

"Test Methods in Spectrophotometry: Stray-Light Determination." Richard E. Poulson; *Applied Optics*; Vol. 3-No. 1; January, 1964.

This paper describes methods and materials used in straylight determination for the purpose of evaluating the performance of UV Spectrophotometers.

reprint 6215

"Use of Infrared Analyzer for Total Carbon Determination" by C. E. Van Hall and V. A. Stenger, *Water and Sewage Works*, June, 1964.

A discussion of the method employed by the Beckman Carbonaceous Analyzer in determining either the total carbon or the total organic carbon content of waste streams. The authors, both with the Special Services Laboratory of the Dow Chemical Company, describe the advantages and the limitations of the instrument and its application in waste treatment plants. Data is presented which indicates a significant correlation of total organic content with the results of both standard BOD and COD tests.

reprint 6216

"Rapid Determination of Percent Reduction of Iron Ore" by E. Price, *Engineering and Mining Journal*, June, 1964.

A comparison of three methods for determining the percent removal of oxygen from iron oxides—the more conventional methods which might be termed gravimetric and metallization, or a new method based on specific gravity. The new method is described in detail by the author, a member of the Research Division of the Allis Chalmers Manufacturing Company. Use of a Beckman Air Comparison Pycnometer helps to make the new procedure faster and easier than conventional methods.

*Teflon is a registered trade mark of Dupont.

	Name			
Please send the literature checked below:I772 (B-1)I8081 (F-13)	Title			
□ 4022-B (F-1) □ 6206 (F-2) □ 4058 (F-4) □ 6215 (F-8)	Company			
□ 7056 (F-13) □ 6216 (F-9)	Division			
	Phone	Extension		
☐ I am not currently on the mailing list	Street	<u> </u>		
Please add my name.	City	State	Zip	
□ Have a sales engineer call.	Company's major products		,,,,,,	
My application is:				
	Vol VL1			
	voi. vi-t			
	Name			
Please send the literature checked below: □ 772 (B-1) □ 8081 (F-13)	Title			
□ 4022-B (F-1) □ 6206 (F-2) □ 4058 (F-4) □ 6215 (F-8)	Сотрапу			
\Box 7056 (F-13) \Box 6216 (F-9)	Division			
	Phone	Extension		
□ I am not currently on the mailing list for your quarterly magazine. THE ANALYZER.	Street			
Please add my name.	City	State	Zip	
	Company's major products			
My application is:				
	Vol. VI-1			
			·····	
Please send the literature checked below:	Name			
□ 772 (B-1) □ 8081 (F-13) □ 4022-B (F-1) □ 6206 (F-2)	Title			
□ 4058 (F-4) □ 6215 (F-8) □ 7056 (F-13) □ 6216 (F-9)	Company			
	Division	T 4 set s		
	rnone			
□ Lam not currently on the mailing list	Street			
I am not currently on the mailing list for your quarterly magazine, THE ANALYZER. Please add my name.	Street City	State	Zip	
 I am not currently on the mailing list for your quarterly magazine, THE ANALYZER. Please add my name. Have a sales engineer call. 	Street City Company's major products	State	Zip	

First Class Permit No. 204 Fullerton California

BUSINESS REPLY MAIL No postage stamp necessary if mailed in the United States

Postage will be paid by

Beckman[®] INSTRUMENTS, INC.

SCIENTIFIC AND PROCESS INSTRUMENTS DIVISION

FULLERTON, CALIFORNIA 92634

	First Class Permit No Fullerton California
BUSINESS REPLY MAIL No postage stamp necessary if mailed in the United States	
Postage will be paid by	
Beckman® INSTRUMENTS, INC.	
SCIENTIFIC AND PROCESS INSTRUMENTS DIVISION	
FULLERTON, CALIFORNIA 92634	

	First Class Permit No. 204 Fullerton California
BUSINESS REPLY MAIL No postage stamp necessary if mailed in the United States	
Postage will be paid by	
Beckman [®] INSTRUMENTS, INC.	
SCIENTIFIC AND PROCESS INSTRUMENTS DIVISION	
FULLERTON, CALIFORNIA 92634	

Deckinali	INSTRUMENTS, INC.		Bulk rate
	SCIENTIFIC AND PROCESS INSTRUMENTS DIVISION FULLERTON, CALIFORNIA • 92634		U. S. Postage PAID Los Angeles
Return requested			Permit no. 14171
		10220	
	HEALTH PHYSICIST	109290	
	GENERAL FLECTRIC X	RAY	
	MILWAUKEF	1 WISC	

Ô

LITHO IN U.S.A.

S65440 165 100GP

 $\hat{\mathbf{\Omega}}$

BROWN RECORDERS

	Model No.	Fange	Inst. No.	Ser. No.
Stack No. 1	153X12V-W7-8K	∋-50 MV	883278	5100-N
Stack No. 2	153X12V-W5-6K	0 -50 M V	883276	5100-N
Hood Foom No. 1	82 <u>9</u> 1	0 -50 MV	936192	5100
Hood Foom No. ?	99 99	0-50 MV	<u>888277</u>	5100-N
Glove Box	KY153X12-VH-II- (V/6)-6-K	0-50 MV	T117981 0001	5100-N
Spare Unit	Y153X12V-W7-(27) A3	0-50 MV	910657	8540

BROWN RECORDERS

	Model No.	Range	Inst. No.	Ser. No.
Stack No. 1	153X12V-W7-6K	0 -5 0 M V	883278	5100-N
Stack No. 2	153X12V-W5-6K	0 -50 M V	883276	5100-N
Hood Room No. 1	11 11	0 -50 MV	986192	5100
Hood Room No. 2	11 11	0-50 M V	883277	5100 - N
Glove Box	KY153X12-VH-II-	0-50 MV	T1179810001	5100-N
Spare Unit	¥153X12V-W7-(27)A3	0-50 MV	910657	6540

J. Jech

Beckman[®]

instructions

ġ.

BULLETIN 338-A BECKMAN MODEL VL LOGARITHMIC MICRO-MICROAMMETER

TABLE OF CONTENTS

SECTION 1.	Pag	je L
	Description	L
SECTION 2.	INSTALLATION	5
	Space Requirements	5 5 5
SECTION 3.	OPERATION)
	Operating Controls9Adjustment and Maintenance Controls9Operating Procedure10Operating Check List11))) 1
SECTION 4.	MAINTENANCE	3
	Equipment13Loop Gain Test13Calibrating Internal Test Currents14Power Supply Test16Checking Chamber Supply Voltage Regulation16Trouble Shooting17	3 3 4 5 5 7
SECTION 5.	THEORY OF THE INSTRUMENT)
	Operation of the Circuit20The Diode Detector22The Feedback Amplifier23The Power Supply24) 2 3 4
SECTION 6.	PARTS LIST	÷

APRIL 1957 PROCESS INSTRUMENTS DIVISION PRINTED IN USA BECKMAN INSTRUMENTS, INC. FULLERTON, CALIFORNIA

LIST OF ILLUSTRATIONS

<u>Fiqure</u>	Title	Page
1.	Characteristics of Model VL Response to Change in Signal Level	2
2.	Characteristics of Reading Errors	3
3.	Mounting Dimensions	6
4.	Controls, Model VL	8
5.	Junction Box for Check of Chamber Supply Voltage	19
6.	Circuit Diagram - Span Adjustment	19
7.	Basic Circuit Diagram	22
8.	Circuit Diagram, Model VL	27

Section 1.

INTRODUCTION

A. DESCRIPTION

The BECKMAN Model VL Logarithmic Micro-Microammeter measures currents in the range of 10-¹³ to 10-⁷ amperes. The entire range of current measurement is presented on a single meter scale for direct reading, with no range-switching requirements. This makes the instrument particularly useful in measuring widely varying signal currents such as those often encountered when working with ion chambers or photocells, or when measuring high-value resistors or insulators over a wide range of applied voltages. Such currents can be measured on a monitoring or continuous basis, with use of a standard recorder unit, but without need of scale-changing mechanisms.

The instrument accommodates any standard 50-millivolt potentiometer-type recorder; six-decade logarithmic chart paper is available for some types of recorders. The instrument also provides a 200-volt regulated voltage supply suitable for use with ion chambers. Full output signal voltage is made available at low impedance for use with auxiliary equipment such as control units.

B. SPECIFICATIONS

1. Power Supply requirements:

a. Electrical - 60 cycle, 100-130 volts, 170 watts.

- 2. Operating conditions:
 - a. Ambient temperature 20 F to 114 F.
 - b. Warm-up time -
 - (1) Constant temperature oven 30 minutes.
 - (2) Amplifier 2 to 4 minutes

3. Instrument characteristics:

a. Indicator - Milliammeter, 4 1/2-inch scale, divided into 6 logarithmic decades.

- b. Detection range 10^{-13} to 10^{-7} amperes.
- c. Drift 5% of decade per 24 hours after equilibration (equivalent to 1% meter displacement).
- d. Accuracy See Figure 2.
- e. Response Time See Figure 1.

-1-

SIGNAL CURRENT LEVEL (AMPERES)



FIGURE 1. Characteristics of Model VL Response to Changes in Signal Level.



INPUT CURRENT

FIGURE 2. CHARACTERISTICS OF READING ERRORS

f. Chamber supply voltage - 200 volts with maximum drift rate of 100 microvolts/second.

4. Output connections:

a. Recorder, 50-millivolt, potentiometer-type.

b. Full output voltage suitable for operation of control equipment.

Section 2.

INSTALLATION

A. SPACE REQUIREMENTS

The rack-mounted Model VL requires a vertical space of 8 3/4 inches in a standard RETMA 19-inch rack. The chassis extends 10 1/4 inches behind the panel; a total of 15 inches is required to accomodate chassis, connectors and cables.

B. TERMINALS

1. Recorder Terminals (S8). These terminals are located on the rear of the chassis. One terminal is connected to chassis ground. The output impedance is 50 ohms.

2. Regulated B+ Terminal (S6). This terminal is located on the rear of the chassis. The potential of the regulated voltage source is about 200 volts; internal resistance is 2 megohms. The voltage supplied through this terminal is used as a polarizing voltage for ion chambers. Since an ion chamber represents a capacitance between the collectrode (connected to the input of the amplifier) and the polarizing electrode (connected to regulated B+), currents can be generated at the input of the amplifier by changes of voltage across the two electrodes. The important criterion of a voltage suitable for ion chamber work is not the dc regulation of the supply, but the rate of change of voltage across the chamber capacitance. The rate of change of voltage of this supply is less than 100 microvolts per second. This will permit operation of chamber units with capacitances of less than 30 micro-microfarads and with induction currents of less than 10-14 amperes.

3. Input Terminal (S1). This is a standard coaxial chassis connector which requires an Amphenol 82-804 cable connector. The choice of input cable is important if the best performance is required. The cable must be well shielded, and the dielectric should be treated to prevent voltage or current generation between the conductor and the shield when the cable is moved, subjected to vibration, or exposed to changes in temperature. Installation should be made so that the cable will not be subjected to motion or rapid temperature change. For operation at low current levels, a cable is available from Amphenol in which the outside of the polyethylene dielectric is coated with colloidal carbon to reduce the generation of voltages and currents. Normally, no earth-to-chassis connection will be required. If, however, with the instrument connected to the current source, the indicating meter needle moves when a part of the system is touched with the hand, a check should be made to ensure adequate shielding of the cable and the current source. If earth grounding is necessary to correct for inadequate shielding, use only a single ground connection - do not connect both the instrument and the current source to ground.

C. MAKING CONNECTIONS

1. Cable for B+ Voltage. To connect the regulated B+ voltage terminal (S6) to an ionization chamber, use RG8/U, RG9/U, or RG10/U cable or equivalent.

-5-



FRONT VIEW



FIGURE 3. Mounting Dimensions

2. Cable for Recorder. To connect the recorder terminals (S8) to a recorder, use a two-conductor cable with at least 150-volt insulation.

3. Input Cable. Use Amphenol 21-467 with maximum length of 240 feet; or use Amphenol 21-537 cable with maximum length of 170 feet.

4. Connectors. The following table indicates the chassis terminal connectors and the mating cable connectors.



FIGURE 4. Controls, Model VL.

Section 3.

OPERATION

A. OPERATING CONTROLS

1. Power Switch. This toggle switch disconnects the instrument from the power line when it is in the OFF position.

NOTE

The switch does not disconnect power from the thermo-regulated oven. This permits maintenance of oven operating temperature while instrument is turned off to keep the high-insulation components clean and moisture free. This reduces warm-up time of instrument.

2. Signal Selector Switch. This switch selects the signal to be connected to the amplifier as follows:

- Position 1 AMPLIFIER BIAS. The input of the amplifier is connected to ground for adjusting d-c amplifier bias. (See Section 3, Paragraph B, 1, Amplifier Bias Adjust).
- Position 2 10^{-8} AMP. An internal standard current of 10^{-8} amperes is connected to the amplifier when adjusting for a correct reading.
- Position 3 10-11 AMP. This is a second internal standard, of 10-11 amperes.

Position 4 - OPERATE. This position is used when a measurement is being made. Selecting this position conects the input to the amplifier.

3. Diode Bias. This control is used in the same manner as the zero control on a conventional amplifier, except that the meter needle is positioned to read a known current instead of zero. Two such test currents are provided. Though only one is necessary to perform zeroing function, two are needed for sensitivity adjustment. (See Section 3, Paragraph B, 2, Span Control).

4. Push to Recover. This push button is used to provide a positive transient at the input of the amplifier to counteract large transients which occur during switching. If a positive transient occurs, the amplifier will recover quickly; but if transient is negative, the amplifier may lose control and recovery will be slow (several minutes under some conditions). The positive transient thus produced by pushing this button compensates for most negative transients. Repeated use of the button will be necessary to compensate for more severe transients.

B. ADJUSTMENT AND MAINTENANCE CONTROLS.

The following controls are used for setting the instrument initially and for maintenance. They are located behind a removable cover on the instrument panel.

1. AMPLIFIER BIAS ADJUST. This control is used to set the bias on the input stage of the d-c amplifier associated with the diode detector. It permits the setting of the input potential within a few millivolts to minimize current leakage across the insulation of the input cable.

This adjustment is delicate, for the amplifier is without feedback and behaves like an amplifier with a gain of about 700. Because of this high amplifier gain, the meter needs only to indicate within three decades of the lower end of the scale.

2. Span Control. This control is essentially the sensitivity control of the instrument. It varies the feedback resistor so that the meter reads correctly for all current values. In an amplifier using a resistor as a current detector, the feedback resistor can be changed in value so that the meter will read full scale for a particular value of current (thus the instrument is made direct reading, and calibration charts are unnecessary). The function of the SPAN control is exactly the same; however, since the logarithmic scale has no zero, a different procedure must be used. Using two known currents, the span control is adjusted until the instrument indicates both currents correctly. The SPAN control interacts to a certain extent with the Diode Bias Control.

3. 10⁻⁸ AMP and 10⁻¹¹ AMP. These controls are used to adjust the test currents. If known currents are applied to the input of the instrument, the instrument can be used as a null device and the internal test currents adjusted to their proper values.

C. OPERATING PROCEDURE

1. Warm-up

a. Set Power Switch to OFF.

b. Connect 115-volt, 60-cycle source to instrument. Wait 30 minutes for thermoregulated oven to reach thermal equilibrium. If the instrument is turned on immediately and operated, the operating procedure must be repeated after the oven has reached equilibrium.

c. Set Signal Selector Switch to 10^{-8} AMP and turn the Power Switch ON. Allow two to four minutes for instrument to warm up.

2. Making Initial Settings

a. Turn DIODE BIAS Control to bring meter to 10^{-8} amperes.

b. Remove cover for sub-panel controls by taking out the two knurled thumbscrews in the 2×7 -inch plate below the indicator.

c. Adjust amplifier bias.

(1) Set Signal Selector Switch to AMPLIFIER BIAS.

(2) Turn AMPLIFIER BIAS ADJUST Control until meter indicates between 10^{-13} and 10^{-10} amperes. When instrument has first been turned on, it may

be difficult to obtain permanent adjustment. Even after the instrument has stabilized, the adjustment will require several seconds and the meter may wander over a range of a half decade.

(3) Readjust the AMPLIFIER BIAS ADJUST again after several hours of instrument operation, and again every few days as experience indicates.

d. Adjust diode bias and span.

(1) Turn Signal Selector Switch to 10-11 AMP.

(2) Turn DIODE BIAS Control to set meter needle at 10^{-13} amperes.

- (3) Turn Signal Selector Switch to 10-8 AMP.
- (4) Turn SPAN Control to set meter needle at 10^{-10} amperes.
- (5) Turn DIODE BIAS Control to set meter needle at 10-8 amperes.

(6) Alternately set Signal Selector Switch at 10^{-11} AMP and 10^{-8} AMP. The meter should read the correct value when either of the two positions is selected.

(7) If there is a discrepancy in the readings it can be corrected by making small changes in the setting of the SPAN Control and subsequent resetting of the DIODE BIAS control to make the 10^{-8} ampere test current read correctly.

3. Operating Instrument

a. Connect the current source to be measured to the Input terminal.

b. Turn Signal Selector Switch to OPERATE.

c. Read measurement of input on meter.

d. If, during input switching operations, transients occur, producing an off-scale reading in the direction of very small current, momentarily depress the PUSH TO RECOVER Control. The meter needle should jump up-scale when the button is released, and then settle at the proper current value. If severe transients occur, more than one operation of the button may be necessary.

If transients occur at times other than when a switching operation is being performed, intermittent connections may be in the current source, or, more likely, shielding may be inadequate. As a final corrective measure, try connecting an earth ground to one part of the equipment.

D. OPERATING CHECK LIST

NOTE

The following procedure for operating the instrument is intended for use only after instrument has been installed, adjusted, and checked out.

1. Turn Power Switch OFF.

2. Connect instrument to 115-volt, 60-cycle power source. Allow 30 minutes for thermoregulated housing to reach constant temperature.

3. Connect current source to INPUT terminal.

4. Turn Signal Selector Switch to 10^{-8} AMP.

5. Turn Power Switch ON. The amplifier warm-up period is three minutes.

6. Adjust DIODE BIAS Control until meter indicates 10-8 amperes.

7. Set Signal Selector Switch to 10^{-11} AMP. Meter should indicate 10^{-11} amperes.

8. Turn Signal Selector Switch to OPERATE to measure input currents.

Section 4.

MAINTENANCE

A. EQUIPMENT

1. Calibrated current source--- delivering 10^{-11} , 10^{-10} , and 10^{-8} amperes. The BECKMAN 1500 Micro-Micro Current Source is recommended.

2. Student's potentiometer.

3. Strip-chart recorder--- with a sensitivity of 50 millivolts and a chart speed on one inch per minute.

4. Oscilloscope--- with vertical sensitivity of 25 millivolts per inch.

5. Voltmeter---with a sensitivity of 20,000 ohms per volt.

6. A small source of radioactivity and an ion chamber--arranged to produce a current of approximately 10-13 amperes.

7. Autotransformer---equipped with a voltmeter such as Variac or Varitran.

B. LOOP GAIN TEST

The loop gain is a measure of the forward gain of the amplifier with the feedback connection broken or otherwise made inoperative. The significance of this test is that it makes possible calculations regarding the effectiveness of the feedback loop in degenerating the input capacitance (cable capacitance, etc.) to reduce the time response of the instrument. Even through the signal current develops a voltage across the diode, the real signal that the amplifier produces is a voltage between the input and ground. The ratio of these two signals is the loop gain of the amplifier.

Since it is very difficult to measure the actual grid-to-ground voltage under signal conditions, use the following method consisting of steps 1, 2, and 3:

1. With the feedback connection inoperative, insert a variable voltage between input and ground as follows:

a. Set the instrument for operation as instructed in Section 3.

b. Leave the Signal Selector Switch in one of the test current positions.

c. Disconnect the yellow lead from the negative meter terminal.

d. Connect the output of a student's potentiometer or similar voltage source between chassis ground and the input terminal, with positive terminal of voltage source to input terminal.

e. Adjust the voltage source to deliver zero voltage.

f. Turn the Signal Selector Switch to OPERATE.

g. Adjust AMPLIFIER BIAS ADJUST until meter needle is off-scale to left.

h. Increase voltage from student's potentiometer until meter indicates 10^{-13} amperes; note the voltage.

i. Increase the voltage until the meter indicates 10^{-7} amperes; again note the voltage.

j. If the instrument drifts appreciably during the measurement, a more accurate measurement can be made by taking readings going upscale, then another pair of readings going downscale. If time consumed in taking each set of readings is about the same, the true reading will be the average of the two values obtained.

2. Measure the voltage necessary to drive the meter from zero to full scale: (It should be 1 to 1 1/2 millivolts.)

a. Disconnect student's potentiometer from input of instrument.

b. Reconnect yellow lead to negative terminal of meter.

c. Connect an accurate voltmeter with full-scale sensitivity of 1 to 3 volts, and an impedance of at least 50,000 ohms from the negative side of the output meter to chassis ground.

d. Connect positive lead from test meter to ground.

e. Set Signal Selector Switch to the 10-11 AMP position.

f. Turn DIODE BIAS Control until meter of VL indicates 10-¹³ amperes. The test meter should read zero.

g. If meter does not, the mechanical zero of the meter does not coincide with the electrical zero. Short the meter and adjust mechanical zero until the meter indicates 10^{-13} amperes. Unshort the meter.

h. With Signal Selector Switch in the 10^{-8} AMP position, turn DIODE BIAS Control until meter indicates 10^{-7} amperes. The test meter now indicates the fullscale output voltage.

3. Calculate the loop gain, using the following expression:

<u>Voltage Output</u> = Loop Gain Voltage Input

The loop gain should be about 700. If it is less than 500, check the voltages through the amplifier against those marked in the circuit diagram. If tube trouble is suspected, the 12AT7 and the 6C4 tubes are most likely to be at fault.

C. CALIBRATING INTERNAL TEST CURRENTS.

1. Using Two Known Currents. The internal test currents can most easily be calibrated by the use of known currents as follows:

a. Connect a known current of 10^{-8} amperes to the input terminal of the instrument.

b. Use the DIODE BIAS control to set the meter to read 10^{-8} amperes.

c. Set Signal Selector Switch to 10^{-8} AMP.

d. Use 10^{-8} AMP Control to set meter to read 10^{-8} amperes. Now the 10^{-8} ampere internal current is calibrated.

e. To calibrate the 10-¹¹ ampere internal current, perform a similar operation using a known current source of 10-¹¹ amperes.

2. Using one Known Current. An alternate method for calibration requires only a single, known-current source. The actual value of the known current is not important other than it must be within the detection limits of the instrument.

To perform this alternate method, you will need a voltage divider with a 100-to-1 ratio and a high-value resistor which is linear with respect to applied voltage over the voltage range of 1 to 100 volts. The voltage divider is used to set the span. Figure 6 shows a suggested circuit. Voltage source should be at least 100 volts.

The voltage divider can be constructed of 1% wire wound resistors which are commercially available. The high resistor value should be such as to produce currents between 10^{-11} and 10^{-8} amperes for best results. A 10^{-10} ohm resistor is indicated.

It can readily be seen that through operation of the switch the current to the instrument will vary by a factor of 100 or two decades. These two currents can be used to set the span of the instrument in the same manner as the internal test currents are used.

With equipment necessary, perform the alternate calibration method as follows:

a. Set the span. (See "Making Initial Settings" in Section 3, Paragraph C).

b. When the span has been set correctly, connect the known current source to the input terminal. (It is suggested that the value of the current be between 10^{-11} and 10^{-8} amperes for best results.)

c. Using the DIODE BIAS Control, set the meter to indicate the exact value of the known current. The instrument is now set to read currents correctly.

d. Set the Signal Selector Switch to 10^{-8} AMP; use the 10^{-8} AMP Control to set the meter to indicate correctly.

e. Set the Signal Selector Switch to 10^{-11} AMP; use the 10^{-11} AMP Control to set the meter to indicate correctly. The internal current sources are not calibrated.

D. POWER SUPPLY TEST

1. Connect an oscilloscope between ground and B+ or B-. Only a small amount of 120 cycle should be observed. The peak-to- peak voltage at these two points should be less than 2 millivolts at the B+ point and less than 75 millivolts at the B- point. No other noise or spikes should be present, except under some conditions when RF is introduced on the a-c supply. Even then, spikes of 25 to 50 millivolts peak-to-peak may not result in poor instrument performance.

2. Note ripple voltages. If they are higher than normal, amplifier loop gain may be low. The 12AT7 or the 35L6 tubes are most likely to be at fault. If the grid voltage is high the 35L6 tubes may have poor emission.

3. If the shunt resistor across the 35L6 tubes changes value materially, the amplifier regulation may be affected. To check the resistor value, turn instrument off, remove a tube (12AT7, for example) and check resistor with an ohmmeter. This is not a complete check, since the resistor trouble might appear only when heated. An open resistor should result in a high bias on the 35L6 tubes.

4. If the shunt resistor is found to be open, replace it, and check the 35L6 tubes on a tube tester. They could have been damaged during overload conditions.

5. If the B+ and B- voltages are noisy, the voltage reference tube (V5) may be noisy. Connect the oscilloscope between B+ and pin 7 of the 5651 tube base. The noise should appear here if the reference tube is at fault.

E. CHECKING CHAMBER SUPPLY VOLTAGE REGULATION

1. Check B+ terminal output. The voltage should be $200 (\pm 10)$ volts. Abnormal reading may indicate a poor 12AX7 tube or a drift in the value of the resistors producing grid bias for the two halves of the tube.

2. To measure regulation of this supply, obtain the following:

a. Polystyrene capacitor of 50 to 100 micro-microfarads with a leakage current of less than 10^{-14} amperes at 200 volts.

b. Signal current of 10^{-13} amperes. This may be obtained from an ion chamber and a small radium source. Use no cables.

c. In practice, it has been found that a 'small junction box (see Figure 5), when connected directly to the instrument INPUT connector, is useful.

3. With capacitor, signal current, and junction box on hand, proceed as follows:

a. Attach junction box to instrument INPUT connector.

b. Attach a small ion chamber to junction box.

c. Polarize chamber with a battery supply returned to instrument chassis ground.

-16-

d. Use small radium source to adjust chamber current to 10^{-13} amperes as read by the instrument.

e. Install the polystyrene capacitor in the junction box, connecting it between INPUT and the insulated pin jack. Close junction box.

f. Connect capacitor pin jack to positive terminal of a 200-volt battery supply.

g. Measure ion chamber current. Change in reading will be due to leakage currents in polystyrene capacitor.

h. Disconnect capacitor connector from 200-volt battery supply and connect it to ground. The instrument should again read the ion chamber current only. Some drift in reading may occur if the capacitor shows any dielectric soak.

i. Allow the soak currents to dissipate before taking a final reading.

j. If the leakage currents from the capacitor are less than 10^{-13} amperes, return capacitor lead to B+ REG connector and set instrument to read ion chamber current.

k. Observe meter for several minutes to note variations in readings. If variations are less than 10^{-14} amperes, the supply is sufficiently well regulated. A recorder would be convenient in making the observations in this step.

F. TROUBLE SHOOTING

1. The manner in which the instrument warms up may be significant in helping to spot trouble. In normal operation the warm-up sequence is as follows:

a. When the instrument is first turned on, only the rectifier tube is supplied with proper voltage; hence it will show red filaments first.

b. The diode activity produces raw B+ voltage across the load capacitor C1. All other tubes in the amplifier, including the series regulator tubes, are supplied with heater current from regulated B+.

c. While the series regulator tubes are cold, the heaters of the tubes are supplied current through the shunt resistor R5.

d. In normal operation, the resistor R5 supplies only half of the regulated current. Moreover, when the heaters of tubes are cold, the resistance of the heaters is lower than during normal operation, so initially the heaters are supplied with only a fraction of their rated current.

e. When the cathodes reach the temperature at which they start to emit, the series regulator tube passes a little plate current which supplies the other tubes with a little more heater current.

f. This process continues for two or three minutes until all the heaters are operating at nominal currents.

g. Normal progress can be followed by observing the meon tube NE2 associated with the plate of tube V6 and the reference tube V5.

h. About one minute after the instrument is turned on, the NE2 will light.

i. After another 3/4 minute, the voltage reference tube will light.

j. In 1/2 minute after the voltage reference tube has lighted, the full voltage will be attained at B+ and B-.

2. If the voltage reference tube lights when the instrument is first turned on, the heater string is broken and can be isolated as follows:

a. Check plug P3 and socket S3 for proper connection.

b. Inspect all tubes for heater continuity.

c. Disconnect P3 from S3 and test socket continuity between pins 5 and 6 with an ohmmeter; a reading of 600 to 800 ohms should be obtained, if resistance is not in this range, the 35L6 tubes are at fault.

d. Connect ohmmeter between socket pins 1 and 4. Reading of 250 to 350 ohms should be obtained. If reading is outside of this range, tube V3, V4, or V10 may be at fault.

e. Connect ohmmeter between pins 1 and 6 of plug P3. Resistance should be 150 to 200 ohms. This step inspects the heaters of tubes V6 and V9.

f. Make checks individually of heaters of the electrometer and diode tubes.

3. If it is necessary to operate the power supply as a separate unit, disconnect plug P3 and insert a 150-ohm, 5-watt resistor across socket pins 1 and 5.



4

ļ





FIGURE 6. Circuit Diagram - Span Adjustment

Section 5.

THEORY OF THE INSTRUMENT

A. OPERATION OF THE CIRCUIT

Figure 7 is the basic circuit diagram which indicates the fundamental circuit configuration and shows the relationship of the controls in the circuit operation. The circuit is essentially a dc amplifier employing 100% negative feedback to stabilize the amplifier gain and to degenerate the input impedance in order that large input capacitances can be tolerated (such capacitances are most often represented by long input cables) without seriously affecting the response time of the device. Conventional micro-microammeters employ high-value input resistors to develop a voltage from the signal current. As a result, practical wide-range instruments use elaborate sensitivity-changing systems to permit satisfactory measurements over a wide range of input current values.

A logarithmic diode is used at the input instead of the high-value resistor. The diode is operated under such conditions that the voltage developed across it is proportional to the logarithm of the current passing through the diode. The diode is operated in this instrument in such a manner that the diode potential changes 0.15 volt for each decade change in the current passing through the diode. Figure 6 shows the diode connected to the input of the amplifier so that the input signal current must pass through the diode. The voltage at the plate of the logarithmic diode detector is impressed on the amplifier. The output voltage is fed back to the cathode of the logarithmic diode detector is comprised of the SPAN resistor and the resistance of the meter. The SPAN resistor is a sensitivity adjustment. The meterscale is divided into six decades; thus, if the increment of voltage change across the diode per decade change of signal current is 0.15 volt, the SPAN control is adjusted to produce 6 x 0.15 = 0.90 volt across the enitre feedback element when the meter reads full scale.

Since a logarithmic element can have no zero, the diode must develop a voltage across it for any particular current level of operation. When the meter current is zero (and indicating 10^{-13} amperes), the diode has a voltage across it corresponding to a current level of 10^{-13} amperes. It is desirable to maintain the input at ground potential to reduce leakage currents across insulators and cable dielectric. It is also advantageous to return the feedback to ground for convenience in the use of auxiliary control and indicating equipment, and to reduce a-c pickup. Both of these objectives may be reached by inserting a potential drop (in the form of Resistor R3) between the feedback point and the cathode of the diode. A variable-current generator comprising the DIODE BIAS control Rl and resistor R2 is connected to R3 in such a manner that the voltage developed across R3 is exactly equal and opposite in polarity to the voltage developed across the diode by a singal current of 10^{-13} amperes. Since the diode does develop a definite voltage drop from a given input signal, zero adjustment of the amplifier itself must be made with the diode removed from the circuit. This is done by grounding the input terminal, which prevents the feedback from having any influence on the input of the amplifier. When the input is grounded, the amplifier is adjusted until the output meter current is zero (needle indicating 10^{-13} amperes). When the input is disconnected from ground, the amplifier maintains the input at ground potential by applying the output signal to the cathode of the detector. When the input is grounded while adjusting the amplifier bias, the

amplifier behaves like a high-gain amplifier with a gain of about 700 and has a full scale sensitivity of about 1.3 millivolts. Consequently, the adjustment is somewhat difficult, but the residual error is small when the input is ungrounded. Suppose that the meter reads 1/3 full scale before ungrounding the input this corresponds to a voltage across the feedback elements of 0.3 volts; then, if the amplifier has a gain of 700, the input will be found at a potential of 0.3/700, or 0.43 millivolt, from ground when the feedback is reconnected.

The use of the screen-voltage adjustment as a bias control, instead of the conventional cathode-type bias control, is possible because of the type of electrometer tube used in this circuit has an extremely low screen m μ ; the use of the control in this position has the further advantage of reducing the current through the slider contact on the control potentiometer to a few microamperes. As this is a feedback amplifier, the input is maintained at constant potential when signal is applied to the input by virtue of the amplifier's ability to do so. In this case the forward gain of the amplifier is about 700, so the change in input potential with respect to ground is equal to the generated output voltage divided by the gain of the amplifier or for full scale meter deflection (0.90 volt), 0.90/700, or 0.0013 volt.

1. Details of the Input Circuit and the Switching. The input circuit is contained within a thermoregulated housing in order to stabilize the temperaturesensitive elements and to provide adequate protection against moisture for the highinsulation components of the circuit. The Signal Selector Switch serves to connect the input of the circuit to ground, to either of two internal current source, or to the input terminal. The switch is designed in such a manner that poles not connected to the input of the circuit are returned to ground. By this means, the internal current generators are maintained under constant conditions and all contacts of the switch are maintained at ground potential, thus materially reducing the leakage currents across the switch insulators, as well as eliminating insulator soak effects.

The internal current standards are produced by impressing suitable voltages across resistors whose values are 3×10^8 and 3×10^{11} ohms. The nominal voltage across these resistors is 3.0 volts; these voltages are made adjustable for exact calibration against external standards.

The transient-control circuit consists of a neon tube connected to the input and returned to ground through a 15-megohm resistor. A 100 micro-microfarad polystyrene capacitor is connected from the junction of these two elements to -35 volts through a 22-megohm resistor. A button on the front panel is a normallyclosed switch which connects the near end of the 22-megohm resistor to B+. Thus, in the non-operating condition, the far end of the neon tube is at ground potential so that current leakage across the base of the tube is kept at a minimum, and the capacitor is discharged through the 15-megohm and 22-megohm resistors and assumes B- voltage (-35 volts). The two resistors act as a voltage divider which maintains the voltage across the neon tube below the firing potential. When the panel control is released, the far end of the capacitor is connected to B+ again. Thus 120 volts appears across the neon tube. Since this is well above the firing potential, the tube fires, transferring a portion of the condenser charge to the input capacitance of the amplifier. This charge tends to make the input of the amplifier positive, thus introducing a transient which is easily handled by the amplifier. A second



FIGURE 7. Basic Circuit Diagram

neon tube serves to illuminate the electrode of the first neon tube to insure proper firing. A one-megohm resistor limits the current to the diode, preventing overloads. The far end of the one megohm resistor is connected to ground through a 100 micromicrofarad polystyrene capacitor. This capacitor is in essence a cable simulator which tends to make the amplifier performance independent of added cable capacitance. The diode is shunted by a 50 micro-microfarad capacitor to reduce noise and to stabilize the circuit against transients.

B. THE DIODE DETECTOR

The didode is a Raytheon CK5704 cathode-type tube with a heater rating of 6.3 volts and 150 milliamperes. In this application, the heater is operated at about 2.0 volts and 80 milliamperes. The important criterion for satisfactory performance is not the heater power but the temperature of the cathode. It can be shown theoretically that the cathode temperature determines the slope of the logarithmic character istic of the device. It has been determined experimentally that the best slope for operation over the range of detection encompassed by this instrument corresponds to a decade increment of 0.15 volt. This means that if the current level is changed by a factor of ten, the voltage across the diode will change by 0.15 volt.

It is important that the cathode be maintained at constant temperature; not only will the decade increment change with cathode temperature changes, but the diode bias for a given current level will change. Indeed, the latter factor is the more important, for a 10% change in the heater current is equivalent to a diode bias change of 40% of a decade. The heater current is obtained from a bleeder across the regulated B voltage supply.

The logarithmic relationship between the current and voltage in the diode detector is due to the fact that the energies of the electrons emitted from the hot cathode surface have Gaussian distribution. If the plate of such a diode is maintained at a fixed voltage negative with respect to the cathode, these electrons emitted from the cathode with energies exceeding the plate-to-cathode voltage will succeed in overcoming the retarding influence of this voltage and reach the plate. Those electrons with less energy will return to the cathode. If the plate-tocathode potential is reduced, more of the electrons emitted from the cathode will have sufficient energy to reach the plate; hence, the plate current will increase. The change in voltage necessary to increase or decrease the plate current tenfold is the decade increment of voltage.

If, instead of controlling the voltage across the diode, we connect the diode to a high-impedance current generator requiring the diode to supply electrons, the diode will adjust the potential between its plate and cathode until the resultant electron flow from cathode to plate meets the demand of the current generator. The voltage across the diode is unique for a given current. If the demand of the current generator changes, then the diode will adjust itself to a new plate-to-cathode potential. Hence, the voltage across the diode is a direct measure of the current flowing through the diode.

If the plate potential were made positive with respect to the cathode, it might be expected that the plate voltage-plate current relation would follow Child's Law. However, this is not the case, for the work function of the cathode surface can be considered as a retarding force on the electrons trying to escape from the cathode surface in just the same manner as a negative potential applied to the plate. This results in the diode behaving in a logarithmic fashion until the plate assumes a positive potential equal to that of the work function of the cathode surface.

The work function of the cathode surface is another variable which will require consideration. It turns out that changes in the work function are responsible for bias drift. The work function is dependent on the material (specifically the surface of the material) and also in part on the history of the treatment of the material. If new diodes are placed in a circuit where the work function of the surface can be measured, the work function can be observed to 'age' in a somewhat predictable manner. During the first week, the change in work function may be several tenths of a volt. (Bear in mind that the full scale voltage sensitivity of the instrument is 0.9 volt.) At the end of three weeks' aging, the rate of change will have leveled off to a value of the order of one to five tenths of a millivolt per day (0.1% of full scale per day). It is important that aged diodes be used to insure low drift characteristics in the amplifier performance.

C. THE FEEDBACK AMPLIFIER

The input tube of the amplifier is a subminiature electrometer tube type CK5886, manufactured by Raytheon. It is operated as a tetrode with a screen voltage of about $4 \frac{1}{2}$ volts and a plate potential of 12 volts. The grid bias is fixed at 3.0 volts and a screen voltage adjustment serves as bias control. In this circuit, the grid current is less than 10-14 amperes. The second stage amplifier is a 12AT7 tube operated as a non-inverting voltage amplifier exhibiting a gain of about 20. Coupling to the output stage is achieved through a neon tube which drops the d-c level of the signal without loss of gain. A 6C4 tube used as the output tube is operated as a cathode follower with a gain of about 0.6. The output signal appears between cathode and ground. The full output signal is applied to the cathode of the diode detector with a polarity opposite that of the signal across the diode. Since the diode impedance is a function of the input signal, the amplifier is required to operate over a wide range of input impedances. In order to stabilize the loop for a broad range of frequencies, a special filter section is used between the 1st and 2nd amplifier stages. The network is designed to produce 45° of phase shift over a wide band of frequencies.

D. THE POWER SUPPLY

The power supply is designed to supply 150 milliamperes of regulated current at about 120 volts. Chassis ground is effected at such a point on the main bleeder that the supply produces a positive 85 volts with respect to the chassis and a negative 35 volts. The high current of the supply is to meet the demands of the diode detector heater, for the heater current must be well regulated to insure good performance. The main bleeder across the supply is comprised of the heaters of all the tubes used in the instrument with the single exception of the rectifier tube. Hence, all the tubes of the instrument are operated under ideal conditions which should result in long service life. Since the power transformer is a regulating transformer, the rectifier tube also is operated under nearly ideal conditions. The power supply regulating amplifier is essentially a three-stage amplifier using a 5651 tube as a voltage reference. The gain of the amplifier, excluding the series regulating tubes, is about 2000. The series regulator is a pair of 35L6 tubes connected in parallel. A shunt resistor across the regulating tubes nominally provides half the current for the supply. The choice of the point on the bleeder for ground provides the best possible regulation for the positive side of the voltage supply; all the critical components are located in this part of the supply. The ripple observed on the B+ voltage is less than 2 millivolts peak-to-peak. The negative supply may have as high as 75 millivolts of ripple.

A second regulated voltage is designed to furnish 200 volts suitable for polarization of ion chambers to be used with the instrument. The amplifier consists of a two-stage amplifier a (12AX7 tube) with a gain of about 2500; regulated B+ of the power amplifier is used as the voltage reference. An output filter on the amplifier further controls the regulation of the output voltage. The internal resistance of the supply is ten megohms. The degree of regulation of the supply can be defined as the rate of change of voltage at the output terminal. Since an ion chamber represents a capacitance coupled from the high-voltage element to the input of the amplifier, and the minimum detectable voltage is known, then the regulation of the polarizing voltage can be determined by application of the relation

 $\frac{I}{C} = \frac{dV}{dT}$

If the chamber is assumed to have an inner-electrode capacitance of 35 micro-microfarads and the maximum induced current to be tolerated is 10^{-14} amperes, then the maximum rate of change of voltage that can be tolerated for the polarizing voltage is 300 microvolts per second. This supply is designed to produce variations of less than 100 microvolts per second. It must be remembered that the dc regulation of supply is not that good.

Section 6.

PARTS LIST

BECKMAN

.

· · · · · · · · · · ·

<u>Symbol</u>

-

Symbol	Description	Part	
R1, 3	51 ohms, ±5%, 1 w, AB GB	7528	
R2, 4	68 ohms, ±5%, 1 w, AB GB	7527	
R5	2,500 ohms, 20 w, Mallory 2HJ2500	7526	
R6, 7	51,000 ohms, ±5%, 1/2 w, AB EB	15067	
R8	36,000 ohms, ±5%, 1/2 w, AB EB	7543	
R9	47,000 ohms, ±5%, 1/2 w, AB EB	15064	
R10, 35	1 megohm, ±5%, 1/2w, AB EB	3180-15	
R11	2 megohms, $\pm 5\%$, $1/2$ w, AB EB	2209	
R12	15,000 ohms, ±5%, 1 w, AB GB	6270 -564- 8	
R13, 21	100,000 ohms, ±5%, 1/2 w, AB EB	15074	
R14	9.1 megohms, $\pm 5\%$, $1/2$ w, AB EB	8856	
R15, 19	5.6 megohms, ±5%, 1/2 w, AB EB	8851	
R16. 37	$0.24 \text{ megohms.} \pm 5\%$. $1/2 \text{ w}$. AB EB	3180-13	
R17	27.000 ohms. ±5%. 1/2 w. AB EB	15085	
R18	8.2 megohms, $\pm 10\%$, $1/2$ w, AB EB	8860	
R20	10 megohms, $\pm 5\%$, $1/2$ w, AB EB	8154	
R22	4.3 megohms, $\pm 5\%$, $1/2$ w, AB EB	8151	
R23	$0.91 \text{ megohm}, \pm 5\%, 1/2 \text{ w}, \text{AB FB}$	8147	
R24	$180.000 \text{ obms. } \pm 5\%. 1/2 \text{ w. AB FB}$	15081	
R25	$50 \text{ obms}, \pm 0.5 \text{ obm}$	70001	
R26	$500 \text{ obms}, \pm 20 \text{ obms})$	/531	
R27	$10 \text{ obms}_{\pm} \pm 1 \text{ obm}$	7550	
R38	$48 \text{ ohms}, \pm 4 \text{ ohms})$	7552	
R28, 39	$22 \text{ megohms, } \pm 5\%$, $1/2 \text{ w}$, AB FB	8852	
R29, 30	$20.000 \text{ obms}_{\pm} \pm 5\% \cdot 1/2 \text{ w}_{\pm} \text{AB FB}$	15057	
R31, 32	$0.33 \text{ mean} = 10\% \cdot 1/2 \text{ w, AB FB}$	8104	
R33	$10.000 \text{ obms.} \pm 5\%. 1/2 \text{ w. AB FB}$	15052	
R34	$24.000 \text{ obms}, \pm 5\%, 1/2 \text{ w}, \text{AB FB}$	7933	
R36	$15 \text{ megohms} = \pm 10\% \cdot 1/2 \text{ w} \text{ AB EB}$	8153	
R40	300 mean = 10% BPC type HBE	5630	
R41	330 000 megohns: ±5% Victoreen Hi-Meg	5623	
R42 43	Potentiometer 5 000 obms Clarostat 43-5000	7520	
R44 45	$100 \text{ obms} \pm 10\% 50 \text{ W} \text{ Obmits}$	08/1	
R46	100 obme + 10% 1/2 w IPC BW - 1/2	9041	
R47	Potentiometer 20 000 obms Helinot 20 000 $A7$	680-4	
R49 R48	100,000,000,000,000,000,000,000,000,000	7628	
R49 50	$2.000 \text{ obm}_{C} \pm 5\%$] w TPC DW=1	7020	
R51	$36\ 000\ \text{ohm}_{2}$ + $5\%\ 1\ \text{w}$ AP CP	6270-556-11	
P50	7500 obm + 5% + 1/2 w AD ED	0270° JJ0 11 5543	
N32 D53	Detentionator 1000 abre Claratet	7500	
n 35 D 5 4	Potentiometer, 1000 ohms, Clarostat	7529	
n 04 D55	Potentioneter, 20,000 onns, Heilpot $\frac{1}{2}$ with rp	7530	
n.J.J D56	12.5 shows = 10%, 1/2 W, AB EB	8098	
D	12.0 Units, $\pm 0/0$, $1/2$ W, IKU BW-1/2 Solootod (to shunt mater W) to 1.00 005	7629	
^N S	Serected (to shunt meter M1 to 1.0±.005 ma)	100/0	
C.7	$20 \times 20 \text{ mrd}$, Mallory FP 234	10063	
υ <u>ζ</u>	$0.0047 \text{ mra}, \pm 20\%, \text{ b} 0.00 \text{ v}, \text{ molded paper},$	15000	
	Sprague /3 P 4/200	15023	
SymbolDescriptionPartC30.047 mfd, $\pm 20\%$, 600 v, molded paper Sprague 73 P 473615029C40.022 mfd, $\pm 20\%$, 600 v, molded paper, Sprague 73 P 2230615026C5, 1520 x 20 mfd, $\pm 20\%$, 600 v, molded paper15021C670.0022 mfd, $\pm 20\%$, 600 v, molded paper, Sprague 73 P 2220615021C81 mfd, 600 v, Tobe Oilmite 0M-6019866C90.068 mfd, $\pm 20\%$, 600 v, molded paper, Sprague 73 P 68306313C100.012 mfd, $\pm 10\%$, 600 v, molded paper, Sprague 73 P 1230613118C110.0025 mfd, $\pm 30\%$ — 10%, 600 v, paper, Sprague 73 P 123066270-509C12500 mmfd, 500 v, mica Sprague 1FM-3510893C13100 mmfd, $\pm 20\%$, 600 v, polystyrene5566C14100 mmfd, $\pm 20\%$, 600 v, mica, Cornell-Dubilier 5R-5118133C14, 18100 mmfd, $\pm 20\%$, 500 v, mica, Cornell-Dubilier 5R-5258132C14100 mmfd, $\pm 20\%$, 500 v, mica, Cornell-Dubilier 5R-5258132C1550 mmfd, 500 v, mica, Cornell-Dubilier 5R-5258132V1, 2Tube, 3516, RCA173V3Tube, 12AT7, RCA15093V4Tube, 12AT7, RCA15093V5Tube, 6C4, RCA15393V11Tube, 5651, RCA1336V2Plug with leads7525V3Tube, 6C4, RCA15093V11Tube, 5674, Raytheon, processed7525V8Tube, 2277, RCA15093S11Transformer (with matched condenser CT)7513 <tr< th=""><th></th><th></th><th>BECKMAN</th></tr<>			BECKMAN
--	----------------	---	------------------
C3 $0.047 \text{ mfd}, \pm 20\%, 600 \text{ v}, molded paper 15029 C4 0.022 \text{ mfd}, 420\%, 600 \text{ v}, molded paper, 15026 C5, 15 20 x 20 mfd, Mallory FP 234 10069 C6, 7 0.022 \text{ mfd}, 4207, 600 \text{ v}, molded paper 15021 Sprague 73 P 2206 15021 9886 C3 1 mfd, 600 v, Tobe Oilmite 0M-601 9886 C3 0.068 \text{ mfd}, \pm 20\%, 600 \text{ v}, molded paper, 313 C10 0.012 \text{ mfd}, \pm 30\%, 600 \text{ v}, molded paper, 313 C11 0.0025 \text{ mfd}, 500 \text{ v}, mica Sprague IFM-35 10893 C12 500 mmfd, 500 v, mica, Cornell-Dubilier 5R-571 8133 C14, 18 100 mmfd, ±20\%, 600 v, polystyrene 5596 C15 500 mmfd, ±20\%, 500 v, mica, Cornell-Dubilier 5R-5025 8132 C14 Tube, 3516, RCA 7523 V1, 2 Tube, 3516, RCA 7524 V3, 6 Tube, 2047, RCA 7523 V4 Tube, 2047, RCA 7523 V5 Tube, 5651, RCA 13361 V10 Tube, 5654, RCA 13361 V10 Tube, 564, RCA 13361 V$	Symbol	Description	<u>Part</u>
Sprague 73 P 473615029C4 $0.022 \text{ mfd}, \pm 20\%, 600 \text{ v}, molded paper, 15026C5, 1520 \times 20 \text{ mfd}, \pm 20\%, 600 \text{ v}, molded paper, 50022 \text{ mfd}, \pm 20\%, 600 \text{ v}, molded paper 50022 \text{ mfd}, \pm 20\%, 600 \text{ v}, molded paper, 50022 \text{ mfd}, \pm 20\%, 600 \text{ v}, molded paper, 5006 \text{ mfd}, \pm 10\%, 600 \text{ v}, molded paper, 5002 \text{ mfd}, \pm 10\%, 600 \text{ v}, molded paper, 5002 \text{ mfd}, \pm 10\%, 600 \text{ v}, molded paper, 5002 \text{ mfd}, \pm 30\%, 600 \text{ v}, molded paper, 5002 \text{ mfd}, \pm 30\%, 500 \text{ v}, mica, Cornell-Dubilier 5R-571 8133C100.025 \text{ mfd}, \pm 30\%, 500 \text{ v}, mica, Cornell-Dubilier 5R-571 8133C110.0025 \text{ mfd}, \pm 20\%, 600 \text{ v}, polystyrene 5596C12500 mmfd, ± 20\%, 600 \text{ v}, polystyrene 5596C13100 mmfd, ± 20\%, 600 \text{ v}, polystyrene 5596C1418C1550 mmfd, ± 20\%, 500 \text{ v}, mica, Cornell-Dubilier 5R-5025 8132C1650 mmfd, ± 50\%, 500 \text{ v}, mica, Cornell-Dubilier 5R-5025 8132V1, 2Tube, 35L6, RCAV3, 6Tube, 12AT7, RCAV4Tube, 26586, RaytheonV7Tube, 5651, RCAV8Tube, CK5704, Raytheon, processedV10Tube, 5684, RaytheonV11Tube, 5686V12Plug with leadsV57549V14Tube, 5680V15Plug, Jenergy Jones P302ABV17Tube, 5680V18Tube, 7593V19Tube, 64, RCAV11Tube, 5680V12Plug with leadsV569S1Socket, 9-pin miniature, Cince 13481V26S21Socket, 9-pin miniature, $	C3	$0.047 \text{ mfd}_{\star} \pm 20\%$, 600 v, molded paper	
C4 $0.022 \text{ mfd}, \pm 20\%, 600 \text{ v}, molded paper, sprague 73 P 2230615026C5, 1520 x 20 mfd, Mallory FP 23410069C6, 70.0022 \text{ mfd}, \pm 20\%, 600 \text{ v}, molded paper15021C81 mfd, 600 v, Tobe Oilmite OM-6019886C90.068 \text{ mfd}, \pm 20\%, 600 \text{ v}, molded paper, sprague 73 P 52306313C100.012 \text{ mfd}, \pm 10\%, 600 \text{ v}, molded paper, sprague 73 P 5230613118C110.0025 \text{ mfd}, \pm 30\% = 10\%, 600 \text{ v}, paper, sprague 73 P 523066270-509C12500 mmfd, ±5%, 500 v, mica Cornell-Dubilier 5R-5118133C13100 mmfd, ±5%, 500 v, mica, Cornell-Dubilier 5R-5188132C14100 mmfd, ±5%, 500 v, mica, Cornell-Dubilier 5R-52258132C1550 mmfd, ±20%, 600 v, polystyrene7560C1650 mmfd, ±20%, 500 v, mica, Cornell-Dubilier 5R-50258132C1950 mmfd, ±20%, 500 v, mica, Cornell-Dubilier 5R-50258132C1950 mmfd, ±20%, 500 v, mica, Cornell-Dubilier 5R-50258132V1, 2Tube, 12A7, RCA7523V4Tube, 12A7, RCA7524V5Tube, CK57C4, Raytheon, processed7525V6Tube, 5257627M1Meter7513T11Transformer (with matched condenser CT)7513P119yith leads7559P119yith leads7559P119yith leads7538P5P119yith leads7538P5P119yith leads7538P5P119yith leads$	00	Sprague 73 P 4736	15029
Solution Sprague 73 P 2206 (structure provided paper) 15026 C5, 15 20 x 20 mfd, Mallory FP 234 10069 C6, 7 0.0022 mfd, 4207, 600 v, molded paper 5prague 73 P 2206 C8 1 mfd, 600 v, Tobe Oilmite OM-601 9886 C9 0.068 mfd, ±20%, 600 v, molded paper, 313 C10 0.012 mfd, ±10%, 600 v, molded paper, 5prague 73 P 12306 13118 C10 0.025 mfd, ±30% — 10%, 600 v, paper, 5prague 2256-A3 1313 C14 0.0025 mfd, ±30% — 10%, 600 v, paper, 5prague 2256-A3 10893 C14 18 100 mmfd, ±20%, 600 v, molded paper 7560 C14 18 100 mmfd, ±20%, 600 v, polystyrene 7560 C17 0.01 mfd, ±20%, 600 v, moles, Cornell-Dubilier 5R-502 8132 9062 C19 50 mmfd, ±20%, 500 v, mica, Cornell-Dubilier 5R-502 8132 9062 C19 50 mmfd, ±20%, 800 v, moles Cornell-Dubilier 5R-502 8132 9062 C19 50 mmfd, ±20%, 500 v, mica, Cornell-Dubilier 5R-502 8132 9062 C19 50 mmfd, ±20%, 800 v, mica, Cornell-Dubilier 5R-502 8132 902070-83734 V10	CA	$0.022 \text{ mfd} \pm 20\%$, 600 v, molded paper.	
C5, 15D0, 20 mfd, Mallory FP 23410069C6, 70.0022 mfd, ±207, 600 v, molded paper15021C70.0022 mfd, ±207, 600 v, molded paper,15021C81 mfd, 600 v, Tobe Oilmite 0M-6019886C90.068 mfd, ±20%, 600 v, molded paper,313C100.012 mfd, ±10%, 600 v, molded paper,13118C110.0025 mfd, ±30% — 10%, 600 v, paper,6270-509C12500 mmfd, 500 v, mica Sprague 1FM-3510893C13100 mmfd, ±50%, 500 v, polystyrene5596C14100 mmfd, ±20%, 600 v, polystyrene7560C170.01 mfd, ±20%, 500 v, mica, Cornell-Dubilier 1DL-3S1 9062C1950 mmfd, ±20%, 600 v, polystyrene7523V1, 2Tube, 1247, RCA7523V3, 6Tube, 1247, RCA7523V4Tube, 12406, RCA13361V5Tube, CK5504, Raytheon, processed7525V9Tube, 1247, RCA13503V11Tube, 5544 CBS7627M1Meter7513T1Transformer (with matched condenser CT)7513P2Plug with leads7559P3Plug, Jones P3022CT9055-0P4Plug, Jones P3022CT9055-0P5Plug, Jones P3022CT7538P5Plug, Jones P3022CT7538P5Plug, Jones P3022CT7539P5Plug, Jones S3022CT7538P5Plug, Jones S3022CT7538P5Socket, Prin miniature, Cince 134817269P5Soc	04	Sprague 73 P 22306	15026
CiteCiteCiteCiteC6, 7 $0.0022 \text{ mfd}, \pm 201, 600 \text{ v}, molded paper15021C81 mdd, 600 v, Tobe Oilmite OM-6019886C90.068 \text{ mfd}, \pm 203, 600 \text{ v}, molded paper,313C100.012 \text{ mfd}, \pm 108, 600 \text{ v}, molded paper,313C110.0025 \text{ mfd}, \pm 308, -10\%, 600 \text{ v}, paper,6270-509C12500 mmfd, 500 v, mica Sprague 1FM-3510893C13100 mmfd, \pm 20\%, 600 \text{ v}, polystyrene5596C14500 mmfd, \pm 20\%, 600 \text{ v}, polystyrene5596C150.01 mfd, \pm 20\%, 500 \text{ v}, mica, Cornell-Dubilier 5R-5T18133C14, 18100 mmfd, \pm 20\%, 500 \text{ v}, mica, Cornell-Dubilier 5B-5Q258132C150.01 mfd, \pm 20\%, 500 \text{ v}, mica, Cornell-Dubilier 5B-5Q258132C1950 mmfd, \pm 5\%, 500 \text{ v}, mica, Cornell-Dubilier 5B-5Q258132C1950 mmfd, \pm 5\%, 500 \text{ v}, mica, Cornell-Dubilier 5B-5Q258132V1, 2Tube, 35L6, RCA7524V3, 6Tube, 12AT7, RCA7523V4Tube, CK3886, Raytheon7553V5Tube, 5551, RCA1733V7Tube, 5651, RCA15033V10Tube, 5AW4 CES7627M1Meter7511Transformer (with matched condenser CT)7513P1P1ug with leads7559P1<9$	C5 15	20×20 mfd Mallory FP 234	10069
Col, 7Col, 0.22 mir, 7.3 P 2220615021Sprague 73 P 220615021C81 mfd, 600 v, Tobe Oilmite 0M-6019886C90.068 mfd, ±20%, 600 v, molded paper,Sprague 73 P 2306313C100.012 mfd, ±10%, 600 v, molded paper,Sprague 73 P 1230613118C110.0025 mfd, 430% — 10%, 600 v, paper,Sprague 2256-AG6270-509C12500 mmfd, 500 v, mica Sprague 1FM-35C14, 18100 mmfd, ±5%, 500 v, mica, Cornel1-Dubilier 5R-5T1C14100 mmfd, ±20%, 600 v, polystyreneC56050 mmfd, ±20%, 600 v, mica, Cornel1-Dubilier 1DL-3S1C1650 mmfd, ±20%, 600 v, mica, Cornel1-Dubilier 1DL-3S1C170.01 mfd, ±20%, 600 v, mica, Cornel1-Dubilier 5R-5Q2C1850 mmfd, ±20%, 600 v, mica, Cornel1-Dubilier 5R-5Q2C1950 mmfd, ±20%, 600 v, mica, Cornel1-Dubilier 5R-5Q2V1, 2Tube, 12AU6, RCAV5Tube, 12AU6, RCAV5Tube, CK5704, Raytheon, processedV7Tube, 6C4, RCAV10Tube, 5AW4 CESV10Tube, 5AW4 CESV11Tube, 5AW4 CESV12Plug with leadsV5Plug, 3-prongV5Socket, 9-pin miniature, Cince 134817269S1Socket, 2-prong Jones S302ABV5Socket, 2-prong Jones S302ABV5Socket, 2-prong Jones S302ABV5Socket, 2-prong Jones S302ABS5Socket, 3-prin miniature, Cinch 56F12865S2Socket, 3-prin Jones S302ABS6 </td <td>C5, 15 C6 7</td> <td>$0.0022 \text{ mfd} \pm 20T - 600 \text{ v}$ molded paper</td> <td></td>	C5, 15 C6 7	$0.0022 \text{ mfd} \pm 20T - 600 \text{ v}$ molded paper	
Sprague 73 P 22209866C81 mfd, 600 v, Tobe Oilmite OM-6019886C90.068 mfd, $\pm 20\%$, 600 v, molded paper,313C100.012 mfd, $\pm 10\%$, 600 v, molded paper,13118C110.0025 mfd, +30%108, 600 v, paper,Sprague 73 P 123066270-509C12500 mmfd, 500 v, mica Sprague 1FM-3510893C13100 mmfd, $\pm 5\%$, 500 v, mica, Cornell-Dubilier 5R-5T18133C14, 18100 mmfd, $\pm 20\%$, 600 v, polystyrene7560C170.01 mfd, $\pm 20\%$, 600 v, mica, Cornell-Dubilier 1DL-3S19062C1950 mmfd, $\pm 20\%$, 500 v, mica, Cornell-Dubilier 5R-5Q258132V1, 2Tube, 35L6, RCA7523V3, 6Tube, 12AUG, RCA7524V5Tube, CK5866, Raytheon, processed7525V4Tube, 12AX7, RCA13361V10Tube, 564, RCA13361V10Tube, 564, RCA13361V11Tube, 564, RCA15093V11Tube, 564, CA7559P10g, Jones P302AB7559P10g, Jones P302AB7559P10g, Jones P302AB7538P5P10g, Jones P302AB7538P10g, 3-prong4536Socket, 9-pin miniature, Cince 134817269Socket, 9-pin miniature, Cinch 56F1286512052S4Socket, 7-ping Jones S302AB5611S5Socket, 7-ping Jones S302AB5611S6Socket, 7-ping Jones S302AB5611S6Socket, 7-ping Jones S302AB5611<	,	5nnaqua 73 P 22206	15021
C3 1 mit, 400 v, 100 0 v, nolded paper, Sprague 73 P 68306 313 C10 0.012 mfd, $\pm 20\%$, 600 v, molded paper, Sprague 73 P 12306 13118 C11 0.0025 mfd, $\pm 10\%$, 600 v, molded paper, Sprague 2256-AG 6270-509 C12 500 mnfd, 500 v, mica Sprague 1FM-35 10893 C13 100 mnfd, $\pm 20\%$, 600 v, polystyrene 5596 C14 100 mnfd, $\pm 20\%$, 600 v, polystyrene 5596 C15 50 mnfd, $\pm 20\%$, 600 v, mica, Cornell-Dubilier 5R-5T1 8133 C14 100 mnfd, $\pm 20\%$, 600 v, mica, Cornell-Dubilier 5R-5Q25 8132 C17 0.01 mfd, $\pm 20\%$, 500 v, mica, Cornell-Dubilier 5R-5Q25 8132 V1, 2 Tube, 12AUG, RCA 7523 V3, 6 Tube, 12AUG, RCA 7524 V5 Tube, CKS704, Raytheon, processed 7525 V9 Tube, 5AW4 CBS 7627 V10 Tube, 5AW4 CBS 7559 P1ug with leads 7559 P1ug with leads 7559 P1ug, 2-prong, Jones P302AB 7538 P1ug, 3-prong 4536 Socket, 9-pin miniature, Cince 13461 7269 Socket, 2-prong Jones S302AB 5614 <td>C0</td> <td>3 prague 13 P 22200</td> <td>9886</td>	C 0	3 prague 13 P 22200	9886
Cy Diversify (1, 200, 000 mills) (2, 00, 000 mills) (2, 00, 000 mills) (2, 0000 mills) (2, 0000 mills) (2, 000 mills) (2, 000 mills		1 (mid, 600 V, 1000 011 mide of 601 001)	,000
Sprague 75 pros00Single 73 product 75	69	$0.008 \text{ mid}, \pm 200, 000 \text{ v}, \text{ morded paper},$	313
CloseDistrict Titler, 1002, 000, molder paper, Sprague 73 P 1230613118Cl10.0025 mfd, +30% — 10%, 600 v, paper, Sprague 2256-AG6270-509Cl2500 mmfd, 500 v, mica Sprague 1FM-3510893Cl3100 mmfd, ±20%, 600 v, polystyrene5596Cl650 mmfd, ±20%, 600 v, polystyrene7560Cl70.01 mfd, ±20%, 600 v, mica, Cornell-Dubilier 5R-5118133Cl4, 18100 mmfd, ±20%, 600 v, mica, Cornell-Dubilier 5R-5228132Cl950 mmfd, ±20%, 600 v, mica, Cornell-Dubilier 5R-5228132V1, 2Tube, 3516, RCA6270-83V3, 6Tube, 12AT7, RCA7523V4Tube, 12AT7, RCA7524V5Tube, 5651, RCA7555V9Tube, 65886, Raytheon7875-0V8Tube, 654, RCA13361V10Tube, 58704, Raytheon, processed7555V9Tube, 624, RCA15093V11Tube, 5844 CBS7669P3Plug with leads7558P4Plug, 2-prong, Jones P302AB7538P5Plug, Jones P302CCT9055-0P9Plug, 3-prong4536S5Socket, 9-pin miniature, Cinch 56F1286512052S6Socket, 4, mphenol 82-8055614S5Socket, 2-prong, Jones S302AB5611S6Socket, 4, mphenol 80C9093S7Socket, 4, mphenol 80C9055-0S8Socket, 4, mphenol 80C9055S8Socket, 4, mphenol 80C9055S8Socket, 4,	C 10	Sprague $73 P 00000$	010
Sprague / 3 P 12306 15110 C11 0.0025 mfd, ±30%	C10	$0.012 \text{ mra}, \pm 10\%, 600 \text{ v}, \text{ molded paper}, 10006$	13118
C110.0025 mrd, +30k 10k, 600 V, paper, Sprague 2256-AG6270-509C12500 mmfd, 500 v, mica Sprague 1FM-3510893C13100 mmfd, ±5%, 500 v, mica, Cornell-Dubilier 5R-5118133C14, 18100 mmfd, ±20%, 600 v, polystyrene5596C1650 mmfd, ±20%, 600 v, polystyrene7560C170.01 mfd, ±20%, 300 v, mica, Cornell-Dubilier 5R-50258132V1, 2Tube, 35L6, RCA6270-8V3, 6Tube, 12AT7, RCA7523V4Tube, 12AT7, RCA7524V5Tube, 5651, RCA737V7Tube, CK5866, Raytheon, processed7525V9Tube, CK5704, Raytheon, processed7525V1Tube, 5AW4 CBS7627M1Meter7511P2Plug with leads7559P4Plug, 2-prong, Jones P302AB7538P5Plug, 3-prong4536Socket, 9-pin miniature, Cince 134817269S3Socket, 7-pin miniature, Cince 134817269S4Socket, 7-pin miniature, Cince 134817269S5Socket, 7-pin miniature, Cince 134817269S6Socket, Amphenol 80-C2F4596S6Socket, Amphenol 80-PC2F4596S8Socket, Amphenol 80-PC2F4596S8Socket, Amphenol 80-PC2F4596S8Socket, Amphenol 80-PC2F4596S8Socket, Amphenol 80-PC2F4596S8Socket, Amphenol 80-PC2F4596S8Socket, Amphenol 80-PC2F4596 <t< td=""><td></td><td>Sprague 73 P 12300 $(100\% - 100\% - 600)$</td><td>15110</td></t<>		Sprague 73 P 12300 $(100\% - 100\% - 600)$	15110
Sprague 2256-AG0270-0093C12500 mmfd, 500 v, mica Sprague 1FM-3510893C13100 mmfd, $\pm 20\%$, 600 v, polystyrene5596C1650 mmfd, $\pm 20\%$, 600 v, polystyrene7560C170.01 mfd, $\pm 20\%$, 600 v, polystyrene7560C1950 mmfd, $\pm 20\%$, 600 v, polystyrene7560C1950 mmfd, $\pm 20\%$, 300 v, mica, Cornell-Dubilier 1DL-3S19062C1950 mmfd, $\pm 20\%$, 500 v, mica, Cornell-Dubilier 5R-50258132V1, 2Tube, 3516, RCA7523V3, 6Tube, 12AT7, RCA7523V4Tube, 5651, RCA173V7Tube, 6551, RCA173V7Tube, 5651, RCA13361V10Tube, 5644, RCA13361V10Tube, 54W4 CBS7627MlMeter7511T1Transformer (with matched condenser CT)7513P2Plug with leads7559P4Plug, 2-prong, Jones P302AB7538P5Plug, Jones P302CT9055-0P9Plug, 3-prong4536Socket, 3-ppin miniature, Cince 134817269S3Socket, 7-pin miniature, Cince 134817269S4Socket, 3-ppin s08 S302AB5611S6Socket, 3-pping Jones S302AB5611S6Socket, 4-mphenol 82-8055614S2Socket, 5-prong, AN 3102A-14S-555629S8Socket, 5-prong, AN 3102A-14S-555629S8Socket, 5-prong, AN 3102A-14S-555629S8	C11	0.0025 mfd, +30% - 10%, 600 V, paper,	6070-500
C12500 mmrd, 500 v, mica Sprague 1rW3310093C13100 mmrd, ±5%, 500 v, mica, Cornell-Dubilier 5R-5118133C14, 18100 mmrd, ±20%, 600 v, polystyrene5596C1650 mmrd, ±20%, 600 v, mica, Cornell-Dubilier 1DL-3519062C170.01 mrd, ±20%, 500 v, mica, Cornell-Dubilier 5R-50258132V1, 2Tube, 3516, RCA6270-8V3, 6Tube, 12AT7, RCA7523V4Tube, 12AU6, RCA7524V5Tube, 5651, RCA173V7Tube, CK5704, Raytheon, processed7525V8Tube, CK5704, Raytheon, processed7525V9Tube, 5AW4 CES7627M1Meter7511T1Transformer (with matched condenser CT)7513P2Plug with leads7559P3Plug with leads7559P4Plug, 3-prong4536S1Socket, Amphenol 82-8055614S2Socket, 7-pin miniature, Cince 134817269S3Socket, Jones S302AB5611S6Socket, Amphenol 80C9093S7Socket, Amphenol 80C9093S7Socket, Amphenol 80C9093S7Socket, Amphenol 80C9093S8Socket, Amphenol 80C9093S7Socket, Amphenol 80C9093S7Socket, Amphenol 80C9093S7Socket, Amphenol 80C9093S7Socket, Amphenol 80C9093S8Socket, Amphenol 80C2313SW1Switch,		Sprague 2255-AG	10902
Cl3100 mmfd, $\pm 20\%$, 500 v, mica, Cornell-Dubilier 5R-5115133Cl4, 18100 mmfd, $\pm 20\%$, 600 v, polystyrene5596Cl70.01 mfd, $\pm 20\%$, 600 v, mica, Cornell-Dubilier IDL-3SI 9062Cl950 mmfd, $\pm 50\%$, 500 v, mica, Cornell-Dubilier 5R-502581328132V1, 2Tube, 35L6, RCA6270-8V3, 6Iube, 12AT7, RCA7523V4Tube, 12AU6, RCA7524V5Tube, 5651, RCA173V7Tube, CK5886, Raytheon7875-0V8Tube, CK5704, Raytheon, processed7525V9Tube, 5641, RCA13361V10Tube, 5484 CES7627M1Meter7511T1Transformer (with matched condenser CT)7513P2Plug with leads7559P4Plug, 2-prong, Jones P302AE7538P5Plug, Jones P302CCT9055-0P9Plug, 3-prong4536Socket, 9-pin miniature, Cince 134817269S3Socket, 9-pin miniature, Cince 134817269S3Socket, 2-prong Jones S302AB5611S6Socket, 2-prong Jones S302AB5611S6Socket, 3000 CCT7539S5Socket, 4mphenol 80C9093S7Socket, 59rong, AN 3102A-145-555629S6Socket, 59rong, AN 3102A-145-555629S8Socket, 59rong, AN 3102A-145-555629S8Socket, 59rong, AN 3102A-145-555629S8Socket, 59rong, AN 3102A-145-555629S	C12	500 mmfd, 500 v, mica Sprague IFM-35	10093 TI 0100
Cl4, 18 100 mmfd, $\pm 20\%$, 600 v, polystyrene 396 Cl6 50 mmfd, $\pm 20\%$, 600 v, polystyrene 7560 Cl7 0.01 mfd, $\pm 20\%$, 300 v, mica, Cornell-Dubilier IDL-3S1 9062 Cl9 50 mmfd, $\pm 5\%$, 500 v, mica, Cornell-Dubilier IDL-3S1 9062 V1, 2 Tube, 3516, RCA 6270-8 V3, 6 Tube, 12AT7, RCA 7523 V4 Tube, 12AT7, RCA 7523 V5 Tube, 5651, RCA 173 V7 Tube, CKS704, Raytheon, processed 7525 V9 Tube, CK5704, Raytheon, processed 7525 V10 Tube, 5AW4 CBS 7627 M1 Meter 7511 T1 Transformer (with matched condenser CT) 7513 V1 Tube, 5AW4 CBS 7559 P1ug with leads 7556 P1ug with leads 7553 S1 Socket, Amphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, 7-pin miniature, Cinch 56F1	C13	100 mmfd, $\pm 5\%$, 500 v, mica, Cornell-Dubilier 5k-5	11 0100 5506
C16 50 mmfd, $\pm 20\%$, 600 v, polystyrene 7500 C17 0.01 mfd, $\pm 20\%$, 300 v, mica, Cornell-Dubilier IDL-3S1 9062 C19 50 mmfd, $\pm 20\%$, 500 v, mica, Cornell-Dubilier 5R-5Q25 8132 V1, 2 Tube, 35L6, RCA 6270-8 V3, 6 Tube, 12AU7, RCA 7523 V4 Tube, 12AU6, RCA 7524 V5 Tube, 5651, RCA 173 V7 Tube, 6K5866, Raytheon 7875-0 V8 Tube, CK5886, Raytheon, processed 7525 V9 Tube, 6C4, RCA 13361 V10 Tube, 12AX7, RCA 15093 V11 Tube, 5844 CBS 7627 M1 Meter 7511 T1 Transformer (with matched condenser CT) 7513 P2 Plug with leads 7559 P3 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P6 Plug, Jones P302CCT 9055-0 P5 Plug, Jones S302CAB 5614 S2 Socket, Amphenol 82-805 5614 S2 Socket, 5-prong, An 3102A-14S-5S 5629	C14, 18	$100 \text{ mmfd}, \pm 20\%, 600 \text{ v}, \text{ polystyrene}$	3590
C17 0.01 mfd, ±20%, 300 v, mica, Cornell-Dubilier IDL-351 9062 C19 50 mmfd, ±5%, 500 v, mica, Cornell-Dubilier SR-5Q25 8132 V1, 2 Tube, 35L6, RCA 6270-8 V3, 6 Tube, 12AT7, RCA 7523 V4 Tube, 12AU6, RCA 7524 V5 Tube, 12AU6, RCA 7524 V5 Tube, CK5866, Raytheon 7875-0 V8 Tube, CK5704, Raytheon, processed 7525 V9 Tube, 6C4, RCA 13361 V10 Tube, 5AW4 CBS 7627 M1 Meter 7513 T1 Transformer (with matched condenser CT) 7513 P2 Plug with leads 7559 P3 Plug with leads 7555 P4 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P9 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 5-prong, AN 3102A-14S-5S 5629	C16	50 mmfd, $\pm 20\%$, 600 v, polystyrene	7000
C19 50 mmfd, ±5%, 500 v, mica, Cornell-Dubilier SH-5Q25 5132 V1, 2 Tube, 35L6, RCA 6270-8 V3, 6 Tube, 12AT7, RCA 7523 V4 Tube, 12AT7, RCA 7524 V5 Tube, 5651, RCA 173 V7 Tube, CK5866, Raytheon 7875-0 V8 Tube, CK5704, Raytheon, processed 7525 V9 Tube, 6C4, RCA 13361 V10 Tube, 5AW4 CBS 7627 M1 Meter 7511 T1 Transformer (with matched condenser CT) 7513 P2 Plug with leads 7559 P4 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CT 9055-0 P9 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 9-pin M 3102A-14S-5S 5614 S4 Socket, 5-prong, AN 3102A-14S-5S 562	C17	0.01 mfd, ±20%, 300 v, mica, Cornell-Dubilier IDL	-351 9002
V1, 2 Tube, 35L6, RCA 6270-8 V3, 6 Tube, 12AT7, RCA 7523 V4 Tube, 12AUG, RCA 7524 V5 Tube, 5651, RCA 173 V7 Tube, CK5704, Raytheon 7875-0 V8 Tube, CK5704, Raytheon, processed 7525 V9 Tube, CK5704, Raytheon, processed 7525 V10 Tube, CK5704, Raytheon, processed 7525 V9 Tube, CK5704, Raytheon, processed 7525 V10 Tube, 5247, RCA 13361 V10 Tube, 5247, RCA 15093 V11 Tube, 5444 CBS 7627 M1 Meter 7513 P2 Plug with leads 7559 P3 Plug with leads 7559 P4 Plug, Jones P302AB 7538 P5 Plug, Jones P302CT 9055-0 P9 Plug, 3-prong 4536 Socket, 9-pin miniatur	C19	50 mmfd, ±5%, 500 v, mica, Cornell-Dubilier 5H-5Q	25 8132
V3, 6 Tube, 12AT7, RCA 7523 V4 Tube, 12AU6, RCA 7524 V5 Tube, 5651, RCA 173 V7 Tube, CK50866, Raytheon 7875-0 V8 Tube, CK5704, Raytheon, processed 7525 V9 Tube, 6C4, RCA 13361 V10 Tube, 5AW4 CBS 7627 M1 Meter 7511 T1 Transformer (with matched condenser CT) 7513 P2 Plug with leads 7559 P3 Plug with leads 7559 P4 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P9 Plug, Jones P302CCT 9055 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 2-prong Jones S302AB 5614 S2 Socket, 5-prong, AN 3102A-14S-5S 5629 S4 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80-PC2F 4596 T81 Thermoregulator, 58°C, normally closed, Edison Type S1-1A <td< td=""><td>V1, 2</td><td>Tube, 35L6, RCA</td><td>0270-8</td></td<>	V1, 2	Tube, 35L6, RCA	0270-8
V4 Tube, 12AU6, RCA /524 V5 Tube, 5651, RCA 173 V7 Tube, CK5866, Raytheon 7875-0 V8 Tube, CK5704, Raytheon, processed 7525 V9 Tube, 6C4, RCA 13361 V10 Tube, 12AX7, RCA 15093 V11 Tube, 5AW4 CBS 7627 M1 Meter 7511 T1 Transformer (with matched condenser CT) 7513 P2 Plug with leads 7559 P3 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P9 Plug, 3-prong 4336 S1 Socket, Amphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, 2-prong Jones S302AB 5611 S6 Socket, 3-prong AN 3102A-14S-5S 5629 S8 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80-PC2F 4536 T81 Thermoregulator, 58°C, normally closed, Edison Type	V3, 6	Tube, 12AT7, RCA	7523
V5 Tube, 5651, RCA 173 V7 Tube, CK5386, Raytheon 7875-0 V8 Tube, CK5704, Raytheon, processed 7525 V9 Tube, 6C4, RCA 13361 V10 Tube, 12AX7, RCA 15093 V11 Tube, 5AW4 CBS 7627 M1 Meter 7511 T1 Transformer (with matched condenser CT) 7513 P2 Plug with leads 7559 P3 Plug, 2-prong, Jones P302AB 7538 P5 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, 7-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, 2-prong Jones S302AB 5611 S6 Socket, 3-prong Jones S302AB 5611 S5 Socket, 7-pin miniature, Cinch 13481 7269 S3 Socket, 2-prong Jones S302AB 5611 S6 Socket, 3-prong AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80C 9093 S7 Socket, Amphenol 80-PC2F 4586 T81 Thermoregulator, 58°C, normally closed, 7515 SW1 Switch, notary, Oak 7515 SW2	V4	Tube, 12AU6, RCA	7524
V7 Tube, CK5886, Raytheon 7875-0 V8 Tube, CK5704, Raytheon, processed 7525 V9 Tube, CK5704, Raytheon, processed 13361 V10 Tube, 12AX7, RCA 15093 V11 Tube, 5AW4 CBS 7627 M1 Meter 7511 T1 Transformer (with matched condenser CT) 7513 P2 Plug with leads 7569 P3 Plug, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P6 Plug, 3-prong 4536 S1 Socket, Maphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, 2-prong Jones S302AB 5611 S6 Socket, Amphenol 80C 9093 S7 Socket, Amphenol 80C-PC2F 4586 T81 Thermoregulator, 58°C, normall	V5	Tube, 5651, RCA	173
V8 Tube, CK5704, Raytheon, processed 7525 V9 Tube, 6C4, RCA 13361 V10 Tube, 12AX7, RCA 15093 V11 Tube, 5AW4 CBS 7627 M1 Meter 7511 T1 Transformer (with matched condenser CT) 7513 P2 Plug with leads 7559 P3 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P9 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cince 13481 7269 S3 Socket, 2-prong Jones S302AB 5611 S6 Socket, 2-prong Jones S302AB 5611 S6 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80C 9093 S7 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80C 9093 S7 Socket, 5-prong, AN 3102A-14S-5S 5	V7	Tube, CK5886, Raytheon	7875-0
V9 Tube, 6C4, RCA 13361 V10 Tube, 12AX7, RCA 15093 V11 Tube, 5AW4 CBS 7627 M1 Meter 7511 Transformer (with matched condenser CT) 7513 P2 Plug with leads 7559 P3 Plug with leads 7555 P4 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P9 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, 7-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, Jones S302CCT 7539 S5 Socket, 2-prong Jones S302AB 5611 S6 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80C 9093 S7 Socket, Amphenol 80-PC2F 4586 TR1 Thermoregulator, 58°C, normally closed, Edison Type S1-1A 1036 SW1 Switch, rotary, Oak 7515 5629 SW2 Switch, DPST, Arrow H&H 80600	V8	Tube, CK5704, Raytheon, processed	7525
V10Tube, 12AX7, RCA15093V11Tube, 5AW4 CBS7627M1Meter7511T1Transformer (with matched condenser CT)7513P2Plug with leads7559P3Plug with leads7555P4Plug, 2-prong, Jones P302AB7538P5Plug, Jones P302CCT9055-0P9Plug, 3-prong4536S1Socket, Amphenol 82-8055614S2Socket, 7-pin miniature, Cince 134817269S3Socket, 7-pin miniature, Cinch 56F1286512052S4Socket, Jones S302CCT7539S5Socket, 2-prong Jones S302AB5611S6Socket, Amphenol 80C9093S7Socket, 5-prong, AN 3102A-14S-5S5629S8Socket, Amphenol 80-PC2F4586TR1Thermoregulator, 58°C, normally closed, Edison Type S1-1A1036SW1Switch, rotary, Oak7515SW2Switch, DPST, Arrow H&H 806002313SW3Switch, normally closed, Grayhill 40027521F1,F2Fuse, 3 amp516NE2Lamp Neon7522	V9	Tube, 6C4, RCA	13361
V11 Tube, 5AW4 CBS 7627 M1 Meter 7511 T1 Transformer (with matched condenser CT) 7513 P2 Plug with leads 7569 P3 Plug with leads 7555 P4 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P9 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, Open miniature, Cince 13481 7269 S3 Socket, Jones S302CCT 7539 S4 Socket, Jones S302CCT 7539 S5 Socket, Jones S302AB 5611 S6 Socket, Jones S302AB 5611 S6 Socket, Amphenol 80C 9093 S7 Socket, S-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80-PC2F 4536 TR1 Thermoregulator, 58°C, normally closed, Edison Type S1-1A 1036 SW1 Switch, rotary, Oak 7515 552 SW2 Switch, normally closed, Grayhill 4002 7521 F1,F2 Fuse, 3 amp	V10	Tube, 12AX7, RCA	15093
M1 Meter 7511 T1 Transformer (with matched condenser CT) 7513 P2 Plug with leads 7569 P3 Plug with leads 7555 P4 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P9 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, Jones S302CCT 7539 S5 Socket, 2-prong Jones S302AB 5611 S6 Socket, S-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80C 9093 S77 Socket, Amphenol 80-PC2F 4536 TR1 Thermoregulator, 58°C, normally closed, Edison Type S1-1A 1036 SW1 Switch, rotary, Oak 7515 313 Sw3 Switch, normally closed, Grayhill 4002 7521 F1,F2 Fuse, 3 amp 516 NE2 Lamp Neon 7522	V 11	Tube, 5AW4 CBS	7627
T1 Transformer (with matched condenser CT) 7513 P2 Plug with leads 7569 P3 Plug with leads 7553 P4 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P9 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, Jones S302CCT 7539 S5 Socket, 2-prong Jones S302AB 5611 S6 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80C 9093 S7 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80-PC2F 4586 TR1 Thermoregulator, 58°C, normally closed, 7515 SW1 Switch, rotary, Oak 7515 SW2 Switch, DPST, Arrow H&H 80600 2313 SW3 Switch, normally closed, Grayhill 4002 7521 F1,F2 Fuse, 3 amp 516 NE2 Lamp Neo	Ml	Meter	7511
P2 Plug with leads 7569 P3 Plug with leads 7555 P4 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P9 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, Jones S302CCT 7539 S5 Socket, 2-prong Jones S302AB 5611 S6 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80-PC2F 4536 TR1 Thermoregulator, 58°C, normally closed, 1036 SW1 Switch, rotary, Oak 7515 SW2 Switch, DPST, Arrow H&H 80600 2313 SW3 Switch, normally closed, Grayhill 4002 7521 F1,F2 Fuse, 3 amp 516 NE2 Lamp Neon 7522	T1	Transformer (with matched condenser CT)	7513
P3 Plug with leads 7555 P4 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P9 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, Jones S302CCT 7539 S5 Socket, 2-prong Jones S302AB 5611 S6 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80-PC2F 4536 TR1 Thermoregulator, 58°C, normally closed, 1036 SW1 Switch, rotary, Oak 7515 SW2 Switch, DPST, Arrow H&H 80600 2313 SW3 Switch, normally closed, Grayhill 4002 7521 F1,F2 Fuse, 3 amp 516 NE2 Lamp Neon 7522	P2	Plug with leads	756 9
P4 Plug, 2-prong, Jones P302AB 7538 P5 Plug, Jones P302CCT 9055-0 P9 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Sccket, Jones S302CCT 7539 S5 Socket, 2-prong Jones S302AB 5611 S6 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80-PC2F 4536 TR1 Thermoregulator, 58°C, normally closed, Edison Type S1-1A 1036 SW1 Switch, rotary, Oak 7515 SW2 Switch, DPST, Arrow H&H 80600 2313 SW3 Switch, normally closed, Grayhill 4002 7521 F1,F2 Fuse, 3 amp 516 NE2 Lamp Neon 7522	P3	Plug with leads	7555
P5 Plug, Jones P302CCT 9055-0 P9 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, Jones S302CCT 7539 S5 Socket, 2-prong Jones S302AB 5611 S6 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80-PC2F 4586 TR1 Thermoregulator, 58°C, normally closed, Edison Type S1-1A 1036 SW1 Switch, rotary, Oak 7515 SW3 Switch, normally closed, Grayhill 4002 7521 F1,F2 Fuse, 3 amp 516 NE2 Lamp Neon 7522	P4	Plug, 2-prong, Jones P302AB	7538
P9 Plug, 3-prong 4536 S1 Socket, Amphenol 82-805 5614 S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, Jones S302CCT 7539 S5 Socket, 2-prong Jones S302AB 5611 S6 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80C-PC2F 4596 TR1 Thermoregulator, 58°C, normally closed, Edison Type S1-1A 1036 SW1 Switch, rotary, Oak 7515 SW2 Switch, normally closed, Grayhill 4002 7521 F1,F2 Fuse, 3 amp 516 NE2 Lamp Neon 7522	P5	Plug, Jones P302CCT	9055-0
S1Socket, Amphenol 82-8055614S2Socket, 9-pin miniature, Cince 134817269S3Socket, 7-pin miniature, Cinch 56F1286512052S4Sccket, Jones S302CCT7539S5Socket, 2-prong Jones S302AB5611S6Sccket, Amphenol 80C9093S7Socket, 5-prong, AN 3102A-14S-5S5629S8Socket, Amphenol 80-PC2F4586TR1Thermoregulator, 58°C, normally closed, Edison Type S1-1A1036SW1Switch, rotary, Oak7515SW2Switch, normally closed, Grayhill 40027521F1,F2Fuse, 3 amp516NE2Lamp Neon7522	P9	Plug, 3-prong	4536
S2 Socket, 9-pin miniature, Cince 13481 7269 S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, Jones S302CCT 7539 S5 Socket, 2-prong Jones S302AB 5611 S6 Socket, Amphenol 80C 9093 S7 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80-PC2F 4586 TR1 Thermoregulator, 58°C, normally closed, Edison Type S1-1A 1036 SW1 Switch, rotary, Oak 7515 SW2 Switch, DPST, Arrow H&H 80600 2313 SW3 Switch, normally closed, Grayhill 4002 7521 F1,F2 Fuse, 3 amp 516 NE2 Lamp Neon 7522	S1	Socket, Amphenol 82-805	561 4
S3 Socket, 7-pin miniature, Cinch 56F12865 12052 S4 Socket, Jones S302CCT 7539 S5 Socket, 2-prong Jones S302AB 5611 S6 Socket, Amphenol 80C 9093 S7 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80-PC2F 4586 TR1 Thermoregulator, 58°C, normally closed, 1036 SW1 Switch, rotary, Oak 7515 SW2 Switch, DPST, Arrow H&H 80600 2313 SW3 Switch, normally closed, Grayhill 4002 7521 F1,F2 Fuse, 3 amp 516 NE2 Lamp Neon 7522	\$2	Socket, 9-pin miniature, Cince 13481	726 9
S4 Sccket, Jones S302CCT 7539 S5 Socket, 2-prong Jones S302AB 5611 S6 Sccket, Amphenol 80C 9093 S7 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80-PC2F 4586 TR1 Thermoregulator, 58°C, normally closed, Edison Type S1-1A 1036 SW1 Switch, rotary, Oak 7515 SW2 Switch, DPST, Arrow H&H 80600 2313 SW3 Switch, normally closed, Grayhill 4002 7521 F1,F2 Fuse, 3 amp 516 NE2 Lamp Neon 7522	S3	Socket, 7-pin miniature, Cinch 56F12865	12052
\$5Socket, 2-prong Jones \$302AB5611\$6Socket, Amphenol 80C9093\$7Socket, 5-prong, AN 3102A-14S-5S5629\$8Socket, Amphenol 80-PC2F4586TR1Thermoregulator, 58°C, normally closed, Edison Type \$1-1A1036\$W1Switch, rotary, Oak7515\$W2Switch, DPST, Arrow H&H 806002313\$W3Switch, normally closed, Grayhill 40027521F1,F2Fuse, 3 amp516NE2Lamp Neon7522	S4	Sccket, Jones S302CCT	753 9
S6 Socket, Amphenol 80C 9093 S7 Socket, 5-prong, AN 3102A-14S-5S 5629 S8 Socket, Amphenol 80-PC2F 4586 TR1 Thermoregulator, 58°C, normally closed, Edison Type S1-1A 1036 SW1 Switch, rotary, Oak 7515 SW2 Switch, DPST, Arrow H&H 80600 2313 SW3 Switch, normally closed, Grayhill 4002 7521 F1,F2 Fuse, 3 amp 516 NE2 Lamp Neon 7522	S 5	Socket, 2-prong Jones S302AB	5611
S7Socket, 5-prong, AN 3102A-14S-5S5629S8Socket, Amphenol 80-PC2F4586TR1Thermoregulator, 58°C, normally closed, Edison Type S1-1A1036SW1Switch, rotary, Oak7515SW2Switch, DPST, Arrow H&H 806002313SW3Switch, normally closed, Grayhill 40027521F1,F2Fuse, 3 amp516NE2Lamp Neon7522	S 6	Socket, Amphenol 80C	9093
S8Socket, Amphenol 80-PC2F4586TR1Thermoregulator, 58°C, normally closed, Edison Type S1-1A1036SW1Switch, rotary, Oak7515SW2Switch, DPST, Arrow H&H 806002313SW3Switch, normally closed, Grayhill 40027521F1,F2Fuse, 3 amp516NE2Lamp Neon7522	S7	Socket, 5-prong, AN 3102A-14S-5S	562 9
TR1Thermoregulator, 58°C, normally closed, Edison Type S1-1A1036SW1Switch, rotary, Oak7515SW2Switch, DPST, Arrow H&H 806002313SW3Switch, normally closed, Grayhill 40027521F1,F2Fuse, 3 amp516NE2Lamp Neon7522	S 8	Socket. Amphenol 80-PC2F	4586
Edison Type S1-1A1036SW1Switch, rotary, Oak7515SW2Switch, DPST, Arrow H&H 806002313SW3Switch, normally closed, Grayhill 40027521F1,F2Fuse, 3 amp516NE2Lamp Neon7522	TR1	Thermoregulator, 58°C, normally closed.	
SW1Switch, rotary, Oak7515SW2Switch, DPST, Arrow H&H 806002313SW3Switch, normally closed, Grayhill 40027521F1,F2Fuse, 3 amp516NE2Lamp Neon7522		Edison Type SI-1A	1036
SW2Switch, DPST, Arrow H&H 806002313SW3Switch, normally closed, Grayhill 40027521F1,F2Fuse, 3 amp516NE2Lamp Neon7522	SW 1	Switch, rotary, Oak	7515
SW3Switch, normally closed, Grayhill 40027521F1,F2Fuse, 3 amp516NE2Lamp Neon7522	SW2	Switch, DPST, Arrow H&H 80600	2313
F1,F2Fuse, 3 amp516NE2Lamp Neon7522	SW3	Switch, normally closed, Gravhill 4002	7521
NE2 Lamp Neon 7522	F1.F2	Fuse, 3 amp	516
	NE2	Lamp Neon	7522

• ">



FIGURE 8. Circuit Diagram, Model VL.

1

Beckman Scientific and Process Instruments Division Beckman Instruments, Inc. j HERMAN L. STEUDEL 312 - 583 - 1020 JUNIFER 3-1020 7360 N. LINCOLN AVE. LINCOLNWOOD 46, ILLINOIS



IN STRUCTION MANUAL

BECKMAN MODEL V MICRO-MICROAMMETER

OPERATING PROCEDURE

Without the 1710 Multiple Switch:

1. Turn Power Switch OFF.

2. Connect instrument to 115-volt, 60 cycle power source. Allow 30 minutes for isothermal shield to reach temperature equilibrium.

3. Turn Operational Switch to MANUAL.

4. Connect Current source to INPUT terminal.

5. Check mechanical zero of meter.

6. Set Selector Switch to ZERO.

7. Turn Power Switch ON. Amplifier warm-up period is 3 to 5 minutes.

8. Adjust meter to zero with ZERO control.

9. Turn Selector Switch to desired range (black figures).

10. Read signal value from meter.

With the 1710 Multiple Switch:

1. Turn Power Switch OFF.

2. Connect instrument to 115-volt, 60 cycle power source. Allow 30 minutes for isothermal shield to reach temperature equilibrium.

3. Make 1710 Switch and B + connections.

4. Turn Operational Switch to AUTOMATIC.

5. Check mechanical zero of meter.

6. Set Selector Switch at ZERO. Turn Power Switch ON. Amplifier warm-up period is 3 to 5 minutes.

7. Adjust meter to zero with ZERO CONTROL.

8. Set Selector Switch to desired Input Resistor position (see paragraph I C 2 and Table I of this manual and Instruction Manual for 1710 Multiple Switch).

9. See Instruction Manual, 1710 Multiple Switch, for measurement of signal currents.

BECKMAN INSTRUCTION MANUAL

MODEL V MICRO-MICROAMMETER

TABLE OF CONTENTS

SECTION I – OPERATION	Pa	ige
A. Description.		1
B. Specifications		2
C. Operating Controls		2
1. Power Switch		2
2. Selector Switch		2
3. Operational Switch		3
4. Zero Control		3
5. Indicating Light		3
6. Meter	• •	3
7. Recorder Terminals	••	3
9. Multiple Switch Socket	••	3
10. Input Terminal	••	3
11. Feedback Terminal	•••	4
12. GND and OSC Terminals		4
D. Operating Procedure		4
SECTION II – MAINTENANCE		
A. Circuit Operation		5
1. Vibrating Reed		5
2. Contact Potential.		6
3. Amplifier 1		7
4. Ampliner 2	••	7
6 Synchronous Rectifier Demodulator	• •	2
7. Feedback Circuit	• •	9
8. Regulated Power Supply Amplifier		9
B. Testing Equipment		9
C. Measuring Amplifier Gain		10
D. Voltage Calibration		11
E. Current Calibration		12
F. Output Noise		12
G. Trouble-Shooting		15
PARTS LIST		17

Beckman Division, BECKMAN INSTRUMENTS, INC. • Fullerton, California

CHICAGO, ILLINOIS . MOUNTAINSIDE, NEW JERSEY

SEPTEMBER

Printed in USA

6	hipped , 20,	1964 / Fri	edit No. 46	216 14	voice Al 410-2	um 4350	
I.R. NO	TYT FOR V	REQ	UEST FOR ISSUANC	E OF	S.N	. NO. 6	00245
Jae na	Reck	T to St	IPPING NOTI	CE	۵ÅT	= 126	161
with the grice Dig	2600 da	where Blad			UNI		
MAN GUN	Eller h	· CIC.		COLLECT PREPAID	DATE SHIP	PED	
	FUILER KOR	2 CHII TORNIC	· · ·				
				SHIP VIA			
<i></i>				SOURCE OF MATERIA	L TO BE SH	IPPED	
CHARGE TO:				OUR ORDER NO.	·····		
		39			07	Y REJID	
	<u> </u>			• •	Y I	T . KLU D	
TEDMS		F.O.B. POINT		DATE REC'D		R. NO	
TENHS							
QUANTITY		DESCRIPTION			UNIT	TOTAL	DELIVERY
		*			PRICE	PRICE	WANTED
1	model VI	- Logna, thimic	Micko 1	1116160			
	History	eter	×				
							×
	100 m						
	Note: t	two Copies to tran Hie D	A Repr Dept.	∀ + R			
		· · · · · · · · · · · · · · · · · · ·			APPROX VAL	LUE, IF G.E	PROPERTY
PRO-RATA CHARG	BE FOR INCOMING	TRANSP (FOR BUYER)			495	- 00	
			REASO	N FOR SHIPMENT			
REPAIR	REPLACE	Round 1	lite Ala	ater 1 al	VI 21	· · · · · ·	
REPAIR OR REPLA	CE	in prove we we	<u>y ya na ka</u> iyazi d itana ku	I A. I	<u> </u>	<u></u>	
RETURN FOR CRED	<u>س</u> ۱۲	HIETER GI	<u>435 1013</u> 1	LECKE T	27 - 16		/
TO BE USED ON P	.0.	ypoint KE	CEIPT N	le Visikle	Ver.	MAGE	te
		Cructors		. <u> </u>			
OTHER	CTEN DV.				1		
SHIPMENI KEQUI			DICATE BY AN	X	CREDIT AC	CT. NO.	
C. ULESE			VENDARIE BROD			AT NO	
MAIERIAL URAW	n I.NVM.	AUURESEL'S PRUP.	TENDUR S PROP.	UUR FRUF.	UNARUE AU	101. HU.	67
		DE IECTED MATI	SALE		CHARCE IN	CONING TRANS	/ @ P TO
Carl.d.	a tak	REJECTED MAT L.	JALL		UNARUE IN	CONTRACTOR	T ~ 1V.
DELIVER RETURNS	TO	TO BE RETURNED	IF NO, WHY?		CHARGE OU	TGOING TRANS	P. TO
Hic	Knad	NOYES _X					
P			SHIPMENT	APPROVED BY:			
⊧736 R3-53					BU	YER	

. .

<u>,</u>

/

LIST OF ILLUSTRATIONS

Figure 1.	Model V Micro-Microammeter	1
Table I.	Full-Scale Current in Amperes	3
Figure 2.	Rear Panel View, Model V	4
Figure 3.	Simplified Schematic Diagram	5
Figure 4.	Principle of a Variable Capacitor	6
Figure 5.	Change in Voltage with Change in Capacitance	6
Figure 6.	Voltage Relationship in a Charged Variable Capacitor	6
Figure 7.	Voltage Relationship in a Charged System	6
Figure 8.	Error-Sensing Circuit	7
Figure 9.	Modulator Signal Phase	7
Figure 10.	Synchronous Demodulator Circuit	8
Figure 11.	Current Flow in Demodulator	8
Figure 12.	Amplifier Gain Measurement Diagram	9
Table II.	Time Constant Data	11
Figure 13.	Voltage Calibration Diagram	11
Figure 14.	Current Calibration Diagram	11
Figure 15.	Capacitor–Discharge Method, Using a Recorder to Measure Time of Discharge	12
Figure 16.	Pattern of Vibrating Capacitor Shorted at Several Places	13
Figure 17.	Pattern of Vibrating Capacitor Shorted Only at the Ends of the Cycle	13
Figure 18.	Overloaded Signal Condition	13
Figure 19.	Normally-Operating Vibrating Capacitor	13
Figure 20.	Schematic Wiring Diagram of Special Calibrating Equipment	15
Figure 21.	Views of Model V Chassis	16
Figure 22.	Circuit Diagram of Model V (Single Input)	19

SECTION I OPERATION



Figure 1. Model V Micro-Microammeter

A. DESCRIPTION

The Beckman Model V Micro-Microammeter is a sensitive current measuring instrument with a vibrating capacitance modulator. Thirteen measuring ranges from 0.3 micro-microamperes to 0.3 microamperes are available.

Many of the highly desirable features of the instrument, including high response speed and accuracy, are due to the degenerative lowering of the input impedance in the more sensitive ranges to a small fraction (1/1000 to 1/6000) of the value of the input resistor used. This is accomplished in such a way that, with proper adjustment, the voltage between the conductor and the shield on the input cable will never exceed 2 mv plus the contact potential in any current range. The use of a vibrating capacitor and an a-c amplifier eliminates the first stage grid current and zero drift inherent in straight d-c amplifiers.

When positive current from an external source flows through the input circuit it tends to increase its potential. The amplifier output to the feedback circuit keeps the feedback end of the input resistors (R_r , Figure 3) sufficiently negative to bring the input circuit quickly to equilibrium potential.

B. SPECIFICATIONS

1. Input Power: 60 cycles, 103-127 volts, 100 watts.

2. Detection Range: 3×10^{-13} to 3×10^{-7} amperes.

3. Indicator: 200 microamperes, 41/2-inch scale.

4. Ambient Temperature Requirements: 20° F. to 115° F.

- 5. Warm-up Time: Amplifier: 3 to 5 minutes Isothermal Shield: 30 minutes
- 6. Zero Shift: Less than 2 millivolts in 24 hours.

7. Chamber Supply Voltage: 210 volts with maximum drift rate of 30 millivolts per second.

- 8. Output Connections:
 - a. Recorder: 50 millivolts potentiometer type b. Multiple Switch: Beckman 1710 Multiple Switch

9. Accuracy: $\pm 5\%$ on all ranges.

10. Time Constant:

	Input to Ground Capacitance (Cable and Detector)			
Range	None	1500 μ μf		
3×10 ⁻¹³	1 second	4 seconds		
10×10^{-13} to 3×10^{-11}	0.9	1.4		
1×10^{-11} to 3×10^{-7}	0.12	0.12		

11. Output Noise: On 3×10^{-13} ampere range; Less than 3% with no input capacitance Less than 8% with $5000\mu\mu$ f input capacitance

C. OPERATING CONTROLS

1. POWER SWITCH

This toggle switch disconnects the instrument, with the exception of the isothermal shield, from line voltage when in the OFF position. CAUTION: This switch must be OFF whenever the outer case is removed because of the danger of the high voltage present in the chassis.

2. SELECTOR SWITCH

This control, marked AMPERES, is used to select the desired range. When using the Model V alone (without a Beckman 1710 Automatic Multiple Switch), the Selector Switch should be in the position of the desired amperes range (black figures).

When using a 1710 Switch without plug-in resistors, the Model V Selector Switch should be on any one of the points in a group (red figures) which provides four ampere ranges, the particular range within the group being selected by the four sensitivity switches on the 1710 Switch. See Table I.

When using a 1710 Switch with plug-in resistors (which shunt the Model V input resistor), the Model V Selector Switch is on a group number (red figures) and the Model V input resistor in the circuit is listed in Table I. The current sensitivity or plug-in resistor value is calculated from the following equation:

$$R_p = \frac{R_v E_f}{I_f R_v - E_f}$$
, in which

- $E_f =$ Feedback or full scale voltage* in volts
- $I_f =$ Full scale current in amperes
- $R_e = Effective \text{ input resistance } (R_v \text{ and } R_p \text{ in parallel})$ in ohms

 $R_v =$ Model V input resistor in ohms

 $R_p =$ Plug-in resistor in ohms

*Selected by 1710 Switch sensitivity switch as indicated in Table I.

If no plug-in resistor is used in one or more channels of the 1710 Switch, the full-scale current I_f may be read directly from Table I, or calculated from

 $I_f = \frac{E_f}{R_v}$, a reduction of the above equation when R_p equals infinity.

Table I applies to all four ion-chamber stations simulaneously. All four ion-chamber currents must be within a 30 to 1 ratio; if the four currents do not fall within this 30 to 1 ratio, then plug-in resistors must be installed in the 1710 Switch to permit measurements on the chambers not within the 30 to 1 limit. The installation of such plug-in resistors can only lower the sensitivity and can never be used to increase it. The installation of a plug-in resistor increases the absolute magnitude of a full scale reading, e.g., from 3×10^{-13} amperes to 3×10^{-12} or 3×10^{-11} amperes.

3. OPERATIONAL SWITCH

This toggle switch, marked MANUAL and AUTO-MATIC, is turned to the latter position when the instrument is used with a 1710 Multiple Switch. When the instrument is used without a 1710 Multiple Switch, the toggle switch must be in the MANUAL position.

4. ZERO CONTROL

Used to bring the meter needle to zero on the meter scale initially and in subsequent corrections for needle drift.

5. INDICATING LIGHT

Lights a few seconds after the instrument is turned on, indicating that the B + supply is operating. If it does not light with the Power Switch in the ON position, the B supply circuit of the instrument is not functioning properly and should be checked.

6. METER

Indicates the magnitude of the current under measurement. Two meter styles are available; one calibrated in percent and the other with two scales of 3 units and 10 units, full scale, respectively. Full scale deflection on both meters represents a current equal to the setting of the Selector Switch. Intermediate readings are fractions of the current indicated by the Selector Switch. For example, a reading of 30 on the percent scale with the Selector Switch set on 3×10^{-11} amps would show 0.9×10^{-11} amps of input current. A reading of 1.5 on the lower scale of the dual scale meter with the Selector Switch set on 3×10^{-11} amps would represent 1.5×10^{-11} amps of input current.

7. RECORDER TERMINALS

These terminals are used to connect a potentiometertype recorder or controller having a range of 50 mv. NO OTHER CONNECTIONS SHOULD BE MADE TO THESE TERMINALS. The recorder should have uniform chart gradations, with the left side corresponding to zero voltage input to the potentiometer. It is important that no temperature sensitive element, such as a cold-junction temperature compensator, be connected to the potentiometer circuit.

8. REGULATED B+ TERMINAL

Marked B + REG, this coaxial fitting has an internal resistance of 1 megohm and is used to supply a precisely-regulated potential of 210 volts, with respect to ground, to polarize an ionization chamber. B +variation, for ionization chambers, shall be less than 0.3 millivolts when the line voltage changes from 103 to 127 volts ac. Use only shielded cable for this connection. Cable shield must be grounded only to the coaxial fitting.

9. MULTIPLE SWITCH SOCKET

This socket, marked MULT SW, permits the use of a 1710 Multiple Switch. The socket provides the feedback voltage circuit and return circuit for the input resistors in the 1710 Switch. See instruction manual supplied with the 1710 Multiple Switch.

10. INPUT TERMINAL

This connection, denoted INPUT 1, is a standard coaxial socket for an Amphenol 82-805 fitting. The input cable is not furnished, but a standard shielded coaxial cable and fitting is required. The shielding of the input cable and jack is of extreme importance and should not be neglected when connecting a current source to the Model V. If a well-shielded cable is used, a chassis-to-earth ground will not normally be required. If touching the instrument causes the meter needle to move, the input or B+ cable shielding is faulty.

A more sensitive test for faulty shielding can be performed with an oscilloscope connected to OSC terminal. If waveform changes when instrument is touched, shielding is faulty. If a chassis-to-earth ground is utilized to correct for faulty shielding, it should be connected only to the GND terminal on the Model V or to the wire provided on the main power cable. Under no circumstances should both the instrument and the ion-chamber be individually connected to earth ground.

11. FEEDBACK TERMINAL

This terminal, marked FB, is used in checking the voltage sensitivity as described under MAINTENANCE.

12 GND AND OSC TERMINALS

These terminals are used for circuit testing. The vacuum tube voltmeter and oscilloscope used in circuit testing are both connected to these terminals.

OPERATING PROCEDURE

Without the 1710 Multiple Switch:

1. Turn Power Switch OFF.

2. Connect instrument to 115-volt, 60 cycle power source. Allow 30 minutes for isothermal shield to reach temperature equilibrium.

- 3. Turn Operational Switch to MANUAL.
- 4. Connect Current source to INPUT terminal.
- 5. Check mechanical zero of meter.

Selector Switch on Model V Group	Input Resistor Model V	Sensiti	tivity Switch Position on 1710 Switch			
(in red)	(ohms)	A (0.1v)**	B (0.3v)**	C (1.0v)**	D (3.0v)**	
I	3.33×1011	0-3×10 ⁻¹³	0-9×10 ⁻¹³	0-3×10 ⁻¹²	0-9×10-12	
II	1.0×1011	0-10×10 ⁻¹³	0-3×10 ⁻¹²	0-10×10 ⁻¹²	0-3×10 ⁻¹¹	
III	1.0×10 ⁹	0-10×10 ⁻¹¹	0-3×10 ⁻¹⁰	0-10×10 ⁻¹⁰	0-3×10 ⁻⁹	
IV	1.0×10^{7}	0-10×10-9	0-3×10 ⁻⁸	0-10×10 ⁻⁸	0-3×10 ⁻⁷	

Table I. Full-Scale Current in Amperes

**Feedback and full-scale voltage

6. Set Selector Switch to ZERO.

7. Turn Power Switch ON. Amplifier warm-up period is 3 to 5 minutes.

8. Adjust meter to zero with ZERO control.

9. Turn Selector Switch to desired range (black figures).

10. Read signal value from meter.

With the 1710 Multiple Switch:

1. Turn Power Switch OFF.

2. Connect instrument to 115-volt, 60 cycle power source. Allow 30 minutes for isothermal shield to reach temperature equilibrium.

3. Make 1710 Switch and B + connections.

4. Turn Operational Switch to AUTOMATIC.

5. Check mechanical zero of meter.

6. Set Selector Switch at ZERO. Turn Power Switch ON. Amplifier warm-up period is 3 to 5 minutes.

7. Adjust meter to zero with ZERO CONTROL.

8. Set Selector Switch to desired Input Resistor position (see paragraph I C 2 and Table I of this manual and Instruction Manual for 1710 Multiple Switch).

9. See Instruction Manual, 1710 Multiple Switch, for measurement of signal currents.



Figure 2. Rear Panel View, Model V



Figure 3. Simplified Schematic Diagram Model V Circuit

A. CIRCUIT OPERATION

The Model V Micro-Microammeter is a very sensitive current-measuring instrument. The signal current flows through a high value resistor (see Figure 3), and the resultant voltage is measured by an electrometer amplifier. In order to achieve rapid response to small signals in the presence of large input capacitance, e.g., cable capacitance, negative voltage feedback is used. The feedback connection is such that the potential input across the signal terminals of the instrument is maintained constant. Thus, the instrument behaves as a classical ammeter. The actual potential level of the input terminals has been chosen as ground to minimize the effect of insulation leakage currents on the measurements.

The circuit is an electronic servo, in which the amplifier strives to maintain the input at ground potential in spite of the magnitude of voltage developed across the input resistor by the signal current. The amplifier unit is designed to compare the input potential with ground potential and to react to any error signal by producing a driving voltage at the bottom end of the high-value resistor with proper sense to return the input to ground level. The ability of the amplifier to counteract the error signal is determined by the magnitude of the open loop gain. The compensating voltage applied to the input also serves to indicate the input signal magnitude.

In the Model V Micro-Microammeter, the dc error voltage signal is first converted into an ac signal by a vibrating capacitor (C, Figure 3). This ac signal is amplified, then rectified synchronously (to preserve dc polarity sense). The resultant dc voltage is applied to the bottom of the input resistor to complete the servo loop.

1. VIBRATING REED

Consider a variable capacitor with no voltage across it; if the capacitance of the unit is changed by moving the plates farther apart, the voltage across the unit will remain at zero. If, however, the capacitor is charged to 1.0 volt in the closed position and the source of charging voltage removed, the voltage across the capacitor increases as the plates are moved apart. (Q = CE), where Q is the charge, C is the capacitance, and E is the voltage.) Q, the charge, is held constant by removing the voltage source; hence, if C charging, then E must change also. See Figure 4. Since the voltage across the capacitor varies inversely with the capacitance of the unit, let the distance between the condenser plates vary in some regular manner and observe the voltage across the unit. See Figure 5. The farther the plates are apart, the higher the voltage; the closer the plates, the lower the voltage.

The actual voltage which can be obtained from such a system is limited by the practical maximum and minimum capacitances that can be achieved. The total maximum and minimum capacitances of the system are the important quantities. In any real system, shunt capacitance and stray capacitance must be considered. For example, suppose a variable capacitor is used with a maximum capacitance of 75 $\mu\mu$ f and a minimum capacitance of $25\mu\mu$ f. Figure 6 shows that if the capacitor is charged to 1.0 volt in the closed position, 2 volts peak-to-peak ac can be obtained from the system. However, if a fixed capacitor of 75 $\mu\mu$ f is connected across the system as shown in Figure 7, only 0.5 volt peak-to-peak will be obtained.





2. CONTACT POTENTIAL

In a real capacitor, there always exists a potential between the plates, even if the leads to the capacitor unit are carefully shorted. At the surface of all metals there exists a characteristic potential barrier which prevents the escape of electrons out of the body of the metal into free space. In some metals this barrier is small; these metals or alloys are used in phototubes. (The barrier can be considered as a battery in series with the surface of the metal.)

If a capacitor is constructed with metals which have identical contact potentials, then the potential difference between the two surfaces would be zero if the leads to the capacitor were shorted. Where it is desired to construct two surfaces with no potential difference, extreme care must be taken. The surfaces should be made of the same kind of metal, and preferably from the same piece of metal. The metal chosen should be thoroughly cleaned and kept clean. Even with these precautions, it is not always possible to reduce the contact potentials to zero; the surface characteristics can be changed drastically by adsorption, abrasion, chemical changes, physical change of state, etc.

The 'residual' contact potential of the Beckman Vibrode is reduced to a few millivolts and will not be greater than 20 millivolts. Change of reed contact potential results in 'zero drift' in the instrument. The



Figure 6. Voltage Relationship in a Charged Variable Capacitor



Figure 7. Voltage Relationship in a Charged System

'drift' of the contact potential will be less than 2 millivolts per day and will have a wandering characteristic, i.e., the drift will not be unidirectional.

The Beckman Vibrode is a variable capacitor which is sealed in a glass envelope and varied by the application of a 120-cycle magnetic field. The capacitor is hence vibrated at a frequency of 120 cycles. The 120cycle magnetic field is derived from a coil driven from a 60-cycle power line.



Figure 8. Error-Sensing Circuit

Figure 8 shows the error-sensing circuit of the Model V. The 300-megohm resistor R1 serves to isolate the reed capacitance from the large cable capacitance so that maximum ac voltage will be developed across the reed. Resistor R2 and Capacitor C1 serve to couple the reed ac-wise to the amplifier and yet isolate the reed from the bias potential of the input tube of the amplifier.

At 120 cycles, the impedance of the reed is about 30 megohms, the impedance of the coupling capacitor C1 about 10 megohms, and the impedance of the cable capacitance about 200,000 ohms; the efficiency of the conversion is limited chiefly by the reed itself. This system differs from the first example used in that the charging voltage is applied to the system continuously. As a result, the maximum ac voltage that can be obtained has a peak-to-peak value of twice the dc voltage. The conversion efficiency of the modulator unit is about 10%.

3. AMPLIFIER 1

The modulated 120-cycle signal is fed into Amplifier 1. This is a three-stage, ac amplifier with negative voltage feedback to stabilize the amplifier gain. The first two stages are subminiature tubes with filament current supplied from regulated B^+ .

A frequency response correction network is employed in the feedback line to the bottom of the reed to improve the transient response of the amplifier.

It is not convenient to measure the gain of the amplifier itself, so in testing the instrument, the gain of the circuit from the input through Amplifier 1 is measured. This measurement includes the modulator unit. Gains of about 3000 are usually observed for Amplifier 1. Gains of less than 750 necessitate correction.

4. AMPLIFIER 2

This amplifier consists of a power amplifier stage which receives the output signal of Amplifier 1 and a synchronous rectifier demodulator to convert the amplified ac to dc. This is followed by a smoothing



filter network. The resultant output is fed back to the input of the instrument to degenerate the input impedance. The current in the dc modulated signal is also used to drive the ouput meter.

5. PHASE RELATIONS

Figure 9 shows the phase and frequency relationships in the reed modulator unit. (A) is the 60-cycle line that is applied to the reed drive coil. This results in a 120-cycle reed motion (B), because the reed plates repel each other each time a magnetic field is generated by the coil current. The reeds repel each other regardless of the direction of the field. The reeds will be farthest apart when the driving coil current is highest, and closest together when the current is zero.

When a charging voltage is applied to the modulator as in the circuit of Figure 8, the highest part of ac voltage generated will be when the reeds are farthest apart. (C) and (D) of Figure 9 show that when a positive signal is applied, the phase of the generated ac signal is in phase with the 120-cycle reed motion, and that when a negative dc signal input is used, the generated ac signal is 180° out of phase with the reed motion. Thus, the polarity of the input can be determined by comparing the phase of the amplified signal with a signal known to be in phase with the reed motion.



Figure 10. Synchronous Demodulator Circuit

6. SYNCHRONOUS RECTIFIER DEMODULATOR

Consider the circuit in Figure 10. At each rectifier, an arrow indicates the direction of electron flow. The amplifier signal from the reed modulator is applied through transformer T1. A 120-cycle reference signal is applied through T2. This signal, obtained from the 120-cycle ripple on the unregulated power supply, is in phase with the 120-cycle magnetic field that drives





Figure 11. Current Flow in Demodulator

the reed. It has been shown that the phase of the signal from the reed modulator depends on the polarity of the dc charging voltage applied to the reed modulator. The ac 120-cycle reference signal must be at least twice as large as any expected dc signal developed across the input resistor detector.

The circuits of Figure 11 show the current flow through the demodulator circuit under the various conditions. Note that current flows in the circuit for only one part of the cycle. The boxes represent the transformer secondaries with polarity indicated. In (B) and (D) the heavy lines indicate the current flow.

In (B) and (D), current also flows through the other rectifiers; however, the actual voltages on the individual rectifiers are such that unequal currents pass through the rectifiers. (The heavy lines indicate the larger current.) Under conditions of no signal, the two rectifiers pass equal currents; hence they produce no net current in the external resistor.

In actual practice, resistors are placed in series with each rectifier to limit the current and also to equalize the rectifier impedances to prevent "one-sidedness." A large capacitor is placed across the external resistor to smooth the rectified voltage suitably for feeding back to the input of the instrument.

7. FEEDBACK CIRCUIT

The filtered voltage signal developed across the feedback resistors (R38, R27, R26, R25) is isolated from circuit potentials by transformer coupling. Thus, both ends of the feedback resistor must be returned to the circuit. The reference end of the resistor is returned through the zero control, which allows adjustment of the feedback potential in relation to ground to compensate for the contact potential of the vibrating reed modulator.

The other end of the feedback resistor is applied to the low impedance end of the input resistors. This feedback voltage must be well filtered to insure that no ac voltages can appear at the input of the instrument without attenuation, at least by gain of the overall amplifier. The alternating current of 60 cycles and the 120-cycle modulator frequency are to be especially discriminated. Each input resistor has a speciallydesigned filter network associated with it to insure sufficient discrimination at the input. For this reason also, any signal source used with the instrument must be carefully designed to eliminate ac pickup voltages from appearing at the input of the instrument.



Figure 12. Amplifier Gain Measurement Diagram

8. REGULATED POWER SUPPLY AMPLIFIER

The amplifier consists of an RCA 5651 Voltage Reference Tube and a three-stage, direct-coupled amplifier which controls a series regulator tube; current feedback is employed through the bleeder string. The first two stages of the amplifier are subminiature tubes with filamentary cathodes heated by the bleeder string current. The input grid of the amplifier is also coupled through capacitor C13 directly to B⁺ to improve the regulation to power line transients. R54, a 35,000ohm resistor, shunts the current regulator tube and carries half the current under nominal line conditions. The output tube of the amplifier uses two neon tubes in the cathode return to avoid cathode degeneration of the signal.

The regulated B^+ is used to provide polarizing voltage for ionization chambers. R51 and C27 provide additional filtering for this purpose.

B. TESTING EQUIPMENT

The following equipment will be required to test the operation of the Model V.

- 1 Model 250 Simpson Analyzer or equivalent
- 1 Oscilloscope, 2" diameter or larger, 1 volt rms/in. minimum vertical sensitivity, 60 cps horizontal

 3×10^{-13} to 3×10^{-7} amperes as shown in column one, Table II. These current values should be ± 2 percent for accurate checks.

Using these instruments the following tests can be made on the Model V.

- 1. Loop gain test: To check the over-all gain in the two amplifiers.
- 2. Gain of Amplifier 1
- 3. Gain of Amplifier 2
- 4. Vibrating capacitor efficiency
- 5. Voltage calibration on all scales
- 6. Current calibration on all scales
- 7. Response time on all scales

8. Power supply output noise level

C. MEASURING AMPLIFIER GAIN

1. Connect the variable voltage source (0-50 mv) in series with the calibrated voltage source. Adjust the calibrated voltage source to deliver 0.0005 volts.

2. Connect the two voltage sources to the instrument through the INPUT jack. Turn the variable source ON and leave the calibrated source turned OFF.

3. Connect the vacuum tube voltmeter and the oscilloscope to test points OSC and GND.

4. Set the Selector Switch on the Model V to the 3×10^{-11} ampere position. Turn the Operational Switch to MANUAL.

5. Adjust the variable voltage source until the meter needle on the Model V reads zero. This voltage is the contact potential and should be 20 my or less.

6. Turn ON the calibrated (0.0005 volts) voltage source.

7. Read E (see Figure 3) on the vacuum tube voltmeter and E on the lower scale of the Model V Meter. With the Selector Switch set at 3×10^{-11} , full scale is 3 volts.

8. Reading these two voltages, the following specifications are given to indicate proper operation of the main circuit components.

- a. Over-all loop gain
 - $Gain = E_2/E_{in} \ge 1500$

If the over-all loop gain is less than 1500 as shown by this test, proceed to check the gain of Amplifiers 1 and 2.

b. Gain of Amplifier 1

 $Gain = E_1 / E_{in} = E / 0.0005$ $= 2000 \text{ E}_{1} \ge 600$

If the gain of Amplifier 1 is less than 600, consult the wiring diagram (Figure 22) and check the circuit for Amplifier 1.

When the over-all loop gain is low and the gain of Amplifier 1 is correct, check the gain of Amplifier 2.

c. Gain of Amplifier 2

 $Gain = E_2/E_1 \ge 2.5$

If the gain of Amplifier 2 as found by this test is less than 2.5, consult the wiring diagram (Figure 22) and check the circuit for Amplifier 2. Occasionally low gain in this stage is caused by the vibrating capacitor being out of phase with the demodulator; if so, replace vibrating capacitor. Also check C15 and C9. See Figure 22.

d. Response Time: Over-all loop gain of less than 1500, as determined in Paragraph above, may result in increased response time, over-damped response, or impaired accuracy. As the over-all loop gain drops, these effects will occur in the above order. The specification of loop gain of 1500 insures a time constant of 4 seconds or less on the 3×10^{-13} scale. A gain of 1500 on 3×10^{-11} scale is approximately equivalent to a gain of 500 on the 3×10^{-13} scale.

$$T = R_{in}C_s + \frac{R_{in}(C_i + C_c)}{A}$$

where T = Time constant in seconds (effective)

 R_r = input resistor (See Table II)

A = Over-all loop gain

- $C_i =$ Internal input to ground capacitance (See Table II)
- $C_c =$ Input cable capacitance
- $C_s = Built-in shunt capacitance (See Table II)$
- $E_{1} =$ Voltage output of Amplifier 1 measured between Oscilloscope Jack J1 and ground with an a-c vacuum-tube voltmeter
- $E_2 = Voltage sensitivity of Model V. See Table II,$ column headed Volts Full Scale
- E_{in} = Voltage impressed on input to Model V (from input to ground, which opens the feedback loop) for calibration of loop gain only.
- For example: Using the 3×10^{-13} scale

$$R_{v} = 3.33 \times 10^{11}$$

$$C_{c} = 5000 \ \mu\mu f$$

$$C_{i} = 30 \ \mu\mu f$$

$$C_{s} = 3 \ \mu\mu f$$

$$A = 500 \text{ (This corresponds to a gain of 1500 on the 3 \times 10^{-13} \text{ scale})}$$

$$\Gamma = 3 \times 10^{11} \left[3 \times 10^{-12} + \frac{5030 \times 10^{-12}}{500} \right]$$

= 0.9 + 3 = 3.9 sec

A





D. VOLTAGE CALIBRATION

The voltage calibration can be checked by applying a known calibrated voltage between the INPUT and FB (Feedback) terminals. Turn the Selector Switch to the current scale corresponding to the desired voltage range as shown in Table II.

- 10×10^{-13} ampere scale corresponds to .1 volt full scale on the Model V meter (M1)
- 3×10^{-12} ampere scale corresponds to .3 volt full scale on the Model V meter (M1)
- 10×10^{-12} ampere scale corresponds to 1 volt full scale on the Model V meter (M1)
- 3×10^{-11} ampere scale corresponds to 3 volts full scale on the Model V meter (M1)

Turn the Operational Switch to MANUAL. Use only those current scales having input resistors of 1×10^{11} ohms. If the voltage registered on the Model V Meter is not within acceptable limits (1.5 percent) proceed to test Amplifiers 1 and 2 as outlined above. If this test has already been made, check condenser C24 $(50 \ \mu f)$ for leakage or a short circuit. The instrument will tend to read high if this condenser is faulty. If these checks are all satisfactory the feedback resistors R25, R26, R27 and R38 should be investigated. The feedback resistors R25, R26, R27 and R38 can be checked by measuring them with a Wheatstone bridge. Before measuring, disconnect the Model V meter (M1) to avoid damaging it with the high currents in the bridge. See the Parts List in this manual for values and tolerances of the feedback resistors. Inspect resistor R35 by measuring with a micro-ammeter the current through the meter and the resistor when the Model V is in operating order and the Power Switch is ON.

BULLETIN 271-D

 $10,000 M_{12}$ 3 MV $T = \frac{3 \cdot 10^{-3}}{10^{10}} 3 \cdot 10^{-13}$ $3 \cdot 10^{-13}$

1. Page 2. Paragraph B 7 should read:

REVISION

SHEET

- 7. Chamber Supply Voltage: 210 volts with maximum drift rate of 0.3 millivolts per second.
- 2. Page 5, Column 1, third paragraph. First sentence should end: "...(C_v. Figure 3)."
- 3. Page 6, Figure 7. Fixed condenser at the top of the illustration should be labeled "75 μμf" instead of "25 μμf."
- 4. Page 7, Figure 9 B. Sine wave should be shifted 90° forward, so that it is in phase with wave in Figure 9 C.
- 5. Page 8, Figure 11 B. Polarity signs should be reversed on voltmeter in lower right corner.
- 6. Page 10, paragraph C 7. First sentence should read:

7. Read E_1 (see Figure 3) on the vacuum tube voltmeter and E_2 on the lower scale of the Model V Meter.

- 7. Page 10, paragraph C 8 a. After "...the gain of Amplifiers 1 and 2.", add: "Usually, the loop gain will be from 3000 to 6000."
- 8. Page 19. See revised circuit diagram, Figure 22.

E. CURRENT CALIBRATION

The current calibration can be determined by supplying a known current to the instrument, from the Beckman 1500 $\mu\mu$ Current Source, for example. In this case the current source must have an internal resistance at least equal to the input resistance R_r (Figure 3) of the Model V. The input resistor used for each current scale is shown in Table II.

The current calibration CANNOT be found by measuring the resistance between the INPUT and GND terminals, even with the instrument disconnected from the line, because the input resistance under these circumstances is shunted by the insulation leakage resistance.

If the voltage calibration is correct but the current calibration is low, the circuit should be checked for moisture or insulation leaks on the high value resistors.

Another method of current calibration follows.

Capacitor Discharge Method for Calculating Values of Input Resistors: The values of the input resistance for the four highest current ranges can be measured on a Wheatstone bridge. The value for the nine lower current ranges can be determined by measuring the time constant of the resistor in conjunction with a known low-leakage capacitor.

1. Connect a condenser of sufficient capacitance between the INPUT and FB terminals to give a product of $R_r \times c=100$ seconds (at least). A condenser with polystyrene dielectric may be used if it has been found to have a leakage resistance over 500 times the value of the input resistor being tested. Use a .1 mfd capacitor for calibrating the 10⁹-ohm resistor (R41). Use a .001 mfd capacitor for measuring the 10^{11} -ohm resistor (R42 and R43 in parallel), and for the 3.33×10^{11} -ohm resistor (R43). Protect the 'capacitor and input leads from external fields by shielding.

2. Make the meter adjustments as usual, with the Selector Switch set to the correct position to select the resistor to be measured (see Table II), and the Operational Switch set at MANUAL.

3. Connect the capacitor and determine if the dielectric absorption currents are negligible. This may be accomplished by discharging the condenser to 0 volts, and observing by the Model V meter the ability of the condenser to remain at 0 volts. Some condensers will act like voltage generators, remaining a few millivolts on either side of 0.

4. Charge the condenser with an external battery to a voltage between .10 and .15 volt on scale 3×10^{-13} , scale 10×10^{-13} , and scale 10×10^{-11}

5. Allow the condenser to discharge through the input resistor, measuring the number of seconds required for the meter needle to go from 10.0 to 3.67.





6. Divide the figure obtained in Step 5 by the condenser capacitance in thousands of micro-microfarads to find the value of the input resistor in thousands of megohms. The use of a recorder will greatly facilitate this test, and will also increase the precision. See Figure 14 for details of connections for using a recorder. The strip chart recorder will produce a chart similar to that shown in Figure 15. Use the time t to solve the equation

 $R = \frac{t}{C}$, where R = value of resistor in thousands of megohms t = time in seconds C = capacitance in micromicrofarads

F. OUTPUT NOISE

The following test should be performed over a period of at least 3 hours with the Model V connected to a recording unit.

1. Disconnect the instrument from line voltage and remove the dust cover by unscrewing the four bolts located above the rear panel.

2. Remove the metal plate located on the underside of the instrument below the INPUT and MULT SW outlets. This provides access to INPUT 1 and B + REG.

3. Connect a 50 $\mu\mu$ f polystyrene dielectric capacitor (Fast No. A9157BA) between INPUT 1 and B+ REG and solder connections. Make the connection at this point to take advantage of the shielding provided by the Model V chassis. Replace metal plate and dust cover. If the condenser is connected between INPUT 1 and B+ REG externally, it must be mounted in a metal box to shield it from stray radiation and other sources of noise. The 50 $\mu\mu$ f condenser simulates a very large ion chamber for test purposes.

4. Set the Selector Switch on the 3×10^{-13} scale and the Operational Switch on MANUAL. Connect the instrument to line voltage and turn the Power Switch ON. The noise level, as indicated by needle fluctuation, should be less than ± 5 percent.

5. If the noise level exceeds ± 5 percent, the power supply regulation may be faulty. It can be checked by suddenly changing the line voltage from 103 volts ac to 127 volts ac. If the needle kick for this shift in line voltage is less than 10 percent on the meter scale, the power regulation is functioning properly. If the noise level is above 5 percent and power supply regulation is not at fault, replacing the 5651 Tube may reduce the noise level.

The waveforms shown in Figures 16 through 19 are typical but their amplitudes will vary from capacitor to capacitor as a function of the uniformity of contact potential distribution on the surfaces. The most uniform surfaces will produce the least residual signal. Typical capacitors at maximum drive will produce .5 to .7 volt rms at the OSC terminal.



Figure 16. Pattern of Vibrating Capacitor Shorted at Several Places

The oscilloscope patterns in Figures 16 and 17 are typical of those of a shorted vibrating capacitor. Although a shorted vibrating capacitor can appear in several ways, it will usually appear as shown in Figure 16, which indicates a short in several places in the cycle. Figure 17 shows the pattern of a vibrating capacitor shorted only at the extreme ends of the cycle.

The pattern in Figure 18 shows an overloaded signal condition. If the load at the input is reduced or elimi-

nated, the pattern should regain its normal appearance; if it does not, the feedback system is inoperative.



Figure 17. Pattern of Vibrating Capacitor Shorted Only at the Ends of the Cycle



Figure 18. Overloaded Signal Condition

The appearance of a good vibrating capacitor is shown in Figure 19. Here the pattern is almost flat with little harmonic content.



Figure 19. Normally-Operating Vibrating Capacitor

Electronic Tubes: If it is necessary to remove electronic tubes, DO NOT handle tubes or wiring until the Power Switch is turned OFF. When handling the CK 533 AX tubes, hold them by the leads or the top of the envelope, never by the pressed end of the envelope. Be sure that not over 15 milliamperes passes through the filaments of the CK533AX tubes. The V3 (12AU7), V4, V7 (6AQ5) and V8 (6X4) tubes can be checked on a standard tube tester; however, such tests are not necessarily conclusive on V1, V2, V6 and V9 (CK533AX) tubes. All tubes are used as supplied by the manufacturer; in replacing tubes, no selection

1

is required. Actual trial replacement should be made first in cases of suspected faulty tubes. A set of spare tubes is located inside the instrument dust cover.

G. TROUBLE-SHOOTING

Following are some symptoms of malfunctions, with indications of probable sources of trouble in each case.

1. Erratic operation of the Meter (ripple over 1 mv on B+ as seen with an oscilloscope).

- a. May be caused by failure of the Regulated Power Supply to control B+ voltage properly, due to:
 - (1) Defective neon lamps (NE51).

(2) Open filament in tube V1, V2, V9, V6, V3B, or V7.

(3) Open resistor in filament string (R3, 10, 9, 33, 32, 34).

(4) R54 open

b. May be caused by noise in the Regulated Power Supply due to:

(1) Intermittent noise in Resistors R3, 10, 9, 33, 32, 34.

(2)Noisy V5 (5651).

- (3) Microphonic V9 or V6 (CK533AX).
- 2. Low Loop Gain may be due to:
- a. Defective Vibrode C_v.
- b. Improper phasing of Demodulator T2. Correct by selecting C15 for maximum voltage across T2.
- c. Defective (open) C9 or C26.
- d. Defective CR1 or CR2 (1N34).
- e. Shorted transformer T1 or T2.

3. Zero offset over 2%. To measure Zero offset, proceed as follows:

With no signal on the Input, and the Selector Switch at ZERO, adjust the meter needle to zero with the Zero control. Now turn the Selector Switch to 3×10^{-13} . The meter needle should return to a point within $\pm 2\%$ of scale zero. If it does not, the cause may be one of the following:

- a. Defective Capacitor C1.
- b. Defective insulation of Switch SW1.
- c. Charge built up on input plug (P3) by friction with mating plug. Charge will dissipate in several hours.

4. Temperature of thermoregulated housing too high or too low. Should be $58 \pm 2^{\circ}$ C. Adjust thermostat with a strong magnet by turning the adjusting bar, or replace the thermostat if contacts are defective.

5. Excessive noise when Selector Switch is in 3×10^{-13} position may be caused by:

- a. Defective Vibrode C.
- b. Microphonic input tube V1.
- 6. Meter will not read full scale on 3×10^{-11} , 3×10^{-9} , or 3×10^{-7} range. Probable causes are:
- a. Improper phasing of Demodulator T2. Correct by selecting C15 for maximum voltage across T2.
- b. Defective 1N34 (CR1, 2).
- c. Shorted turn in Transformer T1 or T2.
- d. Defective Capacitor C12.
- 7. Meter needle does not zero: may be caused by:
- a. Open Zero control R11.
- b. Open filament in tube V1, 2, 9 or 6.
- c. Open resistor in filament circuit R3, 10, 9, 33, 32, 34.



Figure 20. Schematic Wiring Diagram of Special Calibrating Equipment

10 %

0

3 X 10-7

3X10-11

10

.0005

3.0

3

FOR LOOP GAIN MEASUREMENTS



Figure 21. Views of Model V Chassis

PARTS LIST

	Symbol	Description	Beckman Part
	R1	300 Meg, $\pm 10\%$ RPC Type HBF	5630
	R2	$200 \text{ Meg}, \pm 10\%, 1 \text{ watt, S.S. White}$	5270-494
	• R3	100 ohms, $\pm 5\%$, Sprague 5 NIT	4068
	R4, 36	5.1 Meg, $\pm 5\%$, $\frac{1}{2}$ watt, AB EB	10624
	R5, 6, 31	22 Meg, $\pm 10\%$, $\frac{1}{2}$ watt, AB EB	8853
	R7, 3 0	.22 Meg, $\pm 10\%$, $\frac{1}{2}$ watt, AB EB	8102
	R 8	13 Meg, $\pm 5\%$, $\frac{1}{2}$ watt, AB EB	5659
	~ R9	1,500 ohms, $\pm 5\%$, Sprague, Type 5 NIT	4066
	• R 10	200 ohms, $\pm 5\%$, Sprague, Type 5 NIT	4067
-	R11	Potentiometer, 10 ohms, Helipot C.T. No. 10 CZ	5566
	R12, 29	24,000 ohms, $\pm 10\%$, $\frac{1}{2}$ watt, AB EB	7933
	R13	.27 Meg, $\pm 10\%$, $\frac{1}{2}$ watt, AB EB	8103
	R14	91,000 ohms, $\pm 5\%$, $\frac{1}{2}$ watt, AB EB	15073
	R15	10,000 ohms, $\pm 5\%$, $\frac{1}{2}$ watt, AB EB	15052
	R16	3,000 ohms, $\pm 5\%$, $\frac{1}{2}$ watt, AB EB	8845
	R17	51,000 ohms, $\pm 5\%$, $\frac{1}{2}$ watt, AB EB	15067
	R18	20 Meg, $\pm 5\%$, $\frac{1}{2}$ watt, AB EB	2208
	R19, 51, 55	1 Meg, $\pm 10\%$, $\frac{1}{2}$ watt, AB EB	8148
	R20	270 ohms, $\pm 10\%$, $\frac{1}{2}$ watt, AB EB	8070
	R21	$.24 \text{ Meg}, \pm 5\%, \frac{1}{2} \text{ watt}, \text{AB EB}$	3180-13
	R22	.18 Meg, $\pm 5\%$, $\frac{1}{2}$ watt, AB EB	15081
	R23	244 ohms, \pm 5 ohms	5664
	R24	Not used	
	R25	$487.7 \text{ ohms}, \pm 1 \text{ ohm}$	
	R26	$9/4.3 \text{ onms}, \pm 2 \text{ onms}$	5642
	R2/	3415 onms, \pm 7 onms)	
	R28	Not used	1 4 9 0
24 Carlos and a carl	- K32	$5,000 \text{ onms}, \pm 5\%$, 5 watt, Dale Prod, RS-5	1430
	- K33	1,750 onms, $\pm 5\%$, 5 watt, Dale Prod, RS-5	1429
	K34 D25	8,000 onms, $\pm 5\%$, 5 wall, Dale Prod, RS-5	1428
	R33	1 100 abrea (50% 10 wett Mellory Type 111 without areas	1 0040
	R3/ D29	$1,100 \text{ onms}, \pm 5\%, 10 \text{ wall, Mallory, Type THJ vitreous ename0747 shma \pm 25 \text{ shma}$	el 8840
	R38 R20	$9/47$ onnis, ± 23 onnis 51 Mag $\pm 5\%$ (14 watt AD ED	2042
	R39 R40	.51 Meg, $\pm 5\%$, $\frac{1}{2}$ wall, AB EB	9092
	R40 D41	10 Meg, $\pm 1\%$, which Type CP1, 1 wall, w/ vinyite sleeving	5620
	R41 D42	1,000 Meg, $\pm 5\%$ at 5 v and 58° C. Type EBF of victoreen	5623
	R42 D42	145,000 Meg, $\pm 5\%$, 5 v and 58° C. Type HBM of Victoreen	5622
	R43 D44 45	100 abms + 10% 50 wett Obmite	0941
	R44, 45	100 ohms, $\pm 10\%$, 50 watt, Onmite 100 ohms, $\pm 10\%$, 16 watt, IPC. Two PW 16	9841
	R40 D47	$100 \text{ onms}, \pm 10\%, \%$ watt AD ED	9001
	R4/	.50 Meg, $\pm 5\%$, $\frac{3}{2}$ wall, AD ED	1102
	K40, 49, 30	Not used $5.100 \text{ obm} = 5.0\%$ 1/2 worth AD ED	0100
	КJ2 D52	$3,100$ ohms, $\pm 5\%$, $\frac{1}{2}$ wall, AB EB	ð123 15020
	КЈЈ D54	1,000 ohms, $\pm 3\%$, 72 wall, AB EB 25,000 ohms, $\pm 5\%$, 10 worth Spraws	13038
	КJ4 D56	$33,000 \text{ ohms} \pm 5\%, 10 \text{ watt AP CP} $	0 561 14
	NJU D57 50	20,000 onlines, ± 3.70 , 1 watt, AD OD 02/ 220 ohms $\pm 10\%$ 1/2 watt AD ED	0-204-14
	NJ1, JO D50	220 onnis, $\pm 10\%$, 72 waii, AD ED Not used	0009
	RJ7 D60	$75\ 000\ ohms + 5\% \ 16\ watt \ AP EP$	0775
	R00 D61	10 Mag $\pm 10\%$ 1/2 watt AB ED	0//3
	C	Vibrode	5606
	\sim_{v}	101040	2000

Symbol	Description	Beckman Part
01		
	$20001 \text{ mId}, 600 \text{ v}, \text{ Polystyrene}, \pm 20\%$	5596
C2	2 mfd, 400 v, Mallory No. CB-406	5617
C3	.0001 mtd, 400 v, Polystyrene, $\pm 10\%$, Solar No. SDP-40001	12058
C4, 5	.1 + .1 mfd, Mallory No. CBD-602	5619
C6	.5 mfd, Mallory No. CB-404	5618
C7	.01 mfd, 600 v, paper, $\pm 20\%$, Sprague No. 73P 10306	15024
C8	.005 mfd, 600 v, paper, $+30 - 10\%$, Sprague No. 256-AG	6270-453
C9, 10	.1 + .1 mfd, Mallory No. CBD-602	5619
C11, 14	40 mfd, 450 v, electrolytic, Mallory No. FP-238	144
C12	3,000 mfd, 10 v, Mallory No. WP-032	5621
C13	.00025, 600 v, Polystyrene, $\pm 20\%$, Fast A7877FJ	6270-348
C15	.1 mfd, paper, $\pm 20\%$, Sprague No. 73P 10406	13370
C16	Supplied with transformer T4	5615
C17	.01 mfd, 300 v, mica (low loss case), $\pm 20\%$	
	Cornell-Dubilier No. IDL-3S1	9062
C18	.0001 mfd, 500 v, mica, $\pm 20\%$, Cornell-Dubilier No. SW-5T1	2191
C19	.0022 mfd, 600 v, molded paper, $\pm 20\%$,	
	Sprague No. 73P 22206	15028
C20	.033 mfd, 600 v, paper, $\pm 20\%$, Sprague No. 73P 33306	13370
C21	.25 mfd, 400 v, paper, Mallory No. PT 4025	4651
C22	.01 mfd, 600 v, polystyrene, +20%, Fast No, A7878 CJ	9039
C23	.0005 mfd, 600 v, polystyrene, $\pm 20\%$, Fast No. A6290 FJ	9029
C24. 25	Not used	2022
C26	50 mfd, 25 v. Mallory No. TC-29	5559
C27	1 mfd. paper. + 20%. Sprague No. 73P 10406	13370
C28, 29	5 mfd, 600 v. Mallory CB No. 604	5564
C30	.0005 mfd. 500 v. mica. Sprague No. 1 FM-35	10893
C31	1 mfd, 600 y, paper. + 20% Sprague No. 73P 10406	13370
V1. 2. 6. 9	CK 533. Raytheon	12524
V3	12AU7. RCA	12015
V4. 7	6A05. RCA	12015
V5	5651. RCA	173
V8	6X4. RCA	12013
CR1 2	1N34 Sylvania	738
M1	Meter	5536
T1 2	Transformer (demodulator)	5616
T3	Transformer Power	5502
T4	Transformer, Regulating 6.3 Volt Sola No. 301002 Type 12	5615
P1	Connector base flush male Hubbell No. 7556G	4536
P2	Plug Jones No. P-302-CCT	4550
P3	Plug Amphenol No. 82 805 Teflon	5614
S1	Socket Jones No. S 302AB	5611
S1 S2	Not used	5011
52 53	Socket AN 3102A $1/5$ 55	5620
55 S4	Socket Amphenol No. 80C	0003
5- 11 3	Tin-Jack F F Johnson 105-520 Pad	5677
51, 5 12	Tin-Jack F. F. Johnson 105-520 Red	5670
SW1	Switch Selector	JU28 5676
SWI	DDST Toggle Switch 12A 125V Amoun LI ALL No. 00142	30/0
G W Z CW/0	DIST Toggle Switch, 12A, 125V, AITOW H.&H. NO. 82143 DEST Toggle Switch, 6A, 125V, Custon Harrison N. 027047	2313
$F1 \gamma$	Fuse 2 amp SloBio Littelfuse 212002	1918
F1, 4	Fuse, Δ amp. Stobio Entenuse 515002 Fuse, 0.1 amp. Littelfuse 212100	901/
ג'ס דסו	Thermostat adjusted to 500 C	9842
11/1	rnermosiai, aujusieu io 30 °C	1030



Another Beckman Achievement BECKMAN Model W Industrial pH Meter

-A Husky New Addition to America's Most Distinguished Family of



New Savings-New Speed New Production Advantages

Beckman – pioneer of virtually every major advancement in pH equipment – now presents a significant new development in industrial pH instruments – the new Model W Industrial pH Meter.

This new industrial pH Meter takes advantage of modern electronic engineering, using miniature and subminiature tubes and printed circuits. It is only one-quarter the size of conventional pH instruments ...and is designed to save maintenance costs by sturdy, moisture-proof construction. The instrument circuit is of a completely new Unitized design – made in three separate plug-in sections, any of which can be interchanged in the field to avoid downtime.

Whatever your industry, whatever your pH problem, get the fages about Beckman pH Control Instruments.

Beckman division

Many New and Improved Features

(10) x 10 - 10

ADVANCED DESIGN Unusually compact, utilizing subminiature components and printed circuits.

ACCURACY Signal to recorder accurate within ± 0.02 pH.

SIMPLIFIED MAINTENANCE Easily disassembled for quick maintenance ... plug-in units.

DEPENDABILITY Unaffected by line voltage fluctuations, corrosive atmospheres, vibration, or electrical machinery.

SPLASH PROOF Designed for high humidity and outside exposure conditions.

EXPLOSION PROOF Can be used with air purge in explosive atmospheres. Heavy cast iron case available.

BECKMAN INSTRUMENTS, INC., FULLERTON, CALIFORNIA

80271 . PRINTED IN U.S.A.